

ORTHO TECHNOLOGY **Expect Excellence**





Ortho Technology not only offers several aesthetic bracket lines but takes aesthetic treatment to the next level by offering an extensive selection of Tooth Tone coated archwires, ties, springs, and attachments.



Unique and Exclusive Products

Ortho Technology takes pride in offering a multitude of unique and exclusive products. We are constantly researching, developing, and improving our products to ensure the most innovative and highest quality products are available.



Product Line Variety and Completion

Ortho Technology offers all products necessary to complete any patient treatment. From first impressions to debonding pliers and burs as well as patient hygiene accessories, we offer a thorough selection of products.



Superior Customer Service

Ortho Technology strives to provide complete customer satisfaction with every interaction. By doing this, our goal is to ensure a long-lasting and trusted relationship with every customer.



Competitive, Honest Pricing

The goal of Ortho Technology's pricing structure is to supply quality products at a fair price. We understand that a fair product price is just as important as the product itself.



Same Day Shipment

Most orders ship the same day you place your order. You can have your products as soon as the next business day. Place your order before 2:00 P.M. E.S.T. and request overnight shipping. Otherwise we will ship your order via common carrier.



Hassle-Free Returns

We understand that your purchase may not always meet your needs. In an effort to make the return process as hassle-free as possible, we ask that you please call us. Your satisfaction is very important to us, and your input helps us improve our products and service. When you call, a return authorization number will be assigned. This number will help expedite your exchange, credit, or refund.



YOUR FIRST CHOICE

It's Ortho Technology's goal to be "Your First Choice" in orthodontic supplies. Since 1991 Ortho Technology has offered innovative and unique products that orthodontic professionals worldwide have come to trust. Ortho Technology offers a complete line of orthodontic supplies that includes: brackets, bands, buccal tubes, archwires, adhesives, lab supplies, and patient care products. Product selection couldn't be any easier with Ortho Technology's well-organized full-color catalog, featuring thousands of products you use every day.

COMMITMENT **TO QUALITY**

Ortho Technology is committed as a company to provide your practice with the very best products possible, while delivering the utmost in customer service. Ortho Technology is certified to an internationally recognized quality system, ISO 13485. When you choose Ortho Technology, you're choosing a company that has made a tremendous commitment to making products that are of the highest quality.







Our Systems

Ortho Technology Offers a Full Line of Archwire Solutions

TruFlex™, TruForce™, BetaForce™, Universal Form™, and TRIFECTA® lines are available in Stainless Steel, Beta Titanium, Nickel Titanium, Copper Nickel Titanium with and without stops and Thermal Nickel Titanium. Shaped as Pre-Torqued, Full Form, Euro Form, Standard Form, Reverse Curve of Spee and Retraction archwires. Avail-

able in Tooth Tone®, 8 Strand Braided, 3 Strand Twist, 6 Strand Coaxial, Gold Plated, Spooled and Straight lengths. Available with etched

or dimpled midlines.

See pages 67-85

13 Different Bracket Systems Offer Maximum Choice

Crystal clear sapphire, translucent ceramic, metal-lined ceramic, composite, metal-lined composite, stainless steel, gold plated stainless steel, and self-ligating brackets. Patients will love the choices and you will love the performance of Ortho Technology's extensive line of bracket systems.

See pages 8-47

Creatures of the Sea® Orthodontic Elastics

Easy to use identification system manufactured with high quality medical grade materials to deliver consistent and accurate force. Printed tooth diagram aids in proper placement of elastics. Packaged in light-resistant patient packs.

See pages 100-101

TruFit™ 2.0 Molar Bands with Prewelded Buccal Tubes

Each TruFit 2.0 band has precision shaped tooth anatomy for a tight secure fit with

excellent retention. Strong yet malleable proximal walls that require little to no burnishing. Maximum buccal height allows greater placement and positioning flexibility of buccal tubes and lingual attachments.

See pages 49-62



Your Practice

X7° Orthodontic Instruments and PANTERA° Pliers and Cutters°

Our X7 and PANTERA pliers and cutters are hand crafted using the finest grade surgical steel. X7 Orthodontic Instruments are a premium line, precision manufactured in the USA with USA and/or German stainless steel. The X7 line features a shadow-satin surfacing to minimize glare and provide a long-lasting finish, and the PANTERA line offers a smooth, satin finish. Cutting edges on all inserted PANTERA Instruments are made from a custom blended specialty alloy, and both lines of instruments feature smooth contoured edges for both user and patient comfort.

X7 see pages 170-173 Pantera see pages 174-177

Resilience® Bonding Supplies

The light-activated and primer-activated systems provide the optimum properties for bonding ceramic, metal, or plastic brackets to teeth. Both of the composite adhesives have ideal viscosity which delivers predictable handling characteristics that eliminate bracket flotation and allow for easy flash removal. Produces a quick setting, strong mechanical bond which debonds easily. Contains fluoride to help protect against potential decalcification.

See pages 104-109

Galaxy™ Interproximal Reduction

Diamond Discs are perfect for both reduction and contouring. Smoothly remove tooth mass interproximally to create intra-arch space or open contact points between teeth. Great for use with many popular clear aligner systems.

Discs are available in 22.0 mm or 19.0 mm when space is limited.

See page 145

Match

DuraLock® PLUS Impression Trays

Using DuraLock PLUS trays save time. By eliminating the extra steps of coating the tray with adhesive and waiting for it to dry,

the impression process is simplified and accelerated. DuraLock PLUS trays also eliminate inadequate coverage and/or tissue impingement and ensure an accurate fit every time by offering more size choices than competing trays. Reliable strength and stability are the trademarks of DuraLock PLUS. See pages 128-129



DON'T MISS THESE

GREAT **NEW!** PRODUCTS FROM





STAINLESS STEEL BRACKET SYSTEM **Revolutionary Strength and Consistency** in a One-Piece Bracket

See pages 39-41







Nickel Titanium Archwire

See page 74





See page 112















See page 142



Eruption Appliance Curved Eyelet & Angled Button Tweezers

See page 185



ULTRA-GUARD



See pages 152-153

© ORTHO TECHNOLOGY®



2015 INNOVATION, AESTHETICS AND EXCELLENCE

ORTHO TECHNOLOGY®

We understand that as an Orthodontist, you are constantly searching to improve the treatment experience for you and your patients. To support you, we offer innovative solutions through our many bracket lines. From stainless steel to clear sapphire, translucent ceramic and composite aesthetic to passive and active self-ligating, we have the bracket you're looking for. Patients will love the choices, and you will love the performance of Ortho Technology's extensive line of bracket systems. Expect the highest quality available at a fair price with every product we offer.

SELF-LIGATING



Sensation

Active Ceramic Self-Ligating Bracket System

P.10

SensationV

Active Stainless Steel Self-Ligating Bracket System

P.13

Lotus Plus

Passive Self-Ligating Bracket System

P.16

DISCOVER O EXTENSIVE LINE OF BRACK

AESTHETIC



PURE®

Sapphire **Bracket System**

P.22

Reflections

Ceramic **Bracket System**

P.26

Encore!

Metal-Lined Ceramic **Bracket System**

P.29

Avalon®

Metal-Lined Composite **Bracket System**

P.31

OrthoFlex*

Composite Bracket System

P.32

METAL



OTION

Stainless Steel Bracket System

P.36

PINNACLE

NEW!

Stainless Steel Bracket System

P.39

Marquis

Stainless Steel Bracket System

P.42



Bionic*

Stainless Steel **Bracket System**

P.44



Marquis[™] TruGold*

24K Gold Plated **Bracket System**

P.46

Find us on Facebook **f**

Follow us on Twitter



View us on YouTube








The Benefits are Clear

Sensation Active Ceramic Self-Ligating Brackets benefit both patient and practice by providing uncompromised aesthetics and comfort for your patients with faster treatment times for your practice. Sensational results from start to finish. The ceramic material provides strength and durability as well as a translucent appearance that is unmatched by any other bracket system. The rhodium coating on the cobalt-nickel clip increases the clip's strength for extraordinary durability. A strong mechanical bond is created on the bracket base through the use of alumina spheres. Debonding is predictable and easy. The unique guide rail ensures stable and consistent force during clip opening and closing making archwire changes a breeze. With fewer archwire changes, shorter office visits, and reduced chairtime, the benefits are clear.

Identification System



Passive to Active







Transitional **Passive**

Active



© 2015 Ortho Technology, Inc. Sensation is a trademark of Ortho Technology. Empower is a registered trademark of American Orthodontics and In-Ovation is a registered trademark of DENTSPLY GAC

BRACKETSYSTEMS

Rounded Corners

For added patient comfort

Archwire Slot Design

Enhances sliding mechanics

Extremely resistant to slotwear, cracking, and fracturing

High Quality Rhodium Coated Clip

Manufactured using a Cobalt-Nickel alloy to provide resiliency and durability


Non-Removable Clip

Prevents accidental removal

Unique Bonding Base

Alumina spheres provide superior shear bond strength and predictable debonding

> **No Primers Needed** for Bonding



Removable Color Coded ID Dot

For easy bracket selection

For stable and consistent force during clip opening

Compound Contour Base

Ultra-Smooth, Strong Hooks For added patient comfort and

secure engagement of elastics

Vertical Scribe

For vertical repositioning

Generous Contoured Tie-Wing Area

For easy bracket ligation and patient comfort

Translucent

For pleasing aesthetics

Sensation[™] **Work with Efficiency**

Simple Opening







Simple Closing





And Most Importantly

Predictable Debonding







See page 12 for complete Opening, Closing and Debonding Instructions



Scan the QR Code or Go Online www.orthotechnology.com/sensation For complete Opening, Closing, and Debonding Instructions



On Our Website Sensation™ Instructional Video

BRACKETSYSTEM



Single Patient Kit Description Upper/Lower 5x5*

t (20 brackets)	\$299 00 kit
kits (200 brackets)	
kits (400 brackets)	

In	dividual Brackets*	
1	oack (10 brackets)	\$149.50/pk (\$14.95 ea.)
20	packs (200 brackets)	\$139.50/pk (\$13.95 ea.)
40	packs (400 brackets)	\$129.50/pk (\$12.95 ea.)

¹⁰ per pack *Not available in Canada

Sensation™

Ortho Technology's version of the Roth Prescription*

▲ Maxil	lary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1	Central	+11°	+5°	.126"	R	180-12-1100
111	Centrai	+11	+3	.120	L	180-12-2100
2 2	Lateral	+8°	+9°	.116"	R	180-12-1200
22	Laterai	T-0	Τ3	.110	L	180-12-2200
3 3	Cuspid	-2°	+8°	.138"	R	180-12-1310
3 3	Cuspid with hook	-2	Τ0	.136	L	180-12-2310
4/5 4/5	Bicuspid with hook	-7°	0°	.126"	R	180-12-1910
/5//5	with hook	-/	0	.120	L	180-12-2910

Patient Kits Item Numbers .022" Slot Description Upper/Lower 5x5 Hooks 3, 4, & 5 180-12-2011

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lower Anterior	0°	0°	.112"	U	180-12-6800
3 3 Cuspid	-11°	+7°	.138"	R	180-12-3310
3 3 Cuspid with hook	-11	Τ/	.136	L	180-12-4310
Bicuspid with hook	-17°	0°	.126"	R	180-12-3910
⁷⁵ with hook	-17	U	.120	L	180-12-4910

Sensation™

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxil	lary Arch	Torq.	Ang.	Width	R/L	.022" Slot
111	Central	+17°	+4°	.126"	R	180-22-1100
111	Centrai	+ 1 /	+4	.120	L	180-22-2100
2 2	Lateral	+10°	+8°	.116"	R	180-22-1200
22	<u>Z Z</u> Lateral	+10	+0	.110	L	180-22-2200
2 2	Cuspid	+0°	+8°	.138"	R	180-22-1310
3 3	Cuspid with hook	+0	+0	.136	L	180-22-2310
4/5 4/5	Bicuspid with hook	-7°	0°	.126"	R	180-22-1910
/5//5	with hook	-/	0	.120	L	180-22-2910

Patient Kits Item Numbers Description .022" Slot 180-22-2011 Upper/Lower 5x5 Hooks 3, 4, & 5

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lower Anterior	-6°	0°	.112"	U	180-22-6800
3 3 Cuspid	0°	+3°	.138"	R	180-22-3310
3 3 Cuspid with hook	0	Τ.	.136	L	180-22-4310
Bicuspid by With hook	-12°	0°	.126"	R	180-22-3910
^{5 5} with hook	-12	U	0 1.120	L	180-22-4910

^{*}The Ortho Technology versions of the Roth and McLaughlin, Bennett, & Trevisi Rx's are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Sensation™ SLB Opening Instrument **Unique Clip Opening Tool Specifically Designed** to Prevent Overstressing of the Clip



Sensation Brackets must be opened using the OT-7920 Clip Opening Instrument

Item #: OT-7920

3 or more**\$28.50 ea.**

Using the incorrect instrument Heat sterilizable up to 370°F/188°C may result in clip damage.

Important Clip Opening/Closing Tips:

An important part of the closing procedure is making sure the archwire is securely pressed into the base of the archwire slot before pressing on the clip. This can be done with a ligature director or light finger pressure on the archwire. Failure to do so may deform the clip and cause it to be inoperable.

Both Sensation and Sensation M clips open and close in the same manner, however, each require a unique instrument. While both instruments may look similar, each must be used with the appropriate bracket system.



FOR MORE INFORMATION ABOUT...

Sensation Active Self-Ligating Brackets Opening, Closing and Debonding Instructions

See pages 187





Ideal Treatment Control

Sensation M Active Stainless Steel Self-Ligating Brackets provide ideal control in all stages of treatment. The passive stage of the clip during the initial treatment phase allows for more efficient leveling and aligning. Sliding mechanics are also enhanced during the passive and transitional phases. The active stage of the clip engages the archwire securely into the base of the archwire slot, fully expressing the in-out, angulation and torque built-in to the bracket.

Unlike traditional brackets where forces can diminish when ligatures are stretched into place, the self-ligating clip design offers continuous force levels during wire deflection/unloading. This allows practitioners to adjust time between office visits.

Low Profile Design

Smooth rounded facial contours to optimize patient comfort

High Quality Integral Clip

Manufactured using a Cobalt-Nickel alloy to provide resiliency and durability

Ample Under Tie-Wing Area For application of elastomeric chains

Removable Color Coded ID DotFor easy bracket selection

Fixed Clip Design
Prevents accidental removal

Long Axis Scribe Line

Acts as a clear visual cue for placement at the top and bottom of the bonding base



For added patient comfort and secure engagement of elastics


Compound Contour Base

Provides an ideal bracket fit to tooth surface

80 Gauge Micro-Etched Mesh Base

Ensures maximum bond strength with predictable and safe debonding

Identification System




Passive to Active









BRACKETSYSTEM



Sensation[™] M

Ortho Technology's version of the Roth Prescription*

▲ Maxil	lary Arch	Torq.	Ang.	Width	R/L	.022" Slot
111	Central	+11°	+5°	.114"	R	390-12-1100
	Central	trai +11 +5	Τ.)	.114	L	390-12-2100
2 2	Lateral	+8°	+9°	.106"	R	390-12-1200
2 2	Laterai	Τ0	T9	.100	L	390-12-2200
3 3	Cuspid	-2°	+8°	.116"	R	390-12-1310
3 3	with hook	-2	+0	.110	L	390-12-2310
4/5 4/5	Bicuspid	-7°	0°	.116"	R	390-12-1910
<u> </u>	with hook	-7	U	.110	L	390-12-2910

Sensation™ M

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxilla	ry Arch	Torq.	Ang.	Width	R/L	.022" Slot			
1 1	Central	+17°	+4°	.114"	R	390-22-1100			
	Centrai	Τ17	+4 .	.114	L	390-22-2100			
2 2	Lateral	+10°	.10° .0	+8°	.106"	R	390-22-1200		
_Z Z Ld	Laterai		ΤΟ	.100	L	390-22-2200			
3 3	Cuspid	0°	.00	1 0°	n. 10.	0. 10.	.116"	R	390-22-1310
3 3	with hook	0	то	.110	L	390-22-2310			
4/5 4/5	Bicuspid with hook	-7°	0°	٥٠	.116"	R	390-22-1910		
<u> </u>	with hook	-7		.116	L	390-22-2910			

Single Patient Kit Description

Upper/Low	er 5x5*
-----------	---------

1 kit (20 brackets)	\$209.00	kit
10 kits (200 brackets)	\$189.00	kit
20 kits (400 brackets)	\$169.00	kit
Individual Practices		

inuividuai brackets"	
1 pack (10 brackets)	^{\$} 104.50/pk (^{\$} 10.45 ea.)
	\$94.50/pk (\$9.45 ea.)
40 packs (400 brackets)	\$84.50/pk (\$8.45 ea.)
10 per pack	• • • • •

^{*}Not available in Canada

Patient Kits Description

Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers .022" Slot

390-12-2011

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lower Anterior	0°	0°	.098"	U	390-12-6800
Cuspid	-11°	1.7°	.116"	R	390-12-3310
3 3 Cuspid with hook	-11	Τ/	.110	L	390-12-4310
4/4/ Bicuspid with hook	-17°	0°	.116"	R	390-12-3910
^{75 75} with hook	-17	0	.110	L	390-12-4910

Patient Kits Description

Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers .022" Slot 390-22-2011

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lower Anterior	-6°	0°	.098"	U	390-22-6800
3 3 Cuspid	0°	+3°	.116"	R	390-22-3310
3 3 Cuspid with hook				L	390-22-4310
Bicuspid with hook	-12°	0°	.116"	R	390-22-3910
^{75 75} with hook	12			L	390-22-4910

^{*}The Ortho Technology versions of the Roth and McLaughlin, Bennett, & Trevisi Rx's are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Sensation™M SLB Opening Instrument **Unique Clip Opening Tool Specifically Designed** to Prevent Overstressing of the Clip



Sensation M Brackets must be Item #: OT-8920 Opening Instrument. Using the 3 or more\$28.50 ea. incorrect instrument may result Heat sterilizable up to 370°F/188°C in clip damage.

Important Clip Opening/Closing Tips:

An important part of the closing procedure is making sure the archwire is securely pressed into the base of the archwire slot before pressing on the clip. This can be done with a ligature director or light finger pressure on the archwire. Failure to do so may deform the clip and cause it to be inoperable.

Both Sensation and Sensation M clips open and close in the same manner, however, each require a unique instrument. While both instruments may look similar, each must be used with the appropriate bracket system.







Simple Opening - To open the Sensation M Bracket Clip:

- 1. Place the inside hook of the opening instrument into the opening of the door of the bracket interface.
- 2. Tilt the opening instrument gingivally to avoid possible removal of the clip.
- 3. Pull straight down to open using a gentle controlled force.





Simple Closing - To close the Sensation M Bracket Clip

- 1. Place light pressure to seat the archwire in the slot base by using finger pressure or a ligature director.
- 2. Using light finger pressure, press up and in to close the clip.

Sensation

Active Self-Ligating Bracket Systems

For Sensation™ Ceramic and Metal Self-Ligating Bracket Systems

Suggested Order, Duration, and Objective of Basic Archwire Sequence



Passive Phase Round Archwires

.014", .016", or .018"

TruFlex™

Thermal Nickel Titanium Full Form

.014" Upper Item #: 6001-014 .014" Lower Item #: 6001-114

.016" Upper Item #: 6001-016

.016" Lower Item #: 6001-116

.018" Upper Item #: 6001-018

.018" Lower Item #: 6001-118

- Begin Leveling and Alignment
- Correct Rotations
- Begin Archform Development

TruFlex™ Thermal Nickel Titanium

Round

1 pack......\$21.50/pk 10 packs \$17.25/pk 10 per pack



Transitional Phase Square Archwires

.018" x .018"

TRIFECTA®

Thermal Nickel Titanium Full Form

.018" x .018" Upper Item #: 6003-047

.018" x .018" Lower Item #: 6003-147

- or -

.020" x .020"

TRIFECTA®

Thermal Nickel Titanium Full Form

.020" x .020" Upper Item #: 6003-049

.020" x .020" Lower Item #: 6003-149

- Continue Arch Development
- Complete Leveling and Alignment
- Begin Torque Control
- Begin Archform Coordination

TRIFECTA® Thermal Nickel Titanium

■ Rectangle

1 pack.......\$43.50/pk 10 per pack



Active Phase Rectangular Archwires

.018" x .025"

TruForce™ Stainless Steel Full Form

.018" x .025" Upper Item #: 8000-112

.018" x .025" Lower Item #: 8000-212

- or -

.019" x .025"

BetaForce™ Beta Titanium Full Form

.019" x .025" Upper Item #: 6002-111

.019" x .025" Lower Item #: 6002-211

TruForce™ Stainless Steel Full Form

.019" x .025" Upper Item #: 8000-113

.019" x .025" Lower Item #: 8000-213

- Begin Posterior Consolidation Coordination of Archforms
- Complete Posterior Consolidation and **Archform Coordination**
- Complete Torque Expression

TruForce[™] **Stainless Steel**

■ Rectangle

1 pack......**\$7.10/pk** 10 per pack

BetaForce™ Beta Titanium

■ Rectangle

10 per pack

Above chart is suggested order, duration, and objective of basic archwire sequence and by no means a replacement for professional expertise. © 2015 Ortho Technology, Inc. BetaForce, TruFlex, and TruForce are trademarks and TRIFECTA is a registered trademark of Ortho Technology.









BRACKET SYSTEMS

BRACKETSYSTEMS

Lotus Plus®

Passive Self-Ligating Bracket System

COMPARE TO: Damon® by Ormco

Designed by Orthodontists for Total Control and Maximum Patient Comfort


- Low profile mini-twin design with rhomboid shape for easy placement
- Bracket design provides passive engagement to reduce friction and improve sliding mechanics
- No Special Instruments Required. Opens easily with an explorer, closes with utility plier or tweezer, or try our Lotus Plus specialty instruments for additional control
- Mandibular arch clips open gingivally to minimize occlusal interference, maxillary arch clips open occlusally
- Convertible bracket allows for the engagement of ligature ties during the finishing phase if desired

Save Bonding Time!

Lotus Plus ships in the open position!







© 2015 Ortho Technology, Inc. Lotus Plus is a registered trademark of Ortho Technology. Domestic and international patents pending. Damon is a registered trademark of Ormco Corporation.



BRACKETSYST

Open Design Allows for maximum patient hygiene

Unique Clip Design Allows for easy opening and secure closure of the clip

Easy Opening Clip Opens easily with an explorer



Identification System

By quadrant with removable color coded I.D. dot

Ultra Low Profile

For a more aesthetic appearance Maximum patient comfort and satisfaction Avoids occlusal and incisal interference Guides archwire adjacent to the tooth and increases the action of archwire forces

80 Gauge Foil Mesh Base

Ensures maximum bond strength


Mini-Twin Design With rhomboid shape for easy placement

Convertible Bracket Removable clip allows for

ligature ties during finishing phase if desired

Thermal Nickel Titanium For flexibility and long-life

High Quality



Contoured Tie-Wing

For added patient comfort

Generous and Defined Tie-wing Undercuts

Ensures easier and more secure ligation if needed

Compound Contour Base

Provides simplified and accurate bracket placement

Gingivally Offset Bicuspid Bracket Available

Reduces interference while decreasing bond failure

Easy Opening



Opens Easily with No Special Tools

Simply insert the end of an explorer into the circular opening located just below the archwire slot.



Rotate the explorer in a downward 45° motion, gently easing the clip open. Do not pull down on the clip as this will convert the bracket. The clip will "click" once as it disengages from the top lock and "click" a second time when fully opened.



Low Profile



Damon® Q™ In-Ovation® R

Finish with Control

Convertible Bracket


If desired, during the finishing phase, clips can easily be removed and ligatures can secure the archwire into the slot.

To Convert The Bracket

Pull the clip down to remove when you desire to disengage the clip and convert the bracket.



Secure Closure



Gently apply pressure to the plier or tweezer and the clip will begin to move up. The clip will "click" once as it disengages from the bottom lock and "click" a second time when closed.







BRACKET SYSTEMS



Single Patient Kit Description Upper/Lower 5x5

1 kit (20 brackets)	19.00	kit
10 kits (200 brackets)\$2		
25 kits (500 brackets)	89.00	kit

Individual Brackets

1 pack (10 brackets)	\$109.50/pk (\$10.95 ea.)
	s)
	s) \$94.50/pk (\$9.45 ea.)
10 ner nack	

Lotus Plus® Replacement Clips Available

Ask your Sales Representative for Complete Details



Lotus Plus®

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 Central	+12°	+5°	.114"	R	300-12-1100
TIT Central	+ 1 2	+3	.114	L	300-12-2100
2 2 Lateral	+8°	+9°	.110"	R	300-12-1200
ZZ Lateral	+0	+3	1.110	L	300-12-2200
3 3 Cuspid	0°	+11°	.110"	R	300-12-1310
$\frac{3}{3}$ with hook	0	+''		L	300-12-2310
3 3 Cuspid	0.	+11°	.110"	R	300-12-1300
Cuspia	0	T11	.110	L	300-12-2300
4/14/ Bicuspid	٦.	0.	114"	R	300-12-1910
4/5/4/ Bicuspid with hook	-7°	0°	.114"	L	300-12-2910
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.114"	J	300-12-5900

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Patient Kits Item Numbers Description .022" Slot Upper/Lower 5x5 No Hooks 7322-1000

Upper/Lower 5x5 No Hooks 7322-1000
Upper/Lower 5x5 Hooks on Cuspids 7322-1300
Upper/Lower 5x5 Hooks on 3, 4, & 5 7322-1345

▼ Ma	ndibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
1/2 1/2	Lower Anterior	0°	0°	.102"	U	300-12-6800
3 3	Cuspid	-11°	+5°	.114"	R	300-12-3310
3 3	with hook		1 3		L	300-12-4310
3 3	Cuspid	-11°	+5°	.114"	R	300-12-3300
3 3	Cuspiu	-11	Τ.)	.117	L	300-12-4300
4 4	1st Bicuspid	-17°	0°	.114"	R	300-12-3430
4 4	w/hook & offset	-17	U	.114	L	300-12-4430
4 4	1st Bicuspid with offset	-17°	0°	.114"	U	300-12-6420
	2 nd Bicuspid	-22°	0°	.114"	R	300-12-3530
5 5	w/hook & offset	-22	U	.114	L	300-12-4530
5 5	2 nd Bicuspid with offset	-22°	0°	.114"	U	300-12-6520

Lotus Plus®

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+17°	+4°	.114"	R	300-22-1100
TIT Central	+ 1 /	T "	.117	L	300-22-2100
2 2 Lateral	+10°	+8°	.110"	R	300-22-1200
ZIZ Lateral	110	+0	.110	L	300-22-2200
3 3 Cuspid	0°	+8°	.110"	R	300-22-1340
S Cuspiu	0	+0	.110	L	300-22-2340
3 3 Cuspid	0°	+8°	.110"	R	300-22-1350
$\frac{3 3}{}$ with hook	0	Τ0	.110	L	300-22-2350
3 3 Cuspid	-7°	+8°	.110"	R	300-22-1300
S Cuspiu	-/	+0	.110	L	300-22-2300
3 3 Cuspid with hook	-7°	+8°	.110"	R	300-22-1310
with hook	-7	+0	.110	L	300-22-2310
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.114"	U	300-22-5900
4/14/ Bicuspid	-7°	0.	.114"	R	300-22-1910
$\frac{4/5}{5}$ Bicuspid with hook	-7	U	.114	L	300-22-2910

*The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Patient Kits Description Upper/Lower 5x5 No Hooks Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks on 3, 4, & 5** Upper/Lower 5x5 Hooks on 3, 4, & 5** Upper/Lower 5x5 Hooks on 3, 4, & 5**

**Includes brackets upper 300-22-1350, 300-22-2350 and lower 300-22-3350, 300-22-4350. $\dagger 0^*$ torque on cuspids (3s)

▼ Mandibul	ar Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lowe Ante		-6°	0°	.102"	U	300-22-6800
3 3 Cusp	oid	-6°	+3°	.114"	R	300-22-3300
					L	300-22-4300
3 3 Cusp	vid	0°	+3°	.114"	R	300-22-3340
3 3 Cusp	nu	U	+5	.114	L	300-22-4340
3 3 Cusp		-6°	+3°	.114"	R	300-22-3310
^{3 3} with	hook	-0	+υ	.114	L	300-22-4310
3 3 Cusp		0.0	٠.	11.41	R	300-22-3350
ola with	hook	0°	+3°	.114"	L	300-22-4350
	cuspid	-12°	+2°	.114"	R	300-22-3420
with	offset	12	12	.117	L	300-22-4420
1st Bi	cuspid	1.00	. 20	11 411	R	300-22-3430
4 4 w/hoc	k & offset	-12°	+2°	.114"	L	300-22-4430
ElE 2 nd B	icuspid	1 70	. 3.	114"	R	300-22-3520
	offset	-17°	+2°	.114"	L	300-22-4520
5 5 2 nd B	icuspid	1 7.	. 2°	11 411	R	300-22-3530
o o w/hoo	k & offset	-17°	+2°	.114"	L	300-22-4530

BRACKETSYS

Lotus Plus®

Ortho Technology's version of the Damon[®] Standard Torque Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+15°	+5°	.114"	R	300-72-1190
1 Central	713	+3	.114	L	300-72-2190
2 2 1 242 421	+6°	+9°	.110"	R	300-72-1290
22 Lateral	+0	+9	.110	L	300-72-2290
3 3 Cuspid	+7°	+5°	.110"	R	300-72-1390
33 Cuspiu	+/	+3	.110	L	300-72-2390
3 3 Cuspid	+7°	+5°	.110"	R	300-72-1391
$\frac{3 3}{}$ with hook	+/	+3	.110	L	300-72-2391
$\frac{4}{5}$ $\frac{4}{5}$ Bicuspid	-11°	+2°	.114"	R	300-72-1990
4/5 4/5 Bicuspid	-11	+2	.114	L	300-72-2990
4/14/ Bicuspid	-11°	+2°	.114"	R	300-72-1991
4/5 4/5 With hook	-11	+2	.114	L	300-72-2991
4/5 4/5 Bicuspid	-11°	+2°	.114"	R	300-72-1930
$\frac{ -5 ^{5}}{ }$ w/hook & offset		12	.114	L	300-72-2930

Lotus Plus®

Ortho Technology's version of the Damon[®] Low Torque Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot	
1 1 Central	. 2° . 5°	+2° +5° .	.114"	R	300-72-1170	
TIT Cellulai	72	73	.114	L	300-72-2170	
2 2 Lateral	-5°	+9°	.110"	R	300-72-1270	
ZZ Lateral	-5	+9	9 .110	L	300-72-2270	
3 3 Cuspid	-9°	+5° .110"	R	300-72-1370		
S Cuspia	-9		.110	L	300-72-2370	
3 3 Cuspid	-9°	+5°	.110"	R	300-72-1371	
$\frac{3 3}{}$ with hook	-9	+3	.110	L	300-72-2371	
4/14/ Diamaid	-11°	+2°	.114"	R	300-72-1990	
$\frac{4}{5} \frac{4}{5}$ Bicuspid	-11		.114	L	300-72-2990	
4/14/ Bicuspid	-11°	+2°	.114"	R	300-72-1991	
Bicuspid with hook	-11	+2	.114	L	300-72-2991	

Lotus Plus®

Ortho Technology's version of the Damon® Super Torque Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+22°	+5°	.114"	R	300-72-1180
<u>I I Central</u>	+22	+3	.114	L	300-72-2180
2 2 Lateral	+13°	+9°	.110"	R	300-72-1280
ZZ Lateral	713	Τ3	.110	L	300-72-2280
3 3 Cuspid	+11°	+5°	.110"	R	300-72-1380
3 Cuspid	+11	+5		L	300-72-2380
3 3 Cuspid	+11°	+5°	.110"	R	300-72-1381
with hook	+11	+3	.110	L	300-72-2381
4/ 4/ p:i-l	-11°	+2°	114"	R	300-72-1990
$\frac{4}{5}$ Bicuspid	-11	+2	.114"	L	300-72-2990
4/14/ Bicuspid	-11°	+2°	.114"	R	300-72-1991
$\frac{4}{5}$ Bicuspid with hook	-11	TZ	.114	L	300-72-2991

© 2015 Ortho Technology, Inc. Lotus Plus is a registered trademark of Ortho Technology. *The Ortho Technology version of the Damon Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors. Damon is a registered trademark of Ormco Corporation.

Patient Kits Description

Upper/Lower 5x5 No Hooks Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks on 3, 4, & 5

Item Numbers .022" Slot 6322-1000

6322-1300 6322-1345

▼ Mandibular Arc	h Torq.	Ang.	Width	R/L	.022" Slot
1 Central	-3°	+2°	.102"	R	300-72-3190
1 Centrui	-5	+2	.102	L	300-72-4190
2 2 Lateral	-3°	+4°	.102"	R	300-72-3290
Z Lateral	-3	+4	.102	L	300-72-4290
3 3 Cuspid	. 70		11 4"	R	300-72-3390
3 3 Cuspid	+7°	+5°	.114"	L	300-72-4390
Cuspid	+7°	+5°	.114"	R	300-72-3391
3 with hook	(+/	+5	.114	L	300-72-4391
1st Bicuspic	1	4.0	22.41	R	300-72-3492
4 4 with offse		+4°	.114"	L	300-72-4492
4 4 1st Bicuspie		+4°	.114"	R	300-72-3493
4 4 w/hook & off	set	T*	.114	L	300-72-4493
2nd Bicuspi	d 1	-17° +4°	11.41	R	300-72-3592
5 5 with offse			.114"	L	300-72-4592
2nd Bicuspi	d			R	300-72-3593
5 5 2 W/hook & offs		+4°	.114"	L	300-72-4593

Patient Kits

Upper/Lower 5x5 No Hooks Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks on 3, 4, & 5

Item Numbers

6422-1000 6422-1300 6422-1345

▼ Ma	ndibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
11	Central	-11°	+2°	+2° .102"	R	300-72-3170
-	Certifian		12	.102	L	300-72-4170
2 2	Lateral	-11°	+4°	.102"	R	300-72-3270
2 2	Laterai	-11	+4	.102	L	300-72-4270
	Cuspid	0°	+5°	.114"	R	300-72-3370
3 3	Cuspiu	U	+5	.114	L	300-72-4370
3 3	Cuspid	0°		+5° .114"	R	300-72-3371
3 3	with hook	U	+3		L	300-72-4371
4 4	1st Bicuspid	-12°	+4°	.114"	R	300-72-3492
4 4	with offset	-12	T*	.1114	L	300-72-4492
4 4	1st Bicuspid	-12°	+4°	.114"	R	300-72-3493
4 4	w/hook & offset	-12	T -1	.114	L	300-72-4493
5 5	2 nd Bicuspid	-17°	+4°	114"	R	300-72-3592
2 5	with offset	-17	+4	4° .114"	L	300-72-4592
	2 nd Bicuspid	-17°	+4°	114"	R	300-72-3593
5 5	w/hook & offset	-17	+4	.114"	L	300-72-4593

Patient Kits

Upper/Lower 5x5 No Hooks Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks on 3, 4, & 5

Item Numbers

6522-1000 6522-1300 6522-1345

opper, zemer 5x5com cm 5, 1, 4 5						
▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot	
1 1 Central	-3°	+2°	.102"	R	300-72-3190	Z S
1 1 Central	-5	TZ	.102	L	300-72-4190	Standard Torque Rx
2 2 Lateral	-3°	+4°	.102"	R	300-72-3290	ıdaı Je F
Z Laterai	-5	+4	.102	L	300-72-4290	≈ 3
3 3 Cuspid	.120	+5°	114"	R	300-72-3380	
3 3 Cuspid	+13°	+3	-5° .114"		300-72-4380	
3 3 Cuspid with hook	+13°	1 2°	+5° .114" -		300-72-3381	
3 3 with hook	713	+3			300-72-4381	
4 4 1st Bicuspid	-5°	+4°	.114"	R	300-72-3482	
4 4 with offset	-5	+4	.114	L	300-72-4482	
1st Bicuspid	-5°	+4°	.114"	R	300-72-3483	
4 4 w/hook & offset		' '		L	300-72-4483	
5 5 2 nd Bicuspid	-17°	+4°	.114"	R	300-72-3592	Б.
o with offset	-17	T4	.114	L	300-72-4592	tan
2nd Bicuspid	-17°	+4°	.114"	R	300-72-3593	Standard Torque Rx
5 5 w/hook & offset	-17	+4	.114	L	300-72-4593	≥ 0

BRACKET SYSTEMS

Comfort and Confidence Lotus Plus Passive Self-Ugating Bracket System Ask Your Doctor Today About Lotus Plus Braces

Lotus Plus Wall Poster 24" x 30" Full Color Item #: LOTUS-PLUS-POSTER 1 poster \$9.95 ea.

Marketing Tools to GROW Your Practice!





Lotus Plus® Typodonts



On Our Website Lotus Plus® Instructional Video





Scan the QR Code or Go Online www.orthotechnology.com/LotusPlus

Check out all of our online resources - videos, illustrated instructions, color code charts, and FAQs.

BRACKETSYST

Lotus Plus

Passive Self-Ligating Bracket System Suggested Order, Duration, and Objective of Basic Archwire Sequence

Phase I Round Archwires	Duration	Objective
.013" or .014" TruFlex™ Thermal Nickel Titanium Universal Form .013" Item #: 2001-101 .014" Item #: 2001-102	10 Weeks	Cases with severe crowding
.018" TruFlex™ Thermal Nickel Titanium Universal Form	10 Weeks	Level and align Begin arch development Solve 90% of rotations

TruFlex™ Thermal Nickel Titanium

Round 1 pack......\$21.50/pk 10 per pack

Phase II **Rectangular Archwires**

Item #: 2001-104

Item #: 2001-118

.014" x .025" TruFlex™ Thermal Nickel Titanium Universal Form Item #: 2001-114	10-20 Weeks	Complete leveling and alignment Resolve remaining rotations Start control of torque Consolidation of the anterior spaces Continue arch development without expanders or quad-helix
.018" x .025" TruFlex™	4-6 Weeks	Fine tune torque control Continued development of the arch

TruFlex™ **Thermal Nickel Titanium**

■ Rectangle	
1 pack	\$26.25/pk
10 packs	
20 or more	\$18.15/pk
10 per pack	

Phase III **Mechanics at Work**

Thermal Nickel Titanium Universal Form

.019" x .025"	20-30 Weeks	Finish torque control
TruForce™		Consolidation of posterior spaces
Stainless Steel Universal Form		Adjustments to buccal, lingual,
Item #: 2003-119		anterior, and posterior
		Coordination of the arch form

TruForce™ Stainless Steel

■ Rectangle	
1 pack	\$7.10/pk
5 or more 10 per pack	\$5.95/pk
BetaForce™	
Beta Titanium	

■ Rectangle	
1 pack	\$45.99/pk
5 or more	\$42.99/pk
10 per pack	

Phase IV **Finishing**

.019" x .025"
TruForce™ Stainless Steel
Universal Form
Item #: 2003-119
- or -
.019" x .025"
BetaForce™ Beta Titanium
Universal Form

Item #: 2002-119

10 Weeks

Specific detailing of patient arch form


Above chart is suggested order, duration, and objective of basic archwire sequence and by no means a replacement for professional expertise. © 2015 Ortho Technology, Inc. BetaForce, TruFlex, and TruForce are trademarks of Ortho Technology.











Smooth, Strong and Contoured Tie-Wings

For added patient comfort

Archwire Slot Design Diamond and heat polished for superior sliding mechanics

Highest Quality Material Made from high quality monocrystalline sapphire

Unique Bonding Base Zirconia spheres provide superior shear bond strength and predictable debonding

Ultra-Smooth, Strong Hooks

For added patient comfort and secure engagement of elastics

Removable Color Coded ID Dot

For easy bracket selection


Ample Tie-Wing Area

For easy bracket ligation

Contour Base

Provides simplified and accurate bracket placement







© 2015 Ortho Technology, Inc. PURE is a registered trademark of Ortho Technology. U.S. and International patents pending. Inspire ICE is a trademark of Ormco/Sybron Dental Specialties.



PURE® Patient Aesthetics




The macro photograph shows the clarity of the flawless PURE brackets. Made from pure monocrystalline sapphire for a bracket system with the highest transparency available for outstanding aesthetics. Nearly invisible against any tooth shade and stays crystal clear without staining.

PURE® Patient Comfort



Each PURE bracket has been diamond and heat polished to produce rounded, smooth corners for uncompromised patient comfort. The small bracket size and low profile design also provide a desirable bracket both clinicians and patients will appreciate. The PURE brackets are also functional, the tiewings are easy to ligate, exceptionally strong and highly fracture resistant.

PURE® Bond Strength



SEM photographs show the unique Double Fusing Method™ (DFM). The PURE bonding base is coated with zirconia spheres creating millions of undercuts that mechanically lock with the bracket

adhesive. This process allows for secure retention during the length of treatment and provides for predictable debonding at the end of treatment.

PURE® Contoured Base



The anatomically contoured bracket base allows for quick and easy bracket placement. The customized tooth specific design fits the tooth curvatures perfectly to increase bond strength and make bonding appointments a breeze.

PURE® Predictable Debonding



Remove all ligatures and wires. Remove flash from around the brackets with a high speed debonding bur.



Instruct the patient to bite down onto a cotton roll in the area brackets are being removed. Grasp the bracket with a straight or angled debonding plier at the bracket-adhesive interface. The bracket can be held along the occlusal-gingival or the mesial-distal edge. Gently squeeze the plier and gradually increase pressure until the bracket comes free from the tooth. **Do not twist or pull the bracket.**



The bracket is designed to release at the bracket-adhesive interface and will leave a layer of adhesive on the enamel for safety. This means that most of the adhesive will remain on the tooth after debonding.

Remove remaining adhesive with a high speed debonding bur or an adhesive removing plier. After bulk adhesive is removed follow up with a high or low speed enamel polishing bur.

BRACKETSYSTEMS



PURE®

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	ral +12° +5°	+5°	.126"	R	S18001	S22101
TIT Central	T12	+3	.120	L	S18002	S22102
2 2 Lateral	+8°	+9° .116	.116"	R	S18011	S22111
ZZ Lateral	T-0	T 9	+9 .116	L	S18012	S22112
3 3 Cuspid	0°	+11° .126"	126"	R	S18021	S22121
$\frac{3 3}{}$ with hook	0		L	S18022	S22122	
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.126"	U	\$18030	S22130
4/4/ Bicuspid	-7°	0° .126"	.126"	R	S18031	S22131
$\frac{ f ^{2}}{ f ^{5}}$ with hook	'		.120	L	S18032	S22132

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Single Patient Kit Description Upper/Lower 3x3

Opper/Lower 3x3	
1 kit (12 brackets)	^{\$} 134.40 kit
17 kits (204 brackets)	\$128.40 kit
34 kits (408 brackets)	
Upper/Lower 5x5	
1 kit (20 brackets)	\$224.00 kit
10 kits (200 brackets)	\$214.00 kit
20 kits (400 brackets)	\$194.00 kit
Upper 5-5	
1 kit (10 brackets)	\$112.00 kit
20 kits (200 brackets)	\$107.00 kit
40 kits (400 brackets)	
Individual Brackets	
1 pack (10 brackets)	\$112.00/pk (\$11.20 ea.)
20 packs (200 brackets)	\$107.00/pk (\$10.70 ea.)
40 packs (400 brackets) 10 per pack	\$97.00/pk (\$9.70 ea.)

Patient Kits	Item Numbers		
Description	.018" Slot	.022" Slot	
Upper/Lower 3x3 No Hooks	S18100	S22100	
Upper/Lower 3x3 Hooks on Cuspids	S18300	S22300	
Upper/Lower 5x5 No Hooks	S18310	S22310	
Upper/Lower 5x5 Hooks on Cuspids	S18330	S22330	
Upper/Lower 5x5 Hooks 3, 4, & 5	S18340	S22340	
Upper 5-5 No Hooks	S18320	S22320	
Upper 5-5 Hooks on Cuspids	S18360	S22360	
Upper 5-5 Hooks 3, 4, & 5	S18370	S22370	

▼ Ma	ndibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1/2 1/2	Lower Anterior	0°	0°	.112"	U	S18041	S22141
3 3	Cuspid	-11°	+5°	.126"	R	S18051	S22151
3 3	with hook	-11	Τ)	.120	L	S18052	S22152
4 4	1 st Bicuspid	-17°	0°	.126"	U	S18060	S22160
4 4	1st Bicuspid	-17°	0°	.126"	R	S18061	S22161
4 4	with hook	-17	U	.120	L	S18062	S22162
5 5	2 nd Bicuspid	-22°	0°	.126"	U	S18070	S22170
	2 nd Bicuspid	-22°	0°	.126"	R	S18071	S22171
5 5	with hook	-22	U	.120	L	S18072	S22172

PURE®

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+17°	+4°	.126"	R	S22801
Central	T17	T*	.120	L	S22802
2 2 Lateral	+10°	+8°	.116"	R	S22811
ZZ Lateral	+10	+0	1.110	L	S22812
3 3 Cuspid with hook	0°	+8°	.126"	R	S22821
$\frac{313}{}$ with hook	0	+0		L	S22822
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.126"	U	S22833
Bicuspid with hook	-7°	0°	.126"	R	S22831
with hook	-7	U	.120	L	S22832

^{*}The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Patient Kits Description	Item Numbers .022" Slot
Upper/Lower 3x3 No Hooks	S22810
Upper/Lower 3x3 Hooks on Cuspids	S22800
Upper/Lower 5x5 No Hooks	S22820
Upper/Lower 5x5 Hooks on Cuspids	S22830
Upper/Lower 5x5 Hooks 3, 4, & 5	S22840
Upper 5-5 No Hooks	S22850
Upper 5-5 Hooks on Cuspids	S22860
Upper 5-5 Hooks 3, 4, & 5	S22870

▼ Mai	ndibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
1/2 1/2	Lower Anterior	-6°	0°	.112"	U	S22841
3 3	Cuspid	0°	. 2° 126"	R	S22851	
3 3	with hook	0	+3	+3° .126"	L	S22852
111	1 et Di accessi d	pid -12° +2° .	.126"	R	S22863	
4 4 1 st	1 st Bicuspid		+2	.120	L	S22864
4 4	1 st Bicuspid	-12°	+2°	.126"	R	S22861
4 4	with hook	-12	+2	.120	L	S22862
5 5	2nd Disuspid	-17°	+2°	.126"	R	S22873
د اد	2 nd Bicuspid	-17	+2	.126	L	S22874
5 5	2 nd Bicuspid	1.70	·17° +2° .1	120"	R	S22871
5 5	with hook	-17		.126"	L	S22872

Marketing Tools to GROW Your Practice!





PURE Patient Consultation Center with 25 Full Color Brochures Item #: S200

Stand design may vary.



PURE Wall Poster 24" x 30" Full Color Item #: PURE-POSTER





PURE Patient Brochure Refills 1 Refill Pack of 25 Full **Color Trifold Brochures** Item #: S225

PURE®

Standard Edgewise

Patient Kits Item Numbers .022" Slot Description .018" Slot Upper/Lower 5x5 No Hooks S18350 S22350 Upper/Lower 5x5 Hooks on Cuspids S18450 S22450

Standard Edgewise	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Centrals	0°	0°	.126"	U	S18400	S22400
Laterals	0°	0°	.116"	U	S18410	S22410
Cuspids	0°	0°	.126"	U	S18420	S22420
Cuspids	0°	0°	.126"	UR/ _{LL}	S18421	S22421
with hook			.120	UL/ _{LR}	S18422	S22422
Bicuspids	0°	0°	.126"	U	S18430	S22430
Lower Anterior	0°	0°	.112"	U	S18440	S22440





Rounded Corners

For added patient comfort

Archwire Slot Design


Enhances sliding mechanics

Extremely resistant to slot-wear, cracking, and fracturing

Highest Quality Material

Made from 99.9% pure polycrystalline alumina for unsurpassed strength

Compound Contour Base Provides an ideal bracket fit



Removable Color Coded ID Dot

For easy bracket selection

Translucent

For pleasing aesthetics

Dove Tail Grooves


Provide superior mechanical lock

Generous Tie-Wing Area For easy bracket ligation

Contoured Tie-Wing

For added patient comfort

Identification System



© 2015 Ortho Technology, Inc. Reflections is a trademark of Ortho Technology. CLARITY is a trademark of 3M Unitek.



BRACKETSYS


Reflections

Ceramic Bracket System

Proven Clinical Performance

The Reflections Ceramic Bracket System offers a unique combination of aesthetics, function, and strength. The stateof-the-art dovetail base design has generous undercuts which grab any adhesive firmly and ensure maximum mechanical bond strength, yet debond easily. Reflections Brackets are made from 99.9% pure polycrystalline alumina ensuring the translucency patients want. Reflections Brackets are designed to decrease the occurrence of bracket fracture, increase the torque control, increase the sliding mechanics, and provide the clinical functionality doctors need. Reflections Brackets are the clear alternative for aesthetic treatment.





Rounded Corners

For added patient comfort



Archwire Slot Design

Enhances sliding mechanics Extremely resistant to slot-wear, cracking, and fracturing



Compound Contour Base

Provides ideal bracket fit



Dove Tail Grooves

Provide superior mechanical lock



Generous Contoured Tie-Wing Area

For easy bracket ligation and added patient comfort



Translucent Highest Quality Material

99.9% pure polycrystalline alumina for unsurpassed strength and pleasing aesthetics

Reflections™

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.155"	R	REF-001	REF-101
	712	Τ .	.133	L	REF-002	REF-102
2 2 Lateral	+8°	+9°	.133"	R	REF-011	REF-111
ZZ Lateral	Τ0	+9	.133	L	REF-012	REF-112
2/2 Cuspid	. 0° +11°		.131"	R	REF-021	REF-121
$\frac{3 3}{}$ with hook	0	+11	.131	L	REF-022	REF-122
Bicuspid with hook	-7°	0°	.130"	U	REF-031	REF-131
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	REF-030	REF-130

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Patient Kits	Item Numbers				
Description	.018" Slot	.022" Slot			
Upper/Lower 3x3 Hooks on Cuspids	RF18-300	RF22-300			
Upper/Lower 5x5 Hooks on Cuspids	RF18-1300	RF22-1300			
Upper/Lower 5x5 Hooks 3, 4, & 5	RF18-1345	RF22-1345			
Upper 5-5 Hooks on Cuspids	RF18-2300	RF22-2300			
Upper 5-5 Hooks 3, 4, & 5	RF18-2345	RF22-2345			

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Lower Anterior	0°	0°	.110"	U	REF-041	REF-141
Cuspid	-11°	+5°	.133"	R	REF-051	REF-151
3 3 with hook	-11	Τ .	.133	L	REF-052	REF-152
1st Bicuspid with hook	-17°	0°	.133"	U	REF-061	REF-161
4 4 1 st Bicuspid	-17°	0°	.133"	U	REF-060	REF-160
5 5 2nd Bicuspid with hook	-22°	0°	.133"	U	REF-071	REF-171
5 5 2nd Bicuspid	-22°	0°	.133"	U	REF-070	REF-170



BRACKET SYSTEMS



Reflections™

Ortho Technology's version of the Roth Prescription*
"V" Slot (.025" x .025")

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 Central	+12°	+5°	.155"	R	REF-601
111 Central	+12	+3	.133	L	REF-602
2 2 Lateral	+8°	+9°	.133"	R	REF-611
ZZ Lateral	+0	+9	.133	L	REF-612
Cuspid	0°	+11°	.131"	R	REF-621
$\frac{3}{3}$ Cuspid with hook	0			L	REF-622
$\frac{4/5}{5}$ Bicuspid with hook	-7°	0°	.130"	U	REF-631
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	REF-630

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Single Patient Kit Description Upper/Lower 3x3

Opper/Lower 3x3	
1 kit (12 brackets)	\$92.28 kit
17 kits (204 brackets)	\$83.88 kit
34 kits (408 brackets)	\$80.28 kit
Upper/Lower 5x5	
1 kit (20 brackets)	
10 kits (200 brackets)	\$139.80 kit
20 kits (400 brackets)	\$133.80 kit
Upper 5-5	
1 kit (10 brackets)	\$76.90 kit
20 kits (200 brackets)	\$69.90 kit
40 kits (400 brackets)	\$66.90 kit
Individual Brackets	
1 pack (10 brackets)	\$76.90/pk (\$7.69 ea.)
20 packs (200 brackets)	\$69.90/pk (\$6.99 ea.)
40 packs (400 brackets)	\$66.90/pk (\$6.69 ea.)
10 per pack	

Patient Kits Description Upper/Lower 5x5 Hooks on Cuspids

.022" Slot RF322-1300 RF322-1345

Item Numbers

Upper/Lower 5x5 Hooks 3, 4, & 5 RF322-1345

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lower Anterior	0°	0°	.110"	U	REF-640
Cuspid	-11°	+5°	.133"	R	REF-651
3 3 with hook	-11	+3	.133	L	REF-652
4 4 1st Bicuspid with hook	-17°	0°	.133"	J	REF-661
4 4 1st Bicuspid	-17°	0°	.133"	U	REF-660
5 5 2 nd Bicuspid with hook	-22°	0°	.133"	U	REF-671
5 5 2 nd Bicuspid	-22°	0°	.133"	U	REF-670

Reflections™

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+17°	+4°	.155"	R	REF-801
	T17	T4	.133	L	REF-802
2 2 Lateral	+10°	+8°	.133"	R	REF-811
ZZ Lateral	+10	+0	0 .133	L	REF-812
3 3 Cuspid	0°	+8° .1	.131"	R	REF-821
$\frac{3 3}{}$ with hook	0		.131	L	REF-822
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	REF-830
Bicuspid with hook	-7°	0°	.130"	U	REF-831

*The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Patient Kits Description

Upper/Lower 3x3 Hooks on Cuspids Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

RFM18-300 RFM22-300 RFM18-1300 RFM22-1300 RFM18-1345 RFM22-1345

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
Lower Anterior	-6°	0°	.110"	U	REF-841
3 3 Cuspid	0°	+3°	.133"	R	REF-851
with hook	0 +3 .133	.133	L	REF-852	
4 4 Bicuspid	-12°	+2°	.133"	R	REF-861
4 with hook	-12	+2	.133	L	REF-862
4 4 Bicuspid	-12°	+2°	.133"	R	REF-863
4 4 Bicuspiu	-12			L	REF-864
5 5 2nd Bicuspid	-17°	+2°	.133"	R	REF-871
yith hook	-17	+2	.133	L	REF-872
5 5 2 nd Bicuspid	-17°	+2°	.133"	R	REF-873
J J Z Bicuspia	-17	+2	.133	L	REF-874

Reflections™


Standard Edgewise

Patient Kits Description

Upper/Lower 5x5 No Hooks

Item Numbers .022" Slot RF222-1000

Standard Edgewise	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Centrals	0°	0°	.155"	U	REF-200	REF-300
Laterals	0°	0°	.131"	U	REF-210	REF-310
Cuspids	0°	0°	.131"	U	REF-220	REF-320
Bicuspids	0°	0°	.131"	U	REF-230	REF-330
Lower Anterior	0°	0°	.110"	U	REF-240	REF-340



Encore!

Metal-Lined Ceramic Bracket System

The Aesthetics of Ceramics with the Functionality of Metal

Clarity™ Metal-Reinforced by 3M Unitek


Rounded Corners For added patient comfort

Metal-Lined Archwire Slot For improved sliding mechanics and to withstand higher torquing forces

Nickel Free for Sensitive Patients

Highest Quality Material Made from 99.9% pure polycrystalline alumina for unsurpassed strength

> **Compound Contour Base** Provides for an ideal bracket fit



Removable Color Coded ID Dot

For easy bracket selection

Translucent

For pleasing aesthetics

Dove Tail Grooves

Provide superior mechanical lock


Generous Tie-Wing Area

For easy bracket ligation

Contoured Tie-Wing

For added patient comfort





© 2015 Ortho Technology, Inc. Encore! is a trademark of Ortho Technology. Clarity is a trademark of 3M Unitek.



BRACKETSYSTEMS



Single Patient Kit Description Upper/Lower 3x3

1 kit (12 brackets)	\$129.60 kit
17 kits (204 brackets)	\$118.80 kit
Upper/Lower 5x5	
1 kit (20 brackets)	\$216.00 kit
10 kits (200 brackets)	\$198.00 kit
Individual Brackets	
1 pack (10 brackets)	\$108.00/pk (\$10.80 ea.)
20 packs (200 brackets)	\$99.00/pk (\$9.90 ea.)
10 per pack	• • • • • • • • • • • • • • • • • • • •

Encore!™

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 Central	+12°	+12° +5° .	.155"	R	MC-001	MC-101
	T12	+3	.133	L	MC-002	MC-102
2 2 Lateral	+8°	+8° +9° .	.133"	R	MC-011	MC-111
ZZ Lateral	T-0	T 9	.133	L	MC-012	MC-112
2 Cuspid	0°	0° +11° .131"	121"	R	MC-021	MC-121
$\frac{3 3}{}$ with hook	0		L	MC-022	MC-122	
$\frac{4/5}{5}$ Bicuspid with hook	-7°	0°	.130"	U	MC-031	MC-131
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	MC-030	MC-130

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Patient Kits Description

Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers .018" Slot .022" Slot

MC18-1300 MC22-1300 MC18-1345 MC22-1345

▼ Ma	ndibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1/2 1/2	Lower Anterior	0°	0°	.110"	U	MC-041	MC-141
3 3	Cuspid	-11°	+5°	.133"	R	MC-051	MC-151
ا ا	with hook	-11	+3	.133	L	MC-052	MC-152
4 4	1st Bicuspid with hook	-17°	0°	.133"	U	MC-061	MC-161
4 4	1 st Bicuspid	-17°	0°	.133"	U	MC-060	MC-160
5 5	2 nd Bicuspid with hook	-22°	0°	.133"	U	MC-071	MC-171
5 5	2 nd Bicuspid	-22°	0°	.133"	U	MC-070	MC-170

Encore!™

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 Central	+17°	+4°	.155"	R	MC-701	MC-801
111 Central	T17	T*	.133	L	MC-702	MC-802
2 2 Lateral	+10°	+8°	.133"	R	MC-711	MC-811
ZZ Lateral	+10	+8 .1	.133	L	MC-712	MC-812
3 3 Cuspid	0° +	+8°	3° .131"	R	MC-721	MC-821
$\frac{3 3}{}$ with hook	U	+0	.131	L	MC-722	MC-822
$\frac{4/5}{5}$ Bicuspid with hook	-7°	0°	.130"	U	MC-731	MC-831
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	MC-730	MC-830

^{*}The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Patient Kits Description

Upper/Lower 3x3 No Hook Upper/Lower 3x3 Hooks on Cuspids Upper/Lower 5x5 No Hook Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers .018" Slot .022" Slot MCM18-100 MCM22-100 MCM22-300 MCM18-300 MCM18-1000 MCM22-1000

MCM22-1300

MCM22-1345

MC-874

MCM18-1300

MCM18-1345

MC-774

R/L

▼ Mandibular Arch .022" Slot Width .018" Slot Torq. Ang. Lower -6° 0° .110" U MC-741 MC-841 Anterior MC-751 MC-851 Cuspid 3 3 0° +3° .133" with hook MC-752 MC-852 Bicuspid MC-761 MC-861 4 4 -12° +2° .133" with hook MC-762 MC-862 R MC-763 MC-863 4 4 -12° +2° .133" Bicuspid L MC-764 MC-864 2nd Bicuspid R MC-771 MC-871 5 5 -17° +2° .133" with hook Т MC-772 MC-872 R MC-773 MC-873 5 5 2nd Bicuspid -17° +2° .133"

Avalon®

Metal-Lined Composite Bracket System

The Aesthetics Patients Want with the Functionality

Doctors Need™

Single Patient



- Rounded corners for added patient comfort
- Removable color coded ID dot
- Highest quality material made from medical grade polyurethane for superior strength, translucent for pleasing aesthetics
- Metal-lined archwire slot for improved sliding mechanics and to withstand higher torquing forces. Nickel free for sensitive patients
- Compound contour base provides for an ideal bracket fit. Bonds and debonds just like metal
- Generous tie-wing area for easy bracket ligation. Contoured design for added patient comfort
- Dove tail grooves with "micro rock" coating provide superior mechanical lock

Avalon®

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.155"	R	A18-001	A22-001
1 Central	T12	Τ)	.133	L	A18-002	A22-002
2 2 Lateral	+8°	° +9° .133"	R	A18-011	A22-011	
Z Z Lateral	+0	+9	.133	L	A18-012	A22-012
2 2 Cuspid	0°	0° +11°	1° .131"	R	A18-021	A22-021
$\frac{3 3}{}$ with hook	0	+11	.131	L	A18-022	A22-022
$\frac{4/5}{5}$ Bicuspid with hook	-7°	0°	.130"	U	A18-031	A22-031
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	A18-030	A22-030

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

†0° torque on cuspids (3s).

© 2015 Ortho Technology, Inc. Avalon is a registered trademark of Ortho Technology.

Kit Description	
Upper/Lower 3x3	
1 kit (12 brackets)	\$80.28 kit
17 kits (204 brackets)	\$71.88 kit
Upper/Lower 5x5	

Individual Brackets



Patient Kits Description

Upper/Lower 3x3 Hooks on Cuspids Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

ltem	Numbers
.018" Slot	.022" Slot
AV18-300 [†]	AV22-300 ⁺

AV18-1300[†] AV22-1300[†] AV18-1345[†] AV22-1345[†]

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Lower Anterior	0°	0°	.110"	U	A18-041	A22-041
Lower Anterior	-6°	0°	.110"	С	A18-741	A22-841
Cuspid	-11°	+5°	.133"	R	A18-051	A22-051
3 3 with hook	-11	+5	.133	L	A18-052	A22-052
Cuspid	0°	+3°	.133"	R	A18-751	-
3 3 with hook	0	+3	.133	L	A18-752	_
4 4 1st Bicuspid with hook	-17°	0°	.133"	U	A18-061	A22-061
4 4 1st Bicuspid	-17°	0°	.133"	U	A18-060	A22-060
5 5 2nd Bicuspid with hook	-22°	0°	.133"	U	A18-071	A22-071
5 5 2 nd Bicuspid	-22°	0°	.133"	U	A18-070	A22-070





Rounded Corners

For added patient comfort

Archwire Slot Design Enhances sliding mechanics

Limances sharing mechanics


Extremely resistant to slot-wear, cracking, and fracturing

Highest Quality Material Made from polyurethane for

unsurpassed strength

Compound Contour Base Provides an ideal bracket fit

> Bonds and Debonds Just Like Metal



"Micro Rock" MagnificationProvides millions of additional undercuts to grip the adhesive

No Plastic Primers Needed for Bonding

Removable Color Coded ID Dot

For easy bracket selection

Translucent

For pleasing aesthetics

Generous Tie-Wing Area

For easy bracket ligation

Dove Tail Grooves

Provide superior mechanical lock and easy debonding

Contoured Tie-Wing

For added patient comfort





@ 2015 Ortho Technology, Inc. Ortho Flex is a trademark of Ortho Technology.



BRACKETSYS



Single Patient Kit Description Upper/Lower 3x3

	2 bracke			\$47.88	kit
17 kits	(204 bra	ackets	i)	\$43.08	kit

Upper/Lower 5x5

1 kit (20 brackets)	\$79.80	kit
10 kits (200 brackets)	\$71.80	kit

Individual Brackets

1 pack (10 brackets)......\$39.90/pk (\$3.99 ea.) 20 packs (200 brackets)\$35.90/pk (\$3.59 ea.)

10 per pack

OrthoFlex[™]

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.155"	R	COM-001	COM-101
111 Central	+12	+3	.133	L	COM-002	COM-102
2 2 Lateral	+8°	+9°	.133"	R	COM-011	COM-111
22 Lateral	+0	+9	.133	L	COM-012	COM-112
2 Cuspid	0°	+11°	.131"	R	COM-021	COM-121
$\frac{3 3}{}$ with hook	0	+11	.131	L	COM-022	COM-122
$\frac{4/5}{5}$ Bicuspid with hook	-7°	0°	.130"	U	COM-031	COM-131
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	COM-030	COM-130

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Patient Kits	Item Numbers			
Description	.018" Slot	.022" Slot		
Upper/Lower 3x3 Hooks on Cuspids		CM22-300		
Upper/Lower 5x5 No Hooks		CM22-1000		
Upper/Lower 5x5 Hooks on Cuspids	CM18-1300	CM22-1300		
Unner/Lower 5x5 Hooks 3 4 & 5	CM18-1345	CM22-1345		

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Lower Anterior	0°	0°	.110"	U	COM-041	COM-141
Cuspid	-11°	+5°	.133"	R	COM-051	COM-151
3 with hook	-11	+3	.133	L	COM-052	COM-152
4 4 1st Bicuspid with hook	-17°	0°	.133"	U	COM-061	COM-161
4 4 1st Bicuspid	-17°	0°	.133"	U	COM-060	COM-160
5 5 2nd Bicuspid with hook	-22°	0°	.133"	U	COM-071	COM-171
5 5 2 nd Bicuspid	-22°	0°	.133"	U	COM-070	COM-170

OrthoFlex[™]

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+17°	+4°	.155"	R	COM-701	COM-801
TIT Central	+17	+4	.133	L	COM-702	COM-802
2 2 Lateral	+10°	+8°	.133"	R	COM-711	COM-811
Z Z Lateral	+10	+0	.133	L	COM-712	COM-812
3 3 Cuspid	0°	+8°	.131"	R	COM-721	COM-821
$\frac{3 3}{}$ with hook	U	+0	.131	L	COM-722	COM-822
Bicuspid universal	-7°	0°	.130"	U	COM-730	COM-830
Bicuspid with hook	-7°	0°	.130"	U	COM-731	COM-831

The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Patient Kits Description

Upper/Lower 3x3 Hooks on Cuspids Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers .018" Slot .022" Slot CMM18-300 CMM22-300 CMM18-1300 CMM22-1300

CMM18-1345 CMM22-1345

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Lower Anterior	-6°	0°	.110"	U	COM-741	COM-841
3 3 Cuspid	0°	+3°	.133"	R	COM-751	COM-851
with hook	U	Τ.	.133	L	COM-752	COM-852
1st Bicuspid	-12°	+2°	.133"	R	COM-761	COM-861
4 4 with hook	-12	+2	.133	L	COM-762	COM-862
1st Disuspid	-12°	+2°	122"	R	COM-763	COM-863
4 4 1st Bicuspid	-12	+2	.133"	L	COM-764	COM-864
5 5 2 nd Bicuspid	-17°	+2°	.133"	R	COM-771	COM-871
3 3 with hook	-17	ΤΖ	.133	L	COM-772	COM-872
5 5 2 nd Bicuspid	-17°	+2°	.133"	R	COM-773	COM-873
3 / 3 Z Bicuspiu	<u>''</u>	'-	55	L	COM-774	COM-874

OrthoFlex[™]

Standard Edgewise

Standard Edgewise	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Centrals	0°	0°	.155"	U	COM-400	COM-300
Laterals	0°	0°	.131"	U	COM-401	COM-301
Cuspids	0°	0°	.131"	U	COM-402	COM-302
Bicuspids	0°	0°	.131"	U	COM-403	COM-303
Lower Anterior	0°	0°	.110"	U	COM-404	COM-304





BRACKET SYSTEMS



Compare and See the Difference for Yourself!



Tooth Tone® Coated Preformed Shorty Ties

- Durable plastic coating
- Stainless steel

.010" wire (.012" coated) Item #: 907-105 1 pack\$28.29/pk 5 or more\$26.95/pk 100 per pack



Tooth Tone® Coated Kobayashi Shorty Ties

- Durable plastic coating
- Stainless steel

.012" wire (.014" coated) Item #: 907-115 100 per pack

Tooth Tone® Coated Preformed Ligature Ties

.010" wire (.012" coated) Item #: 907-102 1 pack\$26.85/pk 5 or more\$24.15/pk 100 per pack

Tooth Tone[®]





Tooth Tone® Coated Crimpable Bi-Directional Ball Hooks

Item #: 907-002 Fits archwires up to .022" x .025" 3 or more.......559.95/pk 20 per pack



Tooth Tone® **Nickel Titanium** Open Coil **Springs**

7" Lengths	
.012" x .030" (.014" x .027" co	oated)
Item #: 907-003	
1 pack \$58.9	5/pl
5 or more \$56.9	5/nl

3 springs per pack



OrthoFlex™ Composite Buttons

■ Made in the USA

BRACKETSYSTEMS

Makes Aesthetic Treatment Complete

Tooth Tone **Coated Full Form Archwires**

- Superior aesthetics
- Resistant to staining
- Durable coating
- Low friction

Tooth Tone® **Nickel Titanium**

Round	▲ Upper	▼ Lower
.014"	6000-601	6000-701
.016"	6000-602	6000-702
.018"	6000-603	6000-703
.020"	6000-604	6000-704
■ Rectangle		
.016" x .022"	6000-607	6000-707
.018" x .018"	6000-608	6000-708
.018" x .024"	6000-609	6000-709
.019" x .027"	6000-610	6000-710
.020" x .020"	6000-612	6000-712

Note: Size shown is finished size. Tooth Tone coating adds .002" thickness to diameter.

Round

5 or more \$48.99/pk 5 or more \$52.99/pk 10 per pack

■ Rectangle

1 pack..........\$51.99/pk 1 pack..........\$55.99/pk 10 per pack

Tooth Tone® Stainless Steel

Round	▲ Upper	▼ Lower
.014"	8000-802	8000-902
.016"	8000-803	8000-903
.018"	8000-804	8000-904
.020"	8000-805	8000-905
■ Rectangle		
.016" x .022"	8000-807	8000-907
.017" x .025"	8000-809	8000-909
.018" x .018"	8000-810	8000-910
.018" x .024"	8000-812	8000-912
.019" x .027"	8000-813	8000-913
.020" x .020"	8000-817	8000-917

Note: Size shown is finished size. Tooth Tone coating adds .002" thickness to diameter.

5 or more \$30.99/pk 5 or more \$31.99/pk 10 per pack

■ Rectangle

1 pack.......\$32.99/pk 1 pack......\$33.99/pk 10 per pack





Tooth Tone® Crimpable Mini Stops

- Durable aesthetic coating
- For total aesthetic treatment
- Very comfortable to the patient
- Slides easily onto the archwire
 - 2.0 mm in length

Coated Mini Stops - Small Item #: 907-000

1 pack......**\$54.50/pk** 4 or more.........\$51.75/pk Fits round wires up to .018" 25 per pack

Coated Mini Stops - Large Item #: 907-001

Fits square wires up to .020" x .020" rectangular wires up to .019" x .025", 25 per pack



© 2015 Ortho Technology, Inc. Tooth Tone is a registered trademark of Ortho Technology.





Two-Piece Bracket System

Intelligent Hook Design Hooks are angled at 15° to reduce interference with gingival tissue

> **Smooth Facial Contours** For increased patient comfort

Vertical Scribe Line

Runs entire length of bracket base and pad for easy and more accurate placement on the long axis

Compound Contour Base Superior fit and greater contact surface area for improved retention

The Votion Bracket System utilizes Metal Injection Molding (MIM) technology to achieve extremely complex bracket structures and contours that enhance treatment efficiency.

Permanent Identification System With a removable color coded dot

and permanently recessed I.D. mark


Ample Under Tie-Wing Area Provides easy and secure ligation

Rounded Archwire Slot Contours Mesial and distal slot entry points are

rounded to prevent archwire notching when rotations are present and during retraction/protraction procedures



Identification System



© 2015 Ortho Technology, Inc. Votion is a registered trademark of Ortho Technology.



Low Profile Design

Improved aesthetics and patient comfort

Anti-Rotation

Ortho Technology's version of the Roth Prescription features built-in anti-rotation in the cuspid and bicuspid bracket archwire slots to prevent unwanted rotation during sliding mechanics

Intelligent Hook Design

Hooks are designed with a 15° angle to reduce interference with gingival tissue. The ball hook is further rounded to minimize soft tissue irritation and improve the patient experience

Votion®

Ortho Technology's version of the Roth Prescription* with Anti-Rotation

▲ Max	xillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1	Central	+12°	+5°	.154"	R	350-91-1100	350-92-1100
	Central	T12		.134	L	350-91-2100	350-92-2100
2 2	Lateral	+8°	+9°	.125"	R	350-91-1200	350-92-1200
2 2	Laterai	+0	+3	.123	L	350-91-2200	350-92-2200
3 3	Cuspid	-2°	+13°	.138"	R	350-91-1301	350-92-1301
3 3	Cuspia	-2	Τ13	.136	L	350-91-2301	350-92-2301
3 3	Cuspid	-2	+13°	.138"	R	350-91-1311	350-92-1311
3 3	with hook	-2	713		L	350-91-2311	350-92-2311
4/5 4/5	Bicuspid	-7°	٠.	.130"	R	350-91-1902	350-92-1902
<u></u>	ысиѕріи	-/	0°		L	350-91-2902	350-92-2902
4/5 4/5	Bicuspid	-7°	0°	.130"	R	350-91-1912	350-92-1912
<u>/5 /5</u>	with hook	-7	U		L	350-91-2912	350-92-2912
4/5 4/5	Bicuspid	-7°	0°	.130"	R	350-91-1922	350-92-1922
<u> </u>	w/offset pad	-/	U	.130	Ĺ	350-91-2922	350-92-2922
4/5 4/5	Bicuspid	-7°	0°	.130"	R	350-91-1932	350-92-1932
<u> </u>	w/hook & offset	-7	U	.130	L	350-91-2932	350-92-2932

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Maxillary Anti-Rotation Centrals & Laterals = 0° Cuspids = $+4^{\circ}M$ Bicuspids = $+2^{\circ}D$

Mandibular Anti-Rotation Lower Anterior = 0° Cuspids = $+2^{\circ}M$ Bicuspids = $+4^{\circ}D$

Ample Under Tie-Wing Area

Provides easy and secure ligation, and improves chairside efficiencies with fast ligature application and removal

80 Gauge Micro-Etched Mesh Bonding Base

Each base is micro-etched to ensure optimal adhesive properties and prevent bond failures for increased office productivity

Patient Kits	Item N	umbers
Description	.018" Slot	.022" Slot
Upper/Lower 5x5	350-91-2000	350-92-2000
Upper/Lower 5x5 Hooks on Cuspids	350-91-2010	350-92-2010
Upper/Lower 5x5 Hooks 3, 4, & 5	350-91-2011	350-92-2011
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	350-91-2020	350-92-2020
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	350-91-2030	350-92-2030
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	350-91-2031	350-92-2031

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot		
Lower Anterior	-1°	0°	.098"	U	350-91-6800	350-92-6800		
3 3 Cuspid	-11°	+7° .138"	.138"	R	350-91-3301	350-92-3301		
J/J Caspia			50	L	350-91-4301	350-92-4301		
3 3 Cuspid	-11°	+7°	.138"	R	350-91-3311	350-92-3311		
with hook		17	.130	L	350-91-4311	350-92-4311		
4 4 1st Bicuspid	-17°	0°	.130"	R	350-91-3402	350-92-3402		
4 4 1 bicuspiu	-17	U	.130	L	350-91-4402	350-92-4402		
4 4 1st Bicuspid	-17°	0°	.130"	R	350-91-3412	350-92-3412		
4 4 with hook	-17	0	.130	L	350-91-4412	350-92-4412		
1st Bicuspid	-17°	7° 0°	0° .130"	R	350-91-3422	350-92-3422		
4 4 w/offset pad	-17	U	.130	L	350-91-4422	350-92-4422		
4 4 1st Bicuspid	-17°	0° 13	0°	0°	.130"	R	350-91-3432	350-92-3432
4 4 w/hook & offset	-17	U	.130	L	350-91-4432	350-92-4432		
5 5 2 nd Bicuspid	-22°	0°	.130"	R	350-91-3502	350-92-3502		
J J Z bicuspiu		U	.130	L	350-91-4502	350-92-4502		
5 5 2nd Bicuspid	-22°	0°	.130"	R	350-91-3512	350-92-3512		
3 3 with hook	-22	U	.130"	L	350-91-4512	350-92-4512		
5 5 2 nd Bicuspid	-22°	0°	0° .130"	R	350-91-3522	350-92-3522		
3 3 w/offset pad	-22	U		L	350-91-4522	350-92-4522		
5 5 2nd Bicuspid	-22°	0°	.130"	R	350-91-3532	350-92-3532		
w/hook & offset	-22	U	.130	L	350-91-4532	350-92-4532		

BRACKET SYSTEMS



Single Patient Kit Description Upper/Lower 5x5

1 kit (20 brackets)	\$61.80 kit
10 kits (200 brackets)	\$55.80 kit
50 kits (1,000 brackets)	\$49.80 kit
Individual Brackets	
1 pack (10 brackets)	\$30.90/pk (\$3.09 ea.)
20 packs (200 brackets)	
100 packs (1,000 brackets)	
10 per pack	,

Votion®

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

Patient Kits Description	Item Numbers			
With Torque on Cuspids	.018" Slot	.022" Slot		
Upper/Lower 5x5	350-21-2000	350-22-2000		
Upper/Lower 5x5 Hooks on Cuspids	350-21-2010	350-22-2010		
Upper/Lower 5x5 Hooks on 3, 4, & 5	350-21-2011	350-22-2011		
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	350-21-2020	350-22-2020		
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	350-21-2030	350-22-2030		
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	350-21-2031	350-22-2031		

▲ Ma	xillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
111	Central	+17°	ı 10°	+4° .154"	R	350-21-1100	350-22-1100
11	Central	T17	T '1		L	350-21-2100	350-22-2100
2/2	Lateral	+10°	+8°	.125"	R	350-21-1200	350-22-1200
22	Laterai	+10	+0	.123	L	350-21-2200	350-22-2200
2 2	Cuanid	-7°	+8°	.138"	R	350-21-1300	350-22-1300
3 3	Cuspid	-/	+0	.136	L	350-21-2300	350-22-2300
2/2	Cuspid	-7°	+8°	.138"	R	350-21-1310	350-22-1310
3 3	with hook	-/	+0	.136	L	350-21-2310	350-22-2310
2/2	Committee	0°	+8°	.138"	R	350-21-1340	350-22-1340
3 3	Cuspid	U			L	350-21-2340	350-22-2340
2/2	Cuspid	0°	+8°	.138"	R	350-21-1350	350-22-1350
3 3	with hook	U	+0		L	350-21-2350	350-22-2350
4/ 4/	Diamarid	-7°	0°	.130"	R	350-21-1900	350-22-1900
$\frac{4}{5} \frac{4}{5}$	Bicuspid	-/	U	.130	L	350-21-2900	350-22-2900
4/5 4/5	Bicuspid	-7°	0°	120"	R	350-21-1910	350-22-1910
/5 /5	with hook	-/	U	.130"	L	350-21-2910	350-22-2910
4/5 4/5	Bicuspid	-7°	0°	.130"	R	350-21-1920	350-22-1920
<u> 15 15</u>	with offset	-/	U	0 .130	L	350-21-2920	350-22-2920
4/5 4/5	Bicuspid	-7°	0.	1.20"	R	350-21-1930	350-22-1930
<u> </u>	w/hook & offset	-/	0°	.130"	L	350-21-2930	350-22-2930

Patient Kits Description	Item Numbers			
Without Torque on Cuspids	.018" Slot	.022" Slot		
Upper/Lower 5x5	350-21-2040	350-22-2040		
Upper/Lower 5x5 Hooks on Cuspids	350-21-2050	350-22-2050		
Upper/Lower 5x5 Hooks on 3, 4, & 5	350-21-2051	350-22-2051		
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	350-21-2060	350-22-2060		
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	350-21-2061	350-22-2061		
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	350-21-2062	350-22-2062		

▼ Ma	ndibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1/2 1/2	Lower Anterior	-6°	0°	.098"	U	350-21-6800	350-22-6800
3 3	Cuspid	-6°	+3°	.138"	R	350-21-3300	350-22-3300
3 3	Cuspiu	-0	Τ.J	.130	L	350-21-4300	350-22-4300
3 3	Cuspid	-6°	+3°	.138"	R	350-21-3310	350-22-3310
2 2	with hook	-0	Τ.)	.130	L	350-21-4310	350-22-4310
3 3	Cuspid	0°	+3°	.138"	R	350-21-3340	350-22-3340
2 2	Сизріи				L	350-21-4340	350-22-4340
3 3	Cuspid	0°	+3°	.138"	R	350-21-3350	350-22-3350
ا ت	with hook		13	.130	L	350-21-4350	350-22-4350
4 4	1st Bicuspid	-12°	+2°	.130"	R	350-21-3400	350-22-3400
7 7	т ысазыа				L	350-21-4400	350-22-4400
4 4	1st Bicuspid	-12°	+2°	.130"	R	350-21-3410	350-22-3410
7 7	with hook		'-		L	350-21-4410	350-22-4410
4 4	1st Bicuspid	-12°	+2°	.130"	R	350-21-3420	350-22-3420
7 7	with offset		'-	50	L	350-21-4420	350-22-4420
4 4	1st Bicuspid	-12°	+2°	.130"	R	350-21-3430	350-22-3430
.,,	w/hook & offset				L	350-21-4430	350-22-4430
5 5	2 nd Bicuspid	-17°	+2°	.130"	R	350-21-3500	350-22-3500
	<u>'</u>				L	350-21-4500	350-22-4500
5 5	2 nd Bicuspid	-17°	+2°	.130"	R	350-21-3510	350-22-3510
	with hook				L	350-21-4510	350-22-4510
5 5	2 nd Bicuspid	-17°	+2°	.130"	R	350-21-3520	350-22-3520
3 3	with offset			, ,	L	350-21-4520	350-22-4520
5 5	2 nd Bicuspid	-17°	+2°	.130"	R	350-21-3530	350-22-3530
5 5	w/hook & offset	.,	' -	2 .130	L	350-21-4530	350-22-4530

*The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.



Votion Ideal Occlusion Typodont Item #: G50-122 1 typodont.......\$225.00 ea.



DESIGNED FOR USE WITH...

Votion[®] Bondable Buccal Tubes

See page 65



BRACKET SYSTEMS



Intelligent Hook Design

Smooth rounded hooks reduce soft tissue irritation

Cross Scribe Lines

Unique vertical and horizontal scribe lines in the center of the bracket head for superior visual ques. The patented, raised horizontal line on the bracket base further enhances visibility for optimal bonding placement accuracy

Micro-Etched Pylon Base Design

Created from the original Marquis™ base, Pinnacle's redesigned version features angled pylons that offer increased surface area for superior retention and bond strength



Permanent Identification System

With a removable color coded dot and permanently recessed I.D. mark


Rounded Archwire Slot Contours

Mesial and distal slot entry points are rounded to prevent archwire notching when rotations are present and during retraction/protraction procedures

Optimized Compound Contour Base

Brackets are individually contoured for an ideal fit to each tooth for correct positioning





© 2015 Ortho Technology, Inc. Pinnacle is a trademark of Ortho Technology.







BRACKET SYSTEMS



One-Piece Manufacturing Provides High Strength and Ultimate Base-to-Bracket Head Alignment

Pinnacle revolutionizes the original Marquis™ design, and provides increased efficiencies and productivity for doctors demanding the strength and consistency of a one-piece bracket.

Pinnacle's marked improvements stem from its manufacturing process. Through Metal Injection Molding (MIM) technology, more intricate contours and structures are achieved for reduced friction, optimum patient comfort, and enhanced visual placement marks. A specialized finishing process gives Pinnacle a high-shine surface and qualifies it as one of the most updated true twin-wing brackets in the industry.

Pinnacle™

Ortho Technology's version of the Roth Prescription* with Anti-Rotation

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.130"	R	360-91-1100	360-92-1100
TIT Central	T12	Τ.)	.130	L	360-91-2100	360-92-2100
2 2 Lateral	+8°	+9°	.109"	R	360-91-1200	360-92-1200
ZZ Lateral	70	73	.103	L	360-91-2200	360-92-2200
3 3 Cuspid	-2°	+13°	.117"	R	360-91-1301	360-92-1301
2 3 Cuspia	-2	Τ13	.117	L	360-91-2301	360-92-2301
3 3 Cuspid	-2	+13°	.117"	R	360-91-1311	360-92-1311
with hook	-2	Τ13	/	L	360-91-2311	360-92-2311
$\frac{4}{5}$ $\frac{4}{5}$ Bicuspid	-7°	0.	o° .120"	R	360-91-1902	360-92-1902
$\frac{\sqrt{5}}{\sqrt{5}}$ bicuspiu	-/	0°	.120	L	360-91-2902	360-92-2902
4/5 4/5 Bicuspid	-7°	0°	.120"	R	360-91-1912	360-92-1912
$\frac{75 75}{}$ with hook	- /	U	.120	L	360-91-2912	360-92-2912
Bicuspid Woffset pad	-7°	0°	.120"	R	360-91-1922	360-92-1922
$\frac{\sqrt{5}}{\sqrt{5}}$ w/offset pad	- /	U	.120	L	360-91-2922	360-92-2922
Bicuspid 4/5/4/5 Bicuspid w/hook & offset	-7°	0°	.120"	R	360-91-1932	360-92-1932
$\frac{ f }{ f }$ w/hook & offset	-7		.120	L	360-91-2932	360-92-2932

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Maxillary Anti-Rotation Centrals & Laterals = 0° Cuspids = +4°M Bicuspids = +2°D Mandibular Anti-Rotation Lower Anterior = 0° Cuspids = +2°M Bicuspids = +4°D

Low Profile Smooth Facial Contours

For increased patient comfort

Ample Under Tie-Wing Area

Provides easy and secure ligation

Anti-Rotation in the Slot

Ortho Technology's version of the Roth prescription features anti-rotation in the archwire slot to prevent unwanted rotation during sliding mechanics

Patient Kits	Item Numbers			
Description	.018" Slot	.022" Slot		
Upper/Lower 5x5	360-91-2000	360-92-2000		
Upper/Lower 5x5 Hooks on Cuspids	360-91-2010	360-92-2010		
Upper/Lower 5x5 Hooks 3, 4, & 5	360-91-2011	360-92-2011		
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	360-91-2020	360-92-2020		
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	360-91-2030	360-92-2030		
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	360-91-2031	360-92-2031		

▼ Ma	ndibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1/2 1/2	Lower Anterior	-1°	0°	.096"	U	360-91-6800	360-92-6800
3 3	Cuspid	-11°	+7°	.119"	R	360-91-3301	360-92-3301
J J	Сизріи		.,	5	L	360-91-4301	360-92-4301
3 3	Cuspid	-11°	+7°	.119"	R	360-91-3311	360-92-3311
3 3	with hook	-11	Τ/	.113	L	360-91-4311	360-92-4311
4 4	1st Bicuspid	-17°	0°	.119"	R	360-91-3402	360-92-3402
4 4	i Bicuspiu	-17	0	.113	L	360-91-4402	360-92-4402
4 4	1st Bicuspid	-17°	0°	.119"	R	360-91-3412	360-92-3412
4 4	with hook	-17	0	.119	L	360-91-4412	360-92-4412
4 4	1st Bicuspid	-17°	0°	.119"	R	360-91-3422	360-92-3422
4 4	w/offset pad	-17	U		L	360-91-4422	360-92-4422
4 4	1st Bicuspid	-17°	0°	.119"	R	360-91-3432	360-92-3432
4 4	w/hook & offset	-17	U	.119	L	360-91-4432	360-92-4432
5 5	2 nd Bicuspid	-22°	0°	.119"	R	360-91-3502	360-92-3502
ا د	2 bicuspiu	-2.2	U	.113	L	360-91-4502	360-92-4502
5 5	2 nd Bicuspid	-22°	0°	.119"	R	360-91-3512	360-92-3512
2 2	with hook	-22	U	.119	L	360-91-4512	360-92-4512
5 5	2 nd Bicuspid	-22°	0°	.119"	R	360-91-3522	360-92-3522
2 2	w/offset pad	-22	U	.119	L	360-91-4522	360-92-4522
5 5	2 nd Bicuspid	-22°	0°)° .119"	R	360-91-3532	360-92-3532
داد	w/hook & offset	-22	U	.119	L	360-91-4532	360-92-4532

BRACKETSYSTEMS



Single Patient Kit Description Upper/Lower 5x5

Individual Brackets	
50 kits (1,000 brackets)	\$47.80 kit
10 kits (200 brackets)	\$51.80 kit
1 kit (20 brackets)	\$55.80 kit

\$27.90/pk (\$2.79 ea.)
\$25.90/pk (\$2.59 ea.)
\$23.90/pk (\$2.39 ea.)
•

Pinnacle™

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

Patient Kits Description	Item Numbers
With Torque on Cuspids	.022" Slot
Upper/Lower 5x5	360-22-2000
Upper/Lower 5x5 Hooks on Cuspids	360-22-2010
Upper/Lower 5x5 Hooks on 3, 4, 5	360-22-2011
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	360-22-2020
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	360-22-2030
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5	360-22-2031

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+17°	+4°	+4° .130"	R	360-22-1100
1 1 Central	+17	+4	.130	L	360-22-2100
2 2 Lateral	+10°	+8°	.109"	R	360-22-1200
Z Lateral	+10	_ TO	.109	L	360-22-2200
2 2 Cuantid	-7°	+8°	.117"	R	360-22-1300
33 Cuspid	-/	+0	.117	L	360-22-2300
2/2 Cuspid	-7°	+8°	+8° .117"	R	360-22-1310
$\frac{3 3}{}$ with hook	-/	+0	.117	L	360-22-2310
2 2	0°	0° +8°	.117"	R	360-22-1340
33 Cuspid				L	360-22-2340
ala Cuspid	0°	+8°	.117"	R	360-22-1350
$\frac{3 3}{}$ with hook	0 +8	.117	L	360-22-2350	
4/ 4/ Diamaid	-7°	0°	.120"	R	360-22-1900
$\frac{4/5}{5} \frac{4/5}{5}$ Bicuspid		0	.120	L	360-22-2900
4/14/ Bicuspid	-7°	0°	.120"	R	360-22-1910
$\frac{ f ^{2}}{ f }$ with hook	-7	0	.120	L	360-22-2910
4/ 4/ Bicuspid	-7°	0°	.120"	R	360-22-1920
$\frac{ f ^{5}}{ f ^{5}}$ with offset	-7	U	.120	L	360-22-2920
4/14/ Bicuspid	7°	0.	120"	R	360-22-1930
$\frac{ f }{ f }$ w/hook & offset	-7°	0°	.120"	L	360-22-2930


Patient Kits Description Without Torque on Cuspids	Item Numbers .022" Slot
Upper/Lower 5x5	360-22-2040
Upper/Lower 5x5 Hooks on Cuspids	360-22-2050
Upper/Lower 5x5 Hooks on 3, 4, 5	360-22-2051
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	360-22-2060
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	360-22-2061
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	360-22-2062

▼ Mai	ndibular Arch	Torq.	Ang.	Width	R/L	.022" Slot				
1/2 1/2	Lower Anterior	-6°	0°	.096"	U	360-22-6800				
3 3	Cuspid	-6°	+3°	.119"	R	360-22-3300				
2 2	сизріч				L	360-22-4300				
3 3	Cuspid	-6°	+3°	.119"	R	360-22-3310				
2 3	with hook		13	5	L	360-22-4310				
3 3	Cuspid	0°	+3°	.119"	R	360-22-3340				
3 3	Сизріи	_			L	360-22-4340				
3 3	Cuspid	0°	+3°	.119"	R	360-22-3350				
3 3	with hook	0	+3	.113	L	360-22-4350				
4 4	1st Bicuspid	-12°	+2°	.119"	R	360-22-3400				
4 4	тъпсизріи	12 12 111	12	5	L	360-22-4400				
4 4	1st Bicuspid	-12°	+2°	.119"	R	360-22-3410				
4 4	with hook	-12	ΤΖ	.119	L	360-22-4410				
4 4	1st Bicuspid	-12°	+2°	.119"	R	360-22-3420				
4 4	with offset		+2	.119	L	360-22-4420				
4 4	1st Bicuspid	-12°	+2°	.119"	R	360-22-3430				
4 4	w/hook & offset	' -	12	5	L	360-22-4430				
5 5	and Discounted	-17°	+2°	.119"	R	360-22-3500				
3 3	2 nd Bicuspid	' '	TZ	5	L	360-22-4500				
5 5	2 nd Bicuspid	-17°	. 2.	.119"	R	360-22-3510				
د د	with hook	with hook	with hook	-17 +2		+2	" +2°	.119	L	360-22-4510
FF	2 nd Bicuspid	-17°	+2°	.119"	R	360-22-3520				
5 5	with offset	-17	T2	.119	L	360-22-4520				
-	2 nd Bicuspid	-17°	+2°	. 2° 110"	R	360-22-3530				
5 5	w/hook & offset	-17	+2	.119"	L	360-22-4530				

^{*}The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.



Pinnacle Ideal Occlusion Typodont Item #: G50-360 1 typodont.......\$225.00 ea.







One-Piece Bracket System

Identification System

By quadrant with a removable color coded dot and permanently by a recessed I.D. mark

Chamfered Archwire Slot

For easy archwire engagement Archwire slot is parallel to incisal edge of tooth

Highest Quality Material

Cast from the high quality 17-4 stainless steel for unsurpassed strength

Rhomboid Shape

With defined center scribe line for easy and accurate bracket placement

Superior Shear Bond Strength Makes for less bond failure

Ball Hooks

No additional charge

Contoured Tie-Wing

For added patient comfort Tie-wings are parallel to tooth anatomy

Generous Tie-Wing Area For easy bracket ligation

Ultra Low Profile For added patient comfort

Torque-In-Base


Compound-contoured design Base of bracket is parallel to the incisal edge of tooth

State-Of-The-Art Base Design



- Electrical Discharge Machining (EDM)
- Retention superior to mesh bonding pads
- The extended base of the bicuspid provides a better fit and added retention
- The micro-etched pylons are designed at an acute angle generating miniature anchorage points that embed firmly into the adhesive for maximum bond strength

Identification System



© 2015 Ortho Technology, Inc. Marquis is a trademark of Ortho Technology.

BRACKETSYSTEMS

Single Patient Kit Description

x5	
s)\$55 .80 I	kit
kets) ^{\$} 51.80 I	kit
ackets) ^{\$} 47.80 I	kit

Individual Brackets

1 pack (10 brackets)	[§] 27.90/pk (§2.79 ea.)
20 packs (200 brackets)	^{\$} 25.90/pk (^{\$} 2.59 ea.)
100 packs (1,000 brackets)	\$23.90/pk (\$2.39 ea.)
10 per pack	

Marquis™

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot						
1 Central	+12°	+5°	-5° .155" -		M18-001	M22-001						
$\frac{1}{1}$ Central	T12	T 3	.133	L	M18-002	M22-002						
2 2 Lateral	+8°	+9°	.125"	R	M18-011	M22-011						
22 Lateral	70	+3	.123	Г	M18-012	M22-012						
212 Cuspid	0°	+11°	.130"	R	M18-021	M22-021						
$\frac{3 3}{3}$ with hook	0	0	U	T11	TII	T11	TII	T11	.130	L	M18-022	M22-022
3 3 Cuspid	0°	O°	O.	O°	O°	.110	+11°	+11° .130"	R	M18-023	M22-023	
3 3 Cuspid	U	+11	.130	L	M18-024	M22-024						
4/4/ Bicuspid	-7°	0°	.125"	R	M18-031	M22-031						
$\frac{75 75}{}$ with hook	-7	U	.123	L	M18-032	M22-032						
$\frac{4}{5}$ Bicuspid	-7°	0°	.125"	U	M18-030	M22-030						

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Patient Kits Item Numbers Description .018" Slot .022" Slot Upper/Lower 5x5 No Hooks M18-1000 M22-1000 Upper/Lower 5x5 Hooks on Cuspids M18-1300 M22-1300 Upper/Lower 5x5 Hooks 3, 4, & 5 M18-1345 M22-1345

▼ Ma	andibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot		
1/2 1/2	Lower Anterior	0°	0°	.100"	U	M18-040	M22-040		
3 3	Cuspid	-11°	+5°	.125"	R	M18-051	M22-051		
3 3	with hook		+3	.123	L	M18-052	M22-052		
3 3	Cuspid	-11°	+5°	.125"	R	M18-053	M22-053		
3 3	Cuspiu	- +3	-11	-11	ΤΟ	3 1.123	L	M18-054	M22-054
4 4	1st Bicuspid	-17°	0°	.125"	R	M18-061	M22-061		
4 4	with hook	-17	0	.123	L	M18-062	M22-062		
4 4	1st Bicuspid	-17°	0°	.125"	U	M18-060	M22-060		
5 5	2 nd Bicuspid	-22°	0°	.125"	R	M18-071	M22-071		
5 5	with hook	-22	U	.123	L	M18-072	M22-072		
5 5	2 nd Bicuspid	-22°	0°	.125"	U	M18-070	M22-070		

Marquis™

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot				
1 1 Central	+17°	+5°	.155"	R	320-22-1100				
TIT Central	T17	+3	.133	L	320-22-2100				
2 2 Lateral	+10°	+9°	.125"	R	320-22-1200				
ZZ Lateral	+10	+3	+9 .123	L	320-22-2200				
3 3 Cuspid	0°	+8°	.134"	R	320-22-1350				
$\frac{3 3 }{}$ with hook			J	J	0 +0	+0 .	דכו.	L	320-22-2350
3 3 Cuspid	-7°	+8°	.134"	R	320-22-1310				
$\frac{3 3 }{}$ with hook		+6 .134	.134	L	320-22-2310				
4/4/ Bicuspid	-7°	٥٠	.125"	R	320-22-1900				
$\frac{75 75}{}$ with hook	-/	0°	U	-/ 0	0 1.123	L	320-22-2900		
$\frac{4\sqrt{5}}{5}$ Bicuspid	-7°	0°	.125"	U	320-22-5900				

^{*}The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Patient Kits Description

Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks on 3, 4, & 5 Upper/Lower 5x5 Hooks on 3 Upper/Lower 5x5 Hooks on 3, 4, & 5

Item Numbers .022" Slot

320-22-2010 320-22-2011 320-22-2050[†] 320-22-2051[†]

†0° torque on cuspids (3s)

▼ Mandibular Ar	ch Torq.	Ang.	Width	R/L	.022" Slot	
Lower Anterior	-6°	0°	.100"	U	320-22-6800	
3 3 Cuspid	-6°	+3°	.134"	R	320-22-3310	
with hoo	k °	+3	.134	L	320-22-4310	
Cuspid	0°	+3°	.134"	R	320-22-3350	
3 3 3 3 With hoo	k U	+5 .134	0 +3 1.134	.134	L	320-22-4350
1st Bicusp		0°	.125"	R	320-22-3410	
4 4 with hoo	k -12	U	U	.123	L	320-22-4410
4 4 1st Bicusp	oid -12°	0°	.125"	U	320-22-6400	
5 5 2nd Bicusp	oid -17°	0°	.125"	R	320-22-3510	
5 5 with hoo	k -17	U	.123	Ĺ	320-22-4510	
5 5 2 nd Bicus	pid -17°	0°	.125"	U	320-22-6500	

Marquis™

Standard Edgewise

Marquis Standard	Edgewise Sing	le Patient K	(it Description
Upper/Lower 5x5			

1 kit (20 brackets)	.\$41 . 80	kit
10 kits (200 brackets)	.\$35.80	kit

Marquis Standard Edgewise Individual Brackets


1 pack (10 brackets)	
20 packs (200 brackets)	
10 ner nack	

Patient Kits Description

Upper/Lower 5x5

Item Nu	ımbers
.018" Slot	.022" Slot
M18-2000	M22-2000

Standard Edgewise	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Centrals	0°	0°	.155"	U	M18-400	M22-400
Laterals	0°	0°	.125"	U	M18-410	M22-410
Cuspids	0°	0°	.130"	U	M18-420	M22-420
Bicuspids	0°	0°	.125"	U	M18-430	M22-430
Lower Anterior	0°	0°	.100"	U	M18-440	M22-440



Bionic[™]

Stainless Steel Bracket System

Where Form Meets Function for an Intelligent

Bracket Design™





Micro-Cast Bracket Body


Permanent Vertical Scribe Line

Ample Under Tie-Wing Area

80 Gauge Foil Mesh **Bonding Base**

Torque-In-Base Design





Bionic™

Ortho Technology's version of the Roth Prescription* "V" Slot (.025" x .025")

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.154"	R	B18-501	B22-501
TIT Central	T12	+3	.134	L	B18-502	B22-502
2 Lateral	+8°	+9°	.125"	R	B18-511	B22-511
22 Lateral	Τ0	T 9	.123	L	B18-512	B22-512
2/2 Cuspid	0°	+11°	.138"	R	B18-521	B22-521
$\frac{3 3}{}$ with hook	U	+11	.136	L	B18-522	B22-522
3 3 6	0°	+11°	.138"	R	B18-523	B22-523
33 Cuspid	U	+11	.136	L	B18-524	B22-524
4/14/ Bicuspid	-7°	0°	.130"	R	B18-531	B22-531
with hook	-/	0	.130	L	B18-532	B22-532
$\frac{4/5}{4/5}$ Bicuspid	-7°	0°	.130"	U	B18-533	B22-533
4/5 4/5 Bicuspid with offset	-7°	0°	.130"	U	B18-534	B22-534
4/4/ Bicuspid	-7°	0°	.130"	R	B18-535	B22-535
$\frac{ -5 ^{5}}{ }$ w/hook & offset	-/	U	.130	L	B18-536	B22-536

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

© 2015 Ortho Technology, Inc. Bionic is a trademark of Ortho Technology. OmniArch is a registered trademark of DENTSPLY GAC.

Patient Kits Description

Upper/Lower 5x5 No Hooks Upper/Lower 5x5 Hooks on Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers

.018" Slot .022" Slot BV18-1000 BV22-1000 BV22-1300 BV18-1300 BV18-1345 BV22-1345

V Mar	ndibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1/2 1/2	Lower Anterior	0°	0°	.098"	U	B18-540	B22-540
3 3	Cuspid	-11°	+5°	.138"	R	B18-551	B22-551
2 2	with hook	-11	۲	.130	L	B18-552	B22-552
3 3	Cuspid	-11°	+5°	.138"	R	B18-553	B22-553
J J	Сизрій	-11	+3	.136	L	B18-554	B22-554
4 4	1st Bicuspid	-17°	0°	.130"	R	B18-561	B22-561
דןד	with hook	-17	U	.130	L	B18-562	B22-562
4 4	1st Bicuspid	-17°	0°	.130"	U	B18-563	B22-563
4 4	1st Bicuspid with offset	-17°	0°	.130"	U	B18-565	B22-565
4 4	1st Bicuspid	-17°	0°	120"	R	B18-566	B22-566
4 4	w/hook & offset	-17	0	.130"	L	B18-567	B22-567
5 5	2 nd Bicuspid	-22°	0°	.130"	R	B18-571	B22-571
2 2	with hook	-22	U	.130	L	B18-572	B22-572
5 5	2 nd Bicuspid	-22°	0°	.130"	U	B18-573	B22-573
5 5	2 nd Bicuspid with offset	-22°	0°	.130"	U	B18-575	B22-575
5 5	2 nd Bicuspid	-22°	0°	.130"	R	B18-576	B22-576
2 2	w/hook & offset	-22	J	.130	L	B18-577	B22-577

BRACKETSY



Single Patient Kit Description Upper/Lower 5x5

opper/ zone: 5x5		
1 kit (20 brackets)	\$67.00	kit
10 kits (200 brackets)		
50 kits (1,000 brackets)	\$55.00	kit

Individual Brackets

1 pack (10 brackets)	\$33.50/pk (\$3.35 ea.)
20 packs (200 brackets)	
100 packs (1,000 brackets)	
10 per pack	

Bionic[™] **Bio-Progressive**

Ricketts/Hilgers Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+22°	+5°	.154"	R	B18-201	B22-201
TIT Central	+22	+3	.134	L	B18-202	B22-202
2 Lateral	+14°	+9°	.125"	R	B18-211	B22-211
22 Lateral	+14	TE	.123	L	B18-212	B22-212
212 Cuspid	+7°	+5°	.138"	R	B18-221	B22-221
$\frac{3}{3}$ with hook	+/	+3	.136	L	B18-222	B22-222
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.130"	U	B18-230	B22-230
Bicuspid	-7°	0°	.130"	R	B18-231	B22-231
$\frac{75}{5}$ with hook	-/	U	.130	L	B18-232	B22-232
Bicuspid by Bicuspid w/hook & offset	-7°	0°	.130"	R	B18-233	B22-233
$\frac{\sqrt{5}}{\sqrt{5}}$ w/hook & offset	-/	U	.130	L	B18-234	B22-234
$\frac{4/5}{5}$ Bicuspid with offset	-7°	0°	.130"	U	B18-237	B22-237
$\frac{4/5}{5}$ Bicuspid	0°	0°	.130"	U	B18-236	B22-236
4 Bicuspid	0°	0°	.130"	R	B18-238	B22-238
$\frac{\sqrt{5}\sqrt{5}}{\sqrt{5}}$ with hook		0	.130	L	B18-239	B22-239

 $^{{}^*\}mathit{The Ortho Technology version of the Ricketts/Hilgers\ prescription\ is\ not\ implied\ to\ be\ an\ exact}$ version of any other system, nor do we claim any endorsement of their respective doctors.

Item Numbers
.018" Slot .022" Slot
B418-1300 B422-1300 **Patient Kits Ricketts Configuration** Upper/Lower 5x5 Hooks on Cuspids B418-1300 B422-10° torque and 0° angulation on both upper and lower bicuspids

Hilgers Configuration	.018" Slot	.022" Slot
Upper/Lower 5x5 Hooks on Cuspids	B218-1300	B222-1300
Upper/Lower 5x5 Hooks 3, 4, & 5	B218-1345	B222-1345

Hilgers Offsets ConfigurationUpper/Lower 5x5 Hooks on Cuspids
Upper/Lower 5x5 Hooks 3, 4, & 5 **.022" Slot** B322-1300 .018" Slot B318-1300 B318-1345 B322-1345

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
Lower Anterior	0°	0°	.098"	U	B18-240	B22-240
3 3 Cuspid with hook	+7°	+5°	.138"	R L	B18-251 B18-252	B22-251 B22-252
4 4 1st Bicuspid	-12°	0°	.130"	U	B18-260	B22-260
4 4 1st Bicuspid with hook	-12°	0°	.130"	R L	B18-261 B18-262	B22-261 B22-262
4 4 1st Bicuspid with offset	-12°	0°	.130"	U	B18-267	B22-267
4 4 1st Bicuspid w/hook & offset	-12°	0°	.130"	R L	B18-263 B18-264	B22-263 B22-264
5 5 2 nd Bicuspid	-17°	0°	.130"	U	B18-270	B22-270
5 5 2nd Bicuspid with offset	-17°	0°	.130"	U	B18-277	B22-277
5 5 2nd Bicuspid with hook	-17°	0°	.130"	R L	B18-271 B18-272	B22-271 B22-272
5 5 2 nd Bicuspid w/hook & offset	-17°	0°	.130"	R L	B18-273 B18-274	B22-273 B22-274
4/5 4/5 Bicuspid	0°	0°	.130"	U	B18-278	B22-278
Bicuspid with hook	0°	0°	.130"	R L	B18-275 B18-276	B22-275 B22-276





Marquis™ TruGold™

24K Gold Plated Bracket System

Your Complete Source for Adding the Golden Touch™



- Real 24K gold plating makes brackets significantly brighter than competing bracket systems
- Resistant to wear, peeling, and chipping
- One-piece bracket system



TRY OUR...

TruGold 24K Gold Plated Archwire

See page 85



Marquis™ TruGold™

Ortho Technology's version of the Roth Prescription*

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.155"	R	GM18-001	GM22-001
1 1 Celitial	712	Τ.)	.133	L	GM18-002	GM22-002
2 2 Lateral	+8°	+9°	.125"	R	GM18-011	GM22-011
Z Z Lateral	+6	Τ3	.123	L	GM18-012	GM22-012
3 3 Cuspid	0°	+11°	.130"	R	GM18-021	GM22-021
with hook	"	T11	.130	L	GM18-022	GM22-022
4/5 4/5 Bicuspid	-7°	0°	.125"	R	GM18-031	GM22-031
$\frac{75 75}{}$ with hook	-/	0	.123	L	GM18-032	GM22-032
$\frac{4/5}{5}$ Bicuspid	-7°	0°	.125"	U	GM18-030	GM22-030

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Single Patient Kit Description Upper/Lower 5x5

1 kit (20 brackets)	.\$139.80 kit
10 kits (200 brackets)	.\$123.80 kit

Individual Brackets

See pages 42-43 for complete details regarding the Marquis Bracket System



Patient Kits Description

Upper/Lower 5x5 Hooks Cuspids Upper/Lower 5x5 Hooks 3, 4, & 5

Item Numbers .018" Slot .022" Slot

GM18-1300 GM22-1300 GM18-1345 GM22-1345

▼ Ma	ındibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1/2 1/2	Lower Anterior	0°	0°	0° .100"		GM18-040	GM22-040
3 3	Cuspid	-11°	+5°	.125"	R	GM18-051	GM22-051
	with hook		13	23	L	GM18-052	GM22-052
4 4	1st Bicuspid	-17°	U°	0° .125" -	R	GM18-061	GM22-061
.,,	with hook	' '			L	GM18-062	GM22-062
4 4	1st Bicuspid	-17°	0°		U	GM18-060	GM22-060
5 5	2 nd Bicuspid with hook	-22°	0°	.125"	R	GM18-071	GM22-071
	with hook	-22	U	.123	L	GM18-072	GM22-072
5 5	2 nd Bicuspid	-22°	0°	.125"	U	GM18-070	GM22-070

© 2015 Ortho Technology, Inc. TruGold and Marquis are trademarks of Ortho Technology.

MADE USA

ORTHO TECHNOLOGY

Presentation Typodonts

Brackets bonded 7 to 7 with archwires, buccal tubes on upper and lower first and second molars, and clear ligatures in place. Solid tissue colored gingiva base or clear flexibase with tooth tone colored teeth molded separately and inserted in place prior to bonding. The high quality detail makes these hinged typodonts the perfect presentation aid. All typodonts arrive in a hard plastic case with foam insert perfect for safe, protective storage.



Pinnacle™ Item #: G50-360



Votion® Item #: G50-122



PURE® Item #: G50-105



Spider Screw® Item #: G50-801



Sensation™* Item #: G50-109



Marquis™ TruGold 24K Gold Plated Item #: G50-142



Flexibase Ideal Occlusion Item #: G50-154



TruEase™ Mini Item #: G50-000

Avalon® Metal-Lined Composite 1 typodont	G50-203 \$290.00 ea.
Bionic™ Stainless Steel 1 typodont	
Encore!™ Metal-Lined Ceramic 1 typodont	G50-170 \$429.00 ea.
Eruption Appliance 24k Gold Eyelet and Traction Ch Stainless Steel Brackets	G50-205EA
EZ Space Maintainer®l typodont	G50-098 \$270.00 ea.
EZ Space Maintainer® Mini 1 typodont	\$140.00 ea.
EZ Slider™ Mini	
Flexibase Ideal Occlusion with Rooted Tee 1 typodont	

*N/a+	available	in	Canada	
"IVOT	available	ın	Canaaa	

Lotus Plus® Passive Self-Ligating Ideal Occlusion	G50-622
Lotus Plus® Passive Self-Ligating Malocclusion	G50-623
Marquis™ Stainless Steel 1 typodont	G50-140
Marquis™ TruGold™ 1 typodont	G50-142
OrthoFlex™ Composite 1 typodont	
NEW! Pinnacle™ Stainless Steel	G50-360 . \$225.00 ea.
Plain Typodonts Bondable tooth surface Ideal Occlusion	G50-120 G50-152 G50-151 G50-150

PURE®	
Sapphire	
1 typodont	\$395.00 ea.
Reflections™	
Ceramic	
1 typodont	\$310.00 ea.
Sensation™*	
Active Ceramic Self-Ligating	
1 typodont	\$435.00 ea.
Sensation™ M*	
Active Stainless Steel Self-Ligat	ingG50-119
1 typodont	\$425.00 ea.
Spider Screw® System	G50-801
1 typodont	
Spider Screw® Mini	G50-099
1 typodont	
3 or more	\$130.00 ea.
TruEase™ Mini	G50-000
1 typodont	\$135.00 ea.
3 or more	\$125.00 ea.
Votion®	
Stainless Steel	
1 typodont	\$225.00 ea.

BUCCAL TUBES & MOLAR BANDS

BUCCAL TUBES 6 MOLAR BAND

TRUFIT™2.0

TRUEASE

/OTION[®]

i	TruFit™ 2.0	P. 4	49-61
	Bondable Buccal Tubes	_ P.	49-51
	Weldable Buccal Tubes	- P.	52-53
	Molar Bands		P. 54
	Prewelded Band & Tube Assemblies	_ P.	55-58
	Prewelded Band & Tube Assembly Starter Kits		_ P. 59
	Molar Band Conversion Chart	_ P.	6 0-61
	Prewelded & Bondable Lingual Attachmen	ts_	P. 62
	Bicuspid Bands		P. 63

Weldable Bicuspid Brackets

TruEase™ Bondable Buccal Tubes

Votion® Bondable Buccal Tubes

TruCast™ Bondable Buccal Tubes

P. 63

P. 64

P. 65

P. 66



TruFit 2.0

Bondable Maxillary Buccal Tubes

- Available in the most popular prescriptions
- The first molar gingival ball hooks hold elastics and springs securely
- Designed with the latest technology for patient comfort



Maxillary Single Convertible


	Slot Distal Headgear Slot				Bondable		
	Torque	Offset	Tube	Size	Left	Right	
R	-10°	+10°	-	.018"	492-11-2610	492-11-1610	
Т	-10°	+10°	-	.022"	492-12-2610	492-12-1610	
R	-14°	+14°	-	.018"	492-91-2610	492-91-1610	
Г	-14°	+14°	-	.022"	492-92-2610	492-92-1610	
N	-14°	+10°	-	.018"	492-21-2610	492-21-1610	
IV	-14°	+10°	-	.022"	492-22-2610	492-22-1610	
		1	pack (10	tubes)	\$47.00/pk (\$4.70/tube)	
		5+ p	acks (50	tubes)	\$42.50/pk (\$4.25/tube)	



Maxillary Combination Convertible

Gingival Headgear Tube

	Slot	Distal	Headgear	Slot	Bondable		
	Torque	Offset	Tube	Size	Left	Right	
R	-10°	+10°	.045"	.018"	492-11-2671	492-11-1671	
ì	-10°	+10°	.045"	.022"	492-12-2671	492-12-1671	
R	-14°	+14°	.045"	.018"	492-91-2671	492-91-1671	
	-14°	+14°	.045"	.022"	492-92-2671	492-92-1671	
	-14°	+10°	.045"	.018"	492-21-2671	492-21-1671	
V	-14°	+10°	.045"	.022"	492-22-2671	492-22-1671	
		1	pack (10	tubes)	\$50.50/pk (\$5.05/tube)	
		5+ p	acks (50	tubes)	\$46.00/pk (\$4.60/tube)	



M The Ortho Technology version of the

Standard Edgewise Rx

McLaughlin, Bennett, & Trevisi Rx*



Improved Low-Profile Design



.018" x .025" Aı	ıxiliary Tube
------------------	---------------

	Slot Torque	Distal Offset	Headgea Tube	Slot Size	Bond Left	dable Right
R	-10°	+10°	_	.018"	492-11-2631	492-11-1631
	-10°	+10°	_	.022"	492-12-2631	492-12-1631
R	-14°	+14°	_	.018"	492-91-2631	492-91-1631
	-14°	+14°	_	.022"	492-92-2631	492-92-1631
		1 p	ack (10	tubes)	\$50.50/pk (\$5.05/tube)
		5+ pa	acks (50	tubes)	\$46.00/pk (\$4.60/tube)



Maxillary Combination Convertible

Occlusal Headgear Tube

	Slot	Distal	Headgear	Slot	Bondable		
	Torque	Offset	Tube	Size	Left	Right	
R	-10°	+10°	.045"	.018"	492-11-2672	492-11-1672	
	-10°	+10°	.045"	.022"	492-12-2672	492-12-1672	
R	-14°	+14°	.045"	.018"	492-91-2672	492-91-1672	
	-14°	+14°	.045"	.022"	492-92-2672	492-92-1672	
	-14°	+10°	.045"	.018"	492-21-2672	492-21-1672	
N	-14°	+10°	.045"	.022"	492-22-2672	492-22-1672	
S	0°	0°	.045"	.018"	492-51-2672	492-51-1672	
2	0°	0°	.045"	.022"	492-52-2672	492-52-1672	
	1 pack (10 tubes)				\$50.50/pk (\$5.05/tube)		
		5+ p	acks (50	tubes)	\$46.00/pk (\$4.60/tube)	

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

© 2015 Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. Due to our ongoing efforts to update and improve our buccal tube designs, the product you receive may vary slightly from those depicted. †Bases are non-etched.





Maxillary Triple Convertible

Occlusal Headgear Tube
.018" x .025" Auxiliary Tube

	Slot Torque	Distal Offset	Headgea Tube	r Slot Size	Bond Left	dable Right	
l	-10°	+10°	.045"	.018"	492-11-2652	492-11-1652	
١	-10°	+10°	.045"	.022"	492-12-2652	492-12-1652	7
R	-14°	+14°	.045"	.018"	492-91-2652	492-91-1652	
	-14°	+14°	.045"	.022"	492-92-2652	492-92-1652	
	-14°	+10°	.045"	.018"	492-21-2652	492-21-1652	
V	-14°	+10°	.045"	.022"	492-22-2652	492-22-1652	
Ī		1 p	oack (10	tubes)	\$53.50/pk (\$5.35/tube)		
		5+ pa	acks (50	tubes)	\$48.00/pk (\$4.80/tube)	



Maxillary Combination Non-Convertible

Occlusal Headgear Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Bond Left	dable Right	
	-10°	+8°	.045"	.018"	492-11-2682	492-11-1682	
R	-10	+0	.045	.016	492-11-2002		
٦	-10°	+8°	.045"	.022"	492-12-2682	492-12-1682	
•	0°	0°	.045"	.018"	492-51-2682	492-51-1682	
-	0°	0°	.045"	.022"	492-52-2682	492-52-1682	
		1 p	oack (10	tubes)	\$46.00/pk (\$4.60/tube)		
		5+ pa	acks (50	tubes)	\$42.50/pk (\$4.25/tube)	



Maxillary 1st or 2nd Molar Single Non-Convertible

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Bond Left	dable Right
F	-10°	+8°	-	.018"	-	-
	-10°	+8°	_	.022"	_	-
F	-10°	+10°	-	.018"	492-11-2020	492-11-1020
i	-10°	+10°	-	.022"	492-12-2020	492-12-1020
F	-14°	+14°	-	.018"	492-91-2020	492-91-1020
i	-14°	+14°	_	.022"	492-92-2020	492-92-1020
N.	-14°	+10°	-	.018"	492-21-2020	492-21-1020
Λ	-14°	+10°	-	.022"	492-22-2020	492-22-1020
•	0°	0°	-	.018"	492-51-2020	492-51-1020
2	0°	0°	-	.022"	492-52-2020	492-52-1020
		1 p	oack (10	tubes)	\$44.00/pk (\$4.40/tube)
		5+ pa	acks (50	tubes)	\$40.50/pk (\$4.05/tube)

Bondable						
Left	Right	l				
492-11-2023	492-11-1023	L				
492-12-2023	492-12-1023	Ī				
-	-	l				
-	-					
-	_					
-	-					
492-21-2023	492-21-1023					
492-22-2023	492-22-1023	Ī				
492-51-2023	492-51-1023	Ļ				
492-52-2023	492-52-1023	3				
\$44.00/pk (\$4.40/tube)						
\$40.50/pk (\$	4.05/tube)					

- R The Ortho Technology version of the Roth Rx*
- M The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*
- Standard Edgewise Rx



^{*}The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

TruFit 2.0

Bondable Mandibular Buccal Tubes



Mandibular 1st or **2nd Molar Single Non-Convertible**

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Bond Left	dable Right
F	-25°	+5°	-	.018"	-	-
_	-25°	+5°	_	.022"	-	-
F	-25°	+4°	-	.018"	492-91-4020	492-91-3020
_	-25°	+4°	-	.022"	492-92-4020	492-92-3020
	-10°	0°	-	.018"	492-21-4020	492-21-3020
Λ	-10°	0°	_	.022"	492-22-4020	492-22-3020
-	0°	0°	_	.018"	492-51-4020	492-51-3020
2	0°	0°	_	.022"	492-52-4020	492-52-3020
		1	pack (10	tubes)	\$44.00/pk (\$4.40/tube)	
		5+ p	acks (50	tubes)	\$40.50/pk (\$4.05/tube)

Bondable						
Left	Right					
492-11-4023	492-11-3023					
492-12-4023	492-12-3023					
-	-					
-	-					
492-21-4023	492-21-3023					
492-22-4023	492-22-3023					
492-51-4023	492-51-3023					
492-52-4023	492-52-3023					
\$44.00/pk (\$4.40/tube)						
\$40.50/pk (\$4.05/tube)					

TruFit 2.0 Bases







Mandibular **Lip Bumper Convertible Gingival Tube**

	Slot	Distal	Headgear	Slot	Bond	dable
	Torque	Offset	Tube	Size	Left	Right
R	-25°	+4°	.045"	.018"	492-91-4691	492-91-3691
ı	-25°	+4°	.045"	.022"	492-92-4691	492-92-3691
		1 p	oack (10	tubes)	\$50.50/pk (\$5.05/tube)
		5+ pa	acks (50	tubes)	\$46.00/pk (\$4.60/tube)



Mandibular Double Convertible

.018" x .025" Auxiliary Tube

	Slot	Distal	Headgea	r Slot	Bond	dable	
	Torque	Offset	Tube	Size	Left	Right	
ľ	-25°	+4°	-	.018"	492-91-4631	492-91-3631	
F	-25°	+4°	-	.022"	492-92-4631	492-92-3631	7
N	-20°	0°	-	.018"	492-21-4631	492-21-3631	
1	-20°	0°	-	.022"	492-22-4631	492-22-3631	
		1 p	oack (10	tubes)	\$50.50/pk (\$5.05/tube)		
		5+ pa	acks (50	tubes)	\$46.00/pk (\$4.60/tube)	
							•



	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Bond Left	dable Right
R	-25°	+4°	-	.018"	492-91-4610	492-91-3610
٦	-25°	+4°	_	.022"	492-92-4610	492-92-3610
	-20°	0°	-	.018"	492-21-4610	492-21-3610
Ν	-20°	0°	-	.022"	492-22-4610	492-22-3610
		1 p	ack (10	tubes)	\$47.00/pk (\$4.70/tube)
		5+ pa	acks (50	tubes)	\$42.50/pk (\$4.25/tube)

© 2015 Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. Due to our ongoing efforts to update and improve our buccal tube designs, the product you receive may vary slightly from those depicted.







51



Weldable Buccal Tubes



Maxillary Single Convertible

	Slot Torque	Distal Offset	Headgeaı Tube	Slot Size	Welc Left	lable Right	
	-10°	+10°	-	.018"	490-11-2610	490-11-1610	
R	-10°	+10°	-	.022"	490-12-2610	490-12-1610	
R	-14°	+14°	_	.018"	490-91-2610	490-91-1610	
	-14°	+14°	-	.022"	490-92-2610	490-92-1610	
N	-14°	+10°	-	.018"	490-21-2610	490-21-1610	
IV	-14°	+10°	-	.022"	490-22-2610	490-22-1610	
		1 p	ack (10	tubes)	\$36.50/pk (\$3.65/tube)		
		5+ pa	acks (50	tubes)	\$34.00/pk (\$3.40/tube)	



Maxillary Combination Convertible Gingival Headgear Tube

	Slot	Distal	Headgea	Slot	Weld	lable	
	Torque	Offset	Tube	Size	Left	Right	
R	-10°	+10°	.045"	.018"	490-11-2671	490-11-1671	
	-10°	+10°	.045"	.022"	490-12-2671	490-12-1671	
R	-14°	+14°	.045"	.018"	490-91-2671	490-91-1671	
i	-14°	+14°	.045"	.022"	490-92-2671	490-92-1671	
	-14°	+10°	.045"	.018"	490-21-2671	490-21-1671	
IV	-14°	+10°	.045"	.022"	490-22-2671	490-22-1671	
		1 p	oack (10	tubes)	\$40.00/pk (\$4.00/tube)		
		5+ pa	acks (50	tubes)	\$36.50/pk (\$3.65/tube)	



Maxillary Double Convertible

.018" x .025" Auxiliary Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Welc Left	lable Right	
F	-10°	+10°	-	.018"	490-11-2631	490-11-1631	
	-10°	+10°	-	.022"	490-12-2631	490-12-1631	
F	-14°	+14°	-	.018"	490-91-2631	490-91-1631	
П	-14°	+14°	-	.022"	490-92-2631	490-92-1631	
		1 p	oack (10	tubes)	\$40.00/pk (\$4.00/tube) \$36.50/pk (\$3.65/tube)		
		5+ pa	acks (50	tubes)			

Prewelded free of charge on TruFit 2.0 Molar Bands See page 54

- Available loose or prewelded
- 3.5 mm mesial/distal measurement



Maxillary 1st or 2nd Molar Single Non-Convertible

	Slot	Distal	Headgear	Slot	Weldable		
	Torque	Offset	Tube	Size	Left	Right	
R	-10°	+10°	_	.018"	490-11-2020	490-11-1020	
	-10°	+10°	-	.022"	490-12-2020	490-12-1020	
R	-14°	+14°	-	.018"	490-91-2020	490-91-1020	
	-14°	+14°	-	.022"	490-92-2020	490-92-1020	
	-14°	+10°	-	.018"	490-21-2020	490-21-1020	
N	-14°	+10°	_	.022"	490-22-2020	490-22-1020	
c	0°	0°	-	.018"	490-51-2020	490-51-1020	
S	0°	0°	-	.022"	490-52-2020	490-52-1020	
		1 p	ack (10	tubes)	\$35.50/pk (\$3.55/tube)		
		5+ pa	acks (50	tubes)	\$32.50/pk (\$3.25/tube)		



Maxillary Combination Convertible Occlusal Headgear Tube

\$36.50/pk (\$3.65/tube)

Weldable Distal Headgear **Slot** Right **Torque** +10° .045" .018" 490-11-2672 490-11-1672 -10° .045" +10° .022' 490-12-2672 490-12-1672 -14° +14° .045" .018 490-91-2672 490-91-1672 R $+14^{\circ}$.045" .022' 490-92-2672 490-92-1672 .045" +10° .018" 490-21-2672 490-21-1672 -14° +10° .045" .022' 490-22-2672 490-22-1672 .045" .018" 490-51-2672 490-51-1672 S .045" .022" 490-52-2672 490-52-1672 1 pack (10 tubes) \$40.00/pk (\$4.00/tube)

5+ packs (50 tubes)

© 2015 Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. Due to our ongoing efforts to update and improve our buccal tube designs, the product you receive may vary slightly from those depicted.

TruFit 20

Weldable Buccal Tubes



Maxillary Triple Convertible Occlusal Headgear Tube

.018" x .025" Auxiliary Tube


	Slot	Distal Offset	Headgear Tube	Slot Size	Welc Left	lable	
	Torque	Uliset	Tube	Size	Leit	Right	
R	-10°	+10°	.045"	.018"	490-11-2652	490-11-1652	
	-10°	+10°	.045"	.022"	490-12-2652	490-12-1652	
R	-14°	+14°	.045"	.018"	490-91-2652	490-91-1652	
	-14°	+14°	.045"	.022"	490-92-2652	490-92-1652	
	-14°	+10° .045		.018"	490-21-2652	490-21-1652	
N	-14°	+10°	.045"	.022"	490-22-2652	490-22-1652	
		1 p	oack (10	tubes)	\$45.00/pk (\$4.50/tube)		
		5+ pa	acks (50	tubes)	\$41.00/pk (\$4.10/tube)		



Maxillary Combination Non-Convertible

Occlusal Headgear Tube

	Slot	Distal Offset	Headgear Tube	Slot Size	Weldable Left Right			
	Torque	Oliset	Tube	Size	Leit	Right		
F	-10°	+8°	.045"	.018"	490-11-2682	490-11-1682		
	-10°	+8°	.045"	.022"	490-12-2682	490-12-1682		
Ν	-14°	+10°	.045"	.018"	490-21-2682	490-21-1682		
_	-14°	+10°	.045"	.022"	490-22-2682	490-22-1682		
9	0°	0°	.045"	.018"	490-51-2682	490-51-1682		
	0°	0°	.045"	.022"	490-52-2682	490-52-1682		
	0°	+8°	.045"	.018"	490-41-2682	490-41-1682		
2	0°	+8°	.045"	.022"	490-42-2682	490-42-1682		
		1 p	oack (10	tubes)	\$36.50/pk (\$3.65/tube)			
		5+ pa	acks (50	tubes)	\$34.00/pk (\$3.40/tube)			



Standard Edgewise Rx

McLaughlin, Bennett, & Trevisi Rx*

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.



Mandibular Lip Bumper Convertible

	Slot	Distal	Headgea	r Slot	Weldable				
	Torque	Offset	Tube	Size	Left	Right			
R	-25°	+4°	.045"	.018"	490-91-4691	490-91-3691			
	-25°	+4°	.045"	.022"	490-92-4691	490-92-3691			
		1 p	ack (10	tubes)	\$40.00/pk (\$4.00/tube)				
		5+ pa	acks (50	tubes)	\$36.50/pk (\$3.65/tube)				



Mandibular Single Convertible

	Slot Torque	Distal Offset	Headgea Tube	r Slot Size	Welc Left	lable Right	
	-25°	+4°	-	.018"	490-91-4610	490-91-3610	
R	-25°	+4°	-	.022"	490-92-4610	490-92-3610	
N	-20°	0°	-	.018"	490-21-4610	490-21-3610	
IV	-20°	0°	-	.022"	490-22-4610	490-22-3610	
		1 p	ack (10	tubes)	\$36.50/pk (\$3.65/tube)		
		5+ pa	icks (50	tubes)	\$34.00/pk (\$3.40/tube)		



Mandibular 1st or 2nd Molar Single Non-Convertible

	Slot	Distal	Headgear	Slot	Welc	lable	
	Torque	Offset	Tube	Size	Left	Right	
R	-25°	+4°	-	.018"	490-91-4020	490-91-3020	
	-25°	+4°	_	.022"	490-92-4020	490-92-3020	
N	-10°	0°	_	.018"	490-21-4020	490-21-3020	
	-10°	0°	_	.022"	490-22-4020	490-22-3020	
	0°	0°	_	.018"	490-51-4020	490-51-3020	
2	0°	0°	_	.022"	490-52-4020	490-52-3020	
		1	pack (10	tubes)	\$35.50/pk (\$3.55/tube)		
		5+ p	acks (50	tubes)	\$32.50/pk (\$3.25/tube)		



Mandibular Double Convertible

.018" x .025" Auxiliary Tube

	Slot	Dietal	Headgea	r Slot	Weldable		
	Torque	Offset	Tube	Size	Left	Right	
R	-25°	+4°	-	.018"	490-91-4631	490-91-3631	
	-25°	+4°	-	.022"	490-92-4631	490-92-3631	
N	-20°	0°	-	.018"	490-21-4631	490-21-3631	
I N	-20°	0°	-	.022"	490-22-4631	490-22-3631	
		1	pack (10	tubes)	\$40.00/pk (\$4.00/tube)		
		5+ p	acks (50	tubes)	\$36.50/pk (\$3.65/tube)		



Trufit 20 Molar Bands Available Loose or Prewelded with Our Free Provedding County

with Our Free Prewelding Service



Plain Molar Bands

Quadrant	1st Molars	2nd Molars
UR	495-00-1600-000	495-00-1700-000
UL	495-00-2600-000	495-00-2700-000
LR	495-00-3600-000	495-00-3700-000
LL	495-00-4600-000	495-00-4700-000
	1 band	\$2.10 ea.
	100 or more	\$1.70 ea.

TruFit 2.0 Molar Band Features

Shaped for a Precise, Tight Fit on Every Tooth

Each band has precision shaped tooth anatomy for a tight secure fit with excellent retention.

Medium Temper Malleable Stainless Steel Finished to a Brilliant Mirror Luster

TruFit 2.0 stainless steel bands are strong enough to avoid wall collapses while being malleable enough for snap fittings that require little to no burnishing.

Laser Etched Exterior Identification **System with Micro-Etched Interior**

TruFit 2.0 bands display an exterior laser etched identification system for easy product selection. The interior is micro-etched to increase bond strength and retention.

Band Sizes Available From 29.5 to 44

See pages 60-61 for our band conversion chart.

Precision Preweld Position Buccal Attachments to Molar Bands




Entrance to buccal tube will bisect the mesial/buccal cusp



Upper Molar

Lower Molar

Lingual Cleat Attachments Centered on Molar Bands





Maxillary Prewelded Band and Tube Assemblies

Available With or Without Lingual Cleat

MADE USA



Maxillary
Single
Convertible

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Lin	ngual Cleat Right
	-10°	+10°	-	.018"	493-11-2610-000	493-11-1610-000	491-11-2610-000	491-11-1610-000
Н	-10°	+10°	-	.022"	493-12-2610-000	493-12-1610-000	491-12-2610-000	491-12-1610-000
F	-14°	+14°	-	.018"	493-91-2610-000	493-91-1610-000	491-91-2610-000	491-91-1610-000
	-14°	+14°	-	.022"	493-92-2610-000	493-92-1610-000	491-92-2610-000	491-92-1610-000
	-14°	+10°	-	.018"	493-21-2610-000	493-21-1610-000	491-21-2610-000	491-21-1610-000
N	-14°	+10°	-	.022"	493-22-2610-000	493-22-1610-000	491-22-2610-000	491-22-1610-000
		Individual Assemblies			\$ 6.1 (0 ea.	⁵ 6.10 ea.	
	50 or more				\$5.4	5 ea.	\$5.45 ea.	



Maxillary Combination Convertible Occlusal Headgear Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Lin Left	ngual Cleat Right	
R	-10°	+10°	.045"	.018"	493-11-2672-000	493-11-1672-000	491-11-2672-000	491-11-1672-000	
	-10°	+10°	.045"	.022"	493-12-2672-000	493-12-1672-000	491-12-2672-000	491-12-1672-000	
	-14°	+14°	.045"	.018"	493-91-2672-000	493-91-1672-000	491-91-2672-000	491-91-1672-000	
٦	-14°	+14°	.045"	.022"	493-92-2672-000	493-92-1672-000	491-92-2672-000	491-92-1672-000	
	-14°	+10°	.045"	.018"	493-21-2672-000	493-21-1672-000	491-21-2672-000	491-21-1672-000	
Λ	-14°	+10°	.045"	.022"	493-22-2672-000	493-22-1672-000	491-22-2672-000	491-22-1672-000	
9	0°	0°	.045"	.018"	493-51-2672-000	493-51-1672-000	491-51-2672-000	491-51-1672-000	
5	0°	0°	.045"	.022"	493-52-2672-000	493-52-1672-000	491-52-2672-000	491-52-1672-000	
	Individual Assemblies			mblies	\$6.7	0 ea.	§6.70 ea.		
	50 or more				\$6.0	5.00 ea. §6.00 ea.			



Maxillary Combination Convertible Gingival Headgear Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Li	ngual Cleat Right	
R	-10°	+10°	.045"	.018"	493-11-2671-000	493-11-1671-000	491-11-2671-000	491-11-1671-000	
7	-10°	+10°	.045"	.022"	493-12-2671-000	493-12-1671-000	491-12-2671-000	491-12-1671-000	
R	-14°	+14°	.045"	.018"	493-91-2671-000	493-91-1671-000	491-91-2671-000	491-91-1671-000	
1	-14°	+14°	.045"	.022"	493-92-2671-000	493-92-1671-000	491-92-2671-000	491-92-1671-000	
	-14°	+10°	.045"	.018"	493-21-2671-000	493-21-1671-000	491-21-2671-000	491-21-1671-000	
N	-14°	+10°	.045"	.022"	493-22-2671-000	493-22-1671-000	491-22-2671-000	491-22-1671-000	
	Individual Assemblies			mblies	\$6.7°	0 ea.	\$6.70 ea.		
	50 or more				\$6.0	\$6.00 ea. \$6.00 ea.			

R The Ortho Technology version of the Roth Rx*

The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*

S Standard Edgewise Rx

^{*}The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

^{© 2015} Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. Sizes start at 29.5 up to 44. See band size conversion chart on pages 60-61 for specific details. Additional assemblies available. Custom orders also offered. Allow 3-4 weeks for delivery. †Lingual Sheaths are \$2.00 additional cost. Prices subject to change without notice. All bands are supplied with standard weld position unless custom ordered. Custom orders are final and cannot be returned. Due to our ongoing efforts to update and improve our buccal tube designs, the product you receive may vary slightly from those depicted.





Our Most Popular Assemblies are In-Stock and Ready for Same Day Shipment!

Additional Assemblies Available

MADEUSA

Ask your Sales Representative about our Price-Lock Program



Maxillary Prewelded Band and Tube Assemblies

Available With or Without Lingual Cleat



Maxillary 1st Molar Single Non-Convertible

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Lin Left	ngual Cleat Right	
R	-10°	+10°	-	.018"	493-11-2620-000	493-11-1620-000	491-11-2620-000	491-11-1620-000	
H	-10°	+10°	_	.022"	493-12-2620-000	493-12-1620-000	491-12-2620-000	491-12-1620-000	
	-14°	+14°	-	.018"	493-91-2620-000	493-91-1620-000	491-91-2620-000	491-91-1620-000	
Ť	-14°	+14°	_	.022"	493-92-2620-000	493-92-1620-000	491-92-2620-000	491-92-1620-000	
	-14°	+10°	-	.018"	493-21-2620-000	493-21-1620-000	491-21-2620-000	491-21-1620-000	
N	-14°	+10°	-	.022"	493-22-2620-000	493-22-1620-000	491-22-2620-000	491-22-1620-000	
S	0°	0°	-	.018"	493-51-2620-000	493-51-1620-000	491-51-2620-000	491-51-1620-000	
Ī	0°	0°	-	.022"	493-52-2620-000	493-52-1620-000	491-52-2620-000	491-52-1620-000	
	Individual Assemblies				\$6.0	§6.00 ea.		\$6.00 ea.	
	50 or more				\$5.30	5.30 ea. \$5.30 ea.		0 ea.	



Maxillary 2nd Molar Single Non-Convertible

	Slot	Distal	Headgear	Slot	With Ling	jual Cleat	Without Lin	ngual Cleat
	Torque	Offset	Tube	Size	Left	Right	Left	Right
F	-10°	+10°	_	.018"	493-11-2720-000	493-11-1720-000	491-11-2720-000	491-11-1720-000
ľ	-10°	+10°	_	.022"	493-12-2720-000	493-12-1720-000	491-12-2720-000	491-12-1720-000
	-14°	+14°	_	.018"	493-91-2720-000	493-91-1720-000	491-91-2720-000	491-91-1720-000
1	-14°	+14°	_	.022"	493-92-2720-000	493-92-1720-000	491-92-2720-000	491-92-1720-000
	-14°	+10°	-	.018"	493-21-2720-000	493-21-1720-000	491-21-2720-000	491-21-1720-000
N	-14°	+10°	_	.022"	493-22-2720-000	493-22-1720-000	491-22-2720-000	491-22-1720-000
9	0°	0°	-	.018"	493-51-2720-000	493-51-1720-000	491-51-2720-000	491-51-1720-000
	0°	0°	-	.022"	493-52-2720-000	493-52-1720-000	491-52-2720-000	491-52-1720-000
	Individual Assemblies			mblies	\$6.00 ea.		\$6.00 ea.	
	50 or more			r more	\$5.30 ea.		§5.30 ea.	



Maxillary Double Convertible

.018" x .025" Auxiliary Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Li Left	ngual Cleat Right
H	-10°	+10°	-	.018"	493-11-2631-000	493-11-1631-000	491-11-2631-000	491-11-1631-000
ď	-10°	+10°	_	.022"	493-12-2631-000	493-12-1631-000	491-12-2631-000	491-12-1631-000
	-14°	+14°	_	.018"	493-91-2631-000	493-91-1631-000	491-91-2631-000	491-91-1631-000
	-14°	+14°	_	.022"	493-92-2631-000	493-92-1631-000	491-92-2631-000	491-92-1631-000
		Individ	ual Asse	mblies	\$6.7	0 ea.	\$6.7	0 ea.
	50 or more			r more	§6.00 ea.		\$6.0	0 ea.



Our Most Popular Assemblies are In-Stock and Ready for Same Day Shipment!

Additional Assemblies Available

MADE USA

Ask your Sales Representative about our **Price-Lock Program**



Maxillary Prewelded Band and Tube Assemblies Available With or Without Lingual Cleat



Maxillary Triple Convertible **Occlusal Headgear Tube** with Hook

.018" x .025" Auxiliary Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Lin	ngual Cleat Right
F	-10°	+10°	.045"	.018"	493-11-2652-000	493-11-1652-000	491-11-2652-000	491-11-1652-000
Н	-10°	+10°	.045"	.022"	493-12-2652-000	493-12-1652-000	491-12-2652-000	491-12-1652-000
F	-14°	+14°	.045"	.018"	493-91-2652-000	493-91-1652-000	491-91-2652-000	491-91-1652-000
Í	-14°	+14°	.045"	.022"	493-92-2652-000	493-92-1652-000	491-92-2652-000	491-92-1652-000
B	-14°	+10°	.045"	.018"	493-21-2652-000	493-21-1652-000	491-21-2652-000	491-21-1652-000
Λ	-14°	+10°	.045"	.022"	493-22-2652-000	493-22-1652-000	491-22-2652-000	491-22-1652-000
	Cleat Individual Assemblies			mblies	\$6.90 ea.		\$6.90 ea.	
	50 or more			r more	^{\$} 6.10 ea.		\$6.10	O ea.

	-10°	+10°	.045"	.018"	496-11-2652-000	496-11-1652-000	-	-
Ť	-10°	+10°	.045"	.022"	496-12-2652-000	496-12-1652-000	-	-
	With Lingเ	ıal Sheath Iı	ndividual A	ssemblies	\$9.9	O ea.		
			50 c	or more	\$9.00	O ea.		



Maxillary Combination Non-Convertible With Occlusal **Headgear Tube**

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	ual Cleat Right	Without Li	ngual Cleat Right
R	-10°	+8°	.045"	.018"	493-11-2682-000	493-11-1682-000	491-11-2682-000	491-11-1682-000
	-10°	+8°	.045"	.022"	493-12-2682-000	493-12-1682-000	491-12-2682-000	491-12-1682-000
N	-14°	+10°	.045"	.018"	493-21-2682-000	493-21-1682-000	491-21-2682-000	491-21-1682-000
Ĭ	-14°	+10°	.045"	.022"	493-22-2682-000	493-22-1682-000	491-22-2682-000	491-22-1682-000
5	0°	0°	.045"	.018"	493-51-2682-000	493-51-1682-000	491-51-2682-000	491-51-1682-000
	0°	0°	.045"	.022"	493-52-2682-000	493-52-1682-000	491-52-2682-000	491-52-1682-000
5	0°	+8°	.045"	.018"	493-41-2682-000	493-41-1682-000	491-41-2682-000	491-41-1682-000
	0°	+8°	.045"	.022"	493-42-2682-000	493-42-1682-000	491-42-2682-000	491-42-1682-000
	Individual Assemblies			mblies	§6.70 ea.		§6.70 ea.	
	50 or more			r more	⁵ 6.00 ea.		§6.00 ea.	



The Ortho Technology version of the Roth Rx*

M The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*

Standard Edgewise Rx

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

© 2015 Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. Sizes start at 29.5 up to 44. See band size conversion chart on pages 60-61 for specific details. Additional assemblies available. Custom orders also offered. Allow 3-4 weeks for delivery. †Lingual Sheaths are \$2.00 additional cost. Prices subject to change without notice. All bands are supplied with standard weld position unless custom ordered. Custom orders are final and cannot be returned. Due to our ongoing efforts to update and improve our buccal tube designs, the product you receive may vary slightly from those depicted.









Mandibular Prewelded Band and Tube Assemblies

Available With or Without Lingual Cleat





Mandibular 1st Molar Single Non-Convertible

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	gual Cleat Right	Without Lis Left	ngual Cleat Right
R	-25°	+4°	-	.018"	493-91-4620-000	493-91-3620-000	491-91-4620-000	491-91-3620-000
	-25°	+4°	_	.022"	493-92-4620-000	493-92-3620-000	491-92-4620-000	491-92-3620-000
IV	-10°	0°	_	.018"	493-21-4620-000	493-21-3620-000	491-21-4620-000	491-21-3620-000
IV	-10°	0°	_	.022"	493-22-4620-000	493-22-3620-000	491-22-4620-000	491-22-3620-000
G	0°	0°	_	.018"	493-51-4620-000	493-51-3620-000	491-51-4620-000	491-51-3620-000
	0°	0°	_	.022"	493-52-4620-000	493-52-3620-000	491-52-4620-000	491-52-3620-000
		Individ	dual Asse	mblies	\$6.0	0 ea.	\$6.0	0 ea.
			50 o	r more	\$5.3	0 ea.	\$5.3°	0 ea.



Mandibular 2nd Molar Single Non-Convertible

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Lingı Left	ual Cleat Right	Without Lin Left	gual Cleat Right
	-25°	+4°	_	.018"	493-91-4720-000	493-91-3720-000	491-91-4720-000	491-91-3720-000
٦	-25°	+4°	_	.022"	493-92-4720-000	493-92-3720-000	491-92-4720-000	491-92-3720-000
N	-10°	0°	-	.018"	493-21-4720-000	493-21-3720-000	491-21-4720-000	491-21-3720-000
	-10°	0°	_	.022"	493-22-4720-000	493-22-3720-000	491-22-4720-000	491-22-3720-000
	0 °	0°	-	.018"	493-51-4720-000	493-51-3720-000	491-51-4720-000	491-51-3720-000
	0°	0°	-	.022"	493-52-4720-000	493-52-3720-000	491-52-4720-000	491-52-3720-000
		Individ	dual Asse	mblies	\$6.0	0 ea.	\$6.0	0 ea.
	50 or more			r more	\$5.30 ea.		\$5.30 ea.	



Mandibular Single Convertible

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Lir Left	igual Cleat Right
F	-25°	+4°	-	.018"	493-91-4610-000	493-91-3610-000	491-91-4610-000	491-91-3610-000
٦	-25°	+4°	_	.022"	493-92-4610-000	493-92-3610-000	491-92-4610-000	491-92-3610-000
N	-20°	0°	-	.018"	493-21-4610-000	493-21-3610-000	491-21-4610-000	491-21-3610-000
I	-20°	0°	-	.022"	493-22-4610-000	493-22-3610-000	491-22-4610-000	491-22-3610-000
		Individ	lual Asse	mblies	\$6.10	0 ea.	\$6.1¢	O ea.
			50 o	r more	\$5.4	5 ea.	\$5.4	5 ea.



Mandibular Double Convertible .018" x .025" Auxiliary Tube

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	With Ling Left	jual Cleat Right	Without Li Left	ngual Cleat Right
R	-25°	+4°	-	.018"	493-91-4631-000	493-91-3631-000	491-91-4631-000	491-91-3631-000
1	-25°	+4°	_	.022"	493-92-4631-000	493-92-3631-000	491-92-4631-000	491-92-3631-000
N	-20°	0°	-	.018"	493-21-4631-000	493-21-3631-000	491-21-4631-000	491-21-3631-000
IV	-20°	0°	_	.022"	493-22-4631-000	493-22-3631-000	491-22-4631-000	491-22-3631-000
		Individ	lual Asse	mblies	\$6.7	0 ea.	\$6.7	0 ea.
Ī	50 or more			r more	\$6.0	§6.00 ea.		0 ea.



Mandibular Lip Bumper Convertible

	Slot	Distal	Lip	Slot	With Ling	jual Cleat	Without Li	ngual Cleat	
	Torque	Offset	Bumper	Size	Left	Right	Left	Right	
	-25°	+4°	.045"	.018"	493-91-4691-000	493-91-3691-000	491-91-4691-000	491-91-3691-000	
Ï	-25°	+4°	.045"	.022"	493-92-4691-000	493-92-3691-000	491-92-4691-000	491-92-3691-000	
		Individ	lual Asse	mblies	\$6.70	0 ea.	\$6.70	O ea.	
			50 o	r more	\$6.0	§6.00 ea.		⁵ 6.00 ea.	

MOLAR T

TruFit 2.0

Assembly Starter Kits

We've Pre-Selected the Most Popular Sizes with the Most Popular Buccal Tubes to **Get You Started!**



Maxillary
Single Convertible
with Lingual Cleats

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Starter Kit
F	-10°	+10°	-	.018"	493-11-1910-080
	-10°	+10°	-	.022"	493-12-1910-080
N	-14°	+10°	-	.018"	493-21-1910-080
IV	-14°	+10°	-	.022"	493-22-1910-080
		1 kit (8	\$475.00 kit		



Maxillary Combination Convertible with Lingual Cleats

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Starter Kit
-	-10°	+10°	.045"	.018"	493-11-1970-080
R	-10°	+10°	.045"	.022"	493-12-1970-080
N	-14°	+10°	.045"	.018"	493-21-1970-080
11	-14°	+10°	.045"	.022"	493-22-1970-080
		l kit (80	\$475.00 kit		



Mandibular Single Convertible with Lingual Cleats

	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Starter Kit
	-25°	+4°	-	.018"	493-11-1010-080
	-25°	+4°	-	.022"	493-12-1010-080
N	-20°	0°	-	.018"	493-21-1010-080
IV	-20°	0°	-	.022"	493-22-1010-080
		1 kit (8	0 assemb	lies)	\$475.00 kit

R The Ortho Technology version of the Roth Rx*

M The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Ask your Sales Representative for details







40 Left and 40 Right Assemblies

TruFit™ 2.0

Band and Tube Assembly Starter Kits
1 kit...........\$475.00 kit (\$5.94/assembly)

Includes: 40 left and 40 right band & tube assemblies

Includes: 40 left and 40 right band & tube assemblies with lingual cleats, 2x band tray organizers dimensions 11 $^{1}/_4$ " W x 1 $^{1}/_2$ " H x 8 $^{1}/_4$ " D

 2 Stackable Band Organizers included with your purchase

Starter Kits Included Sizes	34.5	32	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40
Left Side	2	2	2	4	4	5	5	5	4	3	2	2
Right Side	2	2	2	4	4	5	5	5	4	3	2	2

Bands Sizes are Available From 29.5 to 44

BAND SIZE CONVERSION TRAYS Included With Your Order





Ask Your Account Representative For Details

© 2015 Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. Due to our ongoing efforts to update and improve our buccal tube designs, the product you receive may vary slightly from those depicted.

MOLAR BANDS - -

Converting to Ortho Technology's TruFit™ 2.0 Molar Bands is Easy!

No Need to Learn New Band Sizes

When you convert to TruFit™ 2.0

Regardless of the brand being currently used in your practice when you convert to TruFit 2.0, we will label our band trays with your current sizing, eliminating the need to learn our sizes and leaving your workflow unchanged.

BAND SIZE CONVERSION TRAYS INCLUDED WITH YOUR ORDER

Ask Your Sales Representative for Details

Upper First and Second Molar Band Sizes

Ortho Technology			29.5		30	30.5			31		31.5		32		32.5		33		33.5		34		34.5		35		35.5		36			36.5	
American MR Full Sizes			1	2	3			5	6		7		8		9		10	11			12	13	14	15				17		18		19	
American MR 1/2 Sizes											1.5		2				2.5	3			3.5		4		4.5		5		5.5		6		
Dentaurum Dentaform		2	3	4		5	6		7		8		9		10		11	12			13		14	15	16		17		18		19	20	
GAC-SnapFit				1		2				3		4		5		6			7		8		9		10		11	12	13			14	
Ormco				1		2	3			4		5		6	7	8	9	10			11		12	13	14	15		16		17		18	19
RM0 Truform	1	2	3	4	5		6			7			8		9		10		11	12		13			14		15		16	17		18	
TP GripTite				00		0			1	2			3		5		6	7	8		9	10	11	12	13		14	15	16	17	18	19	
Unitek Victory Series			29.5		30	30.5			31		31.5		32		32.5		33		33.5		34		34.5		35		35.5		36			36.5	

Lower First Molar Band Sizes

Ortho Technology			29.5			30			30.5		31		31.5		32		32.5		33			33.5		34		34.5		35		35.5			36
American MR Full Sizes					1		2		3	4	5	6	7	8		9	10		11		12	13		14				16		17		18	
American MR 1/2 Sizes											1.5		2		2.5				3			3.5				4		4.5		5			
Dentaurum Dentaform			3	4		5		6	7	8	9	10	11	12		13	14			15	16			17		18		19		20	21		22
GAC-SnapFit				1					2			3		4		5		6				7		8		9	10		11	12		13	
Ormco			1									5		6		7		8	9			10	11	12	13	14	15	16		17		18	
RMO Truform	1	2	3		4	5				6		7			8		9		10			11	12		13		14		15	16		17	
TP GripTite						00			0		1	2		4	5	6	7		8		9	10		11		12		13	14	15	16	17	18
Unitek Victory Series			29.5			30			30.5		31		31.5		32		32.5		33			33.5		34		34.5		35		35.5			36

Lower Second Molar Band Sizes

Ortho Technology			29.5		30		30.5		31		31.5		32		32.5		33		33.5		34		34.5		35		35.5		36		36.5		37
American MR Full Sizes					1	2				4	5		6	7		8	9	10			11	12	13		14	15	16		17	18	19	20	
American MR 1/2 Sizes										1.5		2		2.5			3	3.5				4	4.5			5	5.5		6	6.5	7		7.5
Dentaurum Dentaform					1	2		3	4	5			6	7	8	9	10		11		12	13		14	15		16		17	18	19	20	
GAC-SnapFit	1	2	3	4			5		6		7	8		9		10	11		12	13		14		15	16	17		18		19	20		21
Ormco				1	2		3		4		5	6	7		8		9	10	11		12		13		14	15	16		17	18	19		20
RMO Truform									7		8		9		10		11			12	13			14	15	16	17			18	19	20	
TP GripTite					00		0		1	2	3		4	5	6	7	8		9	10	11		12	13	14	15	16	17	18	19	20	21	
Unitek Victory Series			29.5		30		30.5		31		31.5		32		32.5		33		33.5		34		34.5		35		35.5		36		36.5		37

This band conversion chart is for general reference only. Band sizing varies from company to company and may change periodically. Exact accuracy of this data is not expressed, implied or guaranteed by Ortho Technology.

EMOLARBANDS



See page 54





Upper First and Second Molar Band Sizes

Ortho Technology	37		37.5			38		38.5		39		39.5			40		40.5		41		41.5	42		42.5		43	43.5		44		
American MR Full Sizes		20		21	22	23		24		25			27	28			29		30		31			32							
American MR 1/2 Sizes	6.5		7			7.5		8		9			9.5		10		10.5		11		11.5			12		12.5	13		13.5		
Dentaurum Dentaform	21		22		23	24	25	26		27		28	29	30			31	32													
GAC-SnapFit	15		16		17		18	19	20			21	22		23		24	25			26		27			28		29			
Ormco	20		21				22	23	24		25	26			27		28	29	30			31			32						
RM0 Truform		19	20		21	22	23	24	25	26	27	28		29	30	31	32		33	34	35	36									
TP GripTite	20	21	22		23	24	25	26		27		29	30	31			32		33		34			35		36					
Unitek Victory Series	37		37.5			38		38.5		39		39.5			40		40.5		41		41.5	42		42.5		43	43.5		44		

Lower First Molar Band Sizes

Ortho Technology		36.5		37		37.5		38		38.5		39		39.5		40		40.5		41		41.5		42	42.5	43		43.5	44		
American MR Full Sizes	19		20		21	22	23		24	25		26			27	28	29			30		31			32						
American MR 1/2 Sizes	6	6.5	7		7.5			8				8.5			9	9.5		10.5		11		11.5			12	12.5		13	13.5		
Dentaurum Dentaform	23		24	25	26	27	28	29	30					31	32																
GAC-SnapFit	14	15			16	17	18		19		20	21	22		23		24		25		26		27		28		29		30	31	
Ormco	19		20		21		22	23	24	25	26		27		28	29		30	31	32											
RMO Truform	18		19		20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35		36									
TP GripTite	19	20	21		22		23	24	25	27		28		29	30	31		32		33		34			35	36					
Unitek Victory Series		36.5		37		37.5		38		38.5		39		39.5		40		40.5		41		41.5		42	42.5	43		43.5	44		

Lower Second Molar Band Sizes

Ortho Technology		37.5		38		38.5		39			39.5		40		40.5	41	41.5	42								
American MR Full Sizes		21		22		23	24	25	26	27			28	29	30											
American MR 1/2 Sizes				8					8.5		9	9.5	10		10.5		11									
Dentaurum Dentaform		21		22	23	24		25	26	27			28	29												
GAC-SnapFit	22		23		24		25	26		27		28	29		30											
Ormco		21		22		23	24	25			26	27	28		29	30	31	32								
RMO Truform	21	22		23	24	25	26	27		28	29	30	31	32												
TP GripTite		22		23	24	25	26	27	28		29		30	31	32	33	34	35								
Unitek Victory Series		37.5		38		38.5		39			39.5		40		40.5	41										



Prewelded Lingual Attachments for TruFit 2.0 Molar Bands

MADE USA











Prewelded Cleat

Molar	494-025
Cuspid/Bicuspid	494-024
1 cleat	65¢ ea

Prewelded Lug Molar

Item #: 494-028	
1 lug 35 ¢	ea.
Length 4.8 mm	

Prewelded Indented Lingual Sheath with Hook

UR/LL	494-038
UL/LR	494-039
1 sheath	\$3.80 ea.

Prewelded Indented Lingual Sheath

Item #: 494-032 1 sheath.......^{\$}3.15 ea.

Prewelded Button Oval

Oval Flat494-020 Oval Curved494-021 1 button \$1.50 ea.



Bondable Lingual Cleat

Item #: 400-204 1 pack \$28.30/pk 2 or more \$23.95/pk 10 per pack



Weldable Lingual Cleat

A. Molar494-025-10 **B.** Cuspid/Bicuspid494-024-10 1 pack\$**6.55/pk** 10 per pack



Direct Bond Lingual Cleat

Item #: 494-004-10 1 pack^{\$}20.95/pk Micro-etched 10 per pack



Indented Lingual Sheath with Ball Hook

UR/LL.......494-038-10 UL/LR......494-039-10 1 pack\$61.90/pk 10 per pack



Indented Lingual Sheath

Indent provides additional retention to hold palatal arch bar in sheath.

Item #: 494-032-10
1 pack\$33.25/pk
10 per pack

Bondable Attachments Loose

MADE USA



Direct Bond Eyelets Round Base



Direct Bond Eyelets Rectangle Base



Direct Bond Buttons Round Base

Item #: 494-001-10
Micro-etched
3.5 mm diameter curved one
piece base
1 pack\$20.95/pk
3 or more......\$17.95/pk

10 per pack



Direct Bond Buttons Rectangle Base

A. Oval Flat	400-201
B. Oval Curved	400-202
3.88 mm x 2.48 mr	n one piece base
1 pack	\$20.95/pk
3 or more	\$17.95/pk
10 ner nack	



Ortho Technology's **Bicuspid Bands**

Bicuspid Bands

Item #: 104-300









MAPE USA

Weldable Bicuspid Brackets Roth Prescription* and Standard Edgewise





5 brackets \$3.20 ea. Hooks placed on distal wings of bicuspid bracket

▲ Maxillary Arch	Torq.	Ang.	L/R	.018" Slot	.022" Slot
$\frac{4}{5}$ Bicuspid Twin	-7°	0°	R	W805-259	W805-359
75 75 Bicuspia Twiff			L	W805-260	W805-360
4/14/ Bicuspid Twin	-7°	0.	R	W805-261	W805-361
Bicuspid Twin with Hook			L	W805-262	W805-362
$\frac{4}{5}$ Bicuspid Twin	0°	0°	U	W802-009	W802-109

Lugs

Lingual Cleat

BicuspidW400-010 BicuspidW971-024 1 lug**40**¢ **ea.** 1 cleat......**65**¢ **ea.**

▼ Mandibular Arch	Torq.	Ang.	L/R	.018" Slot	.022" Slot
A A Biomenial Today	-17°	0°	R	W805-289	W805-389
4 4 Bicuspid Twin	-17	U	L	W805-290	W805-390
Bicuspid Twin	-17°	0°	R	W805-291	W805-391
4 4 with Hook		0	L	W805-292	W805-392
5 5 Bicuspid Twin	-22°	0°	R	W805-293	W805-393
5 5 Bicuspid Twin		22 0	L	W805-294	W805-394
5 5 Bicuspid Twin	-22°	0°	R	W805-295	W805-395
^{5 5} with Hook	-22	-22 0	L	W805-296	W805-396
45 45 Bicuspid Twin	0°	0°	U	W802-009	W802-109

^{*}The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Store and Organize Brackets, and Band Assemblies!

Stackable Organizer

- Modular interlocking design
- Clearly marked trays for quick product selection
- Constructed with high quality durable plastic

Band Organizer85180K Bracket Organizer......85181K l organizer......**\$33.50 ea.** 5 or more**\$28.99 ea.** Includes: band/bracket tray Organizer dimensions 11 1/4" W x 1 1/2" H x 8 1/4" D



UR 1st/2nd Molar Band Trav Item #: 4002-05

UL 1st/2nd Molar Band Tray Item #: 4002-06

LR 1st/2nd Molar Band Trav Item #: 4002-07

LL 1st/2nd Molar Band Tray Item #: 4002-08

Tray included with any first time kit orders.









63

BUCCALIUBES

TruEase[™]

Bondable Buccal Tubes




Designed With the Latest Technology for Patient Comfort

Adjustable Comfort Ball Hook

Funneled Entrance
Assists for easy
archwire insertion

Low Profile Design
For maximum
patient comfort



ID Coded by Quadrant

Occlusal/Gingival Positioning Indents

Distal Notch



- Stainless Steel Construction
- Molded Compound Contour Base with 80 Gauge Indented Mesh Pad provides an anatomically correct fit that conforms to the tooth ensuring maximum bond strength
- The Indented Pad assists in easy and accurate placement to conform into the buccal groove



- The Extended Base provides maximum coverage for greater bond strength and control
- 4.2 mm Mesial/Distal Dimension for full rotational control

Single Non-Convertible

	Quadrant	Slot Torque	Distal Offset		Extended Slot S .018"	
	Upper Right	-10°	8°		460-11-1020	460-12-1020
F	Upper Left	-10°	8°		460-11-2020	460-12-2020
	Lower Right	-25°	5°		460-11-3020	460-12-3020
F	Lower Left	-25°	5°		460-11-4020	460-12-4020
Λ	Upper Right	-14°	10°		460-21-1020	460-22-1020
Ï	Upper Left	-14°	10°		460-21-2020	460-22-2020
N	Lower Right	-20°	0°		460-21-3020	460-22-3020
٦	Lower Left	-20°	0°		460-21-4020	460-22-4020
	UR/LL	0°	0°	•	460-51-7020	460-52-7020
Ī	UL/LR	0°	0°	•	460-51-8020	460-52-8020
		1 pac	k (10 t	\$62.00/pk (6.20/tube)	
		5+ pack	ks (50 t	\$58.00/pk (5.80/tube)	

R The Ortho Technology version of the Roth Rx*

The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx* Standard Edgewise Rx

TruEase™ Low Profile Mini Bondable Buccal Tubes MHEUSA

Specifically Designed for Partially Erupted Second Molars

Adjustable

Comfort Ball Hook

Funneled Entrance
Assists for easy archwire insertion

Low Profile Design For maximum patient comfort



ID Coded by Quadrant

Occlusal/Gingival Positioning Indents

Allows for easy placement and minor adjustments



- Stainless Steel Construction
- Mesially Placed Compound Contour Base is designed to fit easily and precisely on the second molar mesial cusp to provide an anatomically correct fit
- Pylon Base with Electrical Discharge Machining (EDM) offers retention superior to mesh bonding pads
- Torque in Base allows easier alignment for optimal functional occlusion and case finishing
- 2.5 mm Mesial/Distal Dimension

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

© 2015 Ortho Technology, Inc. TruEase is a trademark of Ortho Technology.

	Quadrant	Slot Torque	Distal Offset		Slot 9 .018"	Size .022"
Ţ	Upper Right	-10°	0°		470-11-1720	470-12-1720
R	Upper Left	-10°	0°		470-11-2720	470-12-2720
	Lower Right	-25°	0°		470-11-3720	470-12-3720
R	Lower Left	-25°	0°		470-11-4720	470-12-4720
N	UR/LL	-14°	0°		470-21-7723	470-22-7723
Ï	UL/LR	-14°	0°		470-21-8723	470-22-8723
S	UR/LL	0°	0°	•	470-51-7720	470-52-7720
	UL/LR	0°	0°	•	470-51-8720	470-52-8720
	1 pack (10 tubes)				\$62.00/pk (6.20/tube)
	5+ packs (50 tubes)				\$58.00/pk (5.80/tube)

R The Ortho Technology version of the Roth Rx³

The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*

S Standard Edgewise Rx



TruEase Buccal Tubes Mini Typodont

Item #: G50-000

1 mini typodont......**\$135.00 ea.** 3 or more**\$125.00 ea.**

Includes: 1x TruEase Buccal Tube and 1x TruEase Low Profile Mini Buccal Tube

Votion® Bondable Buccal Tubes

VOTION

MADEUSA



Integrated Appliance System Designed to Work Seamlessly with Votion Brackets

The Votion Buccal Tube is a single non-convertible tube designed with modern features for practitioners' and patients' needs. Placement and positioning are made easy with side positioning indents, a defined buccal indent, and a compound contour base. Patients benefit from a smooth, low-profile rounded shape, and lowered hook for arguably the most comfortable treatment experience. Votion Buccal Tubes are the finishing touch to the appliance system that is designed to work seamlessly with Votion Brackets.

Lowered Hook

Funneled Entrance with notch assists in easy archwire insertion into the buccal tube



ID Color-Coded

Occlusal/Gingival Positioning Indents

Defined Buccal Indent

Single Non-Convertible

	Quadrant	Slot Torque	Distal Offset		Extende Slot	
	Upper Right	-14°	14°	0	351-11-1620	351-12-1620
Ŀ	Upper Left	-14°	14°	0	351-11-2620	351-12-2620
	Lower Right	-25°	4°	0	351-11-3620	351-12-3620
ď	Lower Left	-25°	4°	0	351-11-4620	351-12-4620
Ν	Upper Right	-14°	10°	0	351-21-1620	351-22-1620
	Upper Left	-14°	10°	0	351-21-2620	351-22-2620
Ν	Lower Right	-20°	0°	0	351-21-3620	351-22-3620
	Lower Left	-20°	0°	0	351-21-4620	351-22-4620
	1 pack (10 tubes)				\$59.90/pk (\$5.99/tube)
	5+ packs (50 tubes)				\$56.90/pk (\$5.69/tube)

R The Ortho Technology version of the Roth Rx*

The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*

- 80-Gauge Micro-Etched Mesh Pad ensures maximum bond strength
- Metal Injection Molding (MIM) Stainless Steel
- Low Profile for maximum patient comfort
- Occlusal/Gingival Positioning Indents for easy placement and minor adjustments
- Smooth Rounded Surfaces and Lowered Hook position and orientation for improved patient comfort
- Defined Buccal Indent and Compound Contoured Base for ease of placement and anatomical fit that conforms to the tooth



SEAMLESS INTEGRATION WITH OUR...

Votion^o Stainless Steel Bracket System

See pages 36-38



Votion® Mini Bondable Buccal Tubes

MADE USA

Specifically Designed for the Most Comfortable Treatment

Funneled Entrance with notch assists in easy archwire insertion into the buccal tube

Lowered Hook



ID Color-Coded

Occlusal/Gingival Positioning Indents

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

© 2015 Ortho Technology, Inc. Votion is a registered trademark of Ortho Technology.

		Slot		Color		Size
	Quadrant	Torque	Offset	Code	.018"	.022"
	Upper Right	-10°	0°	0	352-11-1620	352-12-1620
	Upper Left	-10°	0°	0	352-11-2620	352-12-2620
	Lower Right	-25°	0°	0	352-11-3620	352-12-3620
	Lower Left	-25°	0°		352-11-4620	352-12-4620
Ν	Upper Right	-14°	10°	0	352-21-1620	352-22-1620
П	Upper Left	-14°	10°	0	352-21-2620	352-22-2620
N	Lower Right	-10°	0°	0	352-21-3620	352-22-3620
П	Lower Left	-10°	0°	0	352-21-4620	352-22-4620
	1 pack (10 tubes)				\$59.90/pk (\$5.99/tube)
	5+ packs (50 tubes)				\$56.90/pk (\$5.69/tube)

R The Ortho Technology version of the Roth Rx*

The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx*

BUCCALIL



Positioning Notch To assist with placement

and avoid interference

Removable Color-**Coded Identification** Marked for quick and accurate placement. Red

upper and yellow lower

Low Profile Hook For patient comfort



Indented Mesh Pad

For a precise contoured fit

Mesial/ Distal **Dimension** 3.6 mm

Chamfered **Entry**

Guides wires into buccal tube



Maximum Bond Strength with a Comfortable Anatomical Fit

- Color coded scribe line for height and angulation reference point, Red = Upper, Yellow = Lower
- Indented 80 gauge mesh pad for an anatomical fit and maximum bond strength
- Occlusal and gingival positioning notches for ease of placement and handling

1 pack (10 tubes)	§45.50/pk (§4.55/tube)
5 packs (50 tubes)	\$41.00/pk (\$4.10/tube)
10+ packs (100 tubes)	\$39.00/pk (\$3.90/tube)

© 2015 Ortho Technology, Inc. TruCast is a trademark of Ortho Technology.

Single Non-Convertible - Bondable

	Slot	Distal	Slot	Size
Quadrant	Torque	Offset	.018"	.022"
Upper Right	-10°	+8°	204B-825	204B-225
Upper Left	-10°	+8°	204B-826	204B-226
Lower Right	-25°	+5°	204B-827	204B-227
Lower Left	-25°	+5°	204B-828	204B-228
Upper Right	-14°	+10°	204B-833	204B-233
Upper Left	-14°	+10°	204B-834	204B-234
Upper Right	-10°	0°	204B-831U	204B-231U
Upper Left	-10°	0°	204B-832U	204B-232U
Lower Right	-10°	0°	204B-831L	204B-231L
Lower Left	-10°	0°	204B-832L	204B-232L
Upper Right	0°	0°	204B-821U	204B-221U
Upper Left	0°	0°	204B-822U	204B-222U
Lower Right	0°	0°	204B-821L	204B-221L
Lower Left	0°	0°	204B-822L	204B-222L

The Ortho Technology version of the Roth Rx* M The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx3

Standard Edgewise Rx

*The Ortho Technology versions of these prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.



2.4 mm

Chamfered Entry

Guides wires into buccal tube

• Smooth, compact design for patient comfort

• Chamfered archwire entry point for easy wire placement

• Contoured, mesially extended base for an accurate, precise fit

Low profile hook is adjustable for patient comfort

Positioning notch for quick alignment

Parallel design for secure tube pick-up and placement

• Bondable and weldable

For patient comfort

1 pack (10 tubes)......\$42.50/pk (\$4.25/tube) 5 packs (50 tubes)\$38.50/pk (\$3.85/tube) 10+ packs (100 tubes)\$36.50/pk (\$3.65/tube)



TruCast™ Mini Second Molar **Bondable Buccal Tubes**

Allows for Earlier Treatment of Partially Erupted Second Molars

Quadrant	Slot Torque	Distal Offset	Slot .018"	Size .022"
Lower Right	-25°	0°	304B-833	304B-233
Lower Left	-25°	0°	304B-834	304B-234
Upper Right	-10°	0°	304B-831	304B-231
Upper Left	-10°	0°	304B-832	304B-232
UL/LR	0°	0°	304B-822	304B-222
UR/LL	0°	0°	304B-821	304B-221

WIREPRODUCTS

ARCHWIRE

An extensive selection of high quality wire types and arch forms

TruFlex™ Nickel Titanium Archwire	P.70-71, P.73, P.78-79
TruFlex™ Thermal Nickel Titanium Archwire	P.72-73
TruFlex™ Copper Nickel Titanium Archwire	P.74-75
TRIFECTA® Thermal Nickel Titanium Archwire	P.76
BetaForce™ Beta Titanium Archwire	P.73, P.77
TruForce™ Stainless Steel Archwire	P.73, P.80-83
Tooth Tone™ Coated Nickel Titanium Archwire	P.84
Tooth Tone™ Coated Stainless Steel Archwire	P.84
TruGold™ 24K Gold Plated Archwire	P.85



WIREPRODUCTS


Archwire Forms






Available in: Etched and Dimpled Etched and Dimpled


Universal Form




Euro Form



Standard Form



Pre-Torqued Full Form



Reverse Curve of Spee



WIREPRODU

Full Form Archwire

The broad, natural archform is the preferred shape for most clinical applications. It is ideal for use with bracket systems that have a built-in prescription, and can be used from start to finish with little to no wire bending.

Available in:

TruFlex™

Nickel Titanium Archwire

TruFlex™ Copper Nickel Titanium Archwire

Tooth Tone® Coated Nickel Titanium Archwire

TruForce™ Stainless Steel Archwire

TruGold™ 24K Gold Plated Archwire

TruFlex™

Thermal Nickel Titanium Archwire

TRIFECTA® Thermal Nickel Titanium Archwire

Tooth Tone® Coated Stainless Steel Archwire

BetaForce™ Beta Titanium Archwire

Euro Form Archwire

The euro arch form has a slightly narrower anterior region when compared to the full form shape, and remains narrow in the posterior region for clinicians who prefer a more oval form.

Available in:

TruFlex™ Nickel Titanium Archwire

> TruFlex™ Copper Nickel Titanium Archwire

TruForce™ Stainless Steel Archwire

TruFlex™ Thermal Nickel Titanium Archwire

BetaForce™ Beta Titanium Archwire

Standard Form Archwire

The standard form, also known as a classic arch form, strictly reflects general arch symmetry. Often it is re-shaped and contoured and has enough length for wire bends, loops, and adjustments.

Available in:

TruForce™ Stainless Steel Archwire

Pre-Torqued Full Form

TruFlex Nickel Titanium Pre-Torqued Archwire provides low force, precise 20° anterior torque while enabling earlier introduction of torque in the upper centrals and laterals. Three anterior segment widths are available with etched center marks for accurate placement. Doctors and patients will appreciate the shorter chairtime provided by eliminating archwire bends.

Available in:

TruFlex™ Nickel Titanium Archwire

Universal Form™

The wide universal arch can be used on the upper or lower and creates the broad smile that displays the first and second bicuspids when viewed from the anterior - perfect for use with self-ligating bracket systems.

Available in:

TruFlex™ Nickel Titanium Archwire

TruFlex™ Copper Nickel Titanium Archwire

TruForce™ Stainless Steel Archwire

TruFlex™ Thermal Nickel Titanium Archwire

BetaForce™ Beta Titanium Archwire

Reverse Curve of Spee & Straight Leg Reverse Curve of Spee

The same great features as regular TruFlex NiTi archwires combining the advantages of super elastic, continuous forces and the benefits of a reverse curve. The Straight Leg RCS TruFlex NiTi is designed with additional length. Both are excellent wires for bite opening and closing due to deep or open bites, lower arch expansion, intrusion of the lower anterior teeth, and also provide excellent leveling and alignment.

Available in:

TruFlex™ Nickel Titanium Archwire

Retraction Archwire

TruForce Stainless Steel Retraction Archwires are available in an assortment of looped or brass post archwires for anterior intrusion or retraction, space closure, or elastic engagement.

Available in:

TruForce™ Stainless Steel Archwire

· Stainless Steel Retraction Archwires

Specialty Archwires

TruForce Stainless Steel is also available straight wire lengths, or three different multi-strand archwires: 3-strand twist, 6-strand coaxial, and 8-strand braided archwire. Each multi-strand wire features light force levels, good flexibility, and can be cut and formed without unraveling.

Available in:

TruFlex™ Nickel Titanium Archwire

· Nickel Titanium Spooled Wire

TruForce™ Stainless Steel Archwire

· 3-Strand Stainless Steel Twist Wire Straight Lengths

· 6-Strand Stainless Steel Coaxial Wire Straight Lengths

· 8-Strand Stainless Steel Braided Archwire Stainless Steel Straight Lengths

BetaForce™ Beta Titanium Archwire

. Beta Titanium Straight Lengths








69

Truflex

Nickel Titanium

Patients will appreciate the comfortable force levels provided during the initial treatment phase. Doctors will love the accurate and consistent results in treatment. TruFlex Nickel Titanium productsguaranteed quality at competitive prices.

- Predictable, light continuous force levels
- Bright finish for optimal sliding mechanics
- Resistant to deformation for high spring back
- Consistent accurate arch shapes
- Satisfaction guaranteed



$\overline{\bigcap}$

TruFlex™ Nickel Titanium Full Form Archwires

The broad, natural archform is the preferred shape for most clinical applications. TruFlex NiTi uses super elastic properties and is provided in the austenitic or activated state. TruFlex NiTi offers mid-range force characteristics, consistent force levels, a smooth bright finish, and excellent rebound and memory properties.





MADE USA

	Etched	Midline	Dimpled Midline		
Round	▲ Upper	▼ Lower	▲ Upper	▼ Lower	
.012"	6000-101	6000-201	D6000-101	D6000-201	
.014"	6000-102	6000-202	D6000-102	D6000-202	
.016"	6000-103	6000-203	D6000-103	D6000-203	
.018"	6000-104	6000-204	D6000-104	D6000-204	
.020"	6000-105	6000-205	D6000-105	D6000-205	
■ Rectangle					
.016" x .016"	6000-106	6000-206	D6000-106	D6000-206	
.016" x .022"	6000-107	6000-207	D6000-107	D6000-207	
.017" x .025"	6000-108	6000-208	D6000-108	D6000-208	
.018" x .018"	6000-109	6000-209	D6000-109	D6000-209	
.018" x .025"	6000-110	6000-210	D6000-110	D6000-210	
.019" x .025"	6000-111	6000-211	D6000-111	D6000-211	
.020" x .020"	6000-113	6000-213	D6000-113	D6000-213	
.021" x .025"	6000-112	6000-212	D6000-112	D6000-212	

Round Etched

1 pack	\$17.75/pk
10 packs	\$13.35/pk
20 or more	
10 per pack	•

■ Rectangle Etched

1 pack	\$21.95/pk
10 packs	
20 or more	\$13.25/pk
10 per pack	•

© 2015 Ortho Technology, Inc. TruFlex is a trademark of Ortho Technology.

Round	Dimpled

I pack	•20.75/pk
10 packs	\$17.25/pk
20 or more	\$12.25/pk
10 per pack	

■ Rectangle Dimpled

= Rectangle Dil	пріси
1 pack	\$24.95/pk
10 packs	
20 or more	\$16.95/pk
10 per pack	•









TruFlex™ Nickel Titanium **Euro Form Archwires**

The euro arch form has a slightly more narrow anterior region when compared to the full form shape. TruFlex NiTi uses super elastic properties and is provided in the austenitic or activated state. TruFlex NiTi offers mid-range force characteristics, consistent force levels, a smooth bright finish, and excellent rebound and memory properties.

Round Etched

1 pack	.\$17.75/pk
10 packs	
20 or more	
10 per pack	•

■ Rectangle Etched

1 pack	\$21.95/pk
10 packs	
20 or more	\$13.25/pk
10 ner nack	

Round Dimpled

1 pack	\$20.75/pk
10 packs	\$17.25/pk
20 or more	\$12.25/pk
10 per pack	•

■ Rectangle Dimpled

1 pack	\$24.95/pk
10 packs	
20 or more	\$16.95/pk
10 per pack	•

	Etched	Midline	Dimpled	Midline
Round	▲ Upper	▼ Lower	▲ Upper	▼ Lower
.012"	5000-101	5000-201	D5000-101	D5000-201
.014"	5000-102	5000-202	D5000-102	D5000-202
.016"	5000-103	5000-203	D5000-103	D5000-203
.018"	5000-104	5000-204	D5000-104	D5000-204
.020"	5000-105	5000-205	D5000-105	D5000-205
■ Rectangle				
.016" x .016"	5000-106	5000-206	D5000-106	D5000-206
.016" x .022"	5000-107	5000-207	D5000-107	D5000-207
.017" x .025"	5000-108	5000-208	D5000-108	D5000-208
.018" x .018"	5000-109	5000-209	D5000-109	D5000-209
.018" x .025"	5000-110	5000-210	D5000-110	D5000-210
.019" x .025"	5000-111	5000-211	D5000-111	D5000-211
.020" x .020"	5000-113	5000-213	D5000-113	D5000-213
.021" x .025"	5000-112	5000-212	D5000-112	D5000-212



TruFlex™ Nickel Titanium **Reverse Curve of Spee**

The same great features as regular TruFlex NiTi archwires combining the advantages of super elastic, continuous force levels and the benefits of a reverse curve. RCS TruFlex NiTi are designed for bite opening and closing due to deep or open bites and also provides excellent leveling and alignment.

Round Etched

1 pack	\$28.95/pk
5 or more	\$25.65/pk
10 per pack	

■ Rectangle Etched

1 pack	§30.95/pk
5 or more	\$27.75/pk
10 ner nack	

Round Dimpled

1 pack	\$29.95/pk
5 or more	\$26.25/pk
10 per pack	

■ Rectangle Dimpled

1 pack	
5 or more	· •
10 ner nack	· •

Etched Midline

Round	▲ Upper	▼ Lower
.012"	6000-311	6000-411
.014"	6000-300	6000-400
.016"	6000-301	6000-401
.018"	6000-302	6000-402
.020"	6000-312	6000-412
■ Rectangle		
.016" x .016"	6000-303	6000-403

.016 X .022	6000-304	6000-404
.017" x .025"	6000-305	6000-405
.018" x .018"	6000-310	6000-410
.018" x .025"	6000-306	6000-406
.019" x .025"	6000-307	6000-407
.020" x .020"	6000-309	6000-409
.021" x .025"	6000-308	6000-408

Dimpled Midline

▲ Upper	▼ Lower
D6000-311	D6000-411
D6000-300	D6000-400
D6000-301	D6000-401
D6000-302	D6000-402
D6000-312	D6000-412

D6000-403
D6000-404
D6000-405
_

D6000-306	D6000-406
D6000-307	D6000-407
D6000-309	D6000-409
D6000-308	D6000-408



TruFlex™ Nickel Titanium **Straight Leg Reverse Curve of Spee**

The same great features as regular TruFlex NiTi archwires combining the advantages of super elastic, continuous forces and the benefits of a reverse curve. Straight Leg RCS TruFlex NiTi are designed with additional length. Excellent wire for bite opening and closing, lower arch expansion, intrusion of the lower anterior teeth and also provides excellent leveling and alignment.

Round Etched

1	pack	.\$2	8.	95,	/pk
5	or more	.\$2	5.	65,	/pk
1	0 per pack				

• Rectangle Etched

1 pack	\$30.95/pk
5 or more	\$27.75/pk
10 per pack	-

Etched Midline

Kouna	▲ Upper	▼ Lower
.014"	5000-300	5000-400
.016"	5000-301	5000-401
.018"	5000-302	5000-402
.020"	5000-312	5000-412
■ Rectangle		
.016" x .016"	5000-303	5000-403
.016" x .022"	5000-304	5000-404
.017" x .025"	5000-305	5000-405
.018" x .018"	5000-310	5000-410
.018" x .025"	5000-306	5000-406
.019" x .025"	5000-307	5000-407
.020" x .020"	5000-309	5000-409
.021" x .025"	5000-308	5000-408

MADE USA

Truflex

Thermal Nickel Titanium



TruFlex™ Thermal Nickel Titanium Full Form Archwires

The broad, natural arch form is the preferred shape for most clinical applications. TruFlex Thermal NiTi is provided in the martensitic or soft and pliable state and becomes super elastic once it reaches the transformation temperature of 90°F/32°C and becomes activated or austenitic. TruFlex Thermal NiTi provides gentle, lower force levels and a smooth, bright finish. The extreme pliability of TruFlex Thermal NiTi wires offers easy ligation and therefore larger size wires can be introduced at earlier phases in treatment without compromising patient comfort.

Round Etched

1 pack	\$21.50/pk
10 packs	
20 or more	\$12.35/pk
10 per pack	

■ Rectangle Etched

1 Pack	•20.23/pk
10 packs	\$22.45/pk
20 or more	\$18.15/pk
10 per pack	-

■ Round Dimpled

1 pack	\$24.25/pk
10 or more	\$19.85/pk
10 per pack	

■ Rectangle Dimpled

1 pack	\$28.65/pk
10 or more	\$25.75/pk
10 per pack	_

MADE USA

	Etched Midline		Dimpled Midline	
Round	▲ Upper	▼ Lower	▲ Upper	▼ Lower
.014"	6001-014	6001-114	D6001-014	D6001-114
.016"	6001-016	6001-116	D6001-016	D6001-116
.018"	6001-018	6001-118	D6001-018	D6001-118
.020"	6001-020	6001-120	D6001-020	D6001-120
■ Rectangle				
.014" x .025"	6001-040	6001-140	D6001-040	D6001-140
.016" x .016"	6001-041	6001-141	D6001-041	D6001-141
.016" x .022"	6001-042	6001-142	D6001-042	D6001-142
.016" x .025"	6001-043	6001-143	D6001-043	D6001-143
.017" x .025"	6001-044	6001-144	D6001-044	D6001-144
.018" x .018"	6001-047	6001-147	D6001-047	D6001-147
.018" x .025"	6001-046	6001-146	D6001-046	D6001-146
.019" x .025"	6001-048	6001-148	D6001-048	D6001-148
.020" x .020"	6001-049	6001-149	D6001-049	D6001-149
.021" x .025"	6001-050	6001-150	D6001-050	D6001-150



TruFlex™ Thermal Nickel Titanium Euro Form Archwires

The euro arch form has a slightly more narrow anterior region when compared to our full form shape. Provides gentle low forces for maximum patient comfort, allowing for the introduction of larger wire sizes earlier in treatment without compromising patient comfort. TruFlex Thermal NiTi is extremely soft and pliable at room temperature for easy ligation. Once ligated, its shape memory becomes active with the patient's natural body temperature (90°F/32°C). Highly smooth, with a bright finish for reduced friction.

Round Etched

1 pack	\$21.50/pk
10 packs	\$17.25/pk
20 or more	\$12.35/pk
10 per pack	•

■ Rectangle Etched

1 pack	^{\$} 26.25/pk
10 packs	\$22.45/pk
20 or more	\$18.15/pk
10 per pack	

Round Dimpled

1 pack	\$24.25/pk
10 or more	
10 per pack	

■ Rectangle Dimpled

1 pack	\$28.65/pk
10 or more	\$25.75/pk
10 per pack	

MAREUSA T

	Etched Midline		e Dimpled Midline	
Round	▲ Upper	▼ Lower	▲ Upper	▼ Lower
.014"	5001-102	5001-202	D5001-102	D5001-202
.016"	5001-103	5001-203	D5001-103	D5001-203
.018"	5001-104	5001-204	D5001-104	D5001-204
.020"	5001-105	5001-205	D5001-105	D5001-205
■ Rectangle				
.016" x .016"	5001-106	5001-206	D5001-106	D5001-206
.016" x .022"	5001-107	5001-207	D5001-107	D5001-207
.017" x .025"	5001-108	5001-208	D5001-108	D5001-208
.018" x .018"	5001-109	5001-209	D5001-109	D5001-209
.018" x .025"	5001-110	5001-210	D5001-110	D5001-210
.019" x .025"	5001-111	5001-211	D5001-111	D5001-211
.020" x .020"	5001-113	5001-213	D5001-113	D5001-213



Universal Form™ Archwires MEUSA

The wide universal arch can be used on the upper or lower and creates the broad smile that displays the first and second bicuspids when viewed from the anterior. Perfect for use with self-ligating bracket systems. The Universal Form archwire provides all the features of our TruForce, TruFlex, and BetaForce wires with the added benefit of our carefully designed arch form.


- Perfect archwire sequencing for self-ligating bracket systems
- Precise and consistent force levels
- Excellent resiliency
- Etched midline

Recommended for use with Self-Ligating Brackets

Actual Size: .014" x .025"

	Size	Item #	Force Level	Archwire
Phase 1	.012"	2000-100		
Initial	.013"	2000-101		TruFlex
Early alignment and	.014"	2000-102	Light	Nickel Titanium
unraveling, start arch	.016"	2000-103		Nickei Titanium
development	.018"	2000-104		
	.013"	2001-101		
	.014"	2001-102	Extra Light	TruFlex
	.016"	2001-103	EXITA LIGHT	Thermal Nickel Titanium
	.018"	2001-104		

	Size	Item #	Force Level	Archwire
Phase 2 Intermediate Transitional.	.014" x .025" .016" x .025" .018" x .025"	2000-114 2000-116 2000-118	Light	TruFlex Nickel Titanium
complete leveling and aligning, correct rotations, and begin	.014" x .025" .016" x .022" .016" x .025"	2001-114 2001-216 2001-116	Extra Light	TruFlex
torque control	.016 x .025	2001-118		Thermal Nickel Titanium







See page 86

TruFlex™ Nickel Titanium

Round	
1 pack	\$17.75/pk
10 packs	\$13.35/pk
20 or more	§8.95/pk
10 per pack	

- Rectaligle	
1 pack	\$21.95/pk
10 packs	\$17.75/pk
20 or more	\$13.35/pk
10 nor nack	

TruFlex™ Thermal Nickel Titanium

Round	
1 pack	\$21.50/pk
10 packs	\$17.25/pk
20 or more	\$12.35/pk
10 per pack	

■ Rectangle	
1 pack	\$26.25/pk
10 packs	\$22.45/pk
20 or more	\$18.15/pk
10 per pack	

TruForce™ Stainless Steel

• Round 1 pack 5 or more 100 per pack	/ -
■ Rectangle	
1 pack	\$7.10/pk
5 or more	\$5.95/pk
10 per pack	•

BetaForce™ Beta Titanium

■ Rectangle	
1 pack	\$45.99/pk
5 or more	
10 per pack	

© 2015 Ortho Technology, Inc. TruFlex, TruForce, BetaForce and Universal Form are trademarks of Ortho Technology.

COPPER Nickel Titanium

The Latest Advancement in Nickel Titanium Archwire Technology

The inclusion of copper to the nickel titanium alloy accentuates the thermal properties and working characteristics of the wire.

The lower hysteresis of TruFlex Copper provides a lighter tactile "feel" of the wire at room temperature which translates to additional flexibility providing for easier placement on cases with crowding conditions.

The unique alloy composition not only gives you easier insertion into the brackets, but also provides a greater resistance to deformation. TruFlex Copper provides controlled loading and unloading forces for a greater working time range, resulting in fewer archwire changes needed.

These controlled force levels translate to gentle, continuous forces providing consistent and predictable tooth movement throughout treatment.



Universal Form Available with Stops and without

Available in 27°C or 35°C Force Levels

- Use 27°C Copper When Higher Force Levels are Needed
- Use 35°C Copper When Medium Force Levels are Needed

TruFlex™ Copper Nickel Titanium Universal Form Archwires





10 INDIVIDUALLY PACKAGED ARCHWIRES

■ Round ■ Rectangle Without Stops 1 pack......**\$63.99/pk** 10 packs**\$55.95/pk** 10 individually packaged archwires

Round	■ Rectangle With Stops
	\$73.99/pk
10 packs.	\$65.95/pk
10 individ	ually packaged archwires

MADE USA

Etched Midline

• Round	Transition Temp	Without Stops	With Stops
.013"	27° C	2005-101	2005-101ST
.014"	27° C	2005-102	2005-102ST
.016"	27° C	2005-103	2005-103ST
.018"	27º C	2005-104	2005-104ST
■ Rectan	gle		
.014" x .02	25" 27° C	2005-105	2005-105ST

.016" x .025" 27° C 2005-106 2005-106ST .017" x .025" 27° C 2005-107 2005-107ST .018" x .025" 27° C 2005-108 2005-108ST .019" x .025" 27° C 2005-109 2005-109ST .017" x .025" 35° C 2005-110 2005-110ST .019" x .025" 35° C 2005-111 2005-111ST











TruFlex™ Copper Nickel Titanium Universal Form Archwires

MADE USA

Stops

Etched Midline

Round	■ Rectangle Without Stops
1 pack	\$53.99/pk
10 packs.	\$45.95/pk
10 per pa	ck loose

Round	■ Rectangle With Stops
1 pack	\$63.99/pk
10 packs.	\$56.95/pk
10 per pa	ck loose

				.014" x .025"	27° C	2004-105	2004-105ST
	Transitio	n Without	With	.016" x .025"	27° C	2004-106	2004-106ST
Round	Temp	Stops	Stops	.017" x .025"	27° C	2004-107	2004-107ST
.013"	27° C	2004-101	2004-101ST	.018" x .025"	27° C	2004-108	2004-108ST
.014"	27º C	2004-102	2004-102ST	.019" x .025"	27º C	2004-109	2004-109ST
.016"	27º C	2004-103	2004-103ST	.017" x .025"	35° C	2004-110	2004-110ST
.018"	27º C	2004-104	2004-104ST	.019" x .025"	35° C	2004-111	2004-111ST



TruFlex™ Copper Nickel Titanium **Full Form Archwires**

Etched Midline

|--|

I pack	.*50.95/pk
10 packs	.\$43.95/pk
10 per pack loose	

■ Rectangle

1 pack	.\$55.95/pk
10 packs	\$49.95/pk
10 ner nack loose	

^{\$} 50.95/pk ^{\$} 43.95/pk		Transitio	n Witho	out Stops
oose	Round	Temp	▲ Upper	▼ Lower
	.014"	27° C	2007-102	2007-202
	.016"	27º C	2007-103	2007-203
\$55.95/pk	.016"	35° C	2007-110	2007-210
\$49.95/pk	.018"	27º C	2007-104	2007-204

.018"

Round

.014"

.016"

.016"

.018"

Transition

27° C

27° C

35° C

27º C

Temp ▲ Upper

2008-102

2008-103

2008-104

MADE USA

MADE USA

Ti	ransition	Withou	it Stops
■ Rectangle	Temp	▲ Upper	▼ Lower
.014" x .025"	27° C	2007-105	2007-205
.016" x .022"	27° C	2007-106	2007-206
.016" x .022"	35° C	2007-112	2007-212
.016" x .025"	27º C	2007-107	2007-207
.017" x .017"	35° C	2007-113	2007-213
.017" x .025"	27º C	2007-108	2007-208
.017" x .025"	35° C	2007-114	2007-214
.018" x .025"	27º C	2007-118	2007-218
.018" x .025"	35° C	2007-115	2007-215
.019" x .025"	27° C	2007-109	2007-209
.019" x .025"	35° C	2007-116	2007-216
.020" x .020"	35° C	2007-117	2007-217

Transition Without

Stops

■ Rectangle Temp



TruFlex™ Copper Nickel Titanium Euro Form Archwires

Etched Midline


Round

1 pack	\$50.95/pk
10 packs	
10 per pack loo.	se

■ Rectangle

1 pack	.\$55.95/pk
10 packs	
10 per pack loose	

10 per pack loose	.018"
© 2015 Ortho Technology Inc. TruFlex is a trademark of Ortho	Гесhnology.



•		***************************************	it stops
■ Rectangle	Temp	▲ Upper	▼ Lower
.014" x .025"	27º C	2008-105	2008-205
.016" x .022"	27° C	2008-106	2008-206
.016" x .022"	35° C	2008-112	2008-212
.016" x .025"	27º C	2008-107	2008-207
.017" x .017"	35° C	2008-113	2008-213
.017" x .025"	27º C	2008-108	2008-208
.017" x .025"	35° C	2008-114	2008-214
.018" x .025"	27º C	2008-118	2008-218
.018" x .025"	35° C	2008-115	2008-215
.019" x .025"	27º C	2008-109	2008-209
.019" x .025"	35° C	2008-116	2008-216
.020" x .020"	35° C	2008-117	2008-217



Without Stops

2008-110 2008-210

2008-111 2008-211

▼ Lower

2008-202

2008-203

2008-204

35°C 2007-111 2007-211




75

TRIFECTA®

Thermal Nickel Titanium

The Right Force on the Right Teeth at the Right Time



Anterior Region = Gentle Forces

The incisor segment has the lowest heat-activated force level specific to moving small, single rooted teeth. The archwire provides gentle forces for patient comfort, reduces unnecessary tooth pressure and decreases the possibility of root resorption.

Bicuspid Region = Moderate Forces

This segment of the archwire delivers a slightly increased force level for moving the larger bicuspids or premolars. The progressively increased force provides quick leveling and alignment to speed up treatment time, without compromising patient comfort levels.

Molar Region = Heavy Forces

The posterior segment of the archwire generates the most force in order to align the larger, more difficult to move first and second molars. Regular archwires have only one force level throughout, which would require 2 or 3 wire changes to achieve the same movement you can get with 1 TRIFECTA wire.

TRIFECTA Force Levels in Grams

Wire Size	Anterior	Bicuspid	Molar
.016" x .016"	65	85	120
.018" x .018"	80	110	155
.020" x .020"	100	150	215
.016" x .022"	105	155	220
.017" x .025"	120	195	275
.018" x .025"	140	200	290
.019" x .025"	155	210	300

■ Rectangle	▲ Upper	▼ Lower
.016" x .016"	6003-041	6003-141
.018" x .018"	6003-047	6003-147
.020" x .020"	6003-049	6003-149
.016" x .022"	6003-042	6003-142
.017" x .025"	6003-044	6003-144
.018" x .025"	6003-046	6003-146
.019" x .025"	6003-048	6003-148

MADE USA

■ Rectangle (Etched	
5 or more 10 per pack	

© 2015 Ortho Technology Inc. TRIFECTA is a registered trademark of Ortho Technology.

3-Point Bend Test performed on a 10.0 mm wire span in the Buccal-Lingual plane to allow for curvature of the archwire. Tooth-moving forces at 1.5 mm deflection.

Beta Force Beta Titanium



Force Levels Between Nickel Titanium and Stainless Steel / Nickel-Free

The BetaForce Beta Titanium archwires are nickel-free with properties that fall between the bendability of Stainless Steel and the flexibility of Nickel Titanium. BetaForce wires have a smooth, polished surface for excellent sliding mechanics and provide moderate wire stiffness levels at 60% less than Stainless Steel. Beta Titanium wires offer excellent formability, making it ideal for intermediate and finishing stages of treatment.

Full Form Archwire

Archwire Size	▲ Upper	▼ Lower
.016"	6002-103	6002-203
.018"	6002-104	6002-204
.016" x .016"	6002-106	6002-206
.016" x .022"	6002-107	6002-207
.017" x .025"	6002-108	6002-208
.018" x .018"	6002-109	6002-209
.018" x .025"	6002-110	6002-210
.019" x .025"	6002-111	6002-211
.021" x .025"	6002-112	6002-212

BetaForce™

1 pack\$45.99/pk 5 or more\$42.99/pk

10 per pack

NEW!

Euro Form Archwire

Archwire Size	▲ Upper	▼ Lower
.016"	6002-113	6002-213
.018"	6002-114	6002-214
.016" x .016"	6002-116	6002-216
.016" x .022"	6002-117	6002-217
.017" x .025"	6002-118	6002-218
.018" x .018"	6002-119	6002-219
.018" x .025"	6002-120	6002-220
.019" x .025"	6002-121	6002-221
.021" x .025"	6002-122	6002-222

BetaForce™

1 pack	\$45.99/pk
5 or more	\$42.99/pk
10 per pack	

BetaForce™ Beta Titanium

14" Straight Lengths

.016"	Х	.016"	6002-406
.016"	х	.022"	6002-407
.017"	Х	.025"	6002-408
.018"	Х	.018"	6002-409
.018"	х	.025"	6002-410
.019"	Х	.025"	6002-411
.021"	Х	.025"	6002-412

1 pack	\$45.95/pk
3 or more	
10 ner nack	-

© 2015 Ortho Technology Inc. BetaForce is trademarks of Ortho Technology.

ORTHOGOP, inc

Disposable Archwire Markers Protect Against Cross Contamination

- · Consistent machine dipped wax tips
- Special thick wax tips will not crack or melt
- . Made with sturdy dowels

Dispenser

Item #: 600-500 1 dispenser \$19.95 ea.

- · The attractive dispenser keeps working area organized
- · Precise marking and adhesion to dry metal surfaces
- Inexpensive and disposable

Markers

Item #: 600-401
1 pack
3 or more
100 markers per pack,
Made in the USA



The only sure way to avoid potential cross-contamination is to use the inexpensive disposable markers..."



"ITU-Flex"

Nickel Titanium

TruFlex™ Nickel Titanium Pre-Torqued Full Form Archwires

Ideal for Use with Passive Self-Ligating Bracket Systems


TruFlex Nickel Titanium Pre-Torqued Archwire provides low force, precise 20° anterior torque while enabling earlier introduction of torque in the upper centrals and laterals. Three anterior segment widths are available with etched centermarks for accurate placement. Doctors and patients will appreciate the shorter chair-time provided by eliminating archwire bends.

	Etched Midline		
■ Rectangle	▲ Upper 34 mm	▲ Upper 38 mm	▼ Lower 28 mm
.016" x .022"	2000-122	2000-124	2000-120
.017" x .025"	2000-128	2000-130	2000-126
.018" x .025"	2000-134	2000-136	2000-132
.019" x .025"	2000-140	2000-142	2000-138
.021" x .025"	2000-146	2000-148	2000-144

Not recommended for use with ceramic brackets.

■ Rectangle Etched

1 pack	§65.95/pk
10 packs	\$63.95/pk
10 per pack	



TruFlex™ Nickel Titanium De-Impacting and Separating Springs

- Designed to effectively separate or de-impact posterior teeth with minimal patient discomfort
- Offset leg design and rectangle nickel titanium wire provides the most efficient method of separating and de-impacting
- Looped over safety ends ensure patient comfort

4 weeks from Start to Finish!











© 2015 Ortho Technology, Inc. TruFlex is a trademark of Ortho Technology.

TruFlex™ Nickel Titanium Spools

- Super elastic and resilient
- Continuous mid-range force levels
- Made in the USA

.012"	6000-000
.014"	6000-001
.016"	6000-002
.018"	6000-003


1 spool	\$22.50 ea.
5 or more	
15' per spool	

TruFlex









TruFlex™ Nickel Titanium **Distalizing Springs**

- Distalize molars without the use of headgear
- Achieve 1.0 mm to 1.5 mm distalization per month
- Fits over any size wire
- Force level decreases as the length of spring increases

_	
065° Dia Cut	

MAREUSA =

.010" x .045"	700-100
1 pack	\$44.50/pk
3 or more	\$39.95/pk
3 per pack. 7" lenat	ths



TruFlex™ Nickel Titanium **Open Coil Springs**

MADE USA

- Designed for maximum results
- Delivers continuous force for opening or holding of spaces between teeth

Lengths

.010" x .030"	700-101
.012" x .030"	700-102
.014" x .036"	700-103
1 pack	\$25.95/pk
4 or more	\$23.50/pl
3 springs per pack	c, 7" lengths

Spools	
.010" x .030"	700-000
.012" x .030"	700-001
1 spool	\$24.50 ea.
4 or more	\$22.50 ea.
15" spool	

Manufallan da de la company de



Unloading Forces ← (Tooth-Moving)

Shown At 200% of Actual Size Recommended Active Spring Length is Two Times Passive Spring Length

TruFlex™ Nickel Titanium **Closed Coil Springs** MADE USA

TruFlex Nickel Titanium Closed Coil Springs deliver the ideal amount of consistent retraction force needed to close or maintain any space.

- Delivers constant force levels
- Active until the spring is closed
- Superelastic Nickel Titanium quality
- Provides comfortable. predictable forces during treatment
- Stainless steel key-end eyelets for easy attachment

Light Force 100 gm

9.0 mm (.010" x .030")700-104L 12.0 mm (.010" x .030")700-105L

Medium Force 150 gm

9.0 mm (.011" x .030")700-104M 12.0 mm (.011" x .030")700-105M

Heavy Force 200 gm

9.0 mm (.012" x .030")......700-104H 12.0 mm (.012" x .030")......700-105H

1 pack\$45.50/pk 5 or more\$43.25/pk 10 per pack



TruFlex™ Nickel Titanium TAD **Closed Coil Springs** MADE USA

Specifically for use with the Spider Screw. Can also be used with other TAD systems with a head size smaller than 3.0 mm.

Light Force 100 gm


9.0 mm (.009" x .030")......900-104L 12.0 mm (.009" x .030").....900-105L

Medium Force 150 gm

9.0 mm (.011" x .030")......900-104M 12.0 mm (.011" x .030")900-105M

Heavy Force 200 gm

9.0 mm (.012" x .030")......900-104H 12.0 mm (.012" x .030")900-105H



10 per pack

U.S. Patent Number D595852.



TruForce archwires are manufactured with the highest quality 304VSS medical grade Stainless Steel. These wires have an extremely high tensile strength and accept intricate bends without breaking.

MADE USA

- High quality Stainless Steel
- Consistent performance
- High luster polish for reduced friction
- Etched permanent center mark

TruForce™ Stainless Steel **Full Form Archwires**

	Round	
1	pack	

\$21.95/pk 100 per pack

■ Rectangle

10 or more\$5.95/pk 10 per pack

Round	▲ Upper	▼ Lower
.012"	8000-101	8000-201
.014"	8000-102	8000-202
.016"	8000-103	8000-203
.018"	8000-104	8000-204
.020"	8000-105	8000-205

■ Rectangle	▲ Upper	▼ Lower
.016" x .016"	8000-106	8000-206
.016" x .022"	8000-107	8000-207
.017" x .017"	8000-115	8000-215
.017" x .022"	8000-108	8000-208
.017" x .025"	8000-109	8000-209
.018" x .018"	8000-110	8000-210
.018" x .022"	8000-111	8000-211
.018" x .025"	8000-112	8000-212
.019" x .025"	8000-113	8000-213
.021" x .025"	8000-114	8000-214

TruForce™ Stainless Steel Euro Form Archwires

Round

1 pack......**\$21.95/pk** 100 per pack

■ Rectangle

10 per pack

Kouna	▲ Opper	* Lower
.012"	9000-101	9000-201
.014"	9000-102	9000-202
.016"	9000-103	9000-203
.018"	9000-104	9000-204
.020"	9000-105	9000-205

.016" x .016"	9000-106	9000-206
.016" x .022"	9000-107	9000-207
.017" x .022"	9000-108	9000-208
.017" x .025"	9000-109	9000-209
.018" x .018"	9000-110	9000-210
.018" x .022"	9000-111	9000-211
.018" x .025"	9000-112	9000-212
.019" x .025"	9000-113	9000-213
.021" x .025"	9000-114	9000-214

▲ Upper

▼ Lower

■ Rectangle



TruForce™ Stainless Steel Standard Form Archwires

Round

1 pack......\$21.95/pk 100 per pack

■ Rectangle

10 per pack

● Kouliu	T obbei	* LOWEI
.012"	7000-101	7000-201
.014"	7000-102	7000-202
.016"	7000-103	7000-203
.018"	7000-104	7000-204
.020"	7000-105	7000-205

■ Rectangle	▲ Upper	▼ Lower
.016" x .016"	7000-106	7000-206
.016" x .022"	7000-107	7000-207
.017" x .022"	7000-108	7000-208
.017" x .025"	7000-109	7000-209
.018" x .018"	7000-110	7000-210
.018" x .022"	7000-111	7000-211
.018" x .025"	7000-112	7000-212
.019" x .025"	7000-113	7000-213
.021" x .025"	7000-114	7000-214

© 2015 Ortho Technology Inc. TruForce is a trademark of Ortho Technology.







TruForce™ Stainless Steel Retraction Archwires

MADE USA

- Anterior retraction
- Anterior intrusion
- Space closure
- Elastic engagement

Post or Loop Measurement Measure along incisal edge from distal of lateral to distal of opposite lateral

- Custom formed to exact measurements
- 304 medical grade Stainless Steel

Looped Archwire

10 per pack

Brass Post Archwire

10 per pack

2 Keyhole Wire





Arch Size	.016" x .022"	.017" x .025"	.019" x .025"
24.0 mm	2KL-224	2KL-424	2KL-624
26.0 mm	2KL-226	2KL-426	2KL-626
28.0 mm	2KL-228	2KL-428	2KL-628
30.0 mm	2KL-230	2KL-430	2KL-630
32.0 mm	2KL-232	2KL-432	2KL-632
34.0 mm	2KL-234	2KL-434	2KL-634
36.0 mm	2KL-236	2KL-436	2KL-636
38.0 mm	2KL-238	2KL-438	2KL-638
40.0 mm	2KL-240	2KL-440	2KL-640
42.0 mm	2KL-242	2KL-442	2KL-642
44.0 mm	2KL-244	2KL-444	2KL-644

4 Keyhole Wire





Arch Size	.016" x .022"	.019" x .025"
24.0 mm	4KL-224	4KL-624
26.0 mm	4KL-226	4KL-626
28.0 mm	4KL-228	4KL-628
30.0 mm	4KL-230	4KL-630
32.0 mm	4KL-232	4KL-632
34.0 mm	4KL-234	4KL-634
36.0 mm	4KL-236	4KL-636
38.0 mm	4KL-238	4KL-638
40.0 mm	4KL-240	4KL-640
42.0 mm	4KL-242	4KL-642
44.0 mm	4KL-244	4KL-644

Stainless Steel Brass Post Wire





Arch Size	.016" x .022"	.017" x .025"	.018" x .025"	.019" x .025"
24.0 mm	BPT-224	BPT-424	BPT-524	BPT-624
26.0 mm	BPT-226	BPT-426	BPT-526	BPT-626
28.0 mm	BPT-228	BPT-428	BPT-528	BPT-628
30.0 mm	BPT-230	BPT-430	BPT-530	BPT-630
32.0 mm	BPT-232	BPT-432	BPT-532	BPT-632
34.0 mm	BPT-234	BPT-434	BPT-534	BPT-634
36.0 mm	BPT-236	BPT-436	BPT-536	BPT-636
38.0 mm	BPT-238	BPT-438	BPT-538	BPT-638
40.0 mm	BPT-240	BPT-440	BPT-540	BPT-640
42.0 mm	BPT-242	BPT-442	BPT-542	BPT-642
44.0 mm	BPT-244	BPT-444	BPT-544	BPT-644

T-Loop Wire Wire





Arch Size	.016" x .022"	.017" x .025"	.019" x .025"
24.0 mm	TLP-224	TLP-424	TLP-624
26.0 mm	TLP-226	TLP-426	TLP-626
28.0 mm	TLP-228	TLP-428	TLP-628
30.0 mm	TLP-230	TLP-430	TLP-630
32.0 mm	TLP-232	TLP-432	TLP-632
34.0 mm	TLP-234	TLP-434	TLP-634
36.0 mm	TLP-236	TLP-436	TLP-636
38.0 mm	TLP-238	TLP-438	TLP-638
40.0 mm	TLP-240	TLP-440	TLP-640
42.0 mm	TLP-242	TLP-442	TLP-642
44.0 mm	TLP-244	TLP-444	TLP-644







81

Truforce

TruForce[™]

Stainless Steel

Straight Lengths



- Light force levels with good flexibility
- 302SS Stainless Steel
- Smooth bright finish
- Made in the USA

TruForce™
Stainless Steel
3-Strand Twist
Wire Straight
Lengths

Size	Archwire	Lengths	Spools
.0155"	705-201	705-101	705-106
.0175"	705-202	705-102	705-107
.0195"	705-203	705-103	705-108

Α	rchwire	
1	pack	\$20.95/pk
5	or more	\$17.95/pk
1	0 per pack	-

Ler	ngths	Spools
1 pa	ack\$14.50/pk	1 spool \$31.25 ea
5 or	more\$13.50/pk	5 or more \$25.50 ea
10 E	per pack, 14" lengths	30' per spool

TruForce

6-Strand

Coaxial Wire

Stainless Steel

Straight Lengths



- Light force levels with good flexibility
- 302SS Stainless Steel
- Smooth bright finish

■ Made in the USA				
Size	Archwire	Lengths	Spools	
.0155"	705-501	705-401	705-301	
.0175"	705-502	705-402	705-302	

Archwire	
1 pack	\$20.95/pk
5 or more	\$18.95/pk
10 per pack	


Lengths	Spools
1 pack\$14.50/pk	1 spool \$31.25 ea.
5 or more \$13.50/pk	5 or more \$25.50 ea.
10 per pack, 14" lengths	30' per spool



TruForce™ Stainless Steel 8-Strand Braided Archwire

- Stable for finishing and holding
- 302SS Stainless Steel
- Gentle force levels with excellent torque control
- Smooth bright finish
- Does not unravel when cut
- Etched midline
- Made in the USA

.016" x .016" .016" x .022" .017" x .025" .018" x .025"	710-101 710-102 710-103
.021" x .025"	
1 pack 5 or more 10 per pack	



- .014".....900-101 .016".....900-102 .018".....900-103 .020".....900-104 .016" x .016"900-114 .016" x .022"900-115 .017" x .025"900-116 .018" x .018"900-117 .018" x .025"900-118 .019" x .025"900-119 .021" x .025"900-120
- 302SS Stainless Steel
- Made in the USA

Round

1 pack\$	13.25/pk
4 or more\$	11.75/pk
25 per pack, 14" leng	ths

■ Rectangle

1 pack	\$17.25/pk
4 or more	\$13.95/pk
25 per pack, 1	4" lenaths



TruForce[™] **Stainless Steel Retainer Wire Straight Lengths**

.022"	900-105
.025"	900-106
.026"	900-107
.028"	900-108
.030"	900-109
.032"	900-110
.036"	900-111
.040"	900-112
045"	900-113

- 302SS Stainless Steel
- High tensile strength
- Smooth bright finish
- Made in the USA

1 pack	\$13.25/pk
4 or more	\$11.75/pk
25 per pack, 14'	' lenaths



Open Coil Spring Spools

Open wound Stainless Steel wire designed to open spaces.

- 304V Stainless Steel
- Made in the USA

.010" x .030"	908-101
.010" x .036"	908-103
1 spool	\$8.95 ea.
5 or more	\$7.95 ea.
3' per spool	



Spring Spools

Closed Coil

Closed wound Stainless Steel wire designed to hold spaces.

- 304V Stainless Steel
- Made in the USA

908-102
\$8.95 ea.
\$7.95 ea.

Truforce[™]



TruForce™ Stainless Steel Preformed Ligature Ties



- High quality dead-soft 304 medical-grade Stainless Steel
- Made in the USA

.009"	906-101
.010"	906-102
.011"	906-103
.012"	906-104
1 pack	\$22.95/pk
10 or more	\$20.95/pk
1000 per pack	

TruForce™ Stainless Steel Kobayashi Ties

- High quality dead-soft 304 medical-grade Stainless Steel
- Quickly place retraction hooks on brackets without hooks
- Made in the USA

.010"	905-100
.012"	905-101
.014"	905-102
1 pack	\$18.50/pk
5 or more	\$16.25/pk
100 per pack	


TruForce™ Stainless Steel Preformed Shorty Ligature Ties



- High quality dead-soft 302 medical-grade Stainless Steel
- Twisted ends allow quick placement over brackets

.010"	906-105
1 pack	\$21.95/pk
5 or more	\$19.95/pk
500 ner nack	

TruForce™ Stainless Steel Kobayashi Shorty Ties



- High quality dead-soft 302 medical-grade Stainless Steel
- Spot welded hooks and twisted ends allow quick placement over brackets

.012"	905-105
.014"	905-103
1 pack	\$26.95/pk
10 or more	
100 per pack	-



TruFlex™ Nickel Titanium Chilling Spray

- Quickly cool and shrink hot retainer material
- Turns NiTi wire soft instantly
- Makes fitting or adjusting easy
- Instantly chills wire down to -62°F/-52°C
- Made in the USA

Item #: 01-2200	1
1 can	\$27.50 ea.
3 or more	\$25.50 ea.
10 oz (284 gm)	per can

Hazardous materials, special



TruForce[™] **Stainless Steel Spooled Ligature** Wire

- High quality dead-soft 304 medical-grade Stainless Steel
- Diamond drawn bright finish
- Made in the USA

.010"	720-	100
1 spool	\$39.50	ea
3 or more		
1 lb/454 am per	spool	



See pages 84-85





ooth one®



Tooth Tone® Coated Nickel Titanium Full Form Archwires

The same great features as our TruFlex Nickel Titanium archwires with a durable tooth colored overcoat which offers superior aesthetics, stain resistance, and low friction.



Round	▲ Upper	▼ Lower
.014"	6000-601	6000-701
.016"	6000-602	6000-702
.018"	6000-603	6000-703
.020"	6000-604	6000-704
■ Rectangle		
.016" x .022"	6000-607	6000-707
.018" x .018"	6000-608	6000-708
.018" x .024"	6000-609	6000-709
.019" x .027"	6000-610	6000-710
.020" x .020"	6000-612	6000-712

Round	
1 pack	^{\$} 51.99/pl
5 packs	548.99/pl
10 per pack	

■ Rectangle 5 packs \$52.99/pk 10 per pack

Note: Size shown is finished size. Tooth Tone coating adds .002" thickness to diameter.



Round	▲ Upper	▼ Lower
.014"	8000-802	8000-902
.016"	8000-803	8000-903
.018"	8000-804	8000-904
.020"	8000-805	8000-905
■ Rectangle		
.016" x .022"	8000-807	8000-907
.017" x .025"	8000-809	8000-909
.018" x .018"	8000-810	8000-910
.018" x .024"	8000-812	8000-912
.019" x .027"	8000-813	8000-913
.020" x .020"	8000-817	8000-917

1 pack 5 packs 10 per pack	
■ Rectangle 1 pack 5 packs	\$33.99/pk \$31.99/pk

Note: Size shown is finished size. Tooth Tone coating adds .002" thickness to diameter.

Tooth Tone® Coated Preformed Shorty Ties Durable plastic coating



- Stainless steel

.010" wire <i>(.012</i> " ltem #	,
1 pack 5 or more 100 per pack	

Tooth Tone® Coated Preformed



.010" wire (.012	2" coatea)
ltem #	907-102
	\$26.85/pk \$24.15/pk

© 2015 Ortho Technology Inc. Tooth Tone is a registered trademark of Ortho Technology.

Tooth Tone® Nickel Titanium Open Coil Springs

10 per pack

■ Pound

Consistent force level



Coated for superior aesthetics

7" Lengths

.012" x .030"......907-003 (.014" x .032" coated) 1 pack**\$58.95/pk**

Tooth Tone® Coated Kobayashi Shorty Ties



Durable plastic coating

3 springs per pack

Stainless steel

.012" wire (.014	" coated)
	907-115
1 pack	§65.04/pk
5 or more	§62.94/pk
100 per pack	

24K GOLD PLATED ARCHWIRE



Real 24K Gold Plating

Get all the aesthetic benefits and elegance of gold plating while also getting the features of our Nickel Titanium and Stainless Steel archwires.

Round	▲ Upper	▼ Lower
.012"	GN-012	GN-112
.014"	GN-014	GN-114
.016"	GN-016	GN-116
.018"	GN-018	GN-118
.020"	GN-020	GN-120
■ Rectangle		
.016" x .016"	GN-041	GN-141
.016" x .022"	GN-042	GN-142
.017" x .025"	GN-044	GN-144
.018" x .018"	GN-045	GN-145
.018" x .025"	GN-046	GN-146
.019" x .025"	GN-048	GN-148
.020" x .020"	GN-049	GN-149
.021" x .025"	GN-050	GN-150



Round	
1 pack	§44.95/pk
5 packs	\$40.95/pk
10 per pack	

■ Rectangle 1 pack......\$48.95/pk 10 per pack



● Kouliu	▲ Opper	* Lower
.014"	GSF-014	GSF-114
.016"	GSF-016	GSF-116
.018"	GSF-018	GSF-118
.020"	GSF-020	GSF-120
■ Rectangle		
.016" x .016"	GSF-041	GSF-141
.016" x .022"	GSF-042	GSF-142
.017" x .025"	GSF-044	GSF-144
.018" x .025"	GSF-046	GSF-146
.019" x .025"	GSF-048	GSF-148
.021" x .025"	GSF-050	GSF-150

MADE USA

	Round	
1	pack	.\$37.95/pk
5	packs	.\$34.95/pk
1	0 per pack	

■ Rectangle	
1 pack	\$38.95/pk
5 packs	\$35.95/pk
10 per pack	

TruGold™ 24K Gold Plated Preformed Ligature Ties



- Real 24K gold plating
- Stainless steel
- Made in the USA

.010"	G906-102
1 pack	\$83.85/pk
3 or more	\$80.95/pk
100 per pack	•

Pound



TRY OUR...

Marquis™ TruGold™ **24K Gold Plated Bracket System**

Your Complete Source for Adding the Golden Touch

See page 46



© 2015 Ortho Technology Inc. TruGold and Marquis are trademarks of Ortho Technology.

CRIMPABLE STOP APPLICATOR

Save Chair Time With No Lost Archwire Stops!





Score Crimpable Stop Applicators

Score Crimpable Stop Applicators were designed for the clinician's convenience, ease of use, and to save valuable chair time. The pre-loaded applicator offers convenience for easy handling and accuracy when placing crimpable tube stops on archwires. Once positioned, the Score applicator is easily removed with just a twist.

Small

Small .019"	020-400
Large .030"	020-500
1 pack	\$22.95/pk
4 or more	\$19.95/pk
10 per pack	

Follow These Easy Steps!


Step 1.

Insert wire through tube and slide into desired position



Step 2.

Twist to remove applicator and to release crimpable stop



Step 3.

Crimp stop with crimping pliers to secure stop on archwire



Reduce Emergency Visits!

No More Poking Archwires! Crimpable Mini Stops

- Extremely low profile and small in size, making them more aesthetically pleasing than typical crimpable stops
- Very comfortable to the patient due to the small size
- Slides easily onto the archwire for customized placement
- Crimps easily into place with Weingart Style Utility Pliers
- 2.0 mm in length
- Stainless steel
- Made in the USA



Standard Packs

1 pack	\$39.95/pk
4 or more	\$36.95/pk
25 per pack	•



Mini Stops - Small

Standard Packs020-200 Large Economy Packs 100003 For wire sizes .012" through .018"

Mini Stops - Large

Standard Packs......020-300 Large Economy Packs 100004 For rectangular wires .014" x .025" through .019" x .025" & square wires up to a maximum of .020" x .020"

Large Economy Packs

1 pack	§62.75/pk
4 or more	
50 per pack	•

Tooth Tone® Crimpable Mini Stops







- For total aesthetic treatment
- Very comfortable to the patient

Coated Mini Stops - Small Item #: 907-000

For round wires up to .018" 25 per pack





- Slides easily onto the archwire
- 2.0 mm in length

Coated Mini Stops - Large Item #: 907-001

For square wires up to .020" x .020" rectangular wires up to .019" x .025" 25 per pack

Crimpable Split Stops

- Metal injected molded
- Stainless steel

Item #: 020-100 1 pack......\$30.95/pk For round and rectangular wires up to .021" x .025", 10 per pack





ATTACHMENTS





Perfect for TAD Use! **Crimpable Archwire Hooks** with Split Archwire Opening

- Excellent for use with miniscrew attachments
- Easily and securely attach closed coil springs or any auxiliary attachments
- Can be placed onto archwires inside or outside of the mouth
- Curved 6.8 mm tall and fits archwires up to .022" x .025"
- Made in the USA



Patient Left Patient Right	
1 pack 3 or more	







Split archwire opening

• For total aesthetic treatment

Very comfortable to the patient

Crimps securely to archwire

■ 3.3 mm tall

Item #: 907-002 Fits archwires up to .022" x .025"

1 pack	\$62.50/pk
3 or more	
20 ner nack	• •

Crimpable Ball Hooks with Slide-On Archwire Opening

Made in the USA

Slide-on Ball Hooks: 3.2 mm tall

Item #: 014-214

Fits archwires up to .018" x .025"

Large

Item #: 014-215 Fits archwires up to .021" x .025"

1 pack**\$58.50/pk** 20 per pack





Crimpable Bi-Directional Ball Hooks with Split Archwire Opening

- Easily place onto any size archwire inside or outside of the mouth
- Secure anchor for any auxiliary attachment
- Bi-directional hook for easy ligation
- Smooth for patient comfort
- High quality stainless steel
- Fits archwires up to .022" x .025"
- Made in the USA

Short: 1.6 mm tall

Item #: 20102

3 or more......\$36.95/pk 20 per pack

Long: 3.3 mm tall Item #: 20103

3 or more \$40.95/pk 20 per pack





Crimpable Ball Hooks with Split Archwire Opening

- Made from high quality stainless steel
- Smooth, rounded ball hooks can be placed precisely without removing the archwire
- Crimps to rectangular archwires up to .022" x .025"
- Made in the USA

Split Ball Hooks: 3.3 mm tall Item #: 20105

3 or more \$43.95/pk 20 per pack

Split Ball Hooks: 3.2 mm tall

Item #: 014-114

Fits archwires up to .018" x .025"

Item #: 014-115

Fits archwires up to .021" x .025"

1 pack	.\$58.50/pk
2 or more	
20 per pack	













87



- Quick order entry
- Favorite items list
- Detailed item history
- Detailed order history
- One click order tracking
- State-of-the-art transaction security
- No credit card required bill on open account



2.0 mm length, .022" slot,

Crimpable Cross Tubes

- Easily connect auxiliary and custom attachments to archwire
- Perfect for all orthodontic anchorage systems
- Diverse TAD technique use
- High quality stainless steel
- Made in the USA





Crimpable Double Tubes

- Easily connect auxiliary and custom attachments to archwire
- Perfect for all orthodontic anchorage systems
- Diverse TAD technique use
- High quality stainless steel
- Made in the USA

Item #: 100300	
1 pack	
5 or more	\$50.95/pk
3.0 mm length,	.022" slot,
10 per pack	

Direct Bond Eyelets Round Base

■ Made in the USA

Item #: 100310

10 per pack

Item #: 025-100 3.5 mm diameter two piece flat base 3 or more\$27.95/pk 10 per pack



Direct Bond Eyelets Rectangle Base

Made in the USA

Item #: 025-000 3.75 mm x 3.0 mm one piece curved base 1 pack\$31.45/pk 10 per pack



Direct Bond Buttons Round Base

■ Made in the USA

Item #: 400-202R 3.5 mm diameter curved one piece base 1 pack......**\$20.95/pk** 10 per pack



Perfect for Aligner Systems!



OrthoFlex™ Composite Buttons

- Aesthetically pleasing appearance
- Durable stain resistant material
- Secure dovetail base design
- Made in the USA

Item #: OT-2344 2.85 mm x 2.05 mm one piece base 1 pack......\$42.50/pk 3 or more.......339.50/pk 10 per pack



Direct Bond Buttons Rectangle Base

Made in the USA

A 0 1 Fl-4	400 201
A. Oval Flat	400-201
B. Oval Curv	ved400-202
3.88 mm x 2.	48 mm one piece base
1 pack	\$20.95/pk
3 or more	\$17.95/pk
10 per pack	-

ATTACHMENTS

Güray Bite Raiser

Developed by Dr. Enis Güray Ask The Expert enisg@orthotechnology.com









Create a Temporary Bite Opening for **Deep** or **Cross-Bite Cases** that call for **Unobstructed Tooth Movement**

Easy Chairside Application


Simply shape and mold the Güray Bite Raiser for a custom fit, insert and ligate the proximal wings into the headgear tube and ligate the "T" Spur onto the lingual cleat.

Secure, Yet Removable

Güray Bite Raisers, unlike Bite Turbos and Bite Builders, can be easily disconnected at the lingual cleat to check the occlusion. Simply resecure if additional treatment time is needed with no debonding concerns, saving the clinician both time and money.

Features & Benefits

- Inexpensive and practical
- Easy and unobstructed tooth movement
- Ensure optimal oral hygiene
- Comfortable for patient
- · Easy insertion and removal insertion takes a couple of minutes, removal takes even less time



Bite Raiser Thickness

		0.8mm Child	1.0mm Adult
gt	4.0 mm	Item # 10030	Item # 10032
	4.0 mm Arm	(for occlusal h	eadgear tubes)
Raiser	6.0 mm 🛖	Item # 10031	Item # 10033
3ite Ra	Arm 🛱	(for gingival h	eadgear tubes)

Eliminates

- Patient compliance requirement
- Removable bite plates
- Laboratory work
- Restorative materials

Resulting In

- Less chairtime
- Less treatment time and expense

Güray Bite Raiser (Child)

0.8 mm x 4.0 m	m10030
0.8 mm x 6.0 m	m10031
1 pack	\$84.75/pk
3 or more	\$80.95/pk
10 per pack	

Güray Bite Raiser (Adult)
1.0 mm x 4.0 mm10032
1.0 mm x 6.0 mm10033
1 pack
3 or more ^{\$} 80.95/pk
10 per pack

© 2015 Ortho Technology, Inc. Güray Bite Raiser is a trademark of Ortho Technology. U.S. Patent Number 6,726,473.



Developed by Dr. Enis Güray



Ask the Expert

enisg@orthotechnology.com orthotechnology.com/EzFAQ

EZ Space Maintainer® Adjustable, Bondable Space Maintainer

No Impressions Necessary! Treatment in One Appointment!

Ortho Technology's EZ Space Maintainer is more cost-effective and less time-consuming than traditional space maintainer appliances. It requires no impressions, no laboratory construction, and can be directly bonded during one in-office visit. The EZ Space Maintainer is more aesthetic and hygienic than traditional space maintainer appliances. Simple and easy to use, the EZ Space Maintainer is available as a secure, adjustable space maintainer providing easy maintenance of the mesio-distal dimension of any lost, deciduous teeth or as an adjustable appliance by using the NiTi coil included to regain some space.

Features and Benefits:

- Simple and easy to use
- Direct bond appliance
- Less chairtime, less expense
- No lab fees

ATTACHMENTS

- More hygienic, more aesthetic
- Secure and adjustable space maintainer options
- Color coded identification marks by quadrant on 1st molar


EZ Space Maintainer Starter Kit

Item #: 10022
1 pack
3 or more \$105.95/pk
4x space maintainers per pack,
1 for each quadrant

EZ Space Maintainer Replacement Packs Individual Quadrants

Lower Right	
1 pack 3 or more	\$96.50/pk

EZ Space Maintainer Reference Diagram





EZ Space Maintainer Typodont Item #: G50-098

1 typodont ... \$270.00 ea. 1x typodont, 1 space maintainer per quadrant



EZ Space Maintainer Mini Typodont Item #: G50-101

1 typodont ... \$140.00 ea. 3 or more \$130.00 ea.





COME TECHNOLOGY





PLUS! Regain Space with **Optional Included Spring**



To regain space, simply activate the optional coil spring (shown above) by squeezing the stopper with any plier.

© 2015 Ortho Technology, Inc. EZ Space Maintainer is a registered trademark of Ortho Technology. Case studies photos courtesy of Dr. Serkan Gülec. Domestic and international patents pending.





Developed by Dr. Enis Güray enisg@orthotechnology.com

EZ Slider™

Easy Protraction or Retraction For All TAD Systems

The EZ Slider Creates Parallel Forces that Prevent Unwanted Tipping and Rotations While Allowing the Movement of Multiple Teeth

EZ Sliders are interchangeable auxiliary items that create distal or mesial forces depending on which quadrant they are used. They can easily be attached to any archwire with their "click-in-click-out" secure clip action arms which can be used in conjunction with any brand TAD coil spring and accompanying TAD. EZ Sliders ensure easy and accurate parallel force application that prevent unwanted tooth movements such as rotations and tipping. EZ Sliders are manufactured with medical grade 304 stainless steel.

Hook on Left

- Distal forces on the UL/LR segments
- Mesial forces on the UR/LL segments
 Hook on Right
- Mesial forces on the UL/LR segments
- Distal forces on the UR/LL segments

EZ Slider

Includes 3 Sizes:
Short 12.5 mm – Qty 4
Use with premolars and canines
Medium 20.0 mm – Qty 4
Use with first molars
Long 30.0 mm – Qty 4
Use with second molars

2x hooks on left, 2x hooks on right of each size




Solid tissue colored gingiva base with tooth tone colored teeth. 2x Brackets, 1x K1 Spider Screw, 20.0 mm EZ Slider, TruFlex NiTi Coil Spring, and 2x TruFit Molar Bands.

EZ Slider Mini Typodont

Item #: G50-097 1 mini typodont....\$235.00 ea.



To determine the size needed: Length for distal forces is determined by the amount of space needed to consolidate. Length for mesial forces is not critical.



Case Studies

PROTRACTION - MESIALIZATION

For Mesial Movement the Length of the Slider is not Critical. The Size is Related with the Space Needed for Protraction – Mesialization



Mesial Movement First Molar and Premolar Long (30.0 mm)



Mesial Movement First Molar and Premolar Short (12.5 mm)

RETRACTION - DISTALIZATION

Consecutive use of the Long - Medium - Small EZ Sliders



Distalizing Second Molar Long (30.0 mm) Aug. 2010



Distalizing First Molar Medium (20.0 mm) Oct. 2010



Distalizing PreMolars Short (12.5 mm) Dec. 2010

© 2015 Ortho Technology, Inc. EZ Slider is a trademark of Ortho Technology.

ORTHODONTIC ANCHORAGE

Spider Screw® Temporary Anchorage

Device System PORTHO TECHNOLOGY EXCLUSIVE





Giuliano Maino DDS, MD

Creator of the Spider Screw System

Dr. Maino maintains a private orthodontic practice in Vicenza, Italy, and is Visiting Professor of Orthodontics at Ferrara and Insubria University. He is an active member of the Angle Society of Europe, the Italian Society of Periodontology, American Association of Orthodontists, Italian Society of Orthodontics, and European Orthodontic Society. Dr. Maino is also a Post Graduate of Clinical Training in Orthodontics at Boston University.



Ask The Expert

Dr. Maino is Available to Answer **Your Questions**

Visit our website for many frequently asked TAD questions and answers. Submit your questions and get the answers you need.

www.orthotechnology.com/TADfaq



Orthodontic Head

- Designed to facilitate appliance (wires, springs, etc.) placement
- Features two intersecting .022" slots
- Under tie-wing area also functions as another .022" x .025" slot featuring two intersecting slots of .027" in diameter
- Chamfered inlets simplify insertion of wires or ligatures
- Small head dimension is designed for patient comfort



Transmucosal Portion

- Variable length allows for optimal adaptation to different intraoral mucosa thicknesses during biomaintenance
- Short length for areas of thinly attached gingiva
- Long for areas with thick or freely moveable tissues
- Polished with a special treatment to help avoid soft tissue irritation and make cleaning easier to accomplish



Intrabony Portion

- Thread shape has an asymmetrical profile for easy placement
- Ensures maximum stability and avoidance of bone stress




- Self-drilling and self-tapping makes pre-drilling before insertion unnecessary dependent upon bone structure
- Easy to place while reducing the risk of root damage









Clinical Cases

- Intrusion Posterior Areas
- Direct Anchorage Uprighting and Molar Intrusion
- Class II Correction
- Lower Molar Protraction Indirect Anchorage
- Molar Uprighting and Intrusion

Visit our website for detailed case studies and before/after photos

www.orthotechnology.com/SpiderScrewClinical



Spider Screw[®]

Temporary Anchorage Device System

Made From High Quality Materials

- Grade 5 titanium construction
- Nickel-free for sensitive patients



STERILE R

Self-Ligating K1 Plus Series

Conical: Self-Drilling Thread

The Spider Screw Self-Ligating TAD - K1 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.1 mm drill provided to penetrate the cortical plate. The locking head is color coded pink for easy identification.



Long Neck - 2.0 mm Length Item # 6.5 mm SXL-1506 8.0 mm SXL-1508 10.0 mm SXL-1510

3.9 mm diameter head 1.5 mm diameter body

Self-Ligating K2 Plus Series

Conical: Self-Drilling Thread

The Spider Screw Self-Ligating TAD - K2 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.2 mm drill provided to penetrate the cortical plate. The locking head is color coded green for easy identification.



Long Neck - 2.0 mm Length Item # 6.0 mm SXL-1906 SXL-1907 7.0 mm 9.0 mm SXL-1909 11.0 mm SXL-1911

3.9 mm diameter head 1.9 mm diameter body

K1 Series

Conical: Self-Drilling Thread

The Spider Screw - K1 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.1 mm drill provided to penetrate the cortical plate.



Short Neck - 1.0 mm				
Length	Item #			
6.5 mm	SCR-1506			
8.0 mm	SCR-1508			
10.0 mm	SCR-1510			

3.4 mm diameter head 1.5 mm diameter body



Long Neck - 2.0 mm Length Item # 6.5 mm SCL-1506 8.0 mm SCL-1508 10.0 mm SCL-1510

3.4 mm diameter head 1.5 mm diameter body

K2 Series

Conical: Self-Drilling Thread

The Spider Screw - K2 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.2 mm drill provided to penetrate the cortical plate.



Short Neck - 1.0 mm		
Length	Item #	
6.0 mm	SCR-1906	
7.0 mm	SCR-1907	
9.0 mm	SCR-1909	
11.0 mm	SCR-1911	
2 4 4:-	unantau laaad	

1.9 mm diameter body



Long Neck - 2.0 mm Lenath Item # 6.0 mm SCL-1906 7.0 mm SCL-1907 9.0 mm SCL-1909 11.0 mm SCL-1911 3.4 mm diameter head 1.9 mm diameter body

NEW! Spider Pin™

Cylindrical: Pre-Drilling Thread

The Spider Pin is self-tapping and requires pre-drilling. It is ideal for areas where a reduced size head is required and the infrabony portion is narrow. The simple design increases patient comfort and makes it easy to attach coil springs or elastic chains.



Long Neck - 2.0 mm Length Item # 8.0mm SCI-1308

2.6 mm diameter head 1.3 mm diameter body



Long Neck - 2.0 mm Lenath Item # 10.0mm SCI-1310

2.6 mm diameter head 1.3 mm diameter body

Individual Self-Ligating K1 & K2 Plus Series Screws

1	screw\$1	2	7	.95	ea.
6	or more\$1	1	6	.50	ea.
_					

Packaged sterile

Individual K1& K2 Screws

1 screw	\$97.50	ea
6 or more	\$89.95	ea
Packaged sterile		

Individual Spider Pins

1 pin	\$97.50	ea.
6 or more	\$89.95	ea.
Packaged sterile		

Short Neck: Reduced neck height for thin tissue (anterior and lateral areas) Long Neck: Oversize neck height for soft thick tissue (posterior and lateral areas)







93

ORTHODONTIC ANCHORAGE

Spider Screw® Self-Ligating TAD





Insert wire and turn Spider Screw Self-Ligating TAD head 90° with Self-Ligating Driver to secure

© 2015 Ortho Technology, Inc. Spider Screw is a registered trademark of Ortho Technology. Spider Screw Self-Ligating TAD is patented in the European Union and United States.

Spider Screw Self-Ligating TAD Starter Kit

- 1x 1.1 mm x 8.0 mm drill K1
- 1x 1.3 mm x 9.0 mm drill K2
- 1x Pick-up driver shaft SXL
- 1x Contra angle pick-up driver short SXL
- 1x Self-ligating TAD key
- 1x Handle driver SXL
- 1x Screw driver body
- 1x Organizer

Entire Kit Can Be Sterilized

- Autoclavable up to 273°F/134°C
- Heat Sterilizable up to 356°F/180°C
- Chemiclavable up to 297°F/132°C

Spider Screw Self-Ligating Starter Kit

Item #: CSS-6008

1 kit**\$945.00 ea.**

3 or more.......\$885.00 ea.

Replacement Organizer

Item #: CSS-6000

1 tray\$336.95 ea.

1 Year Warranty on Instruments

Self-Ligating Spider Screw Organizer

Included in Starter Kit



8.0 mm Drill - K1

1.1 mm x 8.0 mm

Item #: FSC-1108-S 1 drill**\$90.00 ea.**



9.0 mm Drill - K2 Item #: FSC-1309



Pick-Up Driver Shaft - SXL



Contra Angle Pick-Up Driver Short - SXL



Self-Ligating TAD Key Item #: DXL-2820

1 driver....**\$134.50 ea.**



Handle Driver - SXL

Item #: DSQ-3424 1 driver **\$96.50 ea.**



Screw Driver Body

Item #: DSX-1690N-S

Optional Items



Contra Angle Pick-Up Driver Long - SXL

Item #: DPQ-3425

1 driver............**\$134.95 ea.** *25.0 mm*



Great for Larger Hands or for a Firmer Grip

Optional Round Replacement End for Screw Driver Body DSX-1690N-S

Item #: RCX-2545



Wide and Round End Handles

provide easier, more ergonomic grip for larger hands

Screw Driver Body Wide Handle

Wide Handle Item #: DSX-1690RC

95

Spider Screw K1, K2 Starter Kit

- 1x 1.1 mm x 8.0 mm drill K1
- 1x 1.3 mm x 9.0 mm drill K2
- 1x Pick-up driver shaft
- 1x Cross driver shaft
- 1x Contra angle pick-up driver
- 1x Handle driver
- 1x Screw driver body
- 1x Organizer

Entire Kit Can Be Sterilized

- Autoclavable up to 273°F/134°C
- Heat Sterilizable up to 356°F/180°C
- Chemiclavable up to 297°F/132°C

Spider Screw Starter Kit

Item #: CSS-4008

1 kit.....**\$945.00 ea.**

Replacement Organizer - K1, K2

Item #: CSS-4000

1 tray......**\$359.50 ea.**

1 Year Warranty on Instruments



Included in Starter Kit



8.0 mm Drill - K1 Item #: FSC-1108-S

1 drill**\$90.00 ea.** 1.1 mm x 8.0 mm

9.0 mm Drill - K2

Item #: FSC-1309 1 drill**\$90.00 ea.**

1.3 mm x 9.0 mm

Pick-Up Driver Shaft Item #: DSP-5052N-S

52.0 mm

Cross Driver Shaft

Item #: DSX-2852N-S

52.0 mm



Contra Angle Pick-Up Driver

Item #: DPQ-2820-S

20.0 mm



Handle Driver

Item #: DSQ-2824-S



Screw Driver Body

Item #: DSX-1690N-S

1 handle........... \$117.95 ea.

Optional Items



- calibrated torque is reached to avoid overloading
- Torque can be set from 5 to 20 Ncm Note: Pickup driver sold
- Helps eliminate screw breakage

 Automatically releases when the Torque Screw Driver Body Item #: DST-1600-S«

1 screwdriver \$582.50 ea.

separately. Sterilizable, see

instructions for details.

Download instructions for use here www.orthotechnology.com/TorqueDriverIFU



10.0 mm Drill

Item #: FSC-1210-S«

1.2 mm x 10.0 mm



Contra Angle Pick-Up Driver

Item #: DPQ-2825 1 driver \$134.95 ea.

25.0 mm



Pick-Up **Handle Driver**

Item #: DPH-2824-S



Contra Angle Cross Driver

Item #: DPX-2830-S

30.0 mm

«Available by special order. Additional shipping time may be required.



Spider Pin™ Temporary Anchorage Device System

For Areas Where a Reduced Head Size is Required

- Simple and easy to use
- Smooth, rounded design for patient comfort
- Perfect for NiTi closed coil spring attachments
- Infrabony portion diameter is 1.3 mm at the widest point
- Ideal for narrow interproximal spaces
- Requires no patient cooperation and reduces treatment time
- Pre-drilling and self-tapping





Spider Pin Starter Kit

1.3 mm Diameter

- 1x 0.9 mm x 10.0 mm drill
- 1x Contra angle pick-up driver pin
- 1x Short handle driver pin
- 1x Pick-up handle driver shaft pin
- 1x Screw driver body
- 1x Organizer

Entire Kit Can Be Sterilized

- Autoclavable up to 273°F/134°C
- Heat Sterilizable up to 356°F/180°C
- Chemiclavable up to 297°F/132°C

Spider Pin Starter Kit

Item #: CSS-3006

1 kit**\$795.00 ea.**

1 Year Warranty on Instruments



Included in Starter Kit



10.0 mm Drill - Pin Item #: FSC-0910-S«

0.9 mm x 10.0 mm



Contra Angle Pick-Up Driver - Pin Item #: DPQ-2322-S«



Pick-Up Handle Driver Shaft - Pin

52.0 mm

Item #: DSP-2352N-S[«]
1 driver**\$134.95 ea.**

Short Handle Driver - Pin Item #: DSQ-2324-S*
1 driver..........\$135.95 ea.



Screw Driver Body Item #: DSX-1690N-S«

1 handle...........\$117.95 ea. 99.0 mm

© 2015 Ortho Technology, Inc. Spider Pin is a trademark of Ortho Technology. «Available by special order. Additional shipping time may be required.



Maximum™ Power Chain

0000

- Superior color stability, fade resistant
- High elasticity and excellent rebound
- Medical grade polyurethane
- Non-Latex*

		000	
	Closed	Short	Long
Gray	430-101	440-101	450-101
Clear	430-102	440-102	450-102
Tooth	430-103	440-103	450-103
Royal Blue	430-104	440-104	450-104
Light Blue	430-105	440-105	450-105
Purple	430-106	440-106	450-106
Lavender	430-107	440-107	450-107
Red	430-108	440-108	450-108
Maroon	430-109	440-109	450-109
■■ Black	430-110	440-110	450-110
White	430-111	440-111	450-111
Yellow	430-112	440-112	450-112
Pink	430-113	440-113	450-113
Aqua	430-114	440-114	450-114
Neon Green	430-115	440-115	450-115
Kelly Green	430-116	440-116	450-116
Dark Green	430-117	440-117	450-117
Orange	430-118	440-118	450-118
Teal	430-119	440-119	450-119
Fuchsia Fuchsia	430-126	440-126	450-126
Turquoise	430-127	440-127	450-127
Violet	430-129	440-129	450-129
Navy Blue	430-130	440-130	450-130
Cobalt Blue	430-133	440-133	450-133
Sky Blue	430-170	440-170	450-170
Blue Pearl	430-172	440-172	450-172
Pearl	430-171	440-171	450-171
Metallic Silver	430-120	440-120	450-120
Metallic Gold	430-121	440-121	450-121
Fluorescent Pink	430-160	440-160	450-160
Glow "Classic"	430-150	440-150	450-150

1 spool......\$22.49 ea. 15' per spool (4.572 m)

*Not made from natural rubber latex. Colors may vary.

		MELLY GREEN AGUN TURRICHEE LURIT BLUE PRED GREEN TELL CORRECT MARCON PRED GREEN TELL SAY LURIT BLUE	
)	Make Color Selection Fun!	ORTHO TECHNOLOGY WWW.orthotechnology.com PAPAL BLEE PEREL METALLIC SHOPE SERVER STANDARDCERST SERVER SHOPE S	

Elastomeric Color Palette

For Selection of Maximum Power Chain & Power Stick Elastomerics

Item #: 400-000 63/8"W x 43/8"H

Extreme™ Elastomeric Chain

Highest possible tensile strength with less deformation over an extended period of time when compared to regular chains. Extreme™ Chain offers superior rebound qualities, as material is thinner and 25% stiffer than maximum chain.



■ Made in the USA

5 or more\$35.49 ea.

15' per spool (4.572 m)

	Closed	Short	Long
Extreme Clear	430-902	440-902	450-902
Extreme Silver	430-920	440-920	450-920



Power Sticks™ Elastomeric Ligatures

- New secure grip tab design on Mini and Mid-Size Power Sticks
- New 24-tie stick choice
- Superior color stability, fade resistant
- Medical grade polyurethane (Non-Latex*)
- Made in the USA

Power Sticks

1,000 per pack (20 sticks of 50), .115" ligatures

MIW Mid-Size Power Sticks

1,008 per pack (42 mid-size sticks of 24), .115" ligatures

Mini Power Sticks

Each Mini Power Stick holds 10 ligatures to accommodate a single arch, reducing the potential risk of cross contamination. 1,000 per pack (100 mini sticks of 10), .115" ligatures

1 pack	.\$22.49/pk
10 packs	
20 or more	



Elastomeric Color Palette

NEW!

See page 98 for details

1	Power Stick	Mid-Size	Mini
Assorted Colors	404-000	N/A	N/A
Gray	404-001	405-001	407-001
Clear (Stain Resistant)	404-002	405-002	407-002
Tooth	404-003	405-003	407-003
Royal Blue	404-004	405-004	407-004
Light Blue	404-005	405-005	407-005
Purple	404-006	405-006	407-006
Lavender	404-007	405-007	407-007
Red	404-008	405-008	407-008
Maroon	404-009	405-009	407-009
■■ Black	404-010	405-010	407-010
White	404-011	405-011	407-011
Yellow	404-012	405-012	407-012
Pink	404-013	405-013	407-013
Aqua	404-014	405-014	407-014
Neon Green	404-015	405-015	407-015
Kelly Green	404-016	405-016	407-016
Dark Green	404-017	405-017	407-017
Orange	404-018	405-018	407-018
Teal	404-019	405-019	407-019
Fuchsia Fuchsia	404-026	405-026	407-026
Turquoise	404-027	405-027	407-027
Violet	404-029	405-029	407-029
Navy Blue	404-030	405-030	407-030
Cobalt Blue	404-033	405-033	407-033
Sky Blue	404-070	405-070	407-070
Blue Pearl	404-072	405-072	407-072
Pearl	404-071	405-071	407-071
Metallic Silver	404-020	405-020	407-020
Metallic Gold	404-021	405-021	407-021
Fluorescent Pink	404-060	405-060	407-060
Glow "Classic"	404-050	405-050	407-050
© 2015 Ortho Tachnology Inc. Mayi	Frances an	d Baway Chieles aus hu	adamaulia af

© 2015 Ortho Technology, Inc. Maximum, Extreme, and Power Sticks are trademarks of Ortho Technology. Colors may vary. *Not made from natural rubber latex.

Elastomeric Organizers Power Sticks™ **Combination Maximum™ Power Chain** Item #: OT-PSD MADE IN THE Item #: OT-CDW Item #: OT-PCD USA Organizers do not 3 1/8 "W x 4 1/4"H x 9 1/8 "D 101/2"W x 51/2"H x 51/4"D 63/4"W x 51/2"H x 9 "D include product.








99







- High quality medical grade materials
- Delivers consistent and accurate force
- Tooth diagram aids in proper placement
- Light resistant patient packs

MADE USA

Non-Latex* Elastics

1	box	\$29.5	0/pk
5	or more	\$28.5	0/pk

Natural Latex Elastics

1	box	.\$30.95/pk
5	or more	.\$29.95/pk

Neon Latex Elastics

1 box	.\$30.95/pk
5 or more	.\$29.95/pk
All elastics are packaged:	50 patient
packs of 100 elastics per p	pack.
*Not made from natural rubber latex	·.

Elastic Directors

Creatures of the Sea

Creatures of the Sea

Great White Shark

Sea Otter

Made in the USA

Elastic directors make placing and removing elastics easy for patients. Item #: 500-307

1 pack	\$22.95/pk
5 or more	
100 per pack, assorted	l blue, pink,
areen, and alow-in-the-	-dark.


















Light 2.5 oz. (71 gm)

Natural Latex







Seal



500-104



Medium
4.5 oz.
(128 gm)

Natural

(120 g,	
Natural Latex	
Neon Latex	
Non-Latex*	



500-201 500-401 500-201LF



500-202 500-203 500-402 500-403 500-202LF 500-203LF

Dolphin

500-204 500-404 500-204LF



Sea Turtle

500-405 500-205LF

Heavy 6.5 oz. (184 gm)





500-301
500-501
500-301LF



500-302 500-502 500-302LF



500-303	
500-503	
500-303LF	



I William !	,
500-304	
500-504	
500-304LF	



500-305 500-505 500-305LF

Force level is determined by measuring the force generated when the elastic is stretched three times its slack diameter. *Not made from natural rubber latex.




Creatures of the Sea® Elastic Organizer

Our elastics dispenser is constructed of high quality clear acrylic plastic and holds 5 boxes of elastics. Designed for counter top or wall mounting.

MADE USA

Item #: OT-OEDS 1 organizer... \$52.99 ea. 3 or more \$46.99 ea. 171/4" wide x 10 " tall x 41/4" deep

1/4" (6.4 mm)



3/8" (9.5 mm)



1/2" (12.7 mm)



Extraoral Latex Elastics

Delivers consistent and accurate force for extraoral applications. High quality medical grade materials.

Force Level 16 oz (454 gm)

Inner diameter

1/4" (6.4 mm).....500-506 5/16" (7.9 mm)...500-507

3/8" (9.5 mm).....500-508 1/2" (12.7 mm)...500-509

1 pack \$23.50/pk 3 or more .. \$21.95/pk 1,000 elastics per pack

© 2015 Ortho Technology, Inc. Creatures of the Sea is a registered trademark of Ortho Technology.

















101

























Superior stretch and rebound

- Non-Latex*
- Dark blue color
- Made in the USA



1,000 per pack

Elast-o-loop Separators

- Smooth injection molded plastic
- Non-Latex*
- Non-Radiopaque
- Made in the USA

Item #: 480-301	
1 pack	\$24.50/pk
3 or more	\$22.95/pk
960 per pack	

Elast-O-Loop2 is a trademark of Ortho Technology *Not made from natural rubber latex.



Radiopaque Separators

- Smooth injection molded plastic
- Non-Latex*
- Made in the USA

Posterior	
Blue	480-302
White	480-304
Anterior	
White	480-305
	\$19.50/pk
3 or more	\$17.65/pk
1,000 per pack	



Rotation Wedges

Preformed tabs quickly snap over bracket tie wings for rotating individual teeth.

Gray	577-303
Clear	
1 pack	
3 or more	\$19.95/pk
100 ner nack	• •



ELASTOMERICIPROD





Elastic Thread

- Excellent stretch and rebound
- Stain resistant
- Non-Latex*
- Made in the USA

Elastic Thread - Solid	
.020" Gray	490-099
.020" Clear	
1 spool	\$10.95 ea.
3 or more	
25' spools	
*Not made from natural rubber latex.	

Elastic Thread - Solid	
.025" Gray	490-200
.025" Clear	490-201
.030" Gray	490-202
.030" Clear	490-203
Elastic Thread - Hollow	
.025" Gray	490-100
.025" Clear	
.030" Gray	490-102
.030" Clear	490-103
1 spool 3 or more 50' spools	

Bumper Sleeve

- Excellent stretch and rebound
- Stain resistant
- Non-Latex*
- Made in the USA

• Made III the USA	
Lip Bumper Sleeve	
.018" Gray	
.018" Clear	490-401
.027" Gray	490-402
.027" Clear	
1 spool 3 or more 10' spools	



We Give Back

Ortho Technology is a proud Founder Level Partner of the American Association of Orthodontists Foundation.

Ortho Technology is committed to making a positive difference in the future of the orthodontic community.

Ortho Technology has always been deeply committed to the orthodontic industry and is pleased to be able to give back to the profession by partnering with the AAOF.





BONDINGSUPPLIES

Resilience

Orthodontic Bonding Solutions







Light-Cure Orthodontic Adhesive

Resilience® Light-Activated Orthodontic Bonding Systems

Provides the optimum properties for bonding ceramic, metal, or plastic brackets using light-cure technology. The smooth, tacky composite resin delivers ideal handling characteristics which eliminates bracket flotation and allows for easy flash clean up. Light activation provides additional working time to ensure accurate bracket placement. Produces a strong light cured chemical bond which also debonds easily. Archwires can be placed immediately after fully curing the brackets. Contains fluoride to help protect against potential decalcification.

A. Light-Activated Syringe Ki

B. Light-Activated Capsule Kit

Separate Refill Components

Conditioner/Primer

An Independent Clinical Evaluation Concluded that Resilience® Light-Cure Adhesive is...

"Recommended due to minimal bond failures."

"Very easy system to use."

Rated "Excellent" in its clinical evaluation.



Porcelain Conditioner Item #: OT-PO	\$20.95 ea.
Metal Primer Item #: OT-ME	
C. Capsule Gun Item #: FC-SELE-004	.\$45.00 ea.

D. Light-Activated Introductory Kit

© 2015 Ortho Technology, Inc. Resilience is a registered trademark of Ortho Technology.

esilie

Int Resi





Light-Cure Flowable Composite





Resilience Low Viscosity Light-Cure Flowable Composite adhesive features a precision controlled syringe which prevents overflow and run-on. It ensures accuracy as it flows easily into difficult areas, but is reliable not to slump or run until light curing occurs. The disposable tip provides precise application, accessibility, and control resulting in excellent aesthetics. It is dependable, strong, and more resistant to shock and stress. It can be used in many ways such as bonding lingual retainers, Invisalign® attachments, occlusal build-ups for opening a posterior bite, indirect bonding, and much more! Contains fluoride to help protect against potential decalcification.

LV Flowable Composite Kit

Item #: FC-RESI-007
1 kit \$84.95 ea.
3 or more \$78.50 ea.
1x 1.6 gm syringe green etchant, 3x 2 gm
flowable composite syringes, 1x 3 mL
light cure enamel bonding sealant resin,
25x Flowthru™ 22 gauge syringe tips,
10x 2 ½ " honding brushes instructions

LV Flowable Composite Syringe

Separate Refill Components

25x Flowthru™ 22 gauge syringe tips





Braided Retainer Wire

TRY OUR...

See page 143

© 2015 Ortho Technology, Inc. Resilience is a registered trademark of Ortho Technology.

Resilience

Orthodontic Bonding Solutions







No-Mix Orthodontic Adhesive

Resilience® Primer-Activated Orthodontic Bonding Systems

Provides the optimum properties for bonding ceramic, metal, or plastic brackets using a primer activated no-mix bonding system. This is a chemically cured system which does not require pre-mixing. The thick tacky composite resin delivers ideal handling characteristics which eliminates bracket flotation and allows for easy flash clean up. Produces a strong chemical bond which also debonds easily. Archwires can be placed 5 minutes after bonding the last bracket. Contains fluoride to help protect against potential decalcification.

A. Primer-Activated Syringe Kit

Separate Refill Components

9 mL per bottle

© 2015 Ortho Technology, Inc. Resilience is a registered trademark of Ortho Technology.



B. Primer-Activated Capsule Kit

Conditioner/Primer

C. Capsule Gun

Item #: FC-SELE-004 1 gun**\$45.00 ea.**

D. Primer-Activated Syringe Introductory Kit

1x 3.5 gm adhesive syringe, 1x 3 mL primer, 5x 2 ½" bonding brushes, 1x instruction sheet

MADE USA



Orthodontic Bonding Solutions

- Green color provides highly visible contrast
- Gel stays where placed without overflowing
- Non-drying water soluble consistency
- Ideal medium viscosity



Resilience

12 gm Gel Etch Syringe

Item #: 506102 1 syringe......**\$22.95 ea.** 3 or more \$20.95 ea. 1x 12 gm gel syringe, and 20x plastic tips

Hazardous materials, special shipping required.

Gel Etch Needle Tips

Item #: 272-002 . 25x plastic Flowthru™ 22 gauge tips

Item #: 501009

2 or more.......**\$40.95 ea.**

4x 4 gm syringes with 20x plastic Flowthru™ 22 gauge tips

Hazardous materials, special







Item #: 506104 1 syringe............ \$14.50 ea. 5 or more......**\$13.50 ea.** 1x 4 gm gel syringe, and 5x plastic 22 gauge tips

Hazardous materials, special





1 bottle \$62.50 ea. 226 gm (8 oz) per bottle

Gel Etch Medium Bottle

Item #: 506106 1 bottle\$25.95 ea. 3 or more\$22.95 ea. 57 gm (2 oz) per bottle

Hazardous materials, special



Liquid Etch Large Bottle Item #: 506107

1 bottle\$25.95 ea. 3 or more \$22.95 ea. 29.6 mL (1fl oz) per bottle

Liquid Etch Small Bottle Item #: 503204

1 bottle \$14.95 ea. 3 or more \$13.95 ea. 7 mL (.23fl oz) per bottle

Hazardous materials, special shipping required.

BONDINGSUPPLIES

Resilience

Orthodontic Bonding Solutions

Light-Cure Band Cement





MADE USA

Ideal Viscosity for Posterior Molar Pads





Low Viscosity Gel Allows for Controlled Use! Resilience® Light-Cure Band Cement

Light-Cure Single Component Push Syringes

A light-activated band cement in a single paste system which provides a longer working time and a rapid setting time. Can be placed directly into bands. Excess flash is easy to clean up prior to the light curing and the cement is easy to see due to the blue color. Contains fluoride to reduce the chance of decalcification.

Light-Cure Band Cement Kit

Item #: FC-RESI-038	
1 kit	[§] 127.95 ea.
3 or more	\$119 . 95 ea.
6x 5 gm push syringes	

Refill Syringes

Item #. 1 C-KL31-041	
1 syringe \$28.95	ea
5 gm per syringe	

Multi-Cure Band Cement



© 2015 Ortho Technology, Inc. Resilience is a registered trademark of Ortho Technology.

Resilience® Glass Ionomer Paste/ Paste Band Cement

A multi-cure fluoride containing glass ionomer band cement in a two paste system which offers a dual curing capability. May be used with or without a curing light. There is approximately a 3 minute working time and 10 minute set time without use of curing light.

Multi-Cure Economy Kit

Item #: FC-GRAN-061
1 kit \$126.95 ea .
3 or more
4x 5 gm base and 4x 5 gm catalyst, 2x 50
sheet mixing pads $(2 \frac{3}{8}^{"} \times 2 \frac{3}{8}^{"})$, $10x$
spatulas. 1x instructions for use

BONDINGSUPPL

MADE USA

Resilience

Orthodontic Bonding Solutions







Resilience® Glass Ionomer Powder/Liquid Band Cement

A two part glass ionomer band cement system of powder and liquid. It contains fluoride which reduces the incidence of decalcification. Resilience chemically bonds to enamel without the use of etchants. It has a low film thickness and will not wash out under bands. The final setting time is 6 to 7 minutes.

© 2015 Ortho Technology, Inc. Resilience is a registered trademark of Ortho Technology.

A. Powder/Liquid Starter Kit Item #: FC-GRAN-002

1 kit**\$70.95 ea.** 3 or more......**\$65.50 ea.** 1x 50 gm powder and 1x 50 mL liquid, 1x 50 sheet mixing pad, 1x measuring scoop, 1x instructions for use

B. Powder/Liquid Economy Kit Item #: FC-GRAN-059

1x 100 gm powder and 2x 50 mL liquid, 1x 50 sheet mixing pad, 1x measuring scoop, 1x instructions for use

Separate Refill Components

Measuring Scoop Item #: FN-PLAS-087....... \$3.15 ea.






BONDINGSUPPLES

Discover The GC Fuji Difference Bonds Perfectly In A Wet Field!


GC Fuji ORTHO™ LC Bracket Adhesive

- Bond in a Wet Field
- No Etch Required
- Significant Fluoride Release and Recharge





- Clinically proven. Fuji ORTHO LC offers significant advantages over composite bonding systems and greatly reduces chairtime.
- Fuji ORTHO LC features low sensitivity to moisture and can be used in the presence of saliva.
- Fuji ORTHO LC bonds metal and porcelain brackets and metal bands to enamel.



GC Fuji ORT

"GC.

€ 0086

Refill Components

GC Fuji ORTHO LC

Powder Refill ltem #: 000029 1 bottle 40 gm per bottle	\$219.95 ea.
Liquid Refill Item #: 000030 1 bottle	\$94.50 ea.



GC Fuji ORTHO™ LC Capsules

Capsules Refill Package

Conditioner

Fuji Capsule Applier



GC Fuji ORTHO™ BAND Paste Pak

Glass Ionomer Orthodontic Band Cement

The Paste Pak Dispenser saves time and reduces waste. It delivers consistent mixes every time - no matter how little or large a mix is required. A blue shade glass ionomer with a 3 minute working and 3 minute final set time.

GC Fuji ORTHO BAND Paste Pak Starter Kit

Refill Components



GC Fuji ORTHO™ LC Band Adhesive Paste Pak Automix

No Mixing! No Measuring! A Perfect Mix Every Time!

- Light-cure for longer working time
- Moisture tolerant for wet field bonding
- Fluoride release aids in prevention of white spots
- Orange mixing tip aids in avoiding light interference
- Blue color makes the cement cleanup easier

GC Fuji ORTHO™ BAND LC

Paste Pak Automix Refill

GC Fuji ORTHO™ Automix Tips Refill



BONDINGSUPPLIES



Reduce Chair Time!

Reduce Bonding Time

The Cybird XD advanced curing light is designed particularly for the orthodontic professional in need of super-fast light curing and bonding speed to reduce chair time. It features two powerful curing modes to choose from for practice preference.

Wide Application Range

The Cybird XD LED produces a variable wavelength range of 430-490 nm, assuring full curability of a wide range of adhesives.

Functional Aesthetics

The ergonomically shaped curing light is well-balanced for ideal handling. Cybird XD's square body prevents roll-off at chairside tables and is shaped for a more natural grasp.

Clever Construction

Built from industrial grade lightweight aluminum. The turbo 8 mm fiber-optic light guide and sleek body design grants superb oral access for effortless usability.

Intuitive Keypad

Practitioners save time with simplicity: a shot button and a mode selector / timer button for various applications.

1 Year Limited Warranty



Light Shields

(Green Indicator Light)

The Traditional High Power Mode (1,500 mW/cm2) allows you to polymerize light-cured adhesives at 5, 10, 15, or 20 second intervals.

(Orange Indicator Light)

Easily switch to the Plasma Emulation Mode (2,700 mW/cm2) for quicker cure times of 1.5, 2, 2.5, or 3 second intervals and achieve chair time efficiencies and reduced orthodontic bonding times.

light dispersion to the gingiva.

When using the plasma emulation mode, it is recommended to use the light shields to direct light only to the appliance and enamel, and prevent any unnecessary

Composite Curing Test Block

Cybird XD Curing Light Plug Type/Voltage

1 1: 1- 4	1000 OF
United Kingdom/220	CXD-9000-UK
Europe/220	CXD-9000-EU
US/Domestic/110	CXD-9000-US

1x handpiece, 1x light guide (8 mm turbo), 1x eye shield, 2x light shields, 1x charger base, 1x power adapter, 1x power cord, 1x composite curing

Cybird Replacement Parts -

Battery Module

Item #: CXD-1200 141.75 ea.

Light Guide

11-8 mm (Turbo)	CXD-1400
11-11 mm	CXD-1401
11-11 mm, 90° Angle	CXD-1500
1 light guide	<u>•</u> 73.50 ea.

Light Shield

10 per pack

Eye Shield

Item #: CXD-1600......**18.38 ea.**

Cybird XD Charger Base

Item #: CXD-2000 173.25 ea.

Power Supply Adapter

US (DC6V, 2.5A)	CXM-9000-US
EU (DC6V, 2.5A)	CXM-9000-EU
UK (DC6V, 2.5A)	CXM-9000-EU
1 nower supply adapte	er ±83.00 ea.

© 2015 Ortho Technology, Inc.

G SUPP

IN SECONDS!



Inventor Dr. Roger Wolk

D.D.S., M.S.

Dr. Roger Wolk is the inventor and developer of the Mini-Mold[™] Technique. He has been a prominent Orthodontist in the Los Angeles area for over 30 years.

Dr. Wolk is a retired Diplomate of the American Board of Orthodontics and Member of the American Dental Association



3-STEP TECHNIQUE



FILL the selected Mini-Mold Tip with light-cure resin





PLACE the **Mini-Mold Tip** on desired tooth





LIGHT-CURE the Mini-Mold Tip for 5-20 seconds

With the Cybird XD Advanced LED Curing Light from Ortho Technology on page 112



Mini-Mold™ **Starter Kit**

- Safe for sensitive patients
- Saves chair time and money
- Invisible aesthetic attachments
- Fast 5-20 seconds cure time

Item #: 51399........ \$167.95/kit Kit includes: 27x Mini-Mold Tips 3 of each style, 2x tip handles, and detailed step-by-step instructions

Different Mini-Mold™ Treatment Tips to Choose From


Package includes: 8x or 10x Mini-Mold Tips (same type). 1x tip handle. and detailed step-by-step instructions. Cold sterilizable.



Mini-Mold™ **Posterior Molar Pad Tips**

Designed to open the bite from the posterior teeth. Item #: 51390 \$62.95/pk

10 per pack




Mini-Mold™ **Wire Bonder Tips**

Ideal for making direct bond lingual retainers. 4.5 mm #: 51393

3 mm #: 51397 §62.95/pk 10 per pack, wire sizes up to .020"







Mini-Mold™ **Button Tips**

Creates a functional low profile button quickly and easily.

Item #: 51394 \$62.95/pk 10 per pack



round wire up to .020". Item #: 51391 §62.95/pk 10 per pack



Mini-Mold™ **Habit-Control Tips**

The small, conical spike, ends thumb sucking and tongue thrusting.

Item #: 51395 §62.95/pk 10 per pack



Creates edgewise .022" slot brackets, ideal for minor tooth adjustments.

Item #: 51396 §62.95/pk 10 per pack



Mini-Mold™ **Bite-Ramp Tips**

Creates a bite-opening ramp on the lingual of the upper anterior teeth.

3 mm #: 51392......**\$62.95/pk** 8 per pack

5 mm #: 51389...... \$62.95/pk 8 per pack

© 2015 Ortho Technology, Inc. All rights reserved worldwide. Mini-Mold U.S. Patent Numbers 5,820,368 and 5,681,163.

Mold A Lateral

The Perfect Solution to Fix Peg Laterals During the Course of Treatment



From the Creator of the Mini-Mold™ Technique Dr. Roger Wolk

D.D.S., M.S.

Dr. Roger Wolk is the inventor and developer of the Mini-Mold[™] Technique. He has been a prominent Orthodontist in the Los Angeles area for over 30 years.

Dr. Wolk is a retired Diplomate of the American Board of Orthodontics and Member of the American Dental Association.

5 Steps - 5 Minutes Turns A Peg/Small Lateral Into A FULL-SIZE Lateral!



5-Step Technique



PREP

Etch and prime (seal) the tooth



LIGHT-CURE

Light cure 10 to 20 seconds depending on material and light intensity



TEST FIT

Select the tip size, fill with composite and test fit the Mold-A-Lateral until incisal edge aligns with centrals



REMOVE

Remove the tip of the Mold-A-Lateral



STEP 5

FINISH

Use a finishing bur for final contour adjustments



Kit Contains 8 REUSABLE Tips!

Mold-A-Lateral Bonding Kit

Item #: 51388

MADE USA

4x right lateral bonder tips, (2x large, 2x small), 4x left lateral bonder tips, (2x large, 2x small), 2x handles, 1x instructions



TRY OUR...

Resilience LV Flowable Composite Kit

Item #: FC-RESI-007

See page 105

© 2015 Ortho Technology Inc. Mold-A-Lateral is a trademark and Resilience is a registered trademark of Ortho Technology.



BONDINGSUPPLES



Fully Exposes Both Arches to Provide Maximum Access

- Direct and indirect bonding procedures
- Full field access
- Interproximal stripping

Pedo (Amber)......603-229

Adult (Amber)......603-230

Pedo (Clear).....602-229

Adult (Clear)......602-230


Autoclavable

Disposable

- Bite registration
- Photography
- Debonding
- Microetching

Advantages:

- Unique intraoral extensions retract buccinator muscles
- Provides unobstructed view and superior buccal access
- Maintains accessibility during all types of procedures
- Safe for autoclave and dry heat sterilization up to 350°F/177°C






High Heat Cheek Retractor

The easy to use, one size fits all Cheek Retractor that provides superior intraoral access. The unique buccal extensions retract powerful buccinator muscles for maximum access. Made of high quality medical grade materials that withstand repeated sterilization.

- Maximum patient comfort
- Easy tension control
- Boosts visibility to intraoral work areas
- Made in the USA

ltem #: 601-190	
1 retractor	ea
3 retractors	ea










Adult

Pedo

Extraoral Cheek Retractors

- Full view for bonding of anterior and posterior teeth
- Does not reflect light from flash
- Cold sterilizable
- Made in the USA

Pedo	ORT-300-021	1 pack \$26.45/pk
Teen	ORT-300-022	3 or more
Adult	ORT-300-020	2 per pack










115

Adult (Blue)......601-230 10 per pack





Intraoral Cheek Retractors

- The best retractor for bonding
- Unrestricted access to the entire mouth
- Available heat sterilizable, autoclavable and disposable
- Made in the USA

Disposable

Pedo (Green)601-220	1 pack\$20.95/pk
Teen (Blue)601-221	3 or more
Adult (Red)601-222	10 per pack

Heat Sterilizable (Amber) Pedo......601-200 Teen601-201 3 or more......\$24.25 ea. Adult601-202 Sterilizable up to 370°F/188°C

Autoclavable (Clear)	
Pedo601-210	1 retractor
Teen601-211	3 or more \$22.95 ea.
Adult601-212	Sterilizable up to 270°F/132°C

Great for Adults and Children During Bonding





Lip Ring Retractor



Help Protect Lips Interproximal Stri

- Sized for adults and children
- Made of a high quality medical grade plastic
- Works great for bonding and intraoral photography
- Cold sterilizable



During	
ipping	

Item #: 54683	
1 pack	\$21.50/pk
5 or more	\$19.50/pk
2 ner nack	



Disposable Tongue Blocks Cherry Scented!

- Enables patient to relax jaw
- Helps maintain a dry field
- Made in the USA

Heat Sterilizable

Small	600-211
Regular	600-212
1 pack	§41.25/pk
2 or more	\$38.25/pk
50 per pack	-

- Can be easily trimmed
- Increases patient comfort

11	Saliva Ejector	600-216
12	Adapter Tubes	600-217
k	1 pack	\$19.25/pk
k	3 or more	\$17.25/pk
	3 per pack, cold si	terilizable only



TRY OUR...

Photo Cheek Retractors



See pages 167-168




Nola Dry Field System

Achieve Total Bonding Success...

MADE USA

- Full Arch Access for Bonding
- Adjustable to Arch Length
- Maximum Patient Comfort
- Excellent Soft Tissue Retraction
- Fully Autoclavable to 280°F/135°C
- Reduces Bond Failure
- Prevents Treatment Delays





The Nola Dry Field System completely eliminates saliva and fully exposes both arches making it ideal for all bonding applications.

The system is easy to assemble, be inserted, and used by a single operator. The entire system is completely autoclavable making clean-up easy.

Dry Field Systems

	y i icia 3	ystems
Sm	all (Red)	300-400
Lar	ge (White)	300-401
1 s	ystem	\$59.75 ea.
3 c	r more	\$55.75 ea.

System Includes: 1x cheek retractor, 1x flexible connector, 1x tongue guard, 2x bilateral suction tubes, 1x central suction "Y" connector with low-volume adapter

Dry Field System Separate Components

(A) Cheek Retractor

Assembly with fle	xible connector
Small (Red)	300-414
Large (White)	300-406
1 retractor	\$26.10 ea.

Adjustable Cheek Retractor Keeps the cheeks back

Keeps the cheeks back and fully exposes both dental arches

(B) Tongue Guards

Small (Red)	300-413
Large (Translucent).	300-403
1 pack	49.50/pk
2 per pack	_

Tongue Guard

Holds the tongue in place and allows the patient to continue swallowing

© Bilateral Suction Tubes

Bilateral Silicone Tubing

Precisely directs saliva flow and ensures unrestricted evacuation for a totally dry environment

(P) "Y" Connectors

Item #: 300-408 1 pack......^{\$}16.95/pk 4 per pack

Central Suction

"Y" connector adapter (low-volume and high volume also available)

Low-Volume Adapters Item #: 300-409

1 pack......\$27.95/pk 4x Adapters and 4x "Y" connectors

High-Volume Adapters

Item #: 300-411 1 pack......^{\$}32.95/pk 4x Adapters and 4x "Y" connectors

G Connecting Arm Inserts

Used for connecting Tongue Guards to Cheek Retractor. Item #: 300-402

orthogap



Bonding Pad Organizer

MADE USA

6 1/4"W x 6 3/4"H x 4 3/4"D

Orthogap Bonding Pads Item #: 600-370


Disposable Bonding Pads

Combination mixing pad and bonding pad, saves your office time and money!

- Adhesive backing holds brackets firmly in place
- Liquid impermeable surface for mixing
- Easy bracket access speeds selection
- Padded for easy dispensing



U.S. Patent Number 5,542,844





Plastic Bonding Trays

- With or without mixing wells
- Adhesive holds brackets in place
- Holds one full bracket set-up

With Mixing Wells Item #: 600-350 Without Mixing Wells



Mixing Pads

- Liquid impermeable surface
- Smooth finished surface

2.5" x 2.5" Mixing Pads
Item #: 600-600
1 pack...........\$8.55/pk
3 or more.......\$7.75/pk
3x 50 sheet pads per pack

3" x 5" Mixing Pads
Item #: 600-610
1 pack
3 or more
3x 50 sheet pads per pack



Tube Series Regular Tip

iteiii #. 000-433	
1 pack	.\$36.95/pk
3 or more	.\$33.50/pk

4 trays of 100 applicators per pack

Microbrush

- Precise placement of solutions
- Reduces use of excess solutions
- Bends easily into any angle
- Non-linting and non-absorbent tips
- Color coding avoids application mix-up
- Made in the USA



1 pack\$39.95/pk

3 or more\$36.95/pk 4 trays of 100 applicators per pack

Microbrush

- Precise placement of solutions
- Reduces use of excess solutions
- Bends easily into any angle
- Non-linting and non-absorbent tips
- Color coding avoids application mix-up
- Made in the USA

Dispenser

Item #: MPD	
1 dispenser ^{\$} 13.25	ea.
4 or more ^{\$} 10.95	ea



1 pack	\$39.95/pk
3 or more	\$36.95/pk
4 trays of 100 applica	ators per pack

UltraBrus

1 kit.....**\$22.75 ea.**

1 dispenser, 1 refill cartridge of

Dispenser Kit

100 applicators

Item #: U1D

Microbrush

- Reduces use of excess solutions
- Bends easily into any angle
- Non-linting and non-absorbent tips
- Color coding avoids application mix-up
- Made in the USA

Dispenser

Item #: MPD	
1 dispenser \$13.25	ea.
4 or more	ea.

Precise placement of solutions

Bonding Brushes

Injection molded plastic handle

Prophy Paste

Enhances bond strength

Cleans teeth prior to bonding

• Great time saver for rebonds

1x handle, 20x disposable cups

Autoclavable Kit: 600-240

No fluoride

No flavored oils

Made in the USA

Specially Formulated for Orthodontists

Prophy Handle and Cups

Preps single tooth for bondingMade in the USA

Item #: 15486

1x finger ring

1 canister......**\$31.50 ea.**

3 or more......**\$28.95 ea.**

200x 2 gm single-dose cups and

- 1/4" stem with 1/2" tuft
- Use for sealant or primers
- Made in the USA

Bonding Brush Tip Holder

Disposable Kit: 600-250

1 kit (red)......**\$26.95/pk**


25x disposable handles with cups

Bonding Brush Tips 3/4" FN-PLAS-048 \$15.75/pk 100 per pack





JltraBrüsh. 1.0



- Precise placement of solutions
- Reduces use of excess solutions
- Bends easily into any angle
- Non-linting and non-absorbent tips
- Color coding avoids application mix-up
- Made in the USA

UltraBrush 1.0 Refill

Item #: U1R200	
1 pack \$34.	50/pk
3 or more \$31.	50/pk
2 trays of 100 applicators p	er pack

Bonding Brushes

- Large economy pack
- 2" disposable bonding brush
- Made in the USA

Item #: 610-100	
1 pack	.\$31.50/p
5 or more	.\$27.95/p
400 per pack	

BONDINGSUPPLESI


Microetcher™ Intraoral Sandblaster

Benefits

- Increases bond strength up to 400%
- Stainless construction and replaceable carbide tip
- Nozzle rotates 360° for easy access

Uses

- Etch all metals, composites and amalgam
- Perfect for adhesive removal for rebonds or for etching the inside of bands
- Surface abrasion procedures in seconds







AFTER Microetcher II

Microetcher IIA Fully Autoclavable Intraoral Sandblaster

Item #: 22005

1x handpiece, 2x abrasive reservoirs, 1x air line, and 1x user guide

Autoclavable up to 269°F/132°C Plastic jar and tubing must be removed and abrasive jar filter removal recommended prior to sterilization







Microetcher II Intraoral Sandblaster

Item #: 22000

1x 60° autoclavable nozzle, 2x abrasive reservoirs,

1x air line, and 1x user guide



Microetcher II

Microetcher Attachment Components



Quick Disconnect Kit

2x tube connectors



Female Disconnect Kit



Male Disconnect

Item #: 44020**\$26.25 ea.**



Handpiece Adapter

Item #: 44075**\$68.75 ea.**

BONDINGSUPPLES



Microcab Plus Dust Containment Unit

- A compact design with large internal working area and hand guards
- Bright, fluorescent light for ideal illumination
- Easy to clean with a pullout filter
- Snap-in disposable windows
- Integrated handle for portability
- Ultra quiet fan and filter
- Made in the USA

Microcab Plus

United States 110V: 91010 1 microcab\$435.00 ea.

Includes: 4x windows, 1x light bulb, 1x power cord and 1x instructions Size: 6 1/2" x 12" x 6 1/2"

International 220VAC: 91354 1 microcab \$445.00 ea.

Which Window Shield Is Right For Your Unit?





Order 2nd Generation Item #: 93817

Replacement Parts

Disposable Windo	ow Shields
1st generation	91286
2nd generation	93817
1 pack	\$41.95/pk
10 per pack	
Filter	

Item #: 49055.......\$26.25 ea. Hand Guard Kit

110 volt.....109150 220 volt......109130 Wall Mount (Optional) Item #: 91518..........\$49.25 ea. Item #: 91350.........\$38.50 ea.

Item #: 59090.........\$26.50 ea.

Light Bulb - 9 watt



TruEtch[™] **Aluminum Oxide**

50 micron12300 Intraoral applications, white in color

90 micron12400 Extraoral applications, tan in color

Made in the USA

1 bottle**\$18.95 ea.** 11b (454 gm) per bottle

© 2015 Ortho Technology, Inc. TruEtch is a trademark of Ortho Technology.



■ Ends air/water line contamination

- Decreases bond failure
- Autoclavable nozzle
- Designed for orthodontists
- Detachable lingual elbow

Item #: 890-100 1 tooth dryer \$105.75 ea. 2 or more......**\$99.95 ea.** 110 Volt

Replacement Parts


Nozzle w/Elbow Tip Item #: 890-101 \$19.50 ea.







HEADGEARIPRODUCTS



Multi-Adjustable Facemask®

Discover The Ideal Choice to Achieve Maximum **Class III Compliance**

- Free-flow air vents
- Full range of motion option
- Maximum patient comfort
- Fully adjustable
- Light-weight
- Screws designed for increased durability



Exclusive Patient Comfort Design

Quick and Easy Fit


- 1. Position bend in vertical bar over nose
- 2. Center forehead rest on forehead and lock in position
- 3. Center chin cup on chin and lock in position
- 4. Place crossbar approximately 2.0 mm below the upper lip and lock in position

FEATURE





Chin cup design features air vents for maximum comfort and air circulation.



HEADGEARIPRODUCTS



Pride

Item #: 550-406



Canadian

Item #: 550-411

Item #: 550-408

Pride



Green Camo Item #: 550-413

Hawaiian Hibiscus Item #: 550-407





HEADGEAR PRODUCTS

Union Jack Item #: 550-410

Wild Cat Print Item #: 550-409















Multi-Adjustable	Facemask
1 facomack	SEC EO

1 facemask \$56.50 ea
5 facemasks \$52.50 ea
10 or more
Includes: 1x hexagon key and 2
sets of chin and forehead nads

See page 101 for Extraoral Elastics

© 2015 Ortho Technology, Inc. Multi-Adjustable Facemask is a registered trademark of Ortho Technology. U.S. patent pending.

Replacement Pad and Key Kit
American Pride406-RK2
Canadian Pride411-RK2
Green Camo413-RK2
Hawaiian Hibiscus407-RK2
Pink Camo412-RK2
Racin' Flames408-RK2
Union Jack410-RK2
Wild Cat Print409-RK2
Solid Colors550-RK2
1 kit \$12.49 ea.
1x hexagon key and 2 sets of
chin and forehead pads



HEADGEARIPRODUCTS







Large

Small

Reverse Pull Facemask

Fully Adjustable Comfort Design Aids in Patient Compliance

- Maximum patient comfort
- Treats maxillary insufficiencies
- Fully adjustable
- Stainless steel and plastic
- Leather padded contact points
- Hexagon key included
- Storage bag included

Small	REV-PULS
	\$63.95 ea.
	1x storaae baa

See page 101 for Extraoral Elastics



Replacement Pad Kit Item #: 100002 1 kit.....**\$36.75 ea.** 1x forehead pad and 1x chin pad



Lip Bumpers


- Non-toxic plastic coating allows for maximum labial pressure
- Smooth finish for added patient comfort
- Made of .045" wire

A. Lip Bumper without Loops

Item #: M3050-00 1 lip bumper \$19.95 ea.

B. Lip Bumpers with **Loops and Bayonet Stop**

With removable, numbered color identification system. Measurements are from the distal side of the stop to the distal side of the opposite stop. #1 76.0 mmM3051-01 #2 80.0 mmM3051-02 #3 84.0 mmM3051-03 #4 88.0 mmM3051-04 #5 92.0 mmM3051-05 1 lip bumper......**\$20.95 ea.**



Safety Release Modules

- Flexible tear-proof strap
- Durable steel clamp design
- Attaches to any neckpad
- Safety release

A. Steel C-Clamp

600 gm Medium......550-110 750 gm Heavy......550-111

1 pack......**\$41.95/pk** 4 or more...........\$37.95/pk 10 per pack

- Steel C-Clamp has large facebow holes
- Plastic C-Clamp is brightly colored

B. Plastic C-Clamp

600 gm Medium550-100 750 gm Heavy550-101

10 per pack

HEADGEARIPRODUCIS



NEW!

Cervical Neckpads

- Poly/cotton blend
- Sturdy washable design
- One hook-up tab
- Bright, fun and stylish colors
- Increase patient compliance
- Can be used with any style safety module

*Royal Blue, Navy Blue, Royal Blue/Denim, Red, Black







High Pull Head Caps

- Woven poly/cotton blend
- Sturdy washable design
- Two hook-up tabs
- Use with any style module
- *Royal Blue, Navy Blue, Royal Blue/Denim, Red, Black

Royal Blue	550-750
Assorted Colors*	
1 pack	\$49.95/pk
5 per pack	





Ortho Performance® Mesh Bags

- Headgear storage
- One bag with multiple uses
- Can be used to hold headgear items, multi-adjustable facemasks or oral hygiene kits
- Solid blue in color with a secure zipper closure
- Made from premium micro mesh material

Royal Blue	560111
1 pack	.\$45.95/pk
3 packs	\$40.95/pk
10 per pack	-
Measures 8 3/4" x 5"	

© 2015 Ortho Technology, Inc. Ortho Performance is a registered trademark of Ortho Technology.

ORTHO TECHNOLOGY®

Full Form Extraoral Facebows

High Quality Material

High quality stainless steel is used in the manufacturing process of our facebows.

Wide Inner Bow Size Measurement

For accurate selection of the proper size, measurements are based on center of loop to center of opposite loop of the inner bow.


#1	Yellow	83.0 mm431-111
#2	Red	90.0 mm431-112
#3	Pink	97.0 mm431-113
#4	Blue	104.0 mm431-114
#5	Green	111.0 mm431-115

White-Universal (inner bow without loops) Inner Bow Shapes

Full Form431-110 Standard441-110

Strong, Durable **Laser Welded Joint**

Highly resistant to breakage and deformation during length of treatment



for easy attachment of safety release module








Short Outer Arch Facebows

- Laser welded joint for unsurpassed strength
- One piece wire design of inner bow eliminates breakage

1 facebow\$25.95	ea.
10 or more\$23.95	ea.

#1 83.0 mmM4835-01 #2 90.0 mmM4835-02 #3 97.0 mmM4835-03 #4 104.0 mmM4835-04 #5 111.0 mmM4835-05



Colored Facebows Universal Standard Form

- Bright colors help increase patient compliance
- Coordinate with neckpads
- Coated with a durable nylon
- Resistant to peeling, fading and cracking

White	550-310
Blue	550-311
Purple	550-312
Lime Green	550-319

Standard form .045" inner arch

ORTHO TECHNOLOGY®

Stainless Steel Impression Trays Edge Lock Design For Superior Impressions



U1 X-Large Upper Item #: 1170-01



U2 Large Upper Item #: 1170-02



U3 Medium Upper Item #: 1170-03



U4 Small Upper Item #: 1170-04



U5 X-Small Upper Item #: 1170-05

Upper



L1 X-Large Lower Item #: 1170-11



L2 Large Lower Item #: 1170-12



L3 Medium Lower Item #: 1170-13



L4 Small Lower Item #: 1170-14



L5 X-Small Lower Item #: 1170-15

Lower

Perforated Design

- Perforated to ensure secure alginate retention
- Edge-locking
- Resists tarnishing

- Surgical grade stainless steel
- Withstands repeated sterilization
- Heat sterilizable up to 370°F/188°C

Impression Trays

1 tray \$12.95	
10 trays \$10.95	ea
20 or more \$8.95	ea





Laminated Bite Wafers


- Aluminum foil between two layers of soft wax
- Adapts to patient's arch

Item #: 600-024	
1 pack	\$86.95/pk
3 or more	\$83.95/pk
84 per pack	

DuraLock PLUS







MADE USA

Additional Benefit

Ortho Technology manufactures the DuraLock PLUS Impression Tray to be either disposable or cold sterilizable.

DuraLock® PLUS Introductory Kit

DuraLock® PLUS Individual Packs

1 pack	\$43.50/pk
5 packs	\$40.50/pk
10 or more	\$37.95/pk
50 per pack (upper or lower)	

© 2015 Ortho Technology, Inc. DuraLock is a registered trademark of Ortho Technology. U.S. Patent Numbers D565185 and D571473.








DuraLock PLUS
IMPRESSION TRAYS

Not Your Ordinary Impression Tray!





No Adhesive Needed

Using DuraLock PLUS trays save time. By eliminating the extra steps of coating the tray with adhesive and waiting for it to dry, the impression process is simplified and accelerated.

Largest Variety of Tray Sizes

DuraLock PLUS trays eliminate inadequate coverage and/or tissue impingement and ensure an accurate fit every time by offering more size choices than competing trays.

Stability and Durability

The reliable strength and stability of DuraLock PLUS trays eliminate impression distortion and/or impression materials pulling away from the tray sides.



(4) Medium

Color Code: Green Upper......501-004U Lower.....501-004L

Color Code: Yellow

Upper501-005U Lower.....501-005L



(6) Extra Large

Color Code: Purple Upper.....501-006U Lower501-006L





129

PRESSION SUPPLIES

Orthoprint®

Fast Set Alginate

orthoprint

EXTRA FAST SETTING

Zhermack



Orthoprint





Mixing Time **30** Seconds



Working Time 1 Minute **5** Seconds



Time in Mouth 45 Seconds

Total Set Time 1 Minute 50 Seconds

- Extra-rapid working and setting times
- Total setting time: 1 minute 50 seconds
- Vanilla flavor works with alginate flavors too
- Modified atmosphere packaging ensures product stability - 5 year guarantee
- Thixotropic, non-drip, and highly elastic
- Dust free
- Gluten free

1.1 lb (500 gm) Bag


ltem #: 302145	
1 bag	^{\$} 14.50 ea.
12 bags	
24 or more	\$12.50 ea.
12 haas per case	

1.1 lb (500 gm) Canister

Item #: 302160	
1 canister	.\$17.50 ea.
1x 1.1 lb bag of alginate in	canister,
1x scoop and 1x vial	
10 canisters per case	

Measuring Set

Item #: 105182	
1 set	\$3.00 ea.
1x scoop and 1x vial	







Violet Mixing Time **45** Seconds



Pink Working Time **60** Seconds



White Total Set Time **3** Minutes

Kromopan® Chromatic Alginate



Allows casting of plaster up to 100 hours after the setting of the impression. Gluten free.

1 lb (454 gm) Bag		
Item #: KRO-100-P		
1 bag\$	15.25	ea.
10 or more\$	14.25	ea.
20 haas ner case		

ground shipment recommended when ordering 20 or more bags/canisters

1 lb (454 gm) Canister

Item #: KRO-100-C 1 canister......\$22.50 ea. 1x powder scoop and 1x water vial

Scoop and Vial

Item #: KRO-100-SV



Alghamix II Mixing Machine

Two-speed compact mechanical mixing unit for mixing alginates, silicones, plasters and gypsums. Controlled mixing due to constant speed and reduces mixing time by 30%. Improved ergonomic bowl position allows for optimal operating conditions.

2 Year Manufacturer Warranty

Mixing Machine

Item #: DS-C303215	
1 machine	\$499.95 ea.
1x alginate mixer, 1x mixing bow and instructions	l, 1x spatula



Fast Setting PVS Impression Material

Perfect for Clear Aligners



EXAFAST™ NDS Monophase

- Single phase technique with one viscosity
- Easy to syringe with enough compression to obtain accurate impressions

Cartridge Dispenser

■ For all NDS impression material Item #: 152101 «

1 dispenser \$170.95 ea. 3 or more \$166.95 ea. (Shown on page 132)

Monophase Material

Item #: 137407
1 pack\$37.50 ea.
3 or more ^{\$} 35.50 ea.
2x 48 mL cartridges (96 mL total)
6x mixing tips and instructions

Mixing Tips Type 2LL - Purple Item #: 153048«

1 pack	\$66.95/pk
	\$63.00/pk
48 per pack, r	nono/heavy body



EXAFAST™ NDS Heavy Body

- Non-slumping, high compression tray material used with double phase impressions
- Stabilizes the tray and prevents distortion
- Pushes the wash material into the gingival margins for accuracy

Heavy Body Material	Mixing Tips
Item #: 137207	Type 2LL - Purple
1 pack	Item #: 153048« 1 pack

"Available by special order. Additional shipping time may be required.



EXAFAST™ NDS Regular Wash

- Slightly thicker wash material used with double phase impressions
- Syringeable wash that is a little more thixotropic and less runny than injection material
- Stays in place until pressure is applied to capture critical detail

Regular Wash Material

Item #: 13/00/
1 pack \$37.50 ea.
3 or more \$35.50 ea.
2x 48 mL cartridges (96 mL total)
6x mixing tips and instructions

Mixing Tips Type 2L - Aqua Item #: 153248*

1 pack	\$66.95/pl
3 or more	\$63.00/pl
48 per pack, regul	lar material

Intraoral Tips -For Wash Material



EXAFAST™ NDS Injection Wash

- Syringeable wash material for double phase impressions
- Thinner wash that is a less thixotropic and more flowable
- Stays in place as pressure is applied to capture critical detail

Injection Type Material

Item #: 137107

1 pack	\$37.50 ea.
3 or more	
2x 48 mL cartridge	s (96 mL total)
6x mixing tips and	

Mixing Tips Type 2S - Pink Item #: 153148«

Intraoral Tips -For Wash Material

Item #: 150100«

1 pack	\$59.95/pk
3 or more	
100 per pack, in	traoral material
(Shown above)	

IMPRESSION SUPPLIES

PVS Bite Registration Material with Outstanding Accuracy







GC America Cartridge Dispenser Type 2 NDS

Item #: 152101« 1 dispenser..... \$170.95 ea.

Item #: 153248«

3 or more \$166.95 ea.

Mixing Tips Type 2L - Aqua

EXABITE™ II NDS Registration

- Bite registration material is easy to apply and will not slump or run
- Ideal cream-like consistency allows patients to occlude without deviation or distortion
- Sets quickly in only 45 seconds and is easy to trim if necessary

Bite Registration Material

12x mixing tips and instructions

Item #: 132412		
1 pack \$49	.95	ea.
3 or more\$47	.95	ea.
2x 48 mL cartridges (96 mL total)		

«Available by special order, Additional shippina time

3 or more.......563.00/pk 48 per pack

may be required.



EXAFAST™ Putty

EXAFAST Putty is a vinyl polysiloxane material, consisting of a base and a catalyst that are easy to manipulate for accurate, void-free impressions. EXAFAST Putty produces crisp, extraordinarily smooth, precise impressions and can be used with either regular or injection wash material.

- The highest compression tray material used with double phase impressions
- Extremely smooth and easy to hand mix. Avoid contact with latex, oils and moisture when mixing
- Sets after only 21/4 minutes in the mouth
- 45 second working time

Easy to Mix & Use





1:1 Base and Catalyst Ratio





Hand Mix and Place in **Impression Tray**

Putty Standard Pack

Item #: 137307
1 pack \$97.50 ea.
3 or more \$95.50 ea.
1x jar putty base (287 mL), 1x jar putty
catalyst (287 mL), 2x scoops and instructions





Encourages patient cooperation while reducing gagging, 100% Gluten Free

☐ Very Very Strawberry	250-101
☐ Marvelous Mint	250-102
☐ Chillin' Cherry	250-103
☐ Bubble Gum	250-104
☐ Groovy Grape	250-105
☐ Frosty Rootbeer	250-106
☐ Rowdy Raspberry	250-107
☐ Alotta Piña Colada	250-108
☐ Cosmic Orange	251-101
☐ Ballistic Blueberry	251-103

☐ Classic Cola	251-104
☐ Bodacious Banana	251-106
☐ Fantastic Fruit Punch	251-108
☐ Rich Milk Chocolate	251-109
☐ Wet n' Wild Watermelon.	251-110
☐ Cotton Candy	250-113
☐ Cookie Dough	251-111
☐ Orange Dreamsicle	251-114
□ French Vanilla	251-115

Alginate Flavorings

1 bottle \$9.95	ea.
10 bottles \$9.45	ea.
20 or more\$8.95	ea.
2 fl oz per hottle	

Alginate Flavoring Organizer

itterii ii. O i i ti c	
1 organizer	.\$37.99 ea.
3 or more	
10 5/8"W x 6 5/8"H x	3"D

alogy Inc. Tactor Croat is a tradamark of Outho Tachmalog



© 2015 Ortho Technology, Inc. Tastes Great is a trademark of Ortho Technology.

wonderfill

the EZ in, EZ out

tongue filler



Remove a palm full of non-toxic Wonderfill and form a domed bridge in the tongue area.

With Wonderfill in the tongue area, pour plaster/stone.

Easy to Use Time Saver!

No Mixing. No Mess. No Wasted Time.

Saves Time!

Produces a sleeker, more presentation-worthy model!

PLASTER/STONE Residue when no tongue filler used

- Tedious, time consuming grinding out of build up
- Messy
- Coarse model appearance

20 Up to Faster!

STONE/PLASTER
ALGINATE
Wonderfill

Let plaster/ stone fully set. Then pop out non-sticky Wonderfill for a clean model, using only water & brush to remove any residue.

Wonderfill™ Tongue Filler

ground shipment recommended when ordering by the case

ALGINATE As the tongue filler

- Mixing
- Messy
- · Hard to form
- Set-up and clean-up time required
- Often needs carve-out of residue

Ortho Model[™] Archival System



Top Loading StyleUse for Archive Cases

Save Valuable Office Space!

Each Set Holds Up to 320 Models

- Inexpensive storage
- Easy to assemble
- Super strong construction
- Made in the USA



Front Loading Style
Use for Active Cases

Ideal Storage

Ortho Technology's Ortho Model™ Archival System offers offices an innovative and simple solution for archiving patient models, freeing up valuable office space.

4 Boxes Per Set

Each set includes:

4 boxes with 4 trays each, for a total of 16 trays. Each tray has 20 model compartments, for a total of 320 model compartments per set.

Top Loading


Item #: 37810

Dimensions: 15"H x 191/2"W x 131/4"D

Front Loading

Item #: 37809

1 set of 4 large boxes (16 trays) Dimensions: 16"H x 183/4"W x 13"D



© 2015 Ortho Technology, Inc. Ortho Model is a trademark of Ortho Technology.



BIG ADVANTAGES! SMALL FOOTPRINT!

Everything You Need in a Vacuum Forming Machine



Vacuum Forming Machine

The Essix® SelectVac™ is an improvement to standard vacuum forming machines in almost every aspect. It's lighter, has a smaller footprint, and has up to 50% more vacuum power than standard units. This unit saves time and is safer thanks to no pre-heating, an auto-vacuum start that turns on only when you need it which means no unnecessary noise, an auto-off two-sided heating system with rotating plate support, and a digital timer. The carbon fiber heating element evenly distributes heat for a more precise finished product, and the model cup and beads allow for embedding tall or untrimmed models. Unlike many standard vacuum machines, this unit has the ability to be used in both Europe and America.

- Vacuum power: 6 (inches Hg)
- No pre-heat required
- Automatic vacuum start & auto-off heating system
- · Rotating plate for two-sided heating dual laminates or thicker materials
- Accepts round or square plastics
- Includes glass beads for adjusting model height
- 2 year warranty

Vacuum Machine

US/Domestic/110VAC....ESVM01 Europe/220VACESVM02 $11\frac{1}{2}$ "H x 9"D x 7"W, (9.2 lbs) Includes: 1x cast rack (plane/cup), glass beads, 1x power cord, 1x 4 mm Allen wrench and user manual

Replacement Glass Beads

Item #: AGRA0067 8.82 oz (250 gm)

Rotating Heating Head

Rotates 180° for ready access to models, the model cup and retainer materials



Rotating Material Tray Rotates 180° for heating dual laminates and thicker mouthguard materials



Reversible Cast Rack (Plane or Cup) for Versatility Accommodate taller models or use the included beads for

working with untrimmed models



Glass Beads

Reduce excess trimming and allow easier release of models - great for use with untrimmed or taller models



Carbon Fiber Heating Element Even heat distribution and no pre-heating required







LABSUPPLIES

Clear Advantage™ Series I Clear Retainer Material

Thermal Forming Splint/Copolyester

Clear Advantage Series I is brilliantly clear and offers excellent conformity to the model for an exact fit. The material is durable and bonds with acrylic for creating modified appliances. Forms clear for a true "invisible" retainer. Sheets are sealed with a protective peel away film on both sides.

.020" (0.5 mm)7500-100 \$32.95	/pk
.030" (0.7 mm)7500-105 \$32.95	/pk
.040" (1.0 mm)7500-110\$32.95	/pk
5" x 5", 50 sheets ner nack	

.060" (1.5 mm)......7500-035\$31.50/pk .080" (2.0 mm)......7500-040......\$38.50/pk 5" x 5", 25 sheets per pack

• Invisible retainers

by Raintree Essix

Clear Advantage Series 1 is a Moisture-Sensitive Product - Vacuum sealed to maintain product integrity.



Clear Advantage™ Series II Durable Retainer Material


Thermal Forming Coping/Polypropylene

Clear Advantage Series II is our most durable material and it will not shrink, crack or break. This material is slightly opaque in its transparency and offers excellent formability, flexibility and long term wear characteristics. The increased strength of this material is ideal for making retainers for bruxism patients. Not bondable with acrylic.

- Extended life invisible retainers
- Clear aligners
- Splints and bite planes
- Temporary bridges with pontics
- Indirect bonding trays
- Non-Latex*

• Non-Latex

- BPA free
- Made in the USA



Compare To

Essix C+® Plastic by Raintree Essix

Sports Advantage™ Mouthguard Material

Thermal Forming Soft EVA

Our Sports Advantage material is a specially formulated • Soft splints EVA plastic that is soft and flexible. This material heats easily and uniformly and conforms perfectly to the model. Custom fit mouthguards allow for excellent oxygen intake and a perfect fit for increased compliance.

.040" (1.0 mm)7500-140\$31.95/pk .080" (2.0 mm)7500-165 32.95/pk .120" (3.0 mm)7500-170 \$42.95/pk .160" (3.8 mm)7500-180 \$42.95/pk

5" x 5", 25 sheets per pack

mouthguards

• Nightguards

Model duplication

 Non-Latex* BPA free

• Made in

Compare To Essix Ace® Plastic

by Raintree Essix

© 2015 Ortho Technology, Inc. Sports Advantage and Clear Advantage are trademarks of Ortho Technology. All other trademarks are the property of their respective owners. *Not made from natural rubber latex.















See page 149 See page 148

Cleans and Brightens Removable Dental Appliances!



Retainer Brite cleans Essix and Hawley retainers - plus nightguards, mouthguards, splints, and most other removable dental appliances.

1 box	.\$11.95/box
10 boxes	\$11.25/box
24 or more	.\$10.25/box
96 tablets ner hox	

Precise adjustable flame

Up to 2 hour burn time Butane sold separately

 Temperatures up to 2450°F/1343°C

- Helps remove plaque and tarter
- Kills germs and bacteria
- Cleans and brightens
- Keeps appliances minty fresh
- Made in the USA

Item #: 18160	
1 box	.\$6.95/box
10 boxes	.56.25/box
24 or more	.\$5.75/box
36 tablets per box	



Micro Torch Item #: MT-770



LABSUPPLES

Universal All-In-One Solder/Flux

Perfectly premixed ratio of silver solder and flux

Universal All-In-One Solder/Flux is ready to use straight out of the package. The perfectly premixed ratio of silver solder and flux in convenient single lengths ensures less mess and it's easy to use. Additionally, it ensures less breakage because you never risk using too much or too little flux.

- Ensures less breakage
- Easy to work with
- Less mess

Item #: 2000	
1 pack	\$56.95/pk
3 or more	\$53.95/pk
.039" dia., 6x 6"	lengths per pack



Silver Solder

- Cadmium free Dead soft
- Smooth flowing
- Melting point 1220°F/660°C
- Made in the USA

ltem #: 2020	
1 spool	\$36.25 ea.
4 or more	.\$32 . 95 ea.
.025" dia., 5DWT, 7	' per spool



TruFlow™Flux

- Smooth, quick flowing
- Consistent white color
- Made in the USA


Item #: 2010	
1 jar	[§] 16.49 ea.
3 or more	^{\$} 14.95 ea.
2 oz (56.7 gm)	per jar

© 2015 Ortho Technology, Inc. TruFlow is a trademark of Ortho Technology.

LABSUPPLES

Help Prevent Anterior Teeth From Flaring!





Size	Item #
28.5 mm	550-001
31.5 mm	550-002
34.5 mm	550-003
37.5 mm	550-004

Size	Item #
40.5 mm	550-005
43.5 mm	550-006
46.5 mm	550-007
49.5 mm	550-008

TruForce™ Stainless Steel Flat Bow Preformed Retainer Wire

- Greater stability and axial control of anteriors
- Allows for detailed contouring of wire to the labial surfaces
- Maximum tooth control not available with round wire
- Improves retainer aesthetics
- Made in the USA

TruForce Flat Bows from Ortho Technology are a single retainer wire with two distinct shapes, rectangular over the anterior teeth and round in the posterior segments. The flat rectangular segment helps prevent anterior teeth from flaring and prevents the labial bow from sliding toward the gingiva, while the round posterior wire anchors the appliance. The preformed cuspid loops come in eight sizes and are a great time saver for the lab.




TruForce™ Stainless Steel Retainer Wire

Made in the USA

.022"	900-105«
.025"	900-106«
.026"	900-107«
.028"	900-108«
.030"	900-109«
.032"	
.036"	
.040"	
.045"	900-113«
1 pack	
4 or more	\$11.75/pk

«Available by special order. Additional shipping time may be required.

14" lengths, 25 per pack



TruForce™ Ball Clasps

Spring tempered stainless steel provides excellent retention of removable appliances without irritating tissue.

Made in the USA

028"

.032"	3000-032
.036"	3000-036
1 pack	\$29.95/pk
1 pack 5 or more	

3000-038



TruForce™ Adams Clasps

- Strong interproximal hold without slipping
- Consistent uniformity
- Spring tempered stainless steel
- Made in the USA


7.0 mm	(0.28")	10017
8.0 mm	(0.31")	10018
9.0 mm	(0.35")	10019
10.0 mm	(0.39")	10110
11.0 mm	(0.43")	10111

1 pack......\$13.95/pk 10 or more.....\$11.95/pk .028" diameter, 10 per pack



- Great for building up impression trays
- White in color

Item #: 600-020
1 pack
3 or more ^{\$} 17.49/pk
2/ 11 04 -+-!





Base Plate Wax

- Durable all-purpose wax
- Tough, pliable and tacky

Item #: 600-035	
1 pack	.\$57.50/pk
3 or more	.549.95/pk
5 lbs per pack	



- Strong adhesive properties
- Long working time

Item #: 600-030	
1 pack	§62.99/pl
3 or more	\$59.99/pl
120 per pack	-



Safety Swivel Key

- No fingers in the mouth
- No fear of dropping keys down patient's throat
- Made in the USA

Item #: SSK	
1 pack	\$25.95/pk
5 or more	\$21.95/pk
10 per pack	



Universal Screw

Extremely narrow sizes for use both in upper and lower arches.

10.0 mm	A0800-10
12.0 mm	A0800-12
1 pack	\$53.45/pk
10 per pack	





Standard Screw for Uppers

For all types of transversal expansion movements in the upper palate.


11.0 mm	A0805-11
14.0 mm	A0805-14
1 pack	\$53.45/pk
10 per pack	•



Super Micro Sectional Screw

For movement of single teeth. The U-shaped guide pin ensures high stability.

9.0 mm	019-891-09A
10.0 mm	019-891-10A
1 pack	^{\$} 94.50/pk
10 per pack	



Three Directional Screw with 2 Activation Points

Allows for both transversal and sagittal expansion.

16.0 mm	
1 screw	^{\$} 40.95 ea.





8.0 mm RPE

35 turns for maximum expansion Item #: A0620-08



9.0 mm RPE

40 turns for maximum expansion Item #: A0620-09



11.0 mm RPE

50 turns for maximum expansion Item #: A0620-11



13.0 mm RPE

60 turns for maximum expansion Item #: A0620-13



Rapid Palatal Expanders

The body design of Leone's Rapid Palatal Expanders is completely smooth, with no rough areas or intersections for the maximum in complete patient comfort and hygiene.

Medical Grade Stainless Steel

• Laser Engraved Numbering


Laser Welded Arms

Rapid Palatal Expanders



Activation Keys

Patient Safety Leash



ExpaBent Rapida

Palatal Expander Featuring Pre-Bent Arms

Produced entirely of biomechanical stainless steel. Arms are bent with anatomical angulation. Expansion capacity and direction arrow are laser engraved on the expander body.

- Medical grade stainless steel
- Laser engraved numbering

1 expander	^{\$} 32.95 ea.
13.0 mm	019-2000-13
11.0 mm	
9.0 mm	



Exspider Symmetrical Expander

A posterior articulation joint features expansion arms oriented onwards to get an arc of a circle basal expansion. The use of such expander is indicated when anterior transversal deficiencies on the maxillary arch occur, while the posterior diameter is correct.

- Delivers symmetric forces
- Laser engraved numbering
- Laser welded arms

9.0 mm	A019-2001
1 expander	§63.95 ea.
5 or more	\$58.45 ea.
Rody: 3/4" x 5/16" arm: 11/8"	



Super Lingual Half-Screw Expander

The Half-Screw may be fabricated on the lingual surface of either the maxillary or mandibular arch for expansion and distalization. The overall width has been reduced resulting in an improved fit, allowing for a more comfortable appliance in the patient's mouth.

- The perfect solution for fixed mandibular expansion
- Use for individual quadrant expansion or molar distalization
- Compact fit for comfortable lingual expansion

Item #: A0623-08	
1 expander\$31.85	ea.
5 or more \$29.95	ea.

FIXEDAPPLIANCES

Helps to Break the Habit!



Tongue Tamers[®]







- Eliminates detrimental forces
- Discourages over-active muscle activity
- Bonds to upper or lower arches
- 80 gauge foil mesh bonding base
- High quality 17-4 stainless steel
- Curved base perfect for lingual placement

Simple and Immediate Bite Opening





- Opens deep bites with immediate results
- Bonds lingually to upper anteriors
- Eliminates occlusal interference
- Comfortable for the patient
- Composite material will not wear enamel
- Reduces bond failure on lower brackets
- Curved base perfect for lingual placement
- Made in the USA

© 2015 Ortho Technology, Inc. Tongue Tamers and Bite-Builder are registered trademarks of Ortho Technology.



TRY OUR...

Curved Direct Bond Instrument

- Designed to securely hold and place lingual brackets or attachments
- Works great with our Tongue Tamers, Bite-Builders, and placement of molar tubes
- Made in the USA

Item #: OT-201
1 instrument..........\$24.95 ea.
5 or more\$22.95 ea.
Heat sterilizable up to 370°F/188°C











NEW!

Premium Eruption Appliances

Facilitate the Eruption of Impacted Teeth

WITH SOLID 14K GOLD TRACTION CHAINS









Made with Premium Grade Materials 16 gauge 14K gold or surgical steel chain

16 gauge 14K gold or surgical steel chain with reduced or low nickel content

Guaranteed

Our strong and durable chain will not break and can be left in the mouth for extended periods of time without discoloring, deterioration, or loss of integrity

Versatile

Eyelet selection for profile and base options for individual case treatment with round shaped lumen (chain links) easily accommodates up to a .030" elastic thread

Comfortable

1" long chain is slender, smooth, lightweight for added patient comfort

Made in the USA

Eruption Appliance with 14K Gold Chain

Eruption Appliance with Surgical Steel Chain

For use with up to .030" Solid Elastic Thread **D.** Round Mesh Button Swivel Base ...S19634

1	pack	\$107.95/pk
4	or more	\$99.50/pk
3	per pack	



TRY OUR...

NEW! Eruption Appliance Curved Eyelet & Angled Button Tweezers

 Specialized reverse action tweezers designed to fit and hold low profile eyelet and button style eruption appliance pads

See page 185

Eruption Appliances with 24K Gold Plated Eyelet and Traction Chain







Bondable to cuspids impacted palatally or labially. The 1" traction chain swings freely within the eyelet. The small .016" diameter and large .022" diameter chain links can be cut easily to achieve the desired length.

Made in the USA

Eruption Appliance with Small Traction Chain

For use with .020" Solid Elastic Thread

Round Base/Sm. Chain	
Rectangle Base/Sm. Chain	G19645
1 pack	\$64.50/pk
3 or more	\$60.95/pk
5 per pack, Autoclavable up to	275°F/135°C

Eruption Appliance with Large Traction Chain

For use with .020", .025", and .030" Solid Elastic Thread

Round Base/Lg. Chain	G19650
Rectangle Base/Lg. Chain	G19651
1 pack	\$67.50/pk
3 or more	
5 per pack, Autoclavable up to	o 275°F/135°C

Eruption Appliance Typodont

Features eruption appliances on the exposed palatal and labial upper cuspids. Both the round base and rectangular base (small chain) eruption appliances are used. This typodont also features stainless steel brackets in Roth Rx.

Eruption Appliance Typodont

Item #: G50-205EA

1 typodont**\$282.50 ea.**

FIXEDAPPLANCES



Braided Retainer Wire

Perfect for Permanent Fixed Lingual Retainers Does Not Unwind When Cut

Used to make permanent fixed lingual retainers following orthodontic treatment, the Braided Retainer Wire features a very thin, rectangular wire that can be adapted with no occlusal interference while ensuring patient comfort. It prevents inadvertent tooth movement associated with active force wires. This three strand, stainless steel, braided retainer wire is annealed or "dead soft" so that it can be easily formed and holds bends to make bondable lingual retainers.

Can be used as semi-permanent orthodontic splinting of the upper or lower incisors and to maintain difficult extraction sites.

■ Made in the USA

Item #: 68292	
1 pack	\$42.95/pk
3 or more	\$35.95/pk
.010" depth x .028"	width x 6.0"
lengths, 10 per pac	·k

Flat Titanium Dead-Soft Lingual Wire

Offers easy formability, and ideal for retaining final torque positions

Used for fixed lingual retainers, this wire features increased interproximal strength for semi-permanent splinting and maintaining difficult extraction sites on upper or lower arches. Titanium wire offers reduced wear rates and is nickel-free, eliminating potential patient allergy issues. The "dead-soft" wire can be easily formed, and its flat dimensions are ideal for retaining final torque positions.

■ Made in the USA

Item #: 68293

Bondable Lingual Retainers

Size Measurement

Size selection is determined by millimeter measurement from center of one tooth to the next.

- Anatomical form and offset pads
- Foil mesh bonding pads
- Heat treatable wire
- Made in the USA





Lower Cuspid to Cuspid

▼ Lower/Mandibular Cuspid to Cuspid (3x3) .032" Diameter Wire

#2	22 0	mm	 6	2225
# 4	22.0	1111111	 0	0203
#4	24 0	mm	 6	8281
#6	26.0	mm	 6	8282
11.0	20.0		_	0202
#8	28.0	mm	 ხ	8283
#0	30.0	mm	 6	8284
	- 1			- / -
Ιр	аск		 °35.9	5/pk
2 -		_	 122 E	0/
3 0	r mor	e	 *33.5	U/pK
F	er pa	ماء		
ש כ	er ba	L.K		



Foil Mesh Bonding Pad

Cuspid to Cuspid Kit

Item #: 68280	
1 kit	§119.95 ea.
2 or more	[§] 112.95 ea.
20x retainers, 4	per each size



Upper Central to Central

▲ Upper/Maxillary Central to Central (1x1) .028" Diameter Wire

4.0 mm	68286
6.0 mm	68287
9.0 mm	68288
1 pack	\$35.95/pk
3 or more	\$33.50/pk

Central to Central Kit

ltem #: 68290
1 kit \$119.95 ea.
2 or more
20x retainers: 6x 4.0 mm,
6x 6.0 mm, and 8x 9.0 mm per kit

BURS AND DISCS



ORTHO TECHNOLOGY

Interproximal Reduction



Interproximal Diamond Burs

It is recommended to begin with the medium diamond to quickly remove the required amount of enamel. Follow up with the fine diamond to smooth the tooth surface creating a precise, parallel end result.

- Consistent grit size and friction grip
- Excellent diamond particle exposure
- For use with high speed handpieces



High Speed Interproximal Kit

Item #: HS-IPR-KIT
1 kit......**\$149.95 ea.**

Includes: 2x medium grit long, 2x fine grit long, 2x medium grit short, 2x fine grit short, 2x fine grit mosquito, and bur block

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

© 2015 Ortho Technology, Inc. Galaxy is a trademark of Ortho Technology.



Fine Grit Long (Red Color ID) Item #: 8859-31-010



Medium Grit Long (Blue Color ID) Item #: 859-31-010

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



Fine Grit Short (Red Color ID) Item #: 8858-31-014



Medium Grit Short (Blue Color ID) Item #: 858-31-014



Fine Grit Mosquito (Red Color ID) Item #: 8392-31-016

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



Interproximal Reduction

Diamond Discs

Smoothly Remove Tooth Mass

Flexible design for both reduction and contouring. Smoothly remove tooth mass interproximally to create intra-arch space or open contact points between teeth. Great for use with many popular aligner systems including Invisaligne.

- Flexible for ease of use
- Disc diameter 19.0 mm and 22.0 mm
- For use with low speed handpieces
- Includes 2 diamond discs and 1 mandrel



A. Double-Sided Diamond Disc 19.0 mm88617 22.0 mm88613 1 pack......**\$28.75/pk** 3 or more......\$25.75/pk Disc thickness 0.17 mm

B. Double-Sided See-Thru Diamond Disc

22.0 mm	88655
1 pack	\$40.75/pk
3 or more	\$38.50/pk
Disc thickness 0.2	0 mm

C. Double-Sided Mesh **Diamond Disc**

19.0 mm	88657
22.0 mm	88656
1 pack	\$40.75/pk
3 or more	\$38.50/pk
Disc thickness 0 20	

D. Single-Sided Diamond Disc 19.0 mm88616 22.0 mm88614 1 pack......\$28.75/pk 3 or more......\$25.75/pk Disc thickness 0.10 mm, mounted

diamond side down

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

Invisalign is a registered trademark of Align Technology.



Diamond 3 Disc Intro Kit

22.0 mm Intro Kit

Item #: DD-KIT

1 kit......**\$91.50 ea.**

Includes: 1x double-sided, 1x singlesided, 1x double-sided see-thru with mandrels and bur block

19.0 mm Intro Kit Item #: DDKIT2

1 kit......**\$91.50 ea.**

Includes: 1x double-sided, 1x singlesided, 1x double-sided mesh with mandrels and bur block

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



Diamond 4 Disc Intro Kit

22.0 mm Intro Kit Item #: DDKIT1

Includes: 1x double-sided, 1x single-sided, 1x double-sided see-thru, 1x double-sided 1 kit......s119.75 ea. mesh with mandrels and bur block

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

ORTHO TECHNOLOGY®

Interproximal Reduction



Interproximal Diamond Strips

Perforated, Flexible

Perforated Diamond Strips are designed for complete control during interproximal stripping, shaping, and contouring. Strips allow easy access and precise manual enamel reduction resulting in a smooth, natural finish. Perforated design provides improved visibility, total control, and more flexibility than solid lightning strips. Strips are color coded; yellow for extra fine, red for fine, and blue for medium.

- Single sided and color coded
- Stainless steel with diamond coating
- Uncoated center gateway segment
- 2 widths: 2.5 mm and 3.75 mm
- 3 grits: extra fine (15 micron) fine (30 micron) medium (45 micron)
- Thickness: extra fine (.08 mm) fine (0.10 mm) medium (0.12 mm)
- Length: 147.0 mm



Interproximal Strips Narrow

Yellow, Extra Fine.	88603
Red, Fine	88602
Blue, Medium	88601
Strip width 2.5 mm, 10	per pack

Interproximal Strips Wide

Yellow,	Extra Fine	288503
Red, Fir	1e	88502
Blue, M	edium	88501
Strin width 3	3 75 mm	10 ner nack

1	pack	\$96.75/pk
	or more	
Δ	utoclavable un to	270°F/132°C

Interproximal Kits

Narrow - 2.5 mm	88600
Wide - 3.75 mm	88500
1 kit	\$89.95/kit
3 or more	\$83.25/kit
3x extra fine, 3x fine	., 3x medium

© 2015 Ortho Technology, Inc. Galaxy is a trademark of Ortho Technology.



- Disposable
- For complete visibility
- Fits discs up to 22.0 mm
- Use with doriot style handpieces

Item #: 88615	
1 pack	
3 or more	\$42.95/pk
15 per pack	

Metal Safety Guard

Protects Your Patients' Lips During Interproximal Stripping





Lip Ring Retractor

Great for Adults and Children During Bonding

- Sized for adults and children
- Made of a high quality medical grade plastic



Works great for bonding and intraoral photography

Debonding

Debonding Carbide Burs

- Ideal shapes for finishing all facial and lingual surfaces
- Specifically designed with orthodontic debonding and finishing in mind
- Controlled debonding without nicking, scratching or abrading the enamel
- Corrosion resistant finish
- Smooth, friction grip shank

Available in Slow & High Speed



Slow Speed Latch Type Shank 2.35 mm width



High Speed Friction Grip Shank 1.6 mm width



18 Fluted Small

Item #: H118S-21 (Slow Speed) Item #: H118S (High Speed) 5.7 mm head length



18 Fluted Tapered

Item #: H218-21 (Slow Speed) Item #: H218 (High Speed) 7.3 mm head length

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



18 Fluted Long

Item #: H118L-21 (Slow Speed) Item #: H118L (High Speed) 8.3 mm head length

Slow Speed Latch Shank

1 pack	\$54.99/pk
3 or more	\$52.99/pk
5 per pack	

High Speed Friction Grip Shank

1	pack	\$50.25/pk
3	or more	\$48.25/pk
5	ner nack	



12 Fluted Long Taper

Item #: H375R-21-016 (Slow Speed) Item #: H375R-31-016 (High Speed) 8.0 mm head length



12 Fluted Football

Item #: H379-21-023 (Slow Speed) Item #: H379-31-023 (High Speed) 4.4 mm head length

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



12 Fluted Long Flame

Item #: H48L-21-012 (Slow Speed) Item #: H48L-31-012 (High Speed) 8.0 mm head length

Slow Speed Latch Shank

1 pack	³ 50.95/pk
3 or more	§47.25/pk
5 per pack	

High Speed Friction Grip Shank

1 pack	\$44.95/pk
3 or more	
5 per pack	



- Designed for controlled removal of bulk bonding adhesive after brackets are debonded
- Ideal for posterior debonding

10 Fluted Finishing Carbide Item #: H134-21-014 (Slow Speed) Item #: H134-31-014 (High Speed) 5.9 mm head length

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

Slow Speed Latch Shank

1 pack	\$48.99/pk
3 or more	
5 per pack	


High Speed Friction Grip Shank 1 pack\$43.50/pk

3	or more	\$4	1.	.50/pk	
5	per pack				

Finishing Diamond Burs

BURS AND DISCS

 Ideal for finishing and reshaping teeth during and after orthodontic treatment





Extra Fine Grit Long (Yellow Color ID)

Item #: 135EF-31-014



Fine Grit Long

(Red Color ID) Item #: 135F-31-014

1 pack	\$43.50/pk
3 or more	
High speed, 5 per	pack

Ideal for equilibration and

Extra Fine Grit Football

(Yellow Color ID) Item #: 379EF-31-023



Fine Grit Football

(Red Color ID) Item #: 379F-31-023

Dry heat sterilizable up to 340°F/170°C Autoclavable up to 250°F/121°C



Debonding Burs Intro Kit

Item #: DB-KIT

Dry heat sterilizable up to 340°F/170°C Autoclavable up to 250°F/121°C

Includes: 2x 12 fluted long taper, 2x 12 fluted long flame, 2x 12 fluted football, 2x high shine enamel polisher (W16F) and bur block

Premium



Debonding Kit

Item #: PREM-DB-KIT 1 kit......**\$187.95 ea.** Dry heat sterilizable up to 340°F/170°C Autoclavable up to 250°F/121°

Includes: 2x 10 fluted finishing carbide, 2x 18 fluted small, 2x 18 fluted long, 2x 18 fluted tapered, 2x finishing fine grit long, 2x finishing extra fine grit long, 2x fine grit football, 2x extra fine grit football, 2x high shine enamel polisher (W16F), and bur block



High Speed Burs Intro Kit

Item #: HS-KIT Dry heat sterilizable up to 340°F/170°C Autoclavable up to 250°F/121°C

Includes: 2x terminator metal cutter, 2x IPR medium grit long, 2x IPR fine grit long, 2x 12 fluted long flame, 2x 12 fluted long taper, 2x 12 fluted football, 2x high shine enamel polisher (W16F) and bur block







147



ORTHO TECHNOLOGY

Polishers



Debonding Polishers

- · Safely remove residual adhesive after debonding
- Optimal shapes to reach difficult areas
- Protective silicone based grit does not damage enamel



Debonding Polisher 9.0 mm Item #: G1F 9.0 mm head length



Debonding Polisher 6.5 mm Item #: G3F 6.5 mm head length 7.5 mm head length

Item #: G2F

Debonding Polisher 7.5 mm

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

Finishing Polishing Points

- 5 micron diamond grit suspended in resilient silicone
- Used after debonding to create a high polish finish







High Shine Enamel Polisher 2.35 mm width

Item #: W16D	
1 pack	\$37.95/pk
3 or more	\$35.25/pk
Slow speed, Max. I	RPM 4,000 - 8,000
10.4 mm head l	ength,
3 per pack	

© 2015 Ortho Technology, Inc. Galaxy is a trademark of Ortho Technology.




Acrylic Polishers

- Ideal for finishing & polishing of acrylic appliances
- Made of a resilient silicone

Assorted Pack	
Item #: GLOSS-ASST	
1 pack	\$40.50/pk
5 or more	\$37.25/pk
Includes: 2x coarse,	2x medium,
2x fine 6 ner nack	

Coarse Grit	G634
Medium Grit	G635
Fine Grit	G636
1 pack	\$40.50/pk
5 or more	\$37.25/pk
Max. RPM 10,000,	24.0 mm head
length, 6 per pack	

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



Polishing Wheels

- Pre-mounted on a mandrel for perfect balance
- Ideal for a smooth finish of vacuum formed plastics and metal appliances

Coarse	1670016
Medium	
Fine	1670091
1 pack	\$34.50/pk
5 or more	\$30.95/pk
12 per pack	

Laboratory Cutters & Burs

Vacuum Form Cutter Spiral Flute

 Ideal for acrylic vacuum form plastic and metal cutting



Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

Terminator Metal Cutter

- Designed to smoothly cut through metal
- Eliminates bounce and vibration
- Perfect for herbst or crown removal



Item #: T-2	
1 pack	\$24.45/pk
3 or more	\$22.45/pk
High speed, 4.0 m	ım head length,
5 ner nack	

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

Cutting Wheel

Quickly Cut Away Plastic from Model!

- High quality steel
- One piece tool

Item #: 88800
1 instrument
3 or more \$24.95 ea.
Max RPM 10,000, diameter 10.0 mm,
Dry heat sterilizable up to 370°F/188°C





Bur Blocks



Lab Burs• Holds 14 burs

Item #: BB-100 1 block**\$40.95 ea.**

High Speed Burs

Holds 27 bursItem #: BB-200

1 block......**\$40.95 ea.**

@ ORTHO TECHNOLOGY

Dry heat sterilizable up to $340^{\circ}F/170^{\circ}C$ or autoclavable up to $250^{\circ}F/121^{\circ}C$

Laboratory Carbide Burs

- Extremely precise tooth geometry provides our straight handpiece cutters excellent cutting efficiency of all acrylics while leaving a very smooth finish
- 3/32" shank for use with low speed handpieces
- Speed and pressure sensitive

 you determine how much
 and how fast
- Coated with super hard gold Titanium Nitride to greatly extend the life of the bur



Medium Taper Item #: 79E-11-040 Low speed,

13.0 mm head length



Large Pointed Taper Item #: 257RE-11-060 Low speed, 14.0 mm head length



Large Egg Shape

Item #: 251EF-11-060 Low speed, 14.0 mm head length

1	bur	\$39	.50	ea.
3	or more	\$36	.95	ea.

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



Lab Burs Intro Kit

Includes: 1x large egg shape, 1x medium taper, 1x large pointed taper, 1x vacuum form cutter spiral flute and bur block

Item #: LAB-KIT	
1 kit	\$164.50 ea.
Dry heat sterilizable u	p to 340°F/170°C
or autoclavable up to	250°F/121°C

BURS AND DISCS

MiniStripper

Compact Interproximal Stripping System

Safety & Control

at Your Fingertips

Designed by an Orthodontist

- Small size provides control and stability in tight spaces
- Compact size is safer to use than other interproximal stripping devices
- · Variety of strip sizes for every interproximal stripping need
- Easy to use different strips can be switched out in seconds
- Strips are available as single-sided or double-sided providing precise stripping of intra-arch space
- All kit components are autoclavable up to 275°F/135°C

Stainless steel construction





MiniStripper

MiniStripper Kit Itam #: E4201

MiniStripper Handle

ltem #: 54308	
1 handle	\$45.50 ea.
3 or more	^{\$} 41.25 ea.
Autoclavable up to 275°F/1.	35℃

© 2015 Ortho Technology, Inc. MiniStripper is a trademark of Ortho Technology. Patent Pending.

6.0 mm Single-Sided **Stainless Steel Strips**

Item #: 54304	
1 pack	(
3 or more	(
10 ner nack strin thickness is 0 10 mm	

4.0 mm Single-Sided **Stainless Steel Strips**

ltem #: 54305	
1 pack	\$19.25/pk
3 or more	\$18.25/pk
10 per pack, strip thickness	s is 0.10 mm

6.0 mm Single-Sided Diamond Strips

Item #: 54306	
1 pack	§48.95/pk
3 or more	§45.50/pk
5 per pack, strip thickness	s is 0.15 mm



ORTHODONTICS

Featured in the

World Journal

of Orthodontics

Actual Size

6.0 mm Double-Sided **Stainless Steel Strips**

Item #: 54309	
1 pack	k
3 or more\$39.25/p	k
10 per pack, strip thickness is 0.16 mm	

4.0 mm Double-Sided **Stainless Steel Strips**

Item #: 54310	
1 pack	\$42.50/pk
3 or more	\$39.25/pk
10 per pack, strip thickness	is 0.14 mm

6.0 mm Micro Saw Strips

olo mini micro saw strips	
Item #: 54307	
1 pack	\$37.50/pk
3 or more	\$34.50/pk
5 per pack	

"The Ideal Alternative to Relief Wax"



Snaps Easier Into Place and Stays More Secure Than Competing Brands



Ortho Technology's Comfort Cover shields the patient's lips and cheeks from the discomfort often associated with orthodontic treatment or activities where the risk of oral laceration could occur, such as playing musical instruments or light contact sports where mouthguards are excessive.

Comfort Cover™

Individually packaged

- Packaged for professional presentation to your patients
- Convenience makes living with braces easier every day. Each patient set of Comfort Covers protects both the upper and the lower arch and is made from non-latex material*
- Great for carrying in book bags, handbags, or gym bags
- Made in the USA

Item #: 98165

1 pack......\$31.95/pk

5 patient sets per pack









PATIENT CA

© 2015 Ortho Technology, Inc. Comfort Cover is a trademark of Ortho Technology. U.S. Patent Number D510439. *Not made from natural rubber latex.



SPECIALLY DESIGNED FOR ATHLETES WITH BRACES

ULTRA-GUARD SPORTS MOUTHGUARD

\$7,500⁰⁰ Dental Warranty

Complete Warranty Details Available Online at orthotechnology.com/ultraguard

ULTRA-GUARD SPORTS MOUTHGUARD

Allows for Tooth Movement During Orthodontic Treatment and Maintains Bite Tab Thickness Unlike Boil & Bite Mouthquards



SPORTS MOUTHGUARD



Protection for Both Upper and Lower Teeth Plus Protects Lips and Gums from Lacerations

JULTRA-GUARD SPORTS MOUTHGUARD

Thick Shock Absorbing Bite Tabs Lessen Chances of Concussion



© 2015 All Images and illustrations. Ortho Technology, Inc. ULTRA-GUARD is a registered trademark of Ortho Technology. U.S. Patent Number D473976.

THE PERFECT MOUTHGUARD FOR ALL SPORTS

Shock Absorbing Bite Tabs

Helps to protect the brain from concussion by dissipating the shock of a blow to the lower jaw away from the brain

Ergonomically Contoured

Conforms to the inside of the mouth with smooth rounded edges, providing maximum comfort for the athlete





Double Protection Shield

Helps to protect oral tissue from laceration by shielding the lips and cheeks from both upper and lower teeth



Professional Quality Optional Helmet Strap

The heavy-duty tear resistant strap connects easily to the facemask of helmets



NOW AVAILABLE IN PINK!

ULTRA-GUARD Color Mouthguards Meet NCAA Regulations

Individual Packaging Match Team Colors

	w/strap	no strap
Red	24001	24101
Green	24002	24102
Yellow	24003	24103
Blue	24004	24104
Black	24005	24105
Clear	24006	24106
Pink	24007	24107
Assorted**	24000	24100

**Pink not included.



Non-Latex*

Same price with or without helmet strap. One size fits all.

1 pack	\$40.95/box
5 or more	
12 per box	

*Not made from natural rubber latex.







153



Ortho Performance® Mesh Bag

- One bag with multiple uses
- Can be used to hold headgear items, multi-adjustable facemasks or oral hygiene kits
- Solid blue in color with a secure zipper closure
- Made from premium micro mesh material

© 2015 Ortho Technology, Inc. Ortho Performance is a registered trademark of Ortho Technology.

Item #: 560111	
1 pack	.\$44.49/pk
3 or more	\$40.95/pk
10 per pack, measures	8 ³ / ₄ " x 5"

Elastobyte Wafers

Provides Relief For Patient Discomfort

- Chewing on the elastobyte may increase blood circulation
- Helps relieve patient's pain and discomfort
- Great for recently moved or separated teeth
- Made in the USA

Standard Packs	Large Economy Packs
Grape Scented860-229	Grape Scented860-12
Cherry Scented860-230	Cherry Scented860-13
Bubble Gum Scented860-231	Bubble Gum Scented860-13
1 pack	
4 or more	
25 per pack	200 per pack



PATIENT CARE

Retainer Cases - Portable Protection!



Strong, High-Impact Plastic in Your Choice of Classic, Glow-In-The-Dark, and Sparkle Finishes

- Appliance key holder style
- Durable long life hinge
- Write-on I.D. area
- 1 ¼" depth fits most appliances
- Vented holes prevent odors
- 17 exciting color options
- Made in the USA

MAKES A GREAT MOUTHGUARD CASE FOR MOUTHGUARDS WITHOUT STRAPS!

Classic

350-000
350-001
350-002
350-003
350-004
350-005
350-006
350-007
350-008
350-009
350-010

Glow-In-The-Dark

Assorted350-020

Sparkle

Assorted......350-030

1 pack\$23.75/pk 25 per pack



- Use for retainers, elastics, aligners
- Smooth rounded compact design
- Secure snap-tight closure
- Durable plastic high gloss finish

Color Retainer Cases

ltem #: 45205	
1 pack	\$23.99/pk
3 or more	
10 / 20 10 10 2/4	

12 per pack, 3"x 1" with 3/4" inside depth. 2x each color, melon, tangerine, aqua, violet, royal blue, and white

Black & White Retainer Cases

ltem #: 45204	
1 pack	\$23.99/pk
5 or more	\$21.99/pk
12 per pack, 3" x 1" with	3/4" inside

depth, assorted black and white combinations



Ortho Performance® Patient Relief Wax

Made in the USA

Patient Relief Wax 1 pack\$27.25

1 pack \$27.25/pk 5 or more \$24.25/pk 50 per pack

Relief Wax Dispenser

Plain Wax In Scented Cases

□ Grape	600-001	□ Blueberry600-009
☐ Orange	600-002	☐ Coconut600-011
☐ Cherry	600-003	☐ Chocolate600-012
☐ Cinnamon	600-004	☐ Cotton Candy600-013
☐ Mint	600-005	Unscented Wax Cases
☐ Bubble Gum	600-007	□ Silver600-010
☐ Pineapple	600-008	☐ Glow-In-The-Dark600-006





- Great for adults
- A must for musicians
- Added protection for athletes
- Longer lasting than regular wax
- Made in the USA

Item #: 600-050	
1 pack	\$33.75/pk
5 or more	\$31.75/pk
50x assorted sparkle	cases per pack

INFALAB Powerful Brace Relief in a Compact Kit!







Brace Relief™ Pocket Kit

- Mirrored refillable case for on-the-go application
- All the tools to treat and prevent mouth sores
- Dramatically reduces patient complaints and emergency discomfort visits

Item #: 7900	
1 pack	76.95/pk
3 or more	72.95/pk
36 per pack, each kit o	contains
4x silicone rods and 1)	x tube of gel

PATIENT CARE



DID YOU KNOW...

Brace Relief™ Gel is an excellent pain relieving lubricant for when placing elastomeric separators and nickel titanium separating springs.

Brace Relief™ Gel also numbs the periodontal ligament to provide your patients with more comfortable treatment when

placing appliances!



Great For Use With Separators and Aligners!

Brace Relief™ Gel

- Refill the pocket kit, place into hygiene kits or give out individually
- Apply to tissue at the first sign of irritation and as needed throughout treatment
- Made with 9% benzocaine, natural myrrh, natural mint flavor and healing oils
- Antiseptic and soothing



Item #: 7903	
1 box	\$126.50 ea.
3 or more	§121.95 ea.
144x 3 gm tubes	per box



Brace Gard Mini Pak

- Convenient mirrored case
- Medical grade silicone will not crumble or deteriorate
- Reusable without leaving residue • Will not melt or become brittle
- Cover one bracket or an entire arch

Item #: 7906	
1 pack	\$56.50/pl
3 or more	\$51.95/pl
80 per pack, each	kit contains
3x silicone rods	



Brace Gard Silicone

- Refill the pocket kit, place into hygiene kits or give out individually
- Thick, pliable and transparent
- Medical grade silicone will not crumble or deteriorate
- Reusable without leaving residue
- Not temperature sensitive like wax



Item #: 7908	
1 box	\$74.50 ea.
3 or more	§69.50 ea.
200x 2 rod nacks	ner hox

PATIENTICARE

Making Remineralization Better

Delivers Complex Supercharged Calcium and Phosphate Ions

MI Paste™ & MI Paste Plus™



MI Paste and MI Paste Plus with fluoride contain RECALDENT™ (CPP-ACP) and are preventive options used to correct the mineral imbalance by delivering a complex of supercharged calcium and phosphate ions to the enamel to form calcium phosphate crystals and help remineralize the teeth. MI Paste and MI Paste Plus with fluoride helps reverse the effects of decalcification. Apply topically 1-2 times per day for 3-5 minutes.





MI Paste™

- Replaces lost minerals
- Protects teeth from decay and decalcification
- Strengthens tooth enamel
- Reduces sensitivity
- Buffers plaque acid
- Water based and sugar free


Mint Flavor

Item #: 003679 1 case**\$168.95/case** 3 or more \$165.95/case 10 tubes of MI Paste (40 am each) per case, mint flavor

Assorted Flavors

Item #: 002265 3 or more \$165.95/case 10 tubes of MI Paste (40 gm each) per case, two tubes of each flavor (melon, mint, strawberry, tutti frutti and vanilla)

MI Paste Plus™



- Contains 900 ppm fluoride ions (0.2% sodium fluoride)
- Superior to fluoride alone
- Increases fluoride uptake into plaque and subsurface enamel
- Depresses demineralization

Mint Flavor*

Item #: 002621

1 case......**\$168.95/case** 3 or more\$165.95/case 10 tubes of MI Paste Plus (40 gm each) per case, mint flavor

Item #: 002888

1 case......**\$168.95/case** 3 or more \$165.95/case 10x tubes of MI Paste Plus (40 gm each) per case, vanilla flavor

Assorted Flavors*

Item #: 002614

3 or more\$165.95/case 10 tubes of MI Paste Plus (40 gm each) per case, two tubes of each flavor (melon, mint, strawberry, tutti frutti and vanilla)

*Not available in Canada





Advanced Teeth WHITENING SYSTEM

Whether your patient is currently in orthodontic treatment or just had their braces removed, the advanced teeth whitening kit will help enhance the smiles you've created.

The whitening gel achieves noticeable results for patients in or out of treatment. While the reusable trays reduce the need for additional lab costs, the gel works easily with custom trays.

Ortho Performance Whitening is an excellent way to promote your practice!



PATIENTICARE



The Perfect Complement to a STRAIGHT SMILE

- Two concentrations available. 16% and 22% Carbamide Peroxide, for treatment flexibility
- · Contains Potassium Nitrate, water based gel keeps teeth hydrated
- Thick gel prevents run-off on to gums to lessen irritation
- Neutral PH is safe for enamel.
- Mint Flavored
- · Wear time: 2 4 hours

Whitening Starter Kits

16% Carbamide Peroxide	ORT50110
22% Carbamide Peroxide	ORT50100
1 case	\$227.70/case
6 or more	\$203.94/case
611	_

6 kits per case, 4x 3 mL syringes,

2 whitening trays, storage case, and instructions

© 2015 Ortho Technology, Inc. Ortho Performance is a registered trademark of Ortho Technology.





PATIENTICARE



Patient Compliance and Preventative Care Made Simple!



Complete Oral Hygiene In One Convenient Take Along Kit

COMPARE TO:GUM° Orthodontic Dental Kit

Ortho Performance® Premium Ortho Hygiene Kit 5 Piece

A compact accessory kit to help patients adjust to their new orthodontic appliances and hygiene routine

Premium Ortho Hygiene Kit Contents

- Premium bi-level ortho tapered head toothbrush (ORT10789)
- Ortho wax in jewel case
- 10 floss threaders in plastic case (ORT20412)
- Interproximal brush-cone profile
- Premium mint floss (ORT48037)
- Plastic case

Assorted colors:



1 case\$158.40/case (\$4.40/kit) 4 or more\$140.40/case (\$3.90/kit)

36 per case, measures 8" x 3"



Ortho Performance Ortho Hygiene Kit 8 Piece

The versatility of our Ortho Performance 8 Piece Hygiene kit encourages patient compliance. Made of premium micro mesh material, the eye-catching blue bag is breathable, yet very durable. Even with the kit pieces included, the bag has storage room for more. Perfect for patients on the go!

Ortho Hygiene Kit Contents

- Compact v-trim 35 tuft orthodontic toothbrush
- Travel v-trim 35 tuft toothbrush
- Cylindrical intraoral proxy brush (20043)
- 10 pack of floss threaders
- 12 yards of waxed dental floss
- 1 case of non-allergenic wax for patient comfort
- Plastic patient mouth mirror
- Tooth brushing 2 minute timer
- Royal blue micro mesh bag (560111)

Item #: ORT48017

1 case\$179.69/case (\$4.99/kit)
4 or more\$165.24/case (\$4.59/kit)
36 kits per case/complete cases only
measures 8 3/4" x 5"

160



NEW! Our **Premium Line** of **Toothbrushes**





Item #: ORT10788

COMPARE TO: GUM® Technique Toothbrush

30 tufts with extended tapered bristles to clean interproximal areas and sweep debris from gums.

Assorted colors:





Ortho Performance® **Premium Bi-Level Ortho Tapered Head Toothbrush** – Adult Soft

Item #: ORT10789

40 tufts with v-trim bristle cut and a flexible neck for easier access while cleaning.

Assorted colors:

🕽 Blue 🛑 Green 🌑 Purple 🛑 Red

COMPARE TO: GUM® "V" Trim Toothbrush / Colgate® Orthodontic Toothbrush



Ortho Performance® Premium Sensitive Full Head Toothbrush – Adult Extra Soft

Item #: ORT16778

COMPARE TO: GUM[®] Summit Toothbrush

37 tufts with bristled tapered filaments that are ideal for periodontal cleaning.

Assorted colors:

🔵 Light Blue 🔵 Sea-Foam Green 🔵 Lavender 🔵 Pink



Ortho Performance® Premium Multi-Directional Compact Head Toothbrush - Adult Soft

Item #: ORT16996

32 tufts with opposing angled bristles and power tip to reach deep between teeth, and tongue cleaner.

Assorted colors:

Blue 🦲 Green 🌑 Purple 🛑 Red

COMPARE TO: GUM* Super Tip* Toothbrush / Oral B* Cross Action* Toothbrush

Ortho Performance Premium Toothbrushes

4 or more........554.00/pk

72 toothbrushes per pack, assorted colors



Ortho Performance Premium Toothbrushes Individually Packaged for Patient Presentation







COMPARE TO: GUM® Eez-Thru® Dental Floss

Ortho Performance® Premium Dental Floss



Premium Dental Floss - Mint Item #: ORT48037

1 pack\$30.25/pk 3 or more\$27.25/pk 4 m each, 72 per pack

Premium Dental Floss - Plain

Item #: ORT48038

4 m each, 72 per pack



Ortho Performance® Premium Floss Threader

COMPARE TO: GUM® Proxabrush® Trav-Ler®

Item #: ORT20412 3 or more\$45.00/case 10 per plastic holder, 100 holders per case



Ortho Performance® Interproximal Brush Item #: ORT20040

1 pack\$26.75/pk 36 per pack

Features a blue grip, cone profile 0.45 mm wire with finger rest for ease of use.

Cone Profile





COMPARE TO: GUM[®] V Trim Travel Toothbrushes

Ortho Performance® **Premium Ortho V-Trim Travel Toothbrushes**

Compact body v-trim bristle cut with 34 tufts, designed to snap securely into place for travel.

Item #: ORT44003	
1 pack	59.00/pk
3 or more	55.00/pk
72x toothbrushes pe	r pack

Assorted colors:



Ortho Performance® V-Trim Toothbrushes

Soft nylon bristles with rounded ends provide maximum patient comfort. Assorted bright colors make brushing fun for the patient. Each 41 tuft brush is individually wrapped.

Item #: 34596 144x toothbrushes per pack

Assorted colors:



Red Green



© 2015 Ortho Technology, Inc. Ortho Performance is a registered trademark of Ortho Technology. All other trademarks are property of their respective owners.

Patient Compliance and Preventative Care Made Simple!










Ortho Performance® **Dual-Ended Toothbrushes**

Provides gentle and effective cleaning around brackets and wires. The bristles are made from high quality, round-ended, soft nylon and are firmly anchored in place.

- The perfect brush for orthodontic usage
- V-trim toothbrush with an interproximal brush attached to the handle

Item #: 34597	
1 pack	\$72.50/pk
3 or more	\$66.95/pk
72x toothbrushes per	pack.

36x translucent blue and 36x clear



Ortho Performance® Interproximal Brushes

Taper or cylinder shaped brush with soft, high quality, roundended nylon bristles with 0.40 mm wire, flexible handle and travel cap. Encourages patients to clean around the brackets and underneath wires.

- Brush head is adjustable for hard to reach areas
- Perfect for on-the-go orthodontic usage with protective cap

A. Tapered	20044
B. Cylindrical	20043
1 pack	.\$26.75/pk
5 or more	
36x brushes per paci	



TRY OUR NEW TOOTHBRUSH DISPENSER!

The Perfect Brushing Station Organizer!

See page 189



Item #: OT-TB DISPLAY



Mint Flavor

Ortho Performance® Adult Disposable Pre-Pasted Toothbrushes

Soft, polished bristles with rounded ends provide maximum patient comfort. Each 39 tuft brush is pre-pasted with transparent, US formulated, non-fluoride, mint flavored toothpaste. Individually packaged in a single cellophane bag. Ideal for patient use prior to operatory visit.

Item #: OK11091	/A
1 pack	\$24.99/pk
3 or more	\$23.99/pk
144x toothbrushes p	er pack

" OPT10017

Assorted colors: Blue Red Yellow Green



Bubble Gum Flavor

Ortho Performance® Child Disposable Pre-Pasted Toothbrushes

Extra soft, polished bristles with rounded ends provide maximum patient comfort. Each 30 tuft brush is pre-pasted with a US made fluoride solution, bubble gum flavored toothpaste. Individually packaged in a single cellophane bag. Ideal for patient use prior to operatory visit.

Item #: ORT10930	
1 pack	\$24.99/pk
3 or more	\$23.99/pk
144x toothhrushes ner	nack

Assorted Sparkle colors:

Pink
Blue

PurpleGreen

PATIENTICAREI

Cleans and Brightens Removable Dental Appliances!



Retainer Brite®

Retainer Brite cleans Essix and Hawley retainers – plus nightguards, mouthguards, splints, and most other removable dental appliances.

- Helps remove plaque and tartar
- Kills germs and bacteria
- Cleans and brightens
- Keeps appliances minty fresh
- Made in the USA

Item #: 18159	
1 box	\$11.95/box
10 boxes	\$11.25/box
24 or more	\$10.25/box
96 tablets ner hox	12 hoxes ner case

§6.95/box
§6.25/box
\$5.75/box
24 boxes per case




■ PHOTOGRAPHY




Clinical Photography Supplies

Intraoral Rhodium Coated Photographic Mirrors









The Durable Surface Coating Produces a Sharp and Precise Image Every Time

A. Occlusal (Child)

Item #: 32700 $3^{1}/_{3}$ " (84.6 mm) x $2^{2}/_{5}$ " (61.5 mm) x 2" (50.8 mm)

B. Occlusal (Child Extra Long)

Item #: 32701

 $4^{1/2}$ " (114.3 mm) x $2^{1/2}$ " (63.5 mm) x $2^{1/4}$ " (57.2 mm)

C. Occlusal (Adult)

Item #: 32702

 $4" (101.6 \text{ mm}) \times 3" (76.2 \text{ mm}) \times 2^{1/2}" (63.5 \text{ mm})$

D. Occlusal (Adult Extra Large)

Item #: 32703

 $5^{1/2}$ " (139.7 mm) x 3" (76.2 mm) x $2^{2/3}$ " (67.8 mm)

E. Lingual

Item #: 32704

 $6" (152.4 \text{ mm}) \times 1^{2/5}" (35.6 \text{ mm}) \times 1^{1/2}" (38.1 \text{ mm})$

F. Buccal (Narrow)

Item #: 32705

5³/₄" (145.0 mm) x 2" (50.8 mm) x 1" (25.4 mm)

G. Buccal

Item #: 32706

 $6^{1}/4$ " (158.8 mm) x $2^{1}/4$ " (57.2 mm) x $1^{3}/4$ " (44.4 mm)





(E) Buccal



(G) Lingual

- Superior color and contrast
- Rhodium coated glass

(C) Occlusal (Adult)

- Tempered glass coated on both sides
- Optimum reflectivity and adhesion
- Resistant to abrasion

Follow us on Twitter


- Extremely durable coating
- Use with digital or 35 mm
- Withstands repeated sterilization

1 mirror......**\$64.50 ea.** Autoclavable up to 270°F/132°C


The surface of PixSure Perfect® rhodium coated mirrors ensures maximum reflectivity, adhesion, and resistance to abrasion













Intraoral Stainless Steel Photographic Mirrors

These Highly Polished Stainless Steel Double-Ended Mirrors Provide the Perfect Imaging Solution for Every Patient

- No double reflections
- Polished stainless steel
- Double ended versatility
- Extremely durable finish
- Use with digital or 35 mm
- Withstands repeated sterilization
- The highly polished surfaces of PixSure Perfect® stainless steel mirrors ensure high reflectivity and are virtually indestructible

1 mirror	\$73.95 ea.
4 or more	\$71.95 ea.
Autoclavable up to 2	73°F/134°C

- **A. Adult Occlusal/Child Occlusal**....32800 2³/₄" (70.0 mm) x 7" (178.0 mm) x 2³/₁₆" (55.0 mm)
- **B.** Adult Occlusal/Adult Buccal32803 2¹³/₁₆" (72.0 mm) x 7¹/₂" (190.0 mm) x 2¹/₈" (34.0 mm)

ease of use

- **C.** Adult Occlusal/Child Occlusal....32802 $2^{7}/8$ " (73.0 mm) x $7^{1}/4$ " (185.0 mm) x $2^{9}/16$ " (66.0 mm)
- **D. Child Occlusal/Child Buccal**......32806 2%/16" (65.0 mm) x 615/16" (176.0 mm) x 113/16" (47.0 mm)
- E. Lingual Upper/Lingual Lower.....32804 113/16" (47.0 mm) x 711/16" (195.0 mm) x 113/16" (47.0 mm)
- **F.** Adult Buccal/Child Buccal32805 2²/₁₆" (62.0 mm) x 7³/₈" (189.0 mm) x 1¹³/₁₆" (47.0 mm)
- **G. Buccal**......32801 15/s" (42.0 mm) x 6 "/ι_δ" (170.0 mm) x 1 5/s" (42.0 mm)

■PHOTOGRAPHY



PixSure Perfect

Clinical Photography Supplies

PixSure Perfect®

Retract EEZ® Photo Lip Retractor

Innovative Design Provides Maximum Exposure

- Features easy grip handles
- Great anatomical fit
- One size fits most
- Excellent for occlusal photographs

Not recommended for Dentronix

- Angled bend helps retract lips
- Smooth rounded edge for patient comfort
- Made in the USA

Item #: 84516 (Clear)	
1 retractor	\$28.50 ea.
3 or more	\$25.50 ea.
Autoclavable up to 290°F	F/143°C
Cold Sterilizable	
Heat sterilization may co	ause hazing



PixSure Perfect®

Stainless Steel Double Ended Photo Cheek Retractors

Two Sizes in One for a Sure Fit and Full Unobstructed View

- Double ended with a large retractor on one side and a small retractor on the other
- High quality stainless steel will not discolor, yellow or become brittle like plastic
- Brushed stainless steel construction
- Accepts any method of sterilization
- Ultra-slim design for maximum comfort

ltem #: DT-159	
1 pack	§63.50/pk
3 or more	\$55.95/pk
Sterilizable up to 320°F/16	50°C
1 pair per pack	







_PHOTOGRAPHY■

PixSure Perfect®

Photo Cheek Retractors Double Ended

High Density Plastic Withstands Repeated Sterilization





- Sized to comfortably fit most adults and children
- Durable plastic design
- Will not turn cloudy
- Lightweight construction
- Will not dry out or crack
- Suitable for most types of posterior photography



Item #: DT-155 (Clear)	
1 pack	\$21.95/pk
3 or more	\$19.95/pk
Autoclavable up to 250°F/12	1°C
1 pair per pack	

PixSure Perfect®

Photo Cheek Retractors Double Ended

Smooth Edges for Patient Comfort

- Sized for adult or child
- Maximum view of posterior
- Easily holds lip and cheeks
- Designed for comfort
- Made in the USA



Autoclavable	(Gray)
--------------	--------

Item #: 600-301	
1 pack	\$35.95/pk
3 or more	\$33.95/pk
Autoclavable up to 270°F/132	?°C
1 pair per pack	

Heat Sterilizable (Amber) Item #: 600-302

item #. 000-302
1 pack\$39.50/pk
3 or more\$36.95/pk
Heat Sterilizable up to 370°F/188°C
1 nair ner nack



Cold Sterilizable (Clear)


Item #: 600-303	
1 pack	\$23.95/pk
3 or more	\$21.50/pk
1 nair nar nack	

PixSure Perfect®

Photo Cheek Retractors

Long Handles for Easy Grasp








- Durable plastic design
- Adult and child sizes
- Smooth edges for patient comfort
- Easily withstands repeated sterilization
- Provides full exposure of the posterior

7		
	Small (Blue)	
	Item #: WW-605451	
	Large (Blue)	
	Item #: WW-605450	
	1 pack	\$35.50/pk
	3 or more	\$31.50/pk

Autoclavable up to 273°F/134°C

1 pair per pack



Cut Your Film Costs By More Than 50%



TruVision

EXTRAORAL X-RAY FILM

100 sheets per box

- Automatic or manual processing
- Superb archival stability
- Sharp image quality
- Works with all brands of chemistries

TruVision® Green Sensitive

OTOS-0512 ...Green Panoramic5" x 12" (12.7 x 30.5 cm)......\$99.50 ea. OTOS-0612 ...Green Panoramic6" x 12" (15.2 x 30.5 cm).....\$107.50 ea. OTOS-0810Green Cephalometric.....8" x 10" (20.3 x 25.4 cm).....\$107.50 ea. 100 sheets per box

Ortho Technology's TruVision® Conversion Chart


Film Type	Green Sensitive
TruVision®	TruVision® Green Sensitive
Agfa®	ST8 G
Kodak®	T-MAT G EKTAVISION

Cephalometric Tracing Pads

- Lays flat, will not curl
- Superior optical clarity
- Matte finish
- .003" thickness acetate
- Generous 8" x 10"
- Made in the USA

Item #: OT-3000	
1 pad	\$38.50 ea.
4 or more	

100 sheets per pad



© 2015 Ortho Technology, Inc. TruVision is a registered trademark of Ortho Technology. All other trademarks are the property of their respective owners.

OF THODONT OINSTRUMENTS



X7 Orthodontic Instruments Premium Cutters and Pliers

- Precision Manufactured in the USA with USA and/or German Stainless Steel
- Hand-Crafted Through 80% of the Process
- High Quality 400 Series Surgical Stainless Steel Material Blend
- Stringent Quality Control for Consistent Performance and Reliability
- Handle Designed for Optimal Fit and Comfort
- Contoured Corners Around the Joint Area for a Slim-Line Look and Feel
- Shadow Satin Surfacing to Minimize Glare and Provide a Long-Lasting Finish
- Lifetime Warranty Against Materials Defects and Craftsmanship**
- Steam Sterilize for at least 4 Minutes at 270°F/132°C or 30 Minutes at 250°F/121°C



Distal End Cutter Item #: 400-100

Designed to cut all types of wires up to .021" x .025" (.53 mm x .64 mm). The head safely holds the cut portion of wire to be discarded.



Distal End Cutter with Long Handle Item #: 400-101

Extra long handles designed for additional user comfort. Cuts all types of wires up to .021" x .025" (.53 mm x .64 mm) and safely holds the cut portion to be discarded.



Flush Cut Distal End Cutter
Item #: 400-103
Designed to cut and hold archwire distally

Designed to cut and hold archwire distally flush to the buccal tube or bracket while reducing debondings. It is designed for all types of wire up to a maximum of .019" x .025" (.48 mm x .64 mm). This cutter comes with 1 silicone insert already loaded in the working end and two extras included in the package. Replacement silicone inserts can be re-ordered as needed.

© 2015 Ortho Technology, Inc. X7 is a registered trademark of Ortho Technology.
*Multi Use Cutter Item #: 400-106 **For full warranty and instructions for care information visit www.orthotechnology.com/x7

ORTHODONTICINSTRUME





Micro Cutter Item #: 400-099

The tapered tips are designed to cut soft wire pins and ligature wires up to .012" (.30 mm). The slim, short curved tips allow for access into difficult areas with 15° offset tips. Cutting edge length 0.39" (9.87 mm).



Micro Mini Pin & Ligature Cutter, 15° Angled Item #: 400-108

Designed to cut soft wire pins and ligature wires up to .012" (.30 mm). The 15° offset tip with ultra slim profile is a plus in the most difficult to access areas. Cutting edge length 0.54" (13.72 mm).



Micro Mini Pin **& Ligature Cutter** Item #: 400-104

Designed to cut soft wire pins and ligature wires up to .012" (.30 mm). The ultra slim profile, pointed and tapered tips allow access in the most difficult areas. Cutting edge length 0.50" (12.67 mm).



Multi-Use Cutter Item #: 400-106

Designed to cut soft wire pins and ligature wires up to .015" (.38 mm). Cutting edge length 0.50" (12.67 mm).



Mini Pin & Ligature Cutter Item #: 400-105

Designed to cut soft wire pins and ligature wires up to .012" (.30 mm). The slim profile is a plus in difficult to access areas. Cutting edge length 0.50" (12.67 mm).



Hard Wire Cutter Item #: 400-107

Designed to cut hard wires up to .021" x .025" (.53 mm x .64 mm). Cutting edge length 0.46" (11.70 mm).



Long Tapered Bird Beak Plier Item #: 400-109

Designed to bend a variety of smooth loops up to .021" x .025" (.53 mm x .64 mm).



Bird Beak Plier with Cutter Item #: 400-111

Designed for smooth bending surfaces for a range of loop sizes while providing a sharp blade for cutting. Cuts up to .021" x .025" (.53 mm x .64 mm).



Bird Beak Plier Item #: 400-110

Designed for smooth wire bending for a range of loops sizes for round wires up to .030" (.76 mm).



Light Wire Bird Beak with Groove Plier Item #: 400-130

Designed with a longer cone than a standard bird beak and is useful for forming precise loops in light wires up to .020" (.51 mm). The groove allows for







locating and duplicating loops.

171

ORTHODONTIC INSTRUMENTS -





Light Wire Bird Beak Plier Item #: 400-136

Designed with a longer cone than a standard bird beak and is useful for forming precise loops in light wires up to .020" (.51 mm).



Posterior Band Removing Plier, Long Handle Item #: 400-120

Designed for clinician comfort. The alignment of the tip and occlusal pad is well suited for posterior band removal. These pliers come with a $\frac{3}{16}$ " (4.76 mm) pad already loaded in the working end and an extra $\frac{1}{4}$ " (6.35 mm) pad for trial purposes. Replacement pads can be re-ordered as needed.



Small Weingart Plier Item #: 400-114

Designed with slim, tapered 20° tips with serrations that offer greater access and versatility for a secure grip of archwires and auxiliaries.



Angulated Debonding Plier Item #: 400-122

Designed with a long handle for added user comfort and high angle design for effective anterior and posterior removal of all bracket style types–stainless steel, ceramic, and plastic brackets.



Standard How Utility Plier Item #: 400-115

Designed with serrated tips to provide useful placement and removal of archwires, pins and other auxiliaries.



Debonding Plier Item #: 400-123

Designed for effective anterior and posterior removal of all bracket types – stainless steel, ceramic, and plastic brackets.



Rectangular Arch Forming Plier Item #: 400-116

Designed for placing first, second, and third order bends up to .070" (1.78 mm) blade width and on wires up to .021" x .025" (.53 mm x .64 mm).



Adhesive Removing Plier Item #: 400-124

Designed with a reversible, replaceable blade that is used to remove the bulk of composite material after debonding brackets and attachments. These pliers come with a 3/16" (4.76 mm) pad already loaded in the working end and an extra 1/4" (6.35 mm) pad for trial purposes. Replacement pads can be re-ordered as needed.



Lingual Arch Forming Plier Item #: 400-132

Designed to create double back and triple back bends in .030" (.76 mm) and .036" (.91 mm) lingual arch wires, up to .036" (.91 mm).



Three Jaw Plier Item #: 400-117

Designed for contouring and bending wires up to .030" (.76 mm).





NiTi Three Jaw Plier Item #: 400-135

Designed for the forming and contouring of all archwires, especially for nickel titanium wires up to .020" (.51 mm).



Stop (V-Bend) Plier Item #: 400-134

Places a V-Bend in archwires, including nickel titanium, to prevent the archwire from traveling. Wires up to .021" x .025" (.53 mm x .64 mm).



NiTi Distal Cinch Back Plier Item #: 400-127

Designed as a multi-use bender for bending NiTi wire up to .025" (.64 mm). Tie-back NiTi with no heat.



Tweed Loop Forming Plier Item #: 400-131

Designed to form loops with .045" (1.14 mm), .060" (1.52 mm), and .075" (1.91 mm) diameters for creating bull and helical loops on wires up to .021" x .025" (.53 mm x .64 mm).



Jaraback Plier Item #: 400-129

Designed to be a versatile light wire bending plier with grooves that assist with precise bending and closing of loops. The serrated surface grips wires up to .020" (.51 mm).



Replacement Teflon™ Tips item #: 400-124T

1/4" (6.35 mm) Teflon replacement tips for Adhesive Removing Plier (Item #: 400-124).

1 pack\$29.99/pk 6 per pack



Crown Band Contouring Plier Item #: 400-128

Designed for reshaping and contouring



Replacement Blades Item #: 400-124TT

5.0 mm replacement blade for Adhesive Removing Plier (Item #: 400-124).

1 blade**\$29.99 ea.**



steel molar bands for a customized fit.



Replacement Teflon™ Tips Item #: 400-120T

3/16" (4.76 mm) Teflon replacement tips for Posterior Band Removing Plier (Item #: 400-120).

1 pack\$29.99/pk 6 per pack



Hollow Chop Contouring Plier Item #: 400-133

Perfect for subtle contouring of arch wires up to .030" (.76 mm).



Replacement Teflon™ Tips Item #: 400-120TT

1/4" (6.35 mm) Teflon replacement tips for Posterior Band Removing Plier (Item #: 400-120).

1 pack\$29.99/pk 6 per pack



Replacement Silicone Inserts Item #: 400-103T

Silicone replacement inserts for Flush Cut Distal End Cutter (Item #: 400-103).



Quality Craftsmanship

All PANTERA Pliers & Cutters are hand crafted in the USA using the finest grade 420 stainless steel. They feature a smooth, satin finish for reduced glare, an easy to clean surface, and are corrosion resistant. Handle size, shape, and weight create perfect balance. The instruments are all hand polished, leaving smooth non-glare contoured edges for both user and patient comfort.

Box Joint Construction

All PANTERA Pliers & Cutters have a box joint hinged construction for durability and continuously smooth operation. Box joint pliers are known for perfect tip alignment, longevity and exceptional performance each and every time. Handles were designed with your comfort and ergonomics in mind, because we know a good plier goes a long way.

Highest Quality Material

The cutting edges of all inserted PANTERA instruments are manufactured from a custom blended specialty alloy which offers many advantages over tool steel and tungsten carbide. The PANTERA plier inserts are highly corrosion resistant and will remain bright without tarnishing throughout the life of the instrument. All of the cutting edges are exceptionally tough without being brittle like tungsten carbide and will provide superior sharpness and long lasting strength.

Optimum Performance

You will appreciate the optimum performance of our PANTERA Pliers & Cutters and we stand behind our instruments with a lifetime warranty.

Sterilization

Stands up to all methods of heat sterilization, up to 370°F/188°C.

PANTERA Pliers and Cutters include a lifetime warranty against materials defects and craftsmanship. Like all Orthodontic companies, 'lifetime' references the expected life of the product. Cutting instruments and wire benders are guaranteed for seven years and utility pliers for ten years. This covers breakage and joint failure under normal use. Sharpening cutters and/or repairing minor tip damage is not covered by warranty as this is considered to be routine maintenance for normal use.

1 instrument \$112.50 ea. 6 instruments \$106.50 ea.

© 2015 Ortho Technology, Inc. PANTERA is a registered trademark of Ortho Technology. *Mini Ligature Cutter Item #: 000104



Mini Distal End Cutter with Safety Hold

Item #: 000100

This slim, small cutter is designed for easy access to tight areas and offers a safety hold feature for securely catching and holding the distal end of wires after cutting. The inserted blade design will make cuts within 0.5 mm of the buccal tube. For use with round wires up to .020" and rectangular wires up to .018" x .025".



Long Handle Distal End Cutter with Safety Hold Item #: 000101

The long handle cutter allows easy access to second molars and uses less force for cutting. Features a safety hold, and the inserted blade design will make cuts within 0.5 mm of the buccal tube. For use with round wires up to .020" and rectangular wires up to .018" x .025".



Standard Distal End Cutter with Safety Hold Item #: 000102

This universal standard size cutter easily cuts distal ends. The inserted blade design will make cuts within 0.5 mm of the buccal tube, and catches wire fragments with a safety hold feature. For use with round wires up to .020" and rectangular wires up to .022" x .028".

ORTHODONTICINSTRUMENTS





Flush Cut Distal End Cutter with Safety Hold

Item #: 000103

This exceptional standard size cutter easily cuts wires flush to the buccal tube with the precision inserted blades and features a safety hold spring design to securely catch and hold the distal ends of wire. For use with round wires up to .020" and rectangular wires up to .022" x .028".



Light Wire Plier Item #: 000109

This versatile wire bending plier with precision tapered tips is used for intricate detail work without damaging the wire. Inserted tip sizes are .025". For use with round wires up to .020" or rectangular wires up to .016" x .022".



Mini Ligature Cutter Item #: 000104

The small and slim design of this ligature cutter was specifically designed for easy access to difficult areas. Features stay sharp diamond honed inserted tips. Maximum cutting capacity, soft wire up to .014".



Standard Bird Beak Item #: 000110

The standard plier for heavy wire bending and detail work. The tips are precision crafted to avoid nicking the wire. Inserted tip sizes are .035". For use with round wires up to .032" or rectangular wires up to .022" x .028".



Standard Ligature Cutter Item #: 000106

This universal standard size ligature cutter features a concave frontal design for maximum visibility and accessibility and stay sharp diamond honed insert tips. Maximum cutting capacity, soft wire up to .020".



Standard Bird Beak with Cutter

Item #: 000111

This unique plier combines the standard bird beak and a heavy duty wire cutter at the base of the jaw for both wire forming and cutting. Inserted tip sizes are .035". For use with round wires up to .032" or rectangular wires up to .022" x .028".



Hard Wire Cutter Item #: 000107

The standard size wire cutter features specially designed diamond honed inserted cutting edges that withstand the rigors of cutting hard wires and maintain sharpness along with perfect tip alignment. For use with round or rectangular wires up to .022" x .028".



Serrated Bird Beak Item #: 000112

This plier provides maximum grip and increased control without damaging the wire during bending. Inserted tip sizes are .028" round and .035" square. For use with round wire up to .030" or rectangular wires up to .022" x .028".



30° Angled Ligature Cutter Item #: 000108

The angled tips of this ligature cutter provide easy access and increased visibility for hard to reach areas. Also features diamond honed inserted tips that stay sharp after repeated use. Maximum cutting capacity, soft wire up to .020".



Standard Weingart Plier Item #: 000113

This versatile utility plier has an angled and tapered head design. The fully serrated tips allow for a secure grip of archwires and auxiliaries without damage. Inserted tips are curved gently to form a convenient working angle and ensure lip and tissue safety. For use with round or rectangular wires up to .022" x .028".

ORTHODONTIC INSTRUMENTS -





Small Weingart Plier Item #: 000114

A smaller, lighter duty version of our Standard Weingart. The slimmer profile allows for easy access when space and visibility are limited. The fully serrated tips allow for a secure grip of archwires and auxiliaries without damage. Tips are curved gently to form a convenient working angle and ensure lip and tissue safety. For use with round or rectangular wires up to .016" x .022".



Weingart Plier, Slim Tips Item #: 000158

This Weingart plier is designed with curved ultra-slim tips and elongated handles for best access to areas with least visibility and space. The fully serrated pads offer secure grip of archwires and auxiliaries without damage. The tips are also rounded to ensure lip and soft tissue safety. For use with round or rectangular wires up to .016" x .022".



Standard How Utility Plier Item #: 000115

This multi-functional straight utility plier has 3.5 mm serrated inserted tips which provide excellent grip and a secure hold without nicking the wire. The noncontacting bowed tips ensure lip and tissue safety. For use with round or rectangular wires up to .022" x .028".



Rectangular Arch Forming Plier (Tweed) Item #: 000116

This all around wire bending plier is also perfect for archwire torquing. The heavy duty inserted tips will not damage the wire and are designed to be parallel at an opening of .020". For use with round or rectangular wires up to .022" x .028".



Three Prong Plier Item #: 000117

A standard size plier with strong, precision tips designed for bending and adjusting archwires or lab wires. All three beaks are smooth to avoid wire damage. For use with round wires up to .030" or rectangular wires up to .022" x .028".



0.5 mm Step Plier Item #: 000118

Easily detail your archwire by making 0.5 mm bayonet style bends with this specialty inserted tip plier. Also used for stepping the wire up or down, in or out and can be used intraorally or extraorally. For use with round or rectangular wires up to .022" x .028".



1.0 mm Step Plier Item #: 000119

Easily detail your archwire by making 1.0 mm bayonet style bends with this specialty inserted tip plier. Also used for stepping the wire up or down, in or out and can be used intraorally or extraorally. For use with round or rectangular wires up to .022" x .028".



Posterior Band Removing Plier Item #: 000120

This plier makes band removal a snap due to the alignment of the 5.0 mm occlusal pad and the curved inserted tip. Easy access to posterior areas, simply place the pad on the occlusal surface and the tip under the band and squeeze gently. Replacement pads available.



Angled Bracket Removing Plier Item #: 000122

The angled head allows for easy access to difficult posterior areas. Align the tips between the edges of the bracket base and the tooth surface to quickly and safely remove all types of brackets. Diamond honed inserted tips for increased sharpness and long life.



Straight Bracket Removing Plier Item #: 000123

The straight head allows for easy access when debonding all types of brackets. Simply wedge the tips between the edges of the bracket base and the tooth surface to quickly and safely remove brackets. Diamond honed inserted tips for increased sharpness and long life.

I ORTHODONTIC INSTRUMEN





Adhesive Removing Plier Item #: 000124

This specially designed plier allows access to all areas of the mouth to quickly remove adhesive after brackets are debonded. Simply place the pad on the occlusal surface of the tooth and catch the carbide tip on the excess adhesive to remove. The replaceable, double-sided carbide tip offer extra long life.



Hook Crimping Plier Item #: 000126

This inserted tip plier is designed to securely hold all types of crimpable hooks. One gentle squeeze is all you will need to easily attach hooks to the archwire. Can be used inside or outside of the mouth.



Replacement Carbide Tip for #000124

R00124 tip^{\$}33.50/pk 1 per pack, double-sided



Replacement Teflon Tips for #000120 and #000124

3 per pack Heat Sterilizable up to 370°F/188°C









Bend NiTi in a Cinch!



NiTi Cinch Back Plier

Bend NiTi in a cinch! Simple one step cinch back of nickel titanium archwires. Easy to bend distal ends inside of the mouth, no heat treating of the wire required. Slim line design for maximum patient comfort and ease of use by the clinician. Universal plier can be used in any quadrant. For use with wire up to .022" x .028".

Item #: OT-800	
1 instrument	\$112.50 ea.
6 instruments	\$106.50 ea.
12 or more	\$99.50 ea.
Accepts all methods of	f sterilization









177



Lifetime Warranty

Falcon Orthodontic Instruments include a *lifetime warranty against materials defects and craftsmanship. Cutting instruments and wire benders are guaranteed for seven years and utility pliers for ten years. This covers breakage and joint failure under normal use. Sharpening cutters and/or repairing minor tip damage is not covered by warranty as this is considered to be routine maintenance for normal use.

Stainless Steel Inserts

Ensure clean wire cuts every time

Balanced Joint Assembly

Sliding lock joints are factory adjusted for precise long lasting tip alignment

Ergonomic Handle

Creates excellent balance and a sure grip

- Surgical Grade Stainless Steel
- Highly Corrosion Resistant
- Diamond Honed Edges
- Precise Tip Alignment

1 instrument	\$107.50	ea.
6 instruments	\$98.50	ea.
12 or more	\$92.50	ea.
Heat sterilizable up to	370°F/188	°C

© 2015 Ortho Technology, Inc. Falcon Orthodontic Instruments is a registered trademark of Ortho Technology. **Pin and Ligature Cutter Item #: OT-1001



15° Pin and Fine Wire Cutter Item #: OT-1000

The tips are angled 15° for a better field vision and to improve access. The tool steel inserted tips have backing ledges for greater strength. Cuts soft wire up to .015", pins and elastics.



Hard Wire Cutter Item #: OT-1000HWC

Specially designed for hard wires. The cutting angle is 15° and the edge has been diamond honed to allow cutting hard wire up to .020" round.



Pin and Ligature Cutter Item #: OT-1001

This cutter plier is standard sized with sharp cutting edges designed to cut soft wires up to .015", pins and ligatures. The cutting angle is 7°. Compact design.

ORTHODONTICINSTRUMENTS

Falcon

Orthodontic Instruments®



Mini Pin and Ligature Cutter Item #: OT-1002

This cutting plier has finer tips for easy access into difficult areas. The cutting angle is 7°. Cuts soft wire up to .012", pins, ligatures and elastic smoothly and quickly.



Long Handled Distal End Cutter Item #: OT-1016XL

Compact head on this cutter provides the safety hold features of the OT-1016 but has long handles. Cuts wire up to .018" x .025".



45° Ligature Cutter Item #: OT-1006

The angled tips provide easy access to the posterior area, eliminating cheek impingement. Also designed to cut lingual ligatures. Cuts soft wire up to .015".



Micro-Miniature **Light Wire Cutter** Item #: OT-1021

Compact tip size allows for easy access to tight intra bracket areas. The cutting angle is 7°. Cuts soft wire up to .012".



Distal End Cutter, Safety Hold Item #: OT-1016

This cutter shear cuts hard wires close to the buccal tube and holds the loose wire. Cuts wire up to .022" x .028".




Hard Wire Cutter Item #: OT-7001

Heavy duty cutting surface easily cuts wires and pins. Strong enough to bend archwire ends after distal cutting. Cuts wire up to .022" x .028".



Distal End Cutter, Flush Cut With Safety Hold Item #: OT-1016FC

Cuts "flush" to the distal end of the buccal tube. The wide spring steel clip holds and catches the wire throughout cutting. Cuts wire up to .018" x .025".



Straight Utility Plier, **How Style**

Item #: OT-110

A general utility plier with 3/32" diameter tips which are serrated. Designed with rounded tip for maximum patient comfort and safety in mind.



Small Distal End Cutter, Safety Hold Item #: OT-1016S

Compact head design provides easy access to hard to reach areas while minimizing cheek impingement. Cuts wire up to .018" x .025".



Offset Utility Plier, **How Style** Item #: OT-111

A general utility plier similar to the OT-110 with machined tips, 3/32" diameter, but at a 45° offset angle. Matching serrated tips hold wires firmly.

ORTHODONTIC INSTRUMENTS -

Falcon

Orthodontic Instruments®



Band Seating Plier Item #: OT-113

Firmly grips the bracket tie-wing aiding in the seating of the band. Designed with serrated tips for maximum grip that prevent slippage during use.



Light Wire Plier Item #: OT-140

Designed for working round wire up to .016" diameter. The longer tapered beaks make it easier to bend small diameter loops. Prevents wire scoring.



Ball Hook Crimping Plier Item #: OT-126

Designed to crimp auxiliary stops, hooks and posts to archwires. The slotted design fits over the base of the hook while holding it tight during crimping.



Cinch Back Plier Item #: OT-157

Use for cinching back or forming round and rectangular wire up to .028". The serrated tips are tapered to a fine point. Perfect for accessing the posterior region.



Angled Ball Hook Crimping Plier

Item #: OT-129

Designed with an angled tip to easily crimp auxiliaries to archwires in the posterior region. Slotted tip holds auxiliaries secure to allow for tight crimping.



Weingart Utility Plier Item #: OT-158

Matching serrated pads hold wire firmly at a convenient working angle. Tapered tip fits easily between brackets. Rounded for safety.



Detailing Step Plier Item #: OT-133

Special double-sided offset head forms bayonet bends. Forms right hand or left hand step-up or step-down bends in size 1.0 mm.



Utility Plier (Weingart Style, Slim tips)

Item #: OT-158BJ

The ultra slim tips are curved for easier access to difficult areas. Matching serrated pads hold wire firmly at a convenient working angle. Beaks are carefully rounded to make this plier completely lip safe.



Bird Beak Plier Item #: OT-139

An extremely popular and versatile utility plier designed for working round wire up to .030". Tips are specially designed to prevent wire scoring.



Three Jaw Plier Item #: OT-200

Precision-aligned tips ensure consistent bends. Gradually tapered for delicate bending. Diamond honed edges. Will bend wires up to .030".

Falcon

Orthodontic Instruments®



Hollow Chop Arch Forming Plier Item #: OT-203B

The smooth inner surface allows for subtle contouring and forming of archwires. Forms without scoring or torquing. Perfect for accentuated curves.



Posterior Band Remover, Long Item #: OT-347L

The ³/16" pad rests on the occlusal of the molar and the non-padded remover tip catches the edge of the band for upward or downward removal. Replacement pads available. See page 182.



Stop (V-Bend) Plier Item #: OT-220

Makes 1.0 mm V-bends anywhere on the archwire to provide positive stops to keep the archwire from creeping. Edges are carefully diamond honed to prevent wire scoring. Use on NiTi and Stainless Steel round archwires up to .020".



Direct Bond Bracket Remover (Wide, Replaceable Tip and Pad) Item #: OT-349RT

Designed for removal of direct bond brackets. The 3/16" pad rests on the incisal/ occlusal edge of the tooth to provide leverage. While the tip fits between bracket and tooth surface, lifting bracket off with no stress. Also removes adhesive on tooth. Replacement tips and pads available. See page 182.



Nance Loop Forming Plier Item #: OT-230

Four steps allow precision forming of 3.0-6.0 mm loops. One beak is precision serrated to prevent wire slippage. Beveled edges prevent scoring of wire.



Loop Forming Plier, Tweed Style Item #: OT-350

Designed for making bull and helical closing loops in square or rectangular wire. Tips are designed to hold wire firmly up to .022" x .028."



Angled Direct Bond Bracket Remover Item #: OT-344

Designed to quickly and safely remove all bonded attachments. The 60° angled tip allows for easy access to difficult anterior and posterior areas.



Omega Loop Forming Plier, Tweed Style Item #: OT-352

Concave beak is contoured for making more uniform Omega Loops with ease. 3 sections of .045", .060", and .075" diameter and .090" length.



Direct Bond Remover, Wide Item #: OT-346RT

The wide replaceable precision tips wedge between the base of the bracket and the tooth surface, safely and easily lifting off the bracket.



Jarabak Plier Item #: OT-355

Ideal for precise bending of lighter wires. Precision grooves ensure accurate bending and closing of loops. Designed with serrated tips. Wires up to .020".

ORTHODONTICINSTRUMENTS

Falcon

Orthodontic Instruments®



Acrylic Debonding Plier Item #: OT-360

Designed to quickly and easily remove bonded orthodontic acrylic appliances, such as RPE's and Herbst. The plier has a hinged platform on one side and a debonding blade on the other.



Replacement Pads for #OT-347L

OT-2076 ³/16" **\$29.99/pk**OT-2077 ¹/4" **\$29.99/pk**4 pads per pack
Heat sterilizable up to 370°F/188°C



Lingual Arch Forming Plier Item #: OT-410

Designed to form double-back and triple-back bends in either .030" or .036" wire for inserting in lingual sheaths. Maintains a precise .140" width.




Replacement Tips for #OT-346RT

OT-1405^{\$}28.95/pk 2x tips, 2x screws, and 1x hexagon wrench per pack



Rectangular Arch Forming Plier Item #: OT-810S

Tweed Style, ideal for forming square or rectangular archwires. Beveled edges prevent scoring of archwire. Bends wire up to .022" x .028".



Replacement Tips for #OT-349RT

OT-1406.......\$28.95/pk 1x blade, 1x occlusal pad, 1x screw, and 1x hexagon wrench per pack Heat sterilizable up to 370°F/188°C



Loop Forming Replacement Tip for #OT-352 and #OT-350

OT-316**\$29.99/pk**1x tip, 1x screw, and
1x hexagon wrench per pack



Looking for MSDS Answers?

If you have questions about the chemical safety of our products, you can find a complete list of downloadable product Material Safety Data Sheets on our website at www.orthotechnology.com/MSDS.





Advanced Ergonomic Design Now With Satin Finish

- Large-Diameter Handle Reduces muscle fatigue for less risk of carpal tunnel syndrome
- Feather-Light Design Allows for a relaxed grip which provides improved tactile sensitivity
- Textured Grip Better gripping and rotational control
- Corrosion Resistant Made with high quality stainless steel material

Falcon Hand Instruments 5 or more\$17.25 ea. Heat sterilizable up to 370°F/188°C



Lotus Plus® SLB Explorer

Ortho Technology

Double Ended Ligature Director and Explorer, Item #: OT-2106 For opening Lotus Plus bracket clip

• Features an archwire director end to engage the archwire into the

See page 17 for Lotus Plus SLB opening instructions.







Ortho Technology

Direct Bond Adhesive Remover, Item #: OT-2102





The universal standard size double-ended distal bender easily crimps the distal ends of all typical wire sizes. Designed to be quick and comfortable.

- Easy access in difficult areas
- Textured grip for total control
- Ideal for precise bending
- Stainless steel design
- .018" and .030" openings perfect for all wires
- Simple to use
- Made in the USA

Item #: OT-320 1 instrument \$34.95 ea. 3 or more......\$32.95 ea. Heat sterilizable up to 370°F/188°C







183



For elastic modules or stainless steel ligature ties.

Item #: OT-8015 1 instrument........... \$69.25 ea. Heat sterilizable up to 370°F/188°C



Narrow Tip, High Polished Finish

Narrow tip designed for easy placement of elastic modules. Item #: OT-8015N 1 instrument \$69.25 ea. Heat sterilizable up to 370°F/188°C








Serrated Trimming Scissors

- Serrated design for controlled trimming without slippage
- Curved tip cuts easily and precisely, even in hard to reach lingual areas
- Made in the USA
- 3 1/2" curved serrated trimming scissors

Item #: OT-180 1 instrument.........\$20.95 ea. Heat sterilizable up to 370°F/188°C









To open the Lotus Plus bracket clip, place the straight arm in the clip opening, and the curved arm in the under tie-wing area and squeeze. The tool will fully open the Lotus Plus bracket without the worry of accidental bracket conversion.





Secure Closure

To close the Lotus Plus bracket clip, place the straight arm at the base of the clip and the curved arm in the under tie-wing area and squeeze.

■ Made in the USA

Item #: OT-210 5 or more......**\$23.00 ea.** Heat sterilizable up to 370°F/188°C

Lotus Plus® Archwire Directing and Clip Closing **Tweezer PLUS! Easily Seats Archwires in Crowded Cases**







Secure Closure

To close the Lotus Plus bracket clip while simultaneously directing the archwire, place the notched arm on the archwire and the other arm at the base of the clip and squeeze.

■ Made in the USA

Item #: OT-211 1 instrument\$25.00 ea. Heat sterilizable up to 370°F/188°C







Lotus Plus® SLB TweezerClip Closing Tool

Secure Closure

To close the Lotus Plus bracket clip use the Lotus Plus SLB Tweezer. See page 17 for instructions. 

Direct Bond Instrument

Stainless steel construction. Reverse action will not release bracket until desired. Designed with serrated tips. 





Direct Bond Tweezers with Bracket Aligner

Reverse action will not release bracket until desired. Aligns brackets by placing the opposite end into the bracket's slot. Stainless steel construction.

Item #: OT-322 1 instrument..........\$21.95 ea. 5 or more\$19.85 ea. Heat sterilizable up to 370°F/188°C



Direct Bond Tweezers

Stainless steel construction. Reverse action, specifically designed for accurate posterior placement of brackets or tubes.

■ Made in the USA

NEW!

Eruption Appliance Angled Button Tweezers

Designed to fit snugly around button pad type. Special clawtype tweezers to grip the button securely until final placement.



Heat sterilizable up to 370°F/188°C

NEW!

Eruption Appliance Curved Eyelet Tweezers

Designed to fit low profile eyelets and mesh pad styles. Special grooves to grip flat low profile pads with gold loops that are commonly difficult to hold.

Keat Buccal Tube Tweezers



The curved body of the tweezers provides an unobstructed view for easier precise placement of buccal tubes.

A. Keat Buccal Tube Tweezers Item #: OT-209

B. Reverse Action Keat Buccal Tube Tweezers

Item #: OT-209R 1 instrument.....

Curved Direct Bond Instrument





- Designed to securely hold and place lingual brackets or attachments
- Works great with our Tongue Tamers, Bite-Builders, and placement of molar tubes

ORTHODONTIC INSTRUMENTS!



Now With Satin FinishStainless steel with spring back action for easy use.

Item #: OT-101	
1 instrument	\$41.50 ea.
3 or more	\$38.50 ea.
Heat sterilizable up	to 370°F/188°C



Boone Gauge Stainless Steel

Accurately measures height of bracket placement from incisal edge. 3.5, 4.0, 4.5, and 5.0 mm.

Item #: OT-000-SS	
1 instrument	\$29.25 ea.
5 or more	\$25.50 ea.
Heat sterilizable up to	o 370°F/188°C



Mershon Band Pusher

- 5 sided angled serrated tip
- Reduced risk of slippage
- Lightweight hollow handle

Item #: OT-206	
1 instrument	\$24.50 ea
5 or more	\$22.50 ea
Heat sterilizable up to	370°F/188°C

The state of the s

Positive Positioner

Now With High-Polished Finish

Precise placement of direct bond tubes or attachments made simple with the spring-loaded mechanism and small, grooved tips. The Positive Positioner is made of lightweight, anodized aluminum.



- Access first and second molars easily
- Slides open smoothly to release tube
- Made in the USA

Item #: OT-2222	
1 instrument	\$46.50 ea
3 or more	\$42.75 ea.
Autoclavable up to	250°F/121°C



Aluminum Height Gauge

Measures height of bracket placement from incisal edge. 3.5, 4.0, 4.5, and 5.0 mm.

.018"	OT-2118-SS
.022"	OT-2122-SS
1 instrument	\$32.25 ea.
5 or more	\$28.25 ea.
Autoclavable up t	o 250°F/121°C



High Heat Height Gauge

Same features as the aluminum gauge at a fraction of the price. Withstands heat up to 350°F/177°C.

.018"	OT-2118-HH
.022"	OT-2122-HH
1 instrument	^{\$} 18.85 ea.
5 or more	\$16.50 ea.



- Metal serrated triangle tip
- Molded high impact plastic
- Safe for dry heat sterilizers
- Made in the USA


Item #: OT-600-1	99
1 bite stick	^{\$} 15.50 ea.
10 or more	^{\$} 14.50 ea.
Sterilizable up to	350°F/177°C

ORTHODONTIC INSTRUMEN

ORTHODONTICINSTRUMENTS

Sensation™ Self Ligating Bracket System Opening Instruments













Sensation and Sensation M Clips open and close in the same manner, however, each requires a unique instrument. While both instruments may look similar, each must be used with the appropriate bracket system. Using the incorrect instrument may result in damaging the clip.

To open the Sensation or Sensation M Bracket Clip:

- Place the inside hook of the opening instrument into the opening of the door of the bracket interface.
- 2. Tilt the opening instrument gingivally to avoid possible removal of the clip.
- 3. Pull straight down to open using a gentle controlled force.



Sensation™ SLB Opening Instrument



Sensation™ M SLB Opening Instrument

Unique Clip Opening Tool specifically designed to prevent overstressing of the clip. To close; lightly apply finger pressure on the occlusal surface of the clip. Make sure the spring clip fully closes so the archwire is securely engaged.

Sensation Ceramic SLB Opening Instrument Item #: OT-7920

Sensation M SLB Opening Instrument Item #: OT-8920



FOR MORE INFORMATION ABOUT...

Sensation Active Self-Ligating Brackets See pages 10-12

Sensation M Active Self-Ligating Brackets See pages 13-14



Color Coded Bite Sticks

Molded high impact plastic with metal serrated triangle tip

■ Made in the USA

Cold Sterilizable

Teal	600-201T
Purple	600-201P
Blue	600-201B

Autoclavable

rear	600A-2011
Purple	600A-201P
Blue	600A-201B



DSX Instrument Lubricant

- Seals microscopic pores to protect against corrosion
- Reduces friction at the joint of the instrument
- Premium grade fluoropolymer
- Silicone based
- Made in the USA

	item #: MK553		
r	1 bottle	\$27.99	ea
	3 or more	\$24.99	ea
	4oz per bottle		

ORGANIZERS

MADE USA

Elastomeric Organizers

- Keeps contaminates out
- Convenient storage and visual inventory

A. Combination Organizer Item #: OT-CDW


6 3/4" W x 5 1/2"H x 9"D

B. Power Chain Organizer Item #: OT-PCD

3 ⁷/₈"W x 4 ¹/₄"H x 9 ⁷/₈"D C. Power Stick Organizer Item #: OT-PSD

10 ½"W x 5 ½"H x 5 ¼"D





Instrument Organizers

A. Combination Organizer Item #: PDW-001

3 or more\$42.99 ea. 6 1/2 "W x 3 7/8 "H x 9 1/4 "D

B. Covered Plier Organizer

Item #: 56146 3 or more\$57.99 ea. 7 1/4"W x 7 1/2"H x 6"D

C. Upright Plier Organizers

White: OT-IHW . Black: OT-IHB . Clear: OT-IHC 1 organizer......\$26.99 ea. 3 or more\$23.99 ea. Holds up to 18 pliers 4 7/8"W x 5 1/4"H x 9"D

D. Round Plier Organizer

Item #: PDW-002 1 organizer......\$27.99 ea. 3 or more\$25.99 ea. 4 3/32"H x 4"D



Brackets, Bands, & Wires

Archwire Organizer

A. With Lid Item #: OT-AWH 73/4"W x 41/2"H x 4"D

B. No Lid Item #: OT-AWHS 73/4"W x 25/8"H x 4"D

1 organizer......**\$64.99 ea.**

C. Bonding Pad Organizer

Item #: OT-BP370 6 1/4 "W x 6 3/4 "H x 4 3/4 "D


D. Stackable Band Organizer Item #: 85180K


Stackable Bracket Organizer

Item #: 85181K

Includes: band/bracket tray 11 1/4"W x 1 1/2"H x 8 1/4"D









Organizers do not include product.



Ortho Model™ Archive System Each Kit Holds Up to 320 Models

• Saves valuable office space

- Inexpensive storage
- Holds 320 models
- Easy to assemble
- Super strong construction

A. Top Loading

Item #: 37810

3 or more \$95.95/kit Dimensions: 15"H x 19 1/2"W x 13 1/4"D 1 set of 4 large boxes (16 trays) Use for finished cases

B. Front Loading

Item #: 37809

3 or more.......\$105.95/kit

Dimensions: 16"H x 18 3/4"W x 13"D 1 set of 4 large boxes (16 trays)

Use for active cases

ground shipment recommended due to weight







ORGANIZERS

MADE USA

B. Front Loading Style Use for Active Cases

Elastics Organizers

- High quality acrylic plastic
- Countertop or wall mounted
- Makes for easy access

A. Elastics Dispenser Item #: OT-OED

1 organizer......**\$65.99 ea.** 13"W x 10 ½"H x 4"D

B. Creatures of the Sea® **Elastics Organizer**

Item #: OT-OEDS 1 organizer......**\$52.99 ea.**

3 or more**\$46.99 ea.** 17 1/4 "W x 10"H x 41/4"D



Misc. Organizers

• Chairside presentations

A. Relief Wax Dispenser

Item #: OT-RWD 1 organizer......**\$53.99 ea.** 8 1/4"W x 8 1/2"H x 4 1/2"D

B. Alginate Flavoring Organizer

Item #: OT-AFC 1 organizer......**\$37.99 ea.** 10 5/8"W x 6 5/8"H x 3"D

C. Adhesive Syringe Organizer

Item #: OT-ASO

3 or more**\$72.99 ea.** 7 1/4"W x 4 1/2"H x 10 1/8"D

D. Toothbrush Dispenser

Item #: OT-TB DISPLAY 8 ⁷/16"W x 11 ¹/4"H x 7 ⁷/16"D







189

NDEX

Shipping Information

Shipping: Most orders ship the same day you place your order. You can have your products as soon as the next business day. Place your order before 2:00 P.M. E.S.T. and request overnight shipping. Otherwise we will ship your order via common carrier.

Hassle-Free Returns: Customer satisfaction is Ortho Technology's primary objective, and we understand, your purchase may not always meet your needs. In an effort to make the return process as hassle-free as possible, we ask that you please call us. Your satisfaction is very important to us, and your comments help us improve our products and service. When you call, a mandatory Return Merchandised Authorization (RMA) number will be assigned. This number will help expedite your exchange return.

Return Policy: Include the Return Merchandise Authorization (RMA) number which can be obtained by contacting our Customer Service Department. Products may be returned for credit within 30 days from Invoice date. Products under warranty must meet terms within warranty period. After 31 days and up to 60 days from Invoice date, return request will be evaluated by a Customer Service Representative as a courtesy to our customers. Returns may be authorized when these conditions are met: 1. Products are in original and unopened packaging and in re-sellable condition. 2. Returns must meet quality inspection standards, appearance and physical count. 3. A 15% restocking fee may be charged. Custom orders and special orders are not accepted for return.

Due to OSHA regulations regarding possible exposure to bloodborne pathogens, Ortho Technology is prohibited from accepting instruments for evaluation without proof of sterilization. Unless proof of sterilization is evident, Ortho Technology reserves the right to either destroy these items or return them to the customer.

Order Shortage Notification Policy: Any shortage of your order must be reported promptly within 30 calendar days of your order receipt to be considered. Ortho Technology is not responsible for any shortage reported after 30 calendar days.

Terms and Billing: Orders from qualified accounts are accepted on the basis of payment within 30 days from the date of invoice. All invoices are to be paid in U.S. Dollars. Prices listed are F.O.B. point of origin. A packing slip and invoice are enclosed with each shipment. A statement covering all open invoices is mailed monthly. A 1.5% (\$2.00 minimum) per month service charge on overdue balances will be added

Canadian Customers: To provide expedited and economical shipping for our Canadian customers, orders will be shipped into Canada to your door by a common carrier. All duties, brokerage, GST, and other fees are applied during billing. This will ensure easy receipt of shipment. For easy payment, we recommend using a credit card. If payment by check is required, please write "In U.S. Funds" on the check. All catalog prices and invoice totals are listed in U.S. Dollars. When purchasing products from Ortho Technology, the customer surrenders and transfers the rights to receive the actual amount of duties refunded. All duties paid will, in turn, be refunded to the customer directly by Ortho Technology. GST will not be refunded on returns. Canadian businesses that are GST registered may claim this amount through an input tax credit.

Continued Improvement: Due to continuing product improvements, Ortho Technology reserves the right to change prices, discontinue products or change the design of products without advance notice. We are not responsible for typographical or photographic errors, and we reserve the right to limit quantities.

-A-
Adams Clasps138
Alginate Flavoring
Alginate130
Alginate Mixing Machine130
Aluminum Oxide121
Applicators118-119
Archwires
Beta Titanium73, 77
Braided82
NEW! Copper74-75
Gold Plated85
Nickel Titanium70-71, 73, 78 Pre-Torqued78
Retraction
Stainless Steel
Thermal Nickel Titanium
Tooth Tone (Coated)35, 84
TRIFECTA
Universal Form73
Archwire Markers, Dispenser77
—B—
Ball Clasps138
Ball Hooks34, 87
Rand Cement Kits and Refills
Fuji Ortho Products110-111
Light-Cure108
Paste/Paste108
Powder/Liquid109
Base Plate Wax139
Bicuspid Bands63
Bicuspid Weldable Brackets63
Bi-Directional Hooks34, 87
Bite Raiser89
Bite Sticks
Bite Wafers
Bite-Builder141 Bondable Buttons, Cleats, Eyelets, Sheaths34, 62, 88
Bondable Retainers143
Bondable Tips, Mini-Mold
Bonding Brushes
Ronding Kits and Pofills
Fuji Ortho Products110-111
Light-Activated Syringe, Capsule104-105
Primer-Activated Syringe, Capsule106
rinner receivated by mige, capsure minimizer
Bonding Pads, Travs118
Bonding Pads, Trays
Brace Relief, Kits, Refills157
Brace Relief, Kits, Refills
Brace Relief, Kits, Refills
Brace Relief, Kits, Refills
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 8 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66 Mini Bondable 64-66
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes 8 Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes .80 Bondable .49-51, 64-66 Mini Bondable .64-66 Weldable .52-53 Burs .144
Brace Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs Interproximal Diamond Burs 144 Debonding Carbide Burs 147
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66 Meldable 52-53 Burs Interproximal Diamond Burs 144 Debonding Carbide Burs 147 Finishing Diamond Burs 147
Bracker Organizer, Stackable
Bracker Organizer, Stackable
Bracket Organizer, Stackable
Bracker Organizer, Stackable .157 Bracket Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable .49-51, 64-66 Mini Bondable .64-66 Weldable .52-53 Burs .144 Debonding Carbide Burs .147 Interproximal Diamond Burs .147 Finishing Diamond Burs .147 Laboratory Carbide Burs .149 Buttons, Composite, Stainless Steel .34, 62, 88 — C Capsules, Adhesive CementRefills .104, 106
Bracker Organizer, Stackable
Bracker Organizer, Stackable
Bracker Organizer, Stackable
Bracker Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable .64-66 Mini Bondable .64-66 Weldable .52-53 Burs .147 .140 Interproximal Diamond Burs .147 Laboratory Carbide Burs .147 Laboratory Carbide Burs .147 Buttons, Composite, Stainless Steel .34, 62, 88 — C — Capsules, Adhesive CementRefills .104, 106 Capsule Gun .104, 106, 110 Cephalometric Tracing Pads .169 Chain, Elastomeric .98
Bracker Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable Bondable .49-51, 64-66 Mini Bondable .64-66 Weldable .52-53 Burs .144 Debonding Carbide Burs .147 Finishing Diamond Burs .147 Finishing Diamond Burs .147 Buttons, Composite, Stainless Steel .34, 62, 88 — C C Capsules, Adhesive CementRefills .104, 106, 110 Cephalometric Tracing Pads .169 Chain, Elastomeric .98 Cheek Retractors
Bracker Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .16-21 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable .64-66 Weldable .52-53 Burs .144 .145 Interproximal Diamond Burs .147 Finishing Diamond Burs .147 Laboratory Carbide Burs .147 Laboratory Carbide Burs .149 Buttons, Composite, Stainless Steel .34, 62, 88 — C Capsules, Adhesive CementRefills .104, 106 Capsule Gun .104, 106, 110 Cephalometric Tracing Pads .169 Chain, Elastomeric .98 Cheek Retractors Extraoral .115 Extraoral .115 Intraoral .116
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs 147 Interproximal Diamond Burs 147 Finishing Diamond Burs 147 Finishing Diamond Burs 147 Buttons, Composite, Stainless Steel 34, 62, 88 — C — Capsules, Adhesive CementRefills 104, 106, 110 Capsule Gun 104, 106, 110 Cephalometric Tracing Pads 169 Chain, Elastomeric 98 Cheek Retractors Extraoral 115 High Hea
Bracker Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable .49-51, 64-66 Mini Bondable .64-66 Weldable .52-53 Burs Interproximal Diamond Burs .144 Debonding Carbide Burs .147 Finishing Diamond Burs .147 Laboratory Carbide Burs .149 Buttons, Composite, Stainless Steel .34, 62, 88 — C — Capsules, Adhesive CementRefills .104, 106, 110 Cephalometric Tracing Pads .169 Chain, Elastomeric .98 Cheek Retractors Extraoral .115 Extraoral .115 High Heat .115 Intraoral .116 Nola .117 Photo <t< td=""></t<>
Bracker Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable .64-66 Meldable .52-53 Burs Interproximal Diamond Burs .144 Debonding Carbide Burs .147 Laboratory Carbide Burs .147 Laboratory Carbide Burs .149 Buttons, Composite, Stainless Steel .34, 62, 88 — C — Capsules, Adhesive CementRefills .104, 106, 110 Cephalometric Tracing Pads .169 Chain, Elastomeric .98 Cheek Retractors Extraoral .115 High Heat .115 Intraoral .116 Nola .117 Photo .167-168
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs 147 Interproximal Diamond Burs 144 Debonding Carbide Burs 147 Laboratory Carbide Burs 147 Laboratory Carbide Burs 149 Buttons, Composite, Stainless Steel 34, 62, 88 — C— Capsules, Adhesive CementRefills 104, 106 Cephalometric Tracing Pads 169 Chain, Elastomeric 98 Cheek Retractors Extraoral 115 Intraoral </td
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs 147 Interproximal Diamond Burs 147 Enishing Diamond Burs 147 Finishing Diamond Burs 147 Buttons, Composite, Stainless Steel 34, 62, 88 — C — Capsules, Adhesive CementRefills 104, 106 Capsule Gun 104, 106, 110 Cephalometric Tracing Pads 169 Chain, Elastomeric 98 Cheek Retractors Extraoral 115 High Heat
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs 144 Interproximal Diamond Burs 147 Interproximal Diamond Burs 147 Finishing Diamond Burs 147 Laboratory Carbide Burs 149 Buttons, Composite, Stainless Steel 34, 62, 88 — C — Capsules, Adhesive CementRefills 104, 106 Capsule, Adhesive CementRefills 104, 106 Chain, Elastomeric 98 Cheek Retractors Extraoral 115 <
Bracker Organizer, Stackable .63, 188 Bracket Systems .63 Bicuspid Twins, Weldable .63 Ceramic .10-12, 26-30 Composite .31-33 Gold Plated .46 Active Self-Ligating .10-15 Passive Self-Ligating .16-21 Sapphire .22-25 NEW! Stainless Steel .13-21, 36-45 Buccal Tubes Bondable .49-51, 64-66 Meldable .52-53 Burs Interproximal Diamond Burs .144 Debonding Carbide Burs .147 Laboratory Carbide Burs .147 Laboratory Carbide Burs .149 Buttons, Composite, Stainless Steel .34, 62, 88 — C — Capsules, Adhesive CementRefills .104, 106 Capsule Gun .104, 106, 110 Cephalometric Tracing Pads .169 Chain, Elastomeric .98 Cheek Retractors Extraoral .115 High Heat .115 Intraoral .116 Nola
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs Interproximal Diamond Burs 144 Debonding Carbide Burs 147 Finishing Diamond Burs 147 Laboratory Carbide Burs 147 Buttons, Composite, Stainless Steel 34, 62, 88 — C— Capsules, Adhesive CementRefills 104, 106 Capsule Gun 104, 106 Cephalometric Tracing Pads 169 Chain, Elastomeric 98 Cheek Retractors Extraoral 115 Intraoral 116 </td
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable Bondable 49-51, 64-66 Mini Bondable 49-51, 64-66 Weldable 52-53 Burs 144 Interproximal Diamond Burs 144 Debonding Carbide Burs 147 Finishing Diamond Burs 147 Laboratory Carbide Burs 147 Buttons, Composite, Stainless Steel 34, 62, 88 — C — Capsules, Adhesive CementRefills 104, 106 Capsule Gun 104, 106, 110 Cephalometric Tracing Pads 169 Chain, Elastomeric 98 Cheek Retractors
Bracke Relief, Kits, Refills 157 Bracket Organizer, Stackable 63, 188 Bracket Systems 63 Bicuspid Twins, Weldable 63 Ceramic 10-12, 26-30 Composite 31-33 Gold Plated 46 Active Self-Ligating 10-15 Passive Self-Ligating 16-21 Sapphire 22-25 NEW! Stainless Steel 13-21, 36-45 Buccal Tubes Bondable 49-51, 64-66 Mini Bondable 64-66 Weldable 52-53 Burs Interproximal Diamond Burs 144 Debonding Carbide Burs 147 Finishing Diamond Burs 147 Laboratory Carbide Burs 147 Buttons, Composite, Stainless Steel 34, 62, 88 — C— Capsules, Adhesive CementRefills 104, 106 Capsule Gun 104, 106 Cephalometric Tracing Pads 169 Chain, Elastomeric 98 Cheek Retractors Extraoral 115 Intraoral 116 </td

Crimpable Tubes, Cross & Double	88
NEW! Curing Light	112
Cutting Wheel	
— D —	
De-Impacting Springs	78
NEW! Dental Floss	
Diamond Disc Safety Guard	146
Diamond Discs	
Diamond Strips	
Direct Bond Buttons, Cleats, Eyelets34, 62	130
Disconnect Kits, Microetcher	
Dust Containment Units	
	121
—E—	101
Elastics, Intraoral, Extraoral100-	101
Elastic Directors	100
Elastic Organizers101,	189
Elastic Thread, Bumper Sleeve	
Elastobyte Wafers	154
NEW! Elastomeric Chain and Sticks, Ligatures98	
NEW! Eruption Appliance	
Etchant, Gel, Liquid	107
Exafast, Exabite Material	132
Expansion Screws, Expanders139-	140
EZ Space Maintainer	90
EZ Slider	
—F—	
Facebows	126
Facebow, Laser Welded, Extraoral	126
Facemask, Multi-Adjustable	123
Facemask, Reverse Pull	
Film, Extraoral X-Ray Film	109
Flat Bow Preformed Retainer Wire	138
NEW! Floss Threaders	
Flosser, Cordless	
Flux	
Fuji Ortho Products110-	
Paste Pak Automix, Band Cement	111
-G-	
Glass Ionomer Kits and Refills	
Fuji Ortho Band Paste Pak	111
Paste/Paste Band Cement	
Powder/Liquid Band Cement	
— H —	. 05
Iland Instruments	102
Hand Instruments	183
Hand Instruments	125
Hand Instruments	125
Hand Instruments	125 160
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits - I — Impression Material, PVS	125 160
Hand Instruments	125 160 132
Hand Instruments	125 160 132 129
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS Impression Trays DuraLock Plus Stainless Steel	125 160 132 129 127
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Ilmpression Material, PVS	125 160 132 129 127 157
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 186
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 186 186
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 186 186
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 146
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 183
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 188 182 146 185 183 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 183 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 183 185 187 184
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 187 184 184
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 188 185 187 184 186 184
Hand Instruments	125 160 132 129 127 157 187 186 182 146 185 187 184 184 184
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 187 184 184 184 184
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 187 184 186 184 186
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 187 184 186 184 186 177 186
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 187 184 1184 1184 1184 1184 1184 1184 1
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 127 157 187 186 182 146 185 183 185 184 186 184 186 177 186 182
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — Impression Material, PVS	125 160 132 129 127 157 187 188 185 187 184 184 184 184 187 187 188 187 188 188 188 188 188 188
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — IImpression Material, PVS	125 160 132 129 127 157 187 186 183 185 184 184 184 184 186 177 186 182 182 182 182
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — IImpression Material, PVS	125 160 132 129 127 157 187 186 182 146 183 183 184 184 184 184 184 184 185 182 182 182 182 182 182 183
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 187 186 182 146 183 185 187 184 184 184 184 185 187 186 187 186 187 186 187 186 187 186 187 187 187 188 188 188 188 188 188 188
Hand Instruments	125 160 132 129 127 157 187 188 182 146 185 187 184 184 184 184 184 182 182 182 182 182 182 183 185 187 187 188 189 189 189 189 189 189 189 189 189
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 187 188 182 146 185 187 184 184 184 184 184 182 182 182 182 182 182 183 185 187 187 188 189 189 189 189 189 189 189 189 189
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — IImpression Material, PVS	125 160 132 129 127 157 187 186 182 146 185 187 184 184 184 184 184 182 182 182 182 183 185 187 186 187 188 188 188 188 188 188 188 188 188
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 186 185 183 184 184 184 184 184 184 185 182 182 182 182 182 182 183 183 184 185 185 185 185 185 185 185 185 185 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 186 185 183 184 184 184 184 184 184 185 182 182 182 182 182 182 183 183 184 185 185 185 185 185 185 185 185 185 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ——1— Impression Material, PVS	125 160 132 129 127 157 186 185 183 184 184 184 184 184 182 182 182 182 182 182 182 182 183 184 185 185 185 185 185 185 185 185 185 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — I — Impression Material, PVS	125 160 132 129 127 186 182 146 183 185 187 184 186 177 186 2-95 97 185 144 145 150
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — I — Impression Material, PVS	125 160 132 129 127 187 186 182 183 185 184 184 184 184 182 186 2-97 185 185 185 185 185 185 185 185 185 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — I — Impression Material, PVS	125 160 132 129 127 157 187 186 182 183 185 184 184 184 187 186 187 187 187 187 187 187 187 187 187 187
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ————————————————————————————————————	125 160 132 129 127 157 186 186 184 186 183 187 186 184 186 177 186 182 182 182 182 185 185 185 185 185 185 185 185 185 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ————————————————————————————————————	125 160 132 129 1157 187 186 184 184 188 187 186 184 184 184 184 184 184 184 184 185 186 187 186 186 187 186 186 187 186 186 186 186 186 186 186 186 186 186
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits — I — I — Impression Material, PVS	125 160 132 129 127 157 186 185 183 185 184 184 184 184 182 186 2-97 185 185 185 185 185 185 185 185 185 185
Hand Instruments NEW! Head Caps, High Pull Hygiene Kits ————————————————————————————————————	125 160 132 129 127 187 186 182 183 185 184 184 184 184 182 182 182 182 182 183 185 185 185 185 185 185 185 185 185 185

Ligatures, Elastomeric98-99
Light-Cure Adhesive104-105, 110
Light-Cure Band Cement108
Lingual Buttons62, 88
Lingual Cleats62
Lingual Retainers
Lingual Sheaths
Lip Bumpers
Lip Protectors151, 156-157
Lip Ring Retractor116, 146
Lubricant, Instrument187
Lugs62-63
-M-
Measuring Scoop, Glass Ionomer109
Mesh Bags
Metal Primer
Micro Screws
Micro Torch
Microcab Plus
Microbrush Applicators119
Microetcher IIA & Accessories120
Microetcher II & Accessories120
Mini-Mold Kit & Tips113
MiniStripper
MI Paste/MI Paste Plus158 Mirrors
MI Varnish
Mixing Pads
Model Archival Boxes134, 189
Molar Bands
Assemblies55-58
Assembly Kits59
Attachments
Bands
Band Trays63
Band Organizer, Stackable63, 188
Buccal Tubes49-53
Mold-A-Lateral114
Mouthguard Material136
NEW! Mouthguards152-153
Multi-Adjustable Facemask122-123
-N-
— N — NEW! Neckpads
— N —
— N — 125 NEW! Neckpads
— N — NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 — O — 0
— N — NEW! Neckpads
 NEW! Neckpads
 NEW! Neckpads
 N — NEW! Neckpads
NEW! Neckpads
NEW! Neckpads
— N — NEW! Neckpads
— N — NEW! Neckpads
— N — NEW! Neckpads
NEW! Neckpads
— N — NEW! Neckpads
N — NEW! Neckpads
NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 O—O— Oral Pain Reliever 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers Adhesive Syringe Organizer 189 Alginate Flavoring Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & w/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 188 Elastomeric Organizer 188 Elastics Organizer 199, 188 Elastics Organizer 101, 189 Instrument Organizers 188 Power Chain Organizer 99, 188 Rower Chain Organizer 99, 188 Rower Stick Organizer 99, 188 Relief Wax Dispenser 156, 189
NEW! Neckpads
N—NEW! Neckpads
N—NEW! Neckpads
N—NEW! Neckpads
— N — NEW! Neckpads
N—NEW! Neckpads
N—NEW! Neckpads
N—NEW! Neckpads
N — NEW! Neckpads
N — NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 — O — — Oral Pain Reliever 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers Adhesive Syringe Organizer 189 Adlesive Syringe Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & w/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 188 Elastoreric Organizer 99, 188 Elastics Organizer 101, 189 Instrument Organizers 188 Power Chain Organizer 99, 188 Power Stick Organizer 99, 188 Relief Wax Dispenser 156, 189 Stackable Band Organizer 63, 188 Stackable Bracket Organizer 63, 188 NEW! Toothbrush Dispenser 163, 188 NEW! Toothbrush Dispenser 189 — P Palatal Expanders 139-140 NEW! Patient Hygiene Kits
N—NEW! Neckpads
N — NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 — O — 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers 189 Adhesive Syringe Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & w/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 188 Elastomeric Organizer 99, 188 Elastics Organizer 101, 189 Instrument Organizers 188 Power Chain Organizer 99, 188 Relief Wax Dispenser 156, 189 Stackable Band Organizer 99, 188 Relief Wax Dispenser 156, 189 Stackable Bracket Organizer 63, 188 NEW! Toothbrush Dispenser 189 — P — Palatal Expanders 139-140 NEW! Patient Hygiene Kits 160 Photo Cheek and Lip Retractors 167-168 Photographic Mirrors 165-166
NEW! Neckpads
NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 O—O— Oral Pain Reliever 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers Adhesive Syringe Organizer 189 Alginate Flavoring Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & w/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 188 Elastomeric Organizer 199, 188 Elastics Organizer 101, 189 Instrument Organizer 99, 188 Power Chain Organizer 99, 188 Relief Wax Dispenser 199, 188 Relief Wax Dispenser 156, 189 Stackable Band Organizer 99, 188 Stackable Band Organizer 99, 188 NEW! Toothbrush Dispenser 156, 189 Stackable Bracket Organizer 63, 188 NEW! Toothbrush Dispenser 160 Photo Cheek and Lip Retractors 167-168 Photographic Mirrors, Stainless Steel 166 Photographic Mirrors, Stainless Steel 166 Photographic Mirrors, Stainless Steel 166 Pliers, Cutters 170-182 Power Chain 98 Power Sticks 99 Power Chain 98 Power Chain 98 Power Chain 98 Power Chain 98 Power Chain 99 Power Chain 99 Power Chain 98 Power Chain 99 Power Chain 98 Power Chain 98 Power Stick Color Palette 98-99 Power Chain 99 Power Chain 99 Power Conditioner, Sealant 104, 106, 110 Prophy Paste 119
NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 — O —
N − NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 − O − 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers Adhesive Syringe Organizer 189 Alginate Flavoring Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & W/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 188 Elastomeric Organizer 99, 188 Elastics Organizer 101, 189 Instrument Organizer 99, 188 Power Chain Organizer 99, 188 Power Stick Organizer 99, 188 Relief Wax Dispenser 156, 189 Stackable Band Organizer 63, 188 Stackable Bracket Organizer 63, 188 NEW! Toothbrush Dispenser 189 − P Palatal Expanders 139-140 NEW! Patient Hygiene Kits 160 Photo Cheek and Lip Retractors 167-168 Photor Cheek and Lip Retractors
NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 O O 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers Adhesive Syringe Organizer 189 Alginate Flavoring Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & w/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 19, 188 Elastoreric Organizer 99, 188 Elastics Organizer 101, 189 Instrument Organizer 99, 188 Power Chain Organizer 99, 188 Relief Wax Dispenser 156, 189 Stackable Band Organizer 63, 188 Stackable Bracket Organizer 63, 188 NEW! Toothbrush Dispenser 189 P P Palatal Expanders 139-140 NEW! Patient Hygiene Kits 160 Photo Cheek and Lip Retractors 167-168 Photo Cheek and Lip Retractors 167-168 Photo
N − NEW! Neckpads 125 No Mix Adhesive 106 Nola Dry Field System, Components 117 − O − 157 Open Coil Springs, Lengths, Spools 34, 79, 82, 84 Organizers Adhesive Syringe Organizer 189 Alginate Flavoring Organizer 133, 189 Archival Model Storage Organizer 134, 189 Archwire Organizer with & W/o Lid 188 Bonding Pad Organizer 118, 188 Bracket Prep Tray Organizer 188 Elastomeric Organizer 99, 188 Elastics Organizer 101, 189 Instrument Organizer 99, 188 Power Chain Organizer 99, 188 Power Stick Organizer 99, 188 Relief Wax Dispenser 156, 189 Stackable Band Organizer 63, 188 Stackable Bracket Organizer 63, 188 NEW! Toothbrush Dispenser 189 − P Palatal Expanders 139-140 NEW! Patient Hygiene Kits 160 Photo Cheek and Lip Retractors 167-168 Photor Cheek and Lip Retractors

Resilience Bonding Supplies104-109

Retainer Cases155

Retainer Brite137, 1	64
Retainer Thermal Forming Plastic	⊺36 143
Retract EEZ Photo Lip Retractor1	167
Retraction Archwires	.81
Extraoral	15
High Heat	115 116
Photo167-1	
Rotation Wedges	02
— S — Safety Release Modules1	124
Safety Swivel Key1	139
Saliva Ejector for Tongue Block	116 110
Scissors, Serrated1	84
Separators Elastomeric1	102
Niti Spring	
Radiopaque	
Silicone Silver Solder	
Space Maintainer	.90
Spider Screw92 Spider Screw Self-Ligating TAD Starter Kit9	-97 as
Splint Material1	136
Sponge Pellets1	18
Springs Closed Coil79,	82
De-Impacting	
Distalizing Springs	
Comparation	70
TAD Closed Coil	.79 -87
— T —	07
TADS	
Closed Coil Springs, NiTi TAD	.79 91
Spider Screw TAD System92	-97
NEW! Spider Pin93, Thermoforming Plastics	97 36
Ties	
Ligature Ties, Elastomeric98	-99
Kobayashi Ties	-84 -85
Preformed Shorty Ties34, 83	-84
Tongue Blocks	
Tonque Tamers1	141
Tooth Dryer & Replacement Parts	
Tooth Tone (Coated)	
Wires, Ties, Coils, Hooks, Stops34-35, 84, 86	
Torches, Micro	157 169
Tracing Pads, Cephalometric NEW! Typodonts20, 25, 38, 41, 47, 64, 90-91,	142
— U — Universal Solder/Flux1	127
Utility Wax Strips	139
_ V —	
Vacuum Form Cutter	
_ W _	
Waterpik Flosser	
Wax, Patient Relief1	
Wire, Niti Archwires and Accessories70-73, 78	-70
NEW! Copper74	-75
Gold Plated	.85
Straight Lengths, Spools	78. 84
Wire. Stainless Steel	
Archwires and Accessories73, 80 Brass Post	
Coaxial	.82
Gold PlatedLigature Wire Spooled, Preformed	
Retainer Wire82, 138, 1	
Retraction	.81
Straight Lengths, Spools	24
Twist Wire	.82
Wonderfill Tongue Filler1 — X —	34
— X — X-Ray Film1	169

Ordering Information



4 Ways To Order



Phone Toll Free 1-800-999-3161

Monday - Thursday 8:00 A.M. until 6:00 P.M. E.S.T. Friday 8:00 A.M. until 5:00 P.M. E.S.T.



Online

www.orthotechnology.com

e-mail: service@orthotechnology.com Place your order online 24 hours a day, 7 days a week. Get up-to-theminute information on new products and service offered.



E-Mail

service@orthotechnology.com

We accept orders via email worldwide. Be sure to include your phone number, and shipping and billing addresses so we are able to confirm your order.



Fax

1-813-864-9605

24 hours a day, 7 days a week.



Payment Methods

We accept Visa, MasterCard and American Express. We also bill net 30 days on approved accounts and accept wire transfers from international accounts.



Due to the fluctuation in cost of gold and silver, Ortho Technology products containing said metals may incur price increases or decreases accordingly throughout the life of this catalog.

Order by Phone US: 1-800-999-3161

US: **1-813-991-5896**

Order Toll Free Monday - Thursday 8:00 A.M. until 6:00 P.M. Friday 8:00 A.M. until 5:00 P.M. E.S.T.

Order by Fax 1-813-864-9605

Place your next order online at

Orthotechnology.com

User Friendly Redefined



THINK GREEN! If you are receiving too many copies of this catalog, help save valuable resources and let us know by calling **1-800-999-3161**

17401 Commerce Park Blvd. Tampa, Florida 33647 USA PRESORTED STANDARD U.S. POSTAGE PAID Tampa, Florida Permit# 3269





